Revisions

SO 0191-5

Rev from	Rev to	ECO #
0506	1107	4484-07

Notes:

1. BD Cat. Number 271045
2. Blank (Sheet) Size: Length: 11" Width: 25.5"
 Number of Pages: 6 Number of Sheets: 1
 Page Size: Length 11" Width 8.5" Final Folded Size: 2.75" x 4.25"
3. Style (see illustrations below): # 4

4. See Specification Control Number 8801241 for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# 2755 Blue
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

VS Controlled by BD Caribe, LTD

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 8801241		Category and Description Package Insert, CampyPak Plus, Microaerophilic System Envelopes	Sheet: 1 of 7 <hr/> Scale: N/A	A

BD BBL™ CampyPak™ Plus Microaerophilic System Envelopes with Palladium Catalyst

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INTENDED USE

The BBL™ CampyPak™ Plus is a system for generating a reduced oxygen, enriched carbon dioxide environment within BBL GasPak™ jars. When properly activated and incubated, the envelope will rapidly produce conditions conducive to the primary isolation and cultivation of microaerophilic organisms, such as *Campylobacter jejuni*, *C. coli* and *C. fetus* subsp. *fetus*.¹ The actual atmosphere in the GasPak system will vary according to the number of inoculated plates and on the elevation of the laboratory above sea level.

SUMMARY AND EXPLANATION

First recognized as a significant agent of diarrhea and acute gastroenteritis in the 1970s,²⁻⁵ *C. jejuni* is a leading cause of bacterial diarrhea in the United States.^{6,7} *C. coli* causes diseases in humans similar to those diseases caused by *C. jejuni* but at a much lower frequency.¹ Morris and Patton have suggested that all stool specimens and rectal swabs submitted for enteric pathogens should be routinely screened for *Campylobacter* spp.¹ *C. fetus* subsp. *fetus* sometimes causes disease in humans; the symptoms are generally extraintestinal, i.e., septicemias and in patients with preexisting disease, and generally do not involve diarrheal disease.¹

Selective isolation of *C. jejuni* and *C. coli* from normal intestinal flora in stool specimens includes cultivation on a selective medium containing sheep blood and antimicrobial agents with incubation at 42°C; specimens suspected to contain *C. fetus* subsp. *fetus* are incubated at 37°C. A reduced oxygen environment enriched with carbon dioxide is essential for optimal growth of these organisms.¹ Exposure to atmospheric oxygen may interfere with the growth of more oxygen-sensitive strains.⁵ Slow growing strains may require 3 days of incubation for optimal growth and recovery.^{1,8}

The CampyPak Plus envelopes will provide in the GasPak 100 and GasPak 150 jars an appropriate atmosphere for the cultivation of *C. jejuni*.

The catalyst pellets in the CampyPak Plus envelope are in a sachet attached to the envelope. This method eliminates the use of catalyst pellet chambers which were included in the jar lids of the GasPak 100 and GasPak 150 systems. The CampyPak Plus envelope also features a unique peel-open feature which allows easy access to the water channel of the envelope.

PRINCIPLES OF THE PROCEDURE

The CampyPak Plus disposable hydrogen + carbon dioxide generator envelope for the creation of an environment suitable for microaerophiles in a GasPak jar is a modification of the GasPak disposable hydrogen + carbon dioxide generator envelope employed for the creation of an anaerobic atmosphere in the GasPak jar.

Hydrogen, generated from a sodium borohydride tablet following the addition of water, combines in the presence of the palladium catalyst with a portion of the oxygen in the jar to form water. Unlike the GasPak envelope, the amount of hydrogen generated is insufficient to combine with the total volume of oxygen, resulting in a reduced oxygen atmosphere. Carbon dioxide is generated from a sodium bicarbonate plus citric acid tablet. The filter paper wick in the CampyPak Plus envelope slows the introduction of water into the tablet chamber, allowing the lid to be placed on the GasPak jar before significant volumes of gases are released.

Each CampyPak Plus envelope has a fresh catalyst sachet attached to the back surface, avoiding problems caused by inactivation of catalyst by hydrogen sulfide and other gases produced by some anaerobes.

REAGENTS

Each CampyPak Plus with Palladium Catalyst disposable hydrogen + carbon dioxide generator envelope consists of:

- One Sodium Borohydride Tablet,
- One Sodium Bicarbonate Plus Citric Acid Tablet,
- One Filter Paper (piece),
- Palladium Catalyst Sachet.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

May cause cancer by inhalation.

All usual precautions attendant to handling hydrogen gas should be observed when employing CampyPak Plus hydrogen + carbon dioxide generator envelopes.

Hydrogen is a flammable gas. A mixture of hydrogen gas with oxygen or air in a confined area will explode if ignited by a spark, flame or other source of ignition. The historical incidence rate of reported hydrogen gas ignition is approximately two (2) incidents per million envelopes.

Do not use envelopes if there is evidence of premature activation, which may cause the package to appear bloated.

Catalyst

- Important: Do not allow the catalyst to come into contact with water. Careful addition of 10 mL of water is necessary.**
- THE APPROPRIATE GasPak SYSTEM RACK MUST BE UTILIZED AT ALL TIMES.
- CampyPak Plus envelopes with frays or cuts on the catalyst sachet surface covering should not be used.
- After use, the jar should be opened and allowed to aerate for approximately 15 sec prior to removing the rack and its contents so as to achieve atmosphere equilibration.

BBL CampyPak with palladium catalyst hydrogen + carbon dioxide generator envelopes are intended for use only in properly maintained and operated GasPak 100 and GasPak 150 systems. Their use in other systems may yield misleading results and may even be hazardous to laboratory personnel. For satisfactory and safe use, carefully follow the instructions provided with each system. Never mix different envelope systems in the same jar (i.e., envelopes with or without integral catalysts).

ANY UNSATISFACTORY ENVELOPE MUST BE CUT OPEN AND EXAMINED FOR INTACT OR PARTIALLY INTACT TABLETS. ANY REMAINING TABLETS SHOULD BE PLACED IN A CONVENIENT SINK AND DISSOLVED IN FLOWING WATER. DO NOT DISCARD INTACT ENVELOPES.

Storage Instructions: On receipt, store envelopes in a dry environment at 20 to 30°C. Avoid storage in close proximity to volatile laboratory chemicals.

CampyPak Plus envelopes are ready to use and should be activated according to instructions below and on the envelope. The expiration date is for unopened and intact envelopes.

Do not open until ready to use.

Product Deterioration: Do not use if envelopes or catalyst sachets show evidence of having been punctured, crushed, distorted, bloated, exposed to moisture, or exhibit other signs of deterioration.

PROCEDURE

Material Provided: CampyPak Plus Microaerophilic System Envelopes with Palladium Catalyst.

Materials Required But Not Provided:

- Culture media in Petri dishes, tubes, bottles or other suitable containers,
- Incubator (35 to 42°C),
- GasPak 100 system – Cat Nos. 260626, 260627 or GasPak 150 system – Cat Nos. 260628, 260629,
- Pipette, 10 mL, or syringe.

Test Procedure:

- Place inoculated plates or tubes (in GasPak tube holder) in the GasPak 100 or GasPak 150 rack. Peel open the corner of the envelope along the dotted line and place the envelope in the clip of the rack with the catalyst sachet toward the outside. One envelope is to be used in each GasPak 100 system and three envelopes are to be used in each GasPak 150 system. Do not crease, fold or crush the envelope.
- Add 10 mL of tap, distilled or deionized water through the open corner of the envelope. Use a pipette or syringe inserting only the tip into the open corner. Never push any object down into the envelope.
- Close the GasPak jar promptly after envelope is activated and finger tighten the lid clamp in the manner described in the insert directions for the system being used. Incubate at conditions appropriate for the organisms being cultured (i.e., 42°C).
- After use, open the jar and allow it to aerate for approximately 15 sec prior to removing the rack with its contents. Do not handle the CampyPak Plus envelope until after the rack has been removed from the jar.

NOTE: Since the CampyPak Plus envelope has an integral catalyst, the use of palladium catalyst pellets in the catalyst reaction chambers in the lids of GasPak jars is not required.

User Quality Control:

- Each system should be tested periodically for its ability to provide adequate conditions for the growth of *Campylobacter jejuni*.**
- If visible condensation does not occur within 60 min after a CampyPak Plus system is activated, the reaction should be terminated by opening the jar. Check the condition of the lid and O-ring gasket as described below. Activate again with one fresh CampyPak Plus generator envelope (GasPak 100 system) or three CampyPak Plus envelopes (GasPak 150 system).
- The following quality control check list, if followed regularly, will help assure proper performance:
 - The CampyPak Plus system, once activated, should be maintained upright at all times.
 - CampyPak Plus generator envelope – The envelope must be kept upright and peeled open along the dotted line at the upper right hand corner, as indicated on each envelope. The gas-generating tablets must be present and not displaced or crushed, and the foil channels and reaction chambers must not be crushed or distorted and must be able to accept 10 mL of water. Improper positioning of the envelope or use of a damaged envelope may lead to imbalances in reactivity. Damage to the channels, which regulate the reaction rate, may result in spattering and “boiling over” of the contents of the envelope and hence inadequate gas output. The envelope must be positioned in the jar with the catalyst facing outward.
 - Plastic lid and clamp – Inspect for cracks, chips, split O-ring grooves, and other irregularities. Most cracks (craze or stress lines) develop on prolonged use and forced or excessive tightening of clamps on lids. Early cracks may appear as fine black lines near the center of the lid or along the upper edges of the GasPak 100 clamp near the screw. Cracked lids and clamps which exhibit these fine black lines should not be used and should be replaced.
 - O-ring gasket – Examine for cracks, bumps, dirt, and improper seating in the O-ring groove. Examine O-ring groove wall for splitting or cracking. The O-ring must be resilient, clean and snugly retained within the O-ring groove for satisfactory operation of the system. Replace as necessary.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI (formerly NCCLS) guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Conditions suitable for cultivation of microaerophilic, capnophilic organisms are achieved within 60 min in a properly operating GasPak jar.

LIMITATIONS OF THE PROCEDURE

The amount of oxygen remaining after 24 to 48 h of incubation is influenced by the number of organisms growing in the jar. When growth of organisms is heavy, the colony size of *Campylobacter* isolates may be reduced. This effect may be exaggerated at high altitudes.

PERFORMANCE CHARACTERISTICS

Prior to release, BBL CampyPak Plus Microaerophilic System Envelopes with Palladium Catalyst are tested for performance. Samples of each lot are placed in GasPak jars along with sheep blood agar plates streaked with three strains of *Campylobacter fetus* subsp. *jejuni* - ATCC™ 29428, ATCC 33291 and ATCC 33292. Envelopes are activated and the GasPak system is carefully observed for the presence of condensation within 60 min of activation. The jars are then placed in a 42 ± 2°C incubator for 48 h. Following incubation, moderate to heavy growth of all three strains of *C. fetus* subsp. *jejuni* is observed.

AVAILABILITY

Cat. No.	Description
271045	BBL™ CampyPak™ Plus Microaerophilic System Envelopes with Palladium Catalyst, Disposable H ₂ + CO ₂ Generator Envelopes, Pkg. of 10 envelopes

REFERENCES

- Morris, G.K., and C.M. Patton. 1985. *Campylobacter*, p. 302-308. In E.H. Lennette, A. Balows, W.J. Hausler, Jr. and H.J. Shadomy (ed.), Manual of clinical microbiology, 4th ed. American Society for Microbiology, Washington, D.C.
- Dekeyser, P., M. Gossuin-DeTrain, J.P. Butzler, and J. Sternon. 1972. Acute enteritis due to related *Vibrio*: first positive stool cultures. J. Infect. Dis. 125:390-392.
- Blaser, M.J., I.D. Berkowitz, F.M. LaForce, J. Cravens, L.B. Reller, and W.-L.L. Wang. 1979. *Campylobacter* enteritis: clinical and epidemiologic features. Ann. Intern. Med. 91:179-185.
- Smith, J.P., K. Durfee, and J.H. Marymont, Jr. 1980. Incidence of *Campylobacter* enteritis in the midwestern United States. Am. J. Med. Technol. 46:81-84.
- Smibert, R.M. 1984. *Campylobacter*, p. 111-118. In N.R. Krieg, and J.G. Holt, (ed.). Bergey's manual of systematic bacteriology, vol. 1. Williams & Wilkins, Baltimore.
- French, M.J., and L.W. Riley. 1984. *Campylobacter* infections in the United States: results of 11 state surveillance. Arch. Intern. Med. 144:1610-1612.
- “Campylobacter infections widespread”. Medical World News. 1984 (November) 10.
- Janssen, D., and A.G. Helstad. 1982. Isolation of *Campylobacter fetus* subsp. *jejuni* from human fecal specimens by incubation at 35 and 42°C. J. Clin. Microbiol. 16:398-399.

BD BBL CampyPak Plus Microaerophilic System Envelopes with Palladium Catalyst

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BENEX Limited