
### Revisions

Rev from	Rev to	ECO #
0504	0904	3050-04

**Notes:**

1. BD Cat. Number 260662
2. Blank (Sheet) Size : Length: 8"      Width: 24"  
 Number of Pages: 12      Number of Sheets: 1  
 Page Size: Length 8"      Width 4"      Final Folded Size: 4" x 4"
3. Style (see illustrations below): #4

4. See Specification Control Number 8807471 for Material Information
5. Ink Colors: Printed two sides  Yes     No  
 No. of Colors: 1      PMS# 2755 Blue
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	<p style="font-size: small;">COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION</p> <p style="font-weight: bold; font-size: large;">BD</p> <p>Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA</p>	<p>Category and Description</p> <p>Package Insert, GasPak CO<sub>2</sub> Pouch System</p>		Sheet: 1 of 13	A
Proofer	Date				Scale: N/A	
Checked By	Date		<p>Part Number:</p> <p style="font-size: 1.2em;">8807471</p>			

**BD BBL™ GasPak™ CO<sub>2</sub> Pouch System**

English: pages 1 - 3 Italiano: pagine 6 - 7 8807471  
 Français : pages 3 - 4 Español: páginas 8 - 9 2004/09  
 Deutsch: Seiten 4 - 6

See symbol glossary at end of insert. / Víz popis symbolů na konci příbalového letáku. / Se symbolglossaret i slutningen af indlægssedlen. / Zie lijst met symbolen aan het einde van de bijsluiter. / Vaadake sūmbolite seletust infolehe lõpus. / Katso pakkauselosteen lopussa olevaa kuvamerkkien sanastoja. / Voir le glossaire des symboles à la fin de la notice. / Siehe Symbol-Erklärungen am Ende der Packungsbeilage. / Δείτε το γλωσσάριο των συμβόλων στο τέλος του ένθετου. / A jelmagyarázata a használati utasítás végén található. / Vedere il glossario dei simboli alla fine del foglio illustrativo. / Žr. informacinio lapelio pabaigoje pateikiamą simbolių glosarijų. / Se i symbolforklaringen på slutten av produktvedlegget. / Zobacz objaśnienie symboli na końcu ulotki. / Consulte o glossário de símbolos no fim do folheto informativo. / Pozri slovník symbolov na konci letáka. / Consulte el glosario de símbolos al final del prospecto. / Se symbolförteckningen vid slutet av bipacksedeln.

Pokyny vám poskytnú miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar.

#### INTENDED USE

The **BBL™ GasPak™ CO<sub>2</sub>** Pouch system provides an atmosphere enriched with 3 to 12% carbon dioxide. When properly activated and incubated, the pouch will rapidly produce conditions conducive to the isolation and cultivation of capnophilic organisms, such as *Neisseria gonorrhoeae*, *N. meningitidis*, and *Haemophilus influenzae*.

#### SUMMARY AND EXPLANATION

The stimulatory effects of carbon dioxide on the growth and reproduction of bacteria have been studied extensively. Wherry and Oliver were among the first to report that the primary isolation of the gonococcus was greatly enhanced by an increased carbon dioxide atmosphere.<sup>1</sup> An environment enriched with up to 10% carbon dioxide is essential for optimal growth of *N. gonorrhoeae* and *N. meningitidis*.<sup>2-4</sup> Optimum growth of *Mycobacterium* and *Brucella* with a 10% carbon dioxide-enriched environment has been documented.<sup>5,6</sup> The stimulatory effect of 5 to 10% carbon dioxide on the growth of pneumococci has also been demonstrated.<sup>7</sup>

The **BBL GasPak CO<sub>2</sub>** Pouch system provides the appropriate atmosphere for growth of a wide variety of capnophilic bacteria and allows observation of the contents without opening the pouch. The **BBL GasPak CO<sub>2</sub>** Pouch system offers the convenience of a single use disposable system for one to three Petri dishes.

#### PRINCIPLES OF THE PROCEDURE

A carbon dioxide-enriched atmosphere is produced in each sealed pouch. The self-contained system consists of a reagent packet of dilute citric acid solution with coloring and a gas-impermeable plastic pouch with an integral sachet containing sodium carbonate. When the system is activated, the citric acid comes in contact with the sodium carbonate and reacts to produce a standardized amount of carbon dioxide. The final concentration of carbon dioxide in the pouch will vary depending on the number of Petri dishes used and the volume of the pouch at the time it is sealed.

#### REAGENTS

Each individual **GasPak CO<sub>2</sub>** Pouch system consists of:

- 1 Incubation pouch with an integral reagent sachet.
- 1 Liquid activating reagent packet (green).

The reagent sachet contains 2.5 g of the following ingredients:

- Sodium carbonate,
- Inert extender.

The liquid activating packet contains 5 mL of:

- 9% Citric acid solution,
- Green dye consisting of FD&C Yellow No. 5 and FD&C Blue No. 1.

**Warnings and Precautions**

For *in vitro* Diagnostic Use.

**Store Instructions:** On receipt, store at 20 – 30°C in a dry environment.

**Product Deterioration:** Do not use systems which show evidence of a torn or punctured pouch, or reagent sachets which have evidence of lost reagents. Do not use liquid activating reagent packets which show evidence of leakage or lost liquid contents.

**PROCEDURE**

**Materials Provided:** Each **GasPak** CO<sub>2</sub> Pouch system includes: 50 Liquid activating reagent packets, 50 Incubation pouches (with integral reagent sachet).

**Materials Required But Not Provided:** Chocolate II Agar, Modified Thayer Martin (MTM II) Agar, Martin-Lewis Agar, or other suitable medium, Ancillary materials and equipment required to inoculate Petri dishes, incubator, **BBL GasPak** Pouch Sealing Bar or alternate sealing method (i.e., heat-sealing device).

**Test Procedure:**

1. Remove an incubation pouch and a liquid activating reagent packet from the protective outer bag.
2. Hold the liquid activating reagent packet upright, by the neck; grasp and tear off the tab of the packet. Insert the opened nozzle into the reagent channel of the incubation pouch.
3. Hold the incubation pouch upright and squeeze the liquid activating reagent packet until all contents are dispensed through the reagent channel. Discard the empty packet.
4. Place the desired specimen inside the **GasPak** CO<sub>2</sub> Pouch and clamp or heat-seal the pouch shut.

Note: Do not evacuate or otherwise restrict the internal volume of the pouch when sealing, as that may cause artificially high carbon dioxide concentrations.

5. Incubate each **GasPak** CO<sub>2</sub> Pouch system at a temperature appropriate for the organism being cultured.
6. The **GasPak** CO<sub>2</sub> Pouch systems may be stacked to conserve incubator space and may be removed at any time to allow maximum observation of the enclosed specimens.
7. After incubation, open the sealing bar (alternatively cut or tear the heat-sealed pouch), remove the Petri dishes and dispose of the pouch in an appropriate manner.

**User Quality Control:** Periodic use of a capnophilic quality control strain, such as *Neisseria gonorrhoeae* ATCC™ 19424, is recommended.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

**RESULTS**

Carbon dioxide-enriched conditions are achieved rapidly. The final concentration of carbon dioxide in the pouch varies with the volume of the pouch at the time it is sealed. The carbon dioxide concentrations achieved are 8 to 12% with one plate, 5 to 8% with two plates, and 3 to 6% with three plates. The CO<sub>2</sub> level is achieved after 1 h.

**PERFORMANCE CHARACTERISTICS**

Prior to release, **GasPak** CO<sub>2</sub> Pouch Systems are tested for performance. Two wax plates are placed inside the pouch and the pouch is activated according to label directions. The pouches are immediately sealed and placed in a 35 ± 2°C incubator. After approximately 1 h, a 3 cc sample of gas is removed from the pouch using a 3 cc syringe with a 1 1/2" 22 gauge needle. The entire sample is injected into a gas chromatograph. The percent carbon dioxide for each pouch is between 3 to 6%.

**AVAILABILITY**

Cat. No.	Description
260662	<b>BBL™ GasPak™</b> CO <sub>2</sub> Pouch System, Box of 50 systems

**REFERENCES**

1. Wherry, W.B., and W.W. Oliver. 1916. Adaptation to certain tensions of oxygen as shown by gonococcus and other parasitic and saprophytic bacteria. *J. Infect. Dis.* 19:288-298.
2. Kellogg, D.S., K.K. Holmes, and G.A. Hill. 1976. Cumitech 4, Laboratory diagnosis of gonorrhoea. Coordinating ed., S. Marcus and J.C. Sherris. American Society for Microbiology, Washington, D.C.
3. Bovre, K. 1984. *Neisseriaceae*, p. 288-310. In N.R. Krieg, and J.G. Holt (ed.), *Bergey's manual of systematic bacteriology*, volume 1. Williams & Wilkins, Baltimore.
4. Devaux, D.L., G.L. Evans, C.W. Arndt, and W.M. Janda. 1987. Comparison of the Gono-Pak system with the candle extinction jar for recovery of *Neisseria gonorrhoeae*. *J. Clin. Microbiol.* 25:571-572.

5. Beam, E.R., and G.P. Kubica. 1968. Stimulatory effects of carbon dioxide on the primary isolation of tubercle bacilli on agar-containing medium. *Am. J. Clin. Pathol.* 50:395-397.
6. Alton, G.G., L.M. Jones, and D.E. Pietz. 1975. W.H.O. laboratory techniques in brucellosis. W.H.O. Monogr. Ser. 55.
7. Kempner, W., and C. Schlayer. 1942. Effect of CO<sub>2</sub> on the growth rate of the pneumococcus. *J. Bacteriol.* 43:387-396.

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Becton, Dickinson and Company. © 2004 BD.
ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo, BBL and GasPak are trademarks of
BENEX Limited