
**Revisions**


| Rev from | Rev to  | ECO #   |
|----------|---------|---------|
| 0503     | 2010/06 | 5308-10 |

**Notes:**


- BD Cat. Number 226306
- Blank (Sheet) Size : Length: 11"      Width: 20.625"  
 Number of Pages: 10      Number of Sheets: 1  
 Page Size: Length 11"      Width 4 1/8"      Final Folded Size: 5 1/2" x 4 1/8"
- Style (see illustrations below): # 4


#1


#2


#3


#4

#5

#6

#7

- See Specification Control Number VS-226306 for Material Information
- Ink Colors: Printed two sides  Yes       No  
 No. of Colors: 1      PMS# Standard Black
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

|                         |      |   |  |          |
|-------------------------|------|---|--|----------|
| Label Design            | Date | COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION |  Becton, Dickinson and Company<br>7 Loveton Circle<br>Sparks, MD 21152 USA |          |
| Proofer                 | Date |   |  |          |
| Checked By              | Date |   |  |          |
| Part Number:<br>8808061 |      | Category and Description<br>Package Insert,<br>BBL Prompt Inoculation System  | Sheet: 1 of 11<br><hr/> Scale: 1:1   | <b>A</b> |

# **BD BBL™ Prompt™** **Inoculation System**

## **For use with the Disc Diffusion Susceptibility Test**

English: pages 1 – 3  
Français : pages 3 – 5  
Español: páginas 5 – 7

8808061  
2010/06

### **INTENDED USE**

The **BBL™ Prompt™** Inoculation System is used to prepare standardized suspensions of bacteria for the Bauer-Kirby disc diffusion antimicrobial susceptibility test procedure. The system may be used for rapidly growing bacteria such as *Enterobacteriaceae*, *Staphylococcus* spp., *Pseudomonas aeruginosa*, *Acinetobacter* spp., enterococci and some nonenterococcal streptococci. It also may be used for *Haemophilus influenzae*.

### **SUMMARY AND EXPLANATION**

The Bauer Kirby procedure is a standardized method for the determination of antimicrobial susceptibility based on the agar gel disc diffusion principle.<sup>1-3</sup> This method is published as a consensus standard by the National Committee for Clinical Laboratory Standards (NCCLS) and is periodically reviewed and updated.<sup>4</sup>

The inoculum concentration has been shown to be an important factor in obtaining valid results with the Bauer Kirby procedure.<sup>2-5</sup> The original Bauer-Kirby procedure employed inoculum prepared from log phase cultures. This requires incubation in broth for several hours to achieve the desired inoculum density, equivalent to the 0.5 McFarland turbidity standard, which produces approximately  $1 \times 10^8$  CFU/mL.<sup>3-5</sup> Several studies have subsequently shown that direct standardization of the inoculum, i.e., without incubation, is an acceptable alternative.<sup>5-7</sup> The direct adjustment of inoculum preparation is the preferred method for fastidious organisms that may grow slowly in broth. The M2-T4 edition of the NCCLS disc diffusion standard contained a revised procedure for *Haemophilus influenzae*, in which a new medium, Haemophilus Test Medium (HTM), was employed.<sup>8</sup> This procedure recommended matching the inoculum density to the 0.5 McFarland standard using a photometric device.

The **BBL Prompt** Inoculation System is a device that allows direct standardization of inoculum without adjustment of turbidity or preincubation.<sup>9-11</sup> This method is considered acceptable for routine testing purposes.<sup>4</sup> It has been shown that the **Prompt** system is also satisfactory for preparing inoculum for *H. influenzae* in the Bauer-Kirby procedure.<sup>12</sup>

### **PRINCIPLES OF THE PROCEDURE**

The **Prompt** Inoculation System wand is touched to several bacterial colonies on a primary isolation plate and is placed in the tube of saline provided with the system. The bacteria are suspended in the saline by agitation with a vortex mixer. When the **Prompt** Inoculation System is used in this manner, an inoculum containing approximately  $1.5 \times 10^8$  colony forming units per mL (CFU/mL) can be expected for most bacteria. This is equivalent to the inoculum density that is achieved by matching the turbidity to that of a 0.5 McFarland standard.<sup>4</sup> Agreement of susceptibility results is greater than 95% when compared to the standardized method.<sup>11</sup> The **Prompt** inoculum should be used within 6 h after preparation for most bacteria. For *H. influenzae*, the inoculum should be used within 3 h. The **Prompt** Inoculation System facilitates inoculation by eliminating: 1) the incubation period and 2) the need for manual adjustment of inoculum density.

### **PRODUCT DESCRIPTION**

The **Prompt** Inoculation System consists of an inoculation wand and a tube of saline. The wand is a polypropylene rod attached to a stopper. At the tip of the wand are cross-hatched grooves designed to hold a specific number of bacteria. The cell suspending solution, consisting of 1 mL of sterile saline, is provided in a plastic tube with a snap off cap.

#### **Warnings and Precautions:**

For *in vitro* Diagnostic Use.

1. Colonies must be selected from FRESH culture plates ( $\leq 24$  h).
2. DO NOT FLAME the plastic inoculation wand. It will melt.
3. To avoid contamination, always keep fingers above the ridge on the inoculation wand.
4. Keep wands covered when not in use.
5. Do not use if colonies are very small ( $< 0.5$  mm in diameter).

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. After inoculation, the **Prompt** inoculation tube and wand and other contaminated materials should be sterilized by autoclaving.

**Storage Instructions:** **Prompt** units should be stored at room temperature, (2 – 27°C).

**Product Deterioration:** **Prompt** inoculation tubes should not be used if the tube or cap is cracked or if the saline is cloudy.

### **SPECIMEN COLLECTION AND TRANSPORT**

Specimens received in the laboratory should be treated in the usual manner for the preparation of a primary culture plate.

### **PROCEDURE**

**Material Provided:** 62 **Prompt** inoculation wands and 60 **Prompt** inoculation tubes.

**Materials Required But Not Provided:** Ancillary culture media, reagents, quality control organisms and laboratory equipment as required.

**Test Procedure:** The **Prompt** Inoculation System may be used for the preparation of inocula in disc diffusion procedures.

#### A. Preparation of Bacterial Suspension

1. Remove the required number of **Prompt** inoculation tubes from the box and place in a test tube rack.
2. Remove an inoculation wand from the box.
3. Holding the wand tip *perpendicular* to the agar surface, touch five isolated colonies greater than 1 mm in diameter. (As a reference, the tip of the wand is 2 mm in diameter.) Do not penetrate the agar. Do not scrape or drag the tip across the colonies.

**NOTE:** If the colonies are small (0.5 mm to 1 mm in diameter), touch ten instead of five. For very small, pinpoint colonies, continue incubation of the primary plate until they reach a diameter of approximately 0.5 mm to 1 mm. If the colony diameter is not likely to reach 0.5 mm (e.g., some streptococci), an alternative method for inoculum preparation should be used.

4. While holding the inoculation wand with one hand, remove a **Prompt** inoculation tube from the rack.
5. Bend the cap of the tube sideways until it snaps off.
6. Place the inoculation wand into the tube and press down with a twisting motion to assure a tight seal.
7. Vortex the tube vigorously for ten s to release the bacteria from the wand tip. If the organism is not released from the wand, let the solution sit for 5 min and vortex again.
8. Repeat steps 2 – 7 for all tubes.
9. The bacterial suspension should be used within 6 h of preparation. If not used immediately after preparation, shake vigorously to resuspend the bacteria just prior to use. For *H. influenzae*, use within 3 h of preparation.

#### B. Inoculation of Mueller Hinton Agar or HTM Agar Plate

1. Remove the inoculation wand from the tube. Discard in an appropriate container to be sterilized.
2. Dip a sterile cotton swab into the bacterial suspension. To remove excess liquid, rotate the swab several times with a firm pressure on the inside wall of the tube above the fluid level.
3. Streak the plate using conventional techniques and proceed with the susceptibility test.<sup>4</sup>

**User Quality Control:** The **Prompt** Inoculation System is used as part of a multicomponent diagnostic test. When employed for the NCCLS inoculum preparation procedure, it should produce the results given in the NCCLS standard for the control cultures.<sup>4</sup>

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

## RESULTS

When used according to directions, the **Prompt** Inoculation System will provide the proper inoculum for the Bauer-Kirby test for most clinical isolates. The susceptibility test results obtained with the **Prompt** Inoculation System show an overall 97.0% agreement with the standard Bauer-Kirby procedure.

## LIMITATIONS OF THE PROCEDURE

1. The **Prompt** Inoculation System should not be used when the colony size is less than 0.5 mm in diameter. Examples of organisms that may not meet the 0.5 mm size requirement are *Streptococcus* spp. other than enterococci, *S. bovis* and *S. agalactiae* (group B).  
Picking colonies which are too small will result in underinoculation and may cause a resistant organism to appear to be susceptible. An alternate method of inoculum preparation should be used for colonies less than 0.5 mm in diameter.
2. Some mucoid organisms such as "stringy" *Klebsiella* or *Pseudomonas* may not adhere to the wand when attempting to pick up the colony. This will be visually apparent. An alternate method of inoculum preparation should be used for such organisms.

## SPECIFIC PERFORMANCE CHARACTERISTICS

Three hundred (300) susceptibility results obtained from inocula prepared by the **Prompt** system technique were compared to the results obtained from inocula prepared according to NCCLS standard procedure.<sup>11</sup> The relative frequencies of genera tested were *Escherichia* (16.7%), *Klebsiella* (10%), *Proteus* (14.7%), *Providencia* (8%), *Morganella* (4%), *Enterobacter* (8.7%), *Serratia* (4.7%), *Pseudomonas* (6.7%), *Staphylococcus* (12.7%), enterococci (8.7%) and other group D and group B *Streptococcus* (5.3%). The interpretive agreement between disc diffusion susceptibility tests inoculated according to the NCCLS standard procedure and the **Prompt** Inoculation System was 97.0% for 3,578 drug bacteria combinations. When held for 6 h, significant differences were noted in zone diameters; however, overall interpreted agreement was 97.9% for the 1187 drug bacteria combinations tested. The overall mean for 742 gram-negative and gram-positive bacterial counts was  $1.85 \times 10^8$  CFU with a 95% confidence interval of  $4.66 \times 10^7$  to  $7.36 \times 10^8$  CFU.<sup>11</sup> In another study with 14 strains of *H. influenzae*, the mean inoculum density of freshly prepared suspensions was  $1.5 \pm 0.1 \times 10^8$  CFU/mL. After 3 h at room temperature, the density was  $0.7 \pm 0.2 \times 10^8$  CFU/mL.<sup>13</sup> In this study, the overall agreement when compared with a photometric device was 98.4%, and with the visual McFarland standard method, 97.9%.<sup>13</sup>

## AVAILABILITY

| Cat. No. | Description  |
|----------|--|
| 226306   | <b>BBL™ Prompt™</b> Inoculation System, 60 Tubes and 62 Inoculation Wands. |

## REFERENCES

1. Bauer, A.W., W.M.M. Kirby, J.C. Sherris, and M. Turck. 1966. Antibiotic susceptibility testing by a standardized single disk method. *Am. J. Clin. Pathol.* 45:493-496.
2. Ericsson, H.M., and J.C. Sherris. 1971. Antibiotic sensitivity testing. Report of an international collaborative study. *Acta. Pathol. Microbiol. Scand. Sec. B, Suppl.* 217:1-90.
3. Thornsberry, C., T.L. Gavan, E.H. Gerlach. 1977. Cumitech 6, New developments in antimicrobial agent susceptibility testing. Coordinating ed., J.C. Sherris. American Society for Microbiology, Washington, D.C.
4. National Committee for Clinical Laboratory Standards. 2003. Approved standard: M2-A8. Performance standards for antimicrobial disk susceptibility tests, 8th ed. National Committee for Clinical Laboratory Standards, Wayne, Pa.
5. Barry, A.L., and C. Thornsberry. 1985. Susceptibility tests: diffusion test procedures, p. 978-987. In E.H. Lennette, A. Balows, W.J. Hausler, Jr., and H.J. Shadomy (ed.), *Manual of clinical microbiology*, 4th ed. American Society for Microbiology, Washington, D.C.
6. Barry, A.L., L.J. Joyce, A.P. Adams, and E.J. Benner. 1973. Rapid determination of antimicrobial susceptibility for urgent clinical situations. *Am. J. Clin. Pathol.* 59:693-699.
7. D'Amato, R.F. and L. Hochstein. 1982. Evaluation of a rapid inoculum preparation method for agar disk diffusion susceptibility testing. *J. Clin. Microbiol.* 15:282-285.
8. Jorgensen, J.H., J.S. Redding, L.A. Maher, and A.W. Howell. 1987. Improved medium for antimicrobial susceptibility testing of *Haemophilus influenzae*. *J. Clin. Microbiol.* 25:2105-2113.
9. Baker, C.N., C. Thornsberry, and R.W. Hawkinson. 1983. Inoculum standardization in antimicrobial susceptibility testing: evaluation of overnight agar cultures and the rapid inoculum standardization system. *J. Clin. Microbiol.* 17:450-457.
10. Wicks, J.H., R.L. Nelson, and G.E. Krejcarek. 1983. Rapid inoculum standardization system: a novel device for standardization of inocula in antimicrobial susceptibility testing. *J. Clin. Microbiol.* 17:1114-1119.
11. Lund, M.E., and R.W. Hawkinson. 1983. Evaluation of the Prompt inoculation system for preparation of standardized bacterial inocula. *J. Clin. Microbiol.* 18:84-91.
12. Marsik, F., G. Evans, J. Fowler, and L. Thompson. 1989. Comparison of the BBL™ Prompt™ system, Abbott A-JUST™ and visual method for the preparation of *Haemophilus influenzae* inoculum for the Bauer-Kirby procedure. *Abstr. Annu. Meet. Am. Soc. Microbiol.* C-67, p. 404.
13. Data on file at BD Diagnostics.

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© 2010 BD.
Prompt is a trademark of 3M.
BD, BD Logo and BBL are trademarks of Becton, Dickinson and Company.