## Revisions

SO 0191-5

Rev from	Rev to	ECO #
0607	2010/10	5541-10

**Notes:**

1. BD Cat. Number 245000
2. Blank (Sheet) Size : Length: 8.5"      Width: 11"  
 Number of Pages: 40    Number of Sheets: 10  
 Page Size: Length 8.5"    Width 5.5"      Final Folded Size: 5.5" x 8.5"
3. Style (see illustrations below): # 5

4. See Specification Control Number N/A for Material Information
5. Ink Colors: Printed two sides  Yes     No  
 No. of Colors: 1                      Standard Black
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	<b>BD</b> Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 8809241JAA		Category and Description Package Insert BBL Crystal Identification Systems Enteric/ Nonfermenter ID Kit REF 245000	Sheet: 1 of 41 <hr/> Scale: NONE	A

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## INTENDED USE

The **BBL Crystal™** Enteric/Nonfermenter (E/NF) Identification (ID) System is for the identification of aerobic gram-negative bacteria that belong to the family *Enterobacteriaceae* as well as some of the more frequently isolated glucose fermenting and nonfermenting gram-negative bacilli.

## SUMMARY AND EXPLANATION

The **BBL Crystal E/NF ID System** is a miniaturized identification method. Many of the tests used are modifications of classical methods. These include tests for fermentation, oxidation, degradation and hydrolysis of various substrates. In addition, there are chromogen linked substrates to detect enzymes that microbes use to metabolize various substrates.<sup>1-5</sup>

The **BBL Crystal E/NF ID kit** is comprised of (i) **BBL Crystal E/NF panel lids**, (ii) **BBL Crystal bases** and (iii) **BBL Crystal™** Enteric/Stool ID Inoculum Fluid (IF) tubes. The lid contains 30 dehydrated substrates on tips of plastic prongs. The base has 30 reaction wells. Test inoculum is prepared with the inoculum fluid and is used to fill all 30 wells in the base. When the lid is aligned with the base and snapped in place, the test inoculum rehydrates the dried substrates and initiates test reactions.

Following an incubation period, the wells are examined for color changes. Color changes result from metabolic activities of the microorganisms. The resulting pattern of the 30 reactions is converted into a ten digit profile number that is used as the basis for identification.<sup>6</sup> Biochemical and enzymatic reaction patterns for the 30 **BBL Crystal E/NF ID** substrates with a wide variety of microorganisms are stored in the **BBL Crystal E/NF ID** database. Identification is derived from a comparative analysis of the reaction pattern of the test isolate to those held in the database. A complete list of taxa that comprises the current E/NF database is provided in Table 1 (page 36).

## PRINCIPLES OF THE PROCEDURE

The tests used in the **BBL Crystal E/NF ID System** are based on microbial utilization and degradation of specific substrates detected by various indicator systems. Fermentation reactions detect the ability of an isolate to metabolize carbohydrates in the absence of atmospheric oxygen, and oxidation reactions are based on the ability of an organism to metabolize the substrate with oxygen as the final electron acceptor. Both reactions are usually detected by the use of a pH indicator in the test substrate. Chromogenic substrates upon hydrolysis produce color changes that can be detected visually. In addition, there are other tests that detect the ability of an organism to hydrolyze, degrade, reduce or otherwise utilize a substrate in the **BBL Crystal ID System**. Reactions employed by various substrates and a brief explanation of the principles employed in the system are described in the "Reagents" section.

## REAGENTS

The **BBL Crystal E/NF ID panel** contains 30 enzymatic and biochemical substrates as described below. Panel location indicates the row and column where the well is located (example: 1J refers to Row 1 in Column J).

**Precautions:** For *in vitro* Diagnostic Use

## Reagents and Principles of Tests Employed in the BBL Crystal E/NF ID System

Panel Location	Active Ingredient	Code	Approx. Amt. (g/10 mL)	Pos.	Neg.	Principle (Reference)
4A	Arabinose	ARA	3.5	Gold/Yellow	Orange/Red	Utilization of carbohydrate results in lower pH and change in indicator (Phenol red). <sup>7-10</sup>
4B	Mannose	MNS	3.0	Gold/Yellow	Orange/Red	
4C	Sucrose	SUC	2.8	Gold/Yellow	Orange/Red	
4D	Melibiose	MEL	1.0	Gold/Yellow	Orange/Red	
4E	Rhamnose	RHA	3.0	Gold/Yellow	Orange/Red	
4F	Sorbitol	SOR	3.5	Gold/Yellow	Orange/Red	
4G	Mannitol	MNT	1.8	Gold/Yellow	Orange/Red	
4H	Adonitol	ADO	2.5	Gold/Yellow	Orange/Red	
4I	Galactose	GAL	1.5	Gold/Yellow	Orange/Red	
4J	Inositol	INO	1.3	Gold/Yellow	Orange/Red	

## Reagents and Principles of Tests Employed in the BBL Crystal E/NF ID System (continued)

Panel Location	Active Ingredient	Code	Approx. Amt. (g/10 mL)	Pos.	Neg.	Principle (Reference)
2A	p-n-p-phosphate	PHO	0.025	Yellow	Colorless	Enzymatic hydrolysis of the colorless aryl substituted glycoside or phosphate ester releases yellow p-nitrophenol. <sup>1,5</sup>
2B	p-n-p α-β-glucoside	BGL	0.025	Yellow	Colorless	
2C	p-n-p-β-galactoside	NPG	0.06	Yellow	Colorless	
2D	Proline nitroanilide	PRO	0.07	Yellow	Colorless	Enzymatic hydrolysis of colorless amide substrate releases yellow p-nitroaniline. <sup>1,5</sup>
2E	p-n-p bis-phosphate	BPH	0.02	Yellow	Colorless	Enzymatic hydrolysis of the colorless aryl substituted glycoside or phosphate ester releases yellow p-nitrophenol. <sup>1,5</sup>
2F	p-n-p-xyloside	BXY	0.03	Yellow	Colorless	
2G	p-n-p-α-arabinoside	AAR	0.03	Yellow	Colorless	
2H	p-n-p-phosphorylcholine	PHC	0.03	Yellow	Colorless	
2I	p-n-p-β-glucuronide	GLR	0.02	Yellow	Colorless	
2J	p-n-p-N-acetyl glucosaminide	NAG	0.04	Yellow	Colorless	
1A	γ-L-glutamyl p-nitroanilide	GGL	0.03	Yellow	Colorless	Enzymatic hydrolysis of the colorless amide substrate releases yellow p-nitroaniline. <sup>1,5</sup>
1B	Esculin	ESC	0.14	Brown/ Maroon	Clear/Straw	Hydrolysis of esculin results in a black precipitate in the presence of ferric ion. <sup>11</sup>
1C	p-nitro-DL-phenylalanine	PHE	0.1	Gold/ Dk. Orange	Yellow	Oxidative deamination of phenylalanine results in a brown color in the presence of ferric ion. <sup>7,11</sup>
1D	Urea	URE	0.2	Aqua/Blue	Yellow/Green	Hydrolysis of urea and the resulting ammonia change the pH indicator color (Bromthymol blue). <sup>7,11,12</sup>
1E	Glycine	GLY	0.7	Aqua/Blue	Yellow/Green	Degradation of glycine results in alkaline metabolites that change color of the pH indicator (Bromthymol blue). <sup>13</sup>
1F	Citrate	CIT	0.8	Aqua/Blue	Yellow/Green	Utilization of citrate results in alkaline metabolites that change color of the pH indicator (Bromthymol blue). <sup>7,14</sup>
1G	Malonic acid	MLO	1.5	Aqua/Blue	Yellow/Green	Utilization of malonate results in alkaline metabolites that change color of the pH indicator (Bromthymol blue). <sup>11</sup>
1H	Triphenyl Tetrazolium chloride	TTC	0.15	Pink/Red*	Clear	Reduction of the tetrazolium compound results in formation of a red formazan. <sup>13</sup>
1I	Arginine	ARG	1.5	Red/Purple	Yellow/Brown	Anaerobic catabolism results in pH rise and change in the color of the indicator (Bromcresol purple). <sup>7,15</sup>
1J	Lysine	LYS	0.5	Red/Purple	Yellow/Brown	

\*Precipitate may or may not be visible.

After use, all infectious materials including plates, cotton swabs, inoculum tubes, filter papers used for oxidase or indole tests and **BBL Crystal** panels must be autoclaved prior to disposal or incinerated.

### STORAGE AND HANDLING/SHELF LIFE

On receipt, store the **BBL Crystal E/NF** kit at 2 – 25°C. DO NOT FREEZE. If the kit or any components are stored refrigerated, each should be brought to room temperature prior to use.

**Lids:** Lids are individually packaged and must be stored unopened. Visually inspect the package for holes or cracks in the foil package. Do not use if the packaging appears to be damaged. Lids in the original packaging, if stored as recommended, will retain expected reactivity until the date of expiration.

**Bases:** Bases are packaged in two sets of ten, in **BBL Crystal** incubation trays. The bases are stacked facing down to minimize air contamination. Store unused bases in the tray, in plastic bag. Empty trays should be used to incubate panels.

**Inoculum Fluid:** **BBL Crystal** Enteric/Stool ID Inoculum Fluid (IF) is packaged in two sets of ten tubes. Visually inspect the tubes for cracks, leaks, etc. Do not use if there appears to be a leak, tube or cap damage or visual evidence of contamination (i.e., haziness, turbidity). Expiration dating is shown on the tube label. **BBL Crystal** Enteric/Stool ID Inoculum Fluid may be used with either **BBL Crystal E/NF** or **RS/E** panels.

### SPECIMEN COLLECTION AND PROCESSING

**BBL Crystal** ID Systems are **not** for use directly with clinical specimens. Use isolates from a blood agar plate such as **Trypticase™ Soy Agar** with 5% Sheep Blood. Use of a **MacConkey Agar** plate is also acceptable. The test isolate must be a pure culture no more than 24 h old. Only cotton-tipped applicator swabs should be used to prepare inoculum, as some polyester swabs may cause problems with inoculation of the panels. (See "Limitations of the Procedure.") Once lids are removed from the sealed pouches, they must be used within 1 h to ensure adequate performance. The plastic cover should remain on the lid until used.

The incubator used should be humidified to prevent evaporation of fluid from the wells during incubation. The recommended humidity level is 40 – 60%. The usefulness of **BBL Crystal** ID Systems or any other diagnostic procedure performed on clinical specimens is directly influenced by the quality of the specimens themselves. It is strongly recommended that laboratories employ methods discussed in the *Manual of Clinical Microbiology* for specimen collection, transport and placement on primary isolation media.<sup>16</sup>

## TEST PROCEDURE

**Materials Provided:** BBL Crystal Enteric/NF kit:

20 BBL Crystal Enteric/NF Panel Lids,

20 BBL Crystal Bases,

20 BBL Crystal Enteric/Stool ID Inoculum Fluid Tubes. Each tube has approximately  $2.2 \pm 0.1$  mL of Inoculum Fluid containing: NaCl 8.50 g, 3-Morpholinopropanesulfonic acid 0.8372 g, purified water to 1000 mL.

2 incubation trays,

1 BBL Crystal E/NF Report Pad.

**Materials Required But Not Provided:** Sterile cotton swabs (*do not use polyester swabs*); Incubator (35 – 37°C non-CO<sub>2</sub> (40 – 60% humidity); BBL Crystal Light Box/Panel Viewer (includes BBL Crystal Color Reaction Charts) with BBL Crystal ID System Electronic Codebook or BBL E/NF Manual Codebook (see “Availability”), or BBL Crystal AutoReader; nonselective culture plate (e.g., Trypticase Soy Agar with 5% Sheep Blood); BBL DMACA Indole Reagent Droppers; BBL Oxidase Reagent Droppers (see “Availability”).

Also required are the necessary equipment and labware used for preparation, storage and handling of clinical specimens.

**Test Procedure:** BBL Crystal E/NF ID System requires oxidase and indole test results. Prior to BBL Crystal E/NF panel set-up, oxidase and indole tests should be performed from a nonselective isolation plate no more than 24 h old. Perform oxidase and indole tests per instructions provided in the package insert for these reagents.

Refer to Procedural Chart Illustrations, page 38.

1. Remove lids from pouch. Discard desiccant. Once removed from the pouch, covered lids should be used within 1 h. Do not use the panel if there is no desiccant in the pouch. See Fig. A.
2. Take an inoculum tube and label with patient’s specimen number. Using aseptic technique, with the tip of a sterile cotton swab (*do not use a polyester swab*) or a wooden applicator stick or disposable plastic loop, pick one well isolated large (2 – 3 mm or larger in diameter) colony (or 4 – 5 smaller colonies of the same morphology) from a blood plate such as Trypticase Soy Agar with 5% Sheep Blood. Use of a MacConkey Agar plate is also acceptable.
3. Suspend colonies in a tube of BBL Crystal Enteric/Stool Inoculum Fluid.
4. Recap tube and vortex for approximately 10 – 15 s.
5. Take a base, and mark the patient’s specimen number on the side wall.
6. Pour entire contents of inoculum fluid into target area of the base. See Fig. B.
7. Hold base in both hands and roll inoculum gently along the tracks until all of the wells are filled. Roll *back* any excess fluid to the target area and place the base on a bench top. See Fig. C.
8. Align the lid so that the labeled end of the lid is on top of the target area of the base. See Fig. D.
9. Push down until a slight resistance is felt. Place thumb on edge of lid towards middle of panel on each side and push downwards simultaneously until the lid snaps into place (listen for two “clicks”). See Fig. E.

**Purity Plate:** Using a sterile loop, recover a small drop from the inoculum fluid tube either before or after inoculating the base and inoculate an agar slant or plate (any appropriate media) for purity check. Discard inoculum fluid tube and cap in a biohazard disposal container. Incubate the slant or plate for 18 – 24 h at 35 – 37°C in a non-CO<sub>2</sub> incubator. The purity plate or slant may also be used for any supplementary tests or serology, if required.

**Incubation:** Place inoculated panels in incubation trays. Ten panels can fit in one tray (5 rows of 2 panels). All panels should be incubated **face down** (larger windows facing up; label facing down) in a non-CO<sub>2</sub> incubator with 40 – 60% humidity. Trays should not be stacked more than two high during incubation. The incubation time for E/NF panels is **18 – 20 h** at 35 – 37°C. See Fig. F.

**Reading:** After the recommended period of incubation, remove the panels from the incubator. All panels should be read **face down** (larger windows up; label facing down) using the BBL Crystal Light Box or Panel Viewer. See Fig. G. Refer to the color reaction chart and/or the chart in the “Reagents” section for an interpretation of the reactions. Use the BBL Crystal E/NF Report Pad to record reactions. Alternatively, the BBL Crystal AutoReader may be used to read the panels.

**Calculation of BBL Crystal Profile Number:** Each test result that is scored positive is given a value of 4, 2, or 1, corresponding to the row where the test is located. A value of 0 (zero) is given to any negative result. The numbers (values) resulting from each positive reaction in each column are then added together. A 10-digit number is generated; this is the profile number.

**Example:**

	A	B	C	D	E	F	G	H	I	J
4	+	+	+	–	–	+	+	–	+	–
2	–	–	+	–	+	–	–	+	+	–
1	+	–	–	–	–	–	–	+	+	+
Profile	5	4	6	0	2	4	4	3	7	1

The resulting profile number and off-line test results (indole and oxidase) should be entered on a PC in which the BBL Crystal ID System Electronic Codebook has been installed, to obtain the identification. A manual codebook is also available. If a PC is not available contact BD Technical Services for assistance with the identification. If using the BBL Crystal AutoReader, organisms are automatically identified by the PC.

**User Quality Control:** Quality control testing is recommended for each lot of panels as follows –

1. Set up a **BBL Crystal E/NF** panel with *Klebsiella pneumoniae* ATCC™ 33495 per recommended procedure (refer to “Test Procedure”).
2. Incubate panel for 18 – 20 h at 35 – 37°C.
3. Read panel with **BBL Crystal** Light Box or Panel Viewer and **BBL Crystal E/NF Color Reaction Chart**; record reactions using the **BBL Crystal E/NF Report Pad**. Alternatively, read the panel on the **BBL Crystal AutoReader**.
4. Compare recorded reactions with those listed in Table 2 (page 37). If discrepant results are obtained, confirm purity of quality control strain before contacting **BD Technical Services**.

Expected test results for additional quality control test strains are also listed in Table 2 (page 37).

#### LIMITATIONS OF THE PROCEDURE

The **BBL Crystal E/NF ID System** is designed for the E/NF taxa provided. Taxa other than those listed in Table 1 are not intended for use in this system.

**BBL Crystal** Identification Systems use a modified microenvironment; therefore, expected values for its individual tests may differ from information previously established with conventional test reactions. The accuracy of the **BBL Crystal E/NF Identification System** is based on statistical use of specially designed tests and an exclusive database.

When antisera are available, the biochemical identification of selected organisms, such as *Salmonella*, *Salmonella* subgroup 3, *Shigella*, enteropathogenic *Escherichia coli* A-D, and *Vibrio cholerae*, should be extended by antigenic analysis.<sup>9,16</sup>

Only cotton-tipped applicator swabs should be used to prepare the inoculum suspension as some polyester swabs may cause the inoculum fluid to become viscous. This may result in insufficient inoculum fluid to fill the wells. Once lids are removed from the sealed pouches they must be used within 1 h to ensure adequate performance. The plastic cover should remain on the lid until used.

The incubator where panels are placed should be humidified to prevent evaporation of inoculum fluid from the wells during incubation. The recommended humidity level is 40 – 60%.

The panels, after inoculation, should only be incubated **face down** (larger windows facing up; label facing down) to maximize the effectiveness of substrates.

Colonies should be taken from a blood agar plate such as **Trypticase Soy Agar with 5% Sheep Blood**. Use of a **MacConkey Agar** plate is also acceptable.

**BBL Crystal** Identification Systems are NOT for use directly with clinical specimens.

#### PERFORMANCE CHARACTERISTICS

**Reproducibility:** In an external study involving three (3) clinical laboratories, the reproducibility of E/NF substrates’ (30) reactions was studied by replicate testing. The reproducibility of individual substrate reactions ranged from 96.3 – 100%. The overall reproducibility of **BBL Crystal E/NF** panel was determined to be 99.6%.

**Accuracy of Identification:** The performance of **BBL Crystal E/NF ID System** was compared to currently available commercial systems using clinical isolates and stock cultures.

In an internal study, the performance of the **BBL Crystal E/NF** was evaluated. Results from 169 enteric and nonenteric isolates (representing 45 species) tested were analyzed. Discrepant identifications were resolved by the use of other commercial systems. These results are shown below:

N = 169	ID Without Supplemental Testing	ID With Supplemental Testing	NO ID or Misidentified
<b>BBL Crystal E/NF</b>	163 (96.4%)	167 (98.8%)	2 (1.2%)

The performance of the **BBL Crystal Enteric/Nonfermenter ID** test was evaluated in three independent clinical laboratories.<sup>13</sup> Both routine isolates arriving in the clinical laboratory as well as previously identified isolates of the clinical trial sites’ choice were utilized to establish performance characteristics.

Out of the 299 fresh clinical isolates tested by the laboratories’ current identification methods, the **BBL Crystal ID System** correctly reported 96.7% (289) including 16 instances where two or three organisms were reported and required supplemental testing to resolve.

Out of the 291 previously identified challenge strains confirmed by the laboratories’ current identification methods, the **BBL Crystal ID System** correctly reported 96.9% (282) including 8 instances where two or three organisms were reported and required supplemental testing to resolve.<sup>13</sup>

## AVAILABILITY


Cat. No.	Description	Cat. No.	Description
245000	<b>BBL Crystal™</b> Enteric/Nonfermenter ID Kit, containing 20 each: <b>BBL Crystal™</b> Enteric/NF Panel Lids, <b>BBL Crystal™</b> Bases, <b>BBL Crystal™</b> Enteric/Stool ID Inoculum Fluid tubes.	245002	<b>BBL Crystal™</b> Identification Systems Enteric/Nonfermenter Manual Codebook.
245031	<b>BBL Crystal™</b> Panel Viewer, Domestic model, 110 V, 60 Hz.	245029	<b>BBL Crystal™</b> Enteric/Stool ID Inoculum Fluid, ctn of 10.
245032	<b>BBL Crystal™</b> Panel Viewer, European model, 220 V, 50 Hz.	245300	<b>BBL Crystal™</b> AutoReader.
245033	<b>BBL Crystal™</b> Panel Viewer, Japanese model, 100 V, 50/60 Hz.	221239	<b>Trypticase™</b> Soy Agar with 5% Sheep Blood, pkg of 20 plates.
245034	<b>BBL Crystal™</b> Panel Viewer Longwave UV Tube.	221261	<b>Trypticase™</b> Soy Agar with 5% Sheep Blood, ctn of 100 plates.
245036	<b>BBL Crystal™</b> Panel Viewer White Light Tube.	261187	<b>BBL™</b> DMACA Indole Reagent Droppers, 50s.
		261181	<b>BBL™</b> Oxidase Reagent Droppers, 50s.

## REFERENCES

1. Edberg, S.C., and C.M. Kontnick. 1986. Comparison of  $\beta$ -glucuronidase-based substrate systems for identification of *Escherichia coli*. *J. Clin. Microbiol.* 24:368-371.
2. Kampfer, P., O. Rauhoff, and W. Dott. 1991. Glycosidase profiles of members of the family *Enterobacteriaceae*. *J. Clin. Microbiol.* 29:2877-2879.
3. Kilian, M., and P. Bulow. 1976. Rapid diagnosis of *Enterobacteriaceae* 1: detection of bacterial glycosidases. *Acta Pathol. Microbiol. Scand. Sect. B.* 84:245-251.
4. Manafi, M., W. Kneifel, and S. Bascomb. 1991. Fluorogenic and chromogenic substrates used in bacterial diagnostics. *Microbiol. Rev.* 55:335-348.
5. Muytjens, H. L., J. van der Ros-van de Repe, and H. A. M. van Druten. 1984. Enzymatic profiles of *Enterobacter sakazakii* and related species with special reference to the  $\alpha$ -glucosidase reactions and reproducibility of the test system. *J. Clin. Microbiol.* 20:684-686.
6. Sneath, P.H.A. 1957. The application of computers to taxonomy. *J. Gen. Microbiol.* 17:201-221.
7. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 1998. *Bailey and Scott's diagnostic microbiology*, 10th ed. Mosby, Inc., St. Louis.
8. Cowan, S.T., and K.J. Steel. 1974. *Manual for the identification of medical bacteria*. 2nd ed. Cambridge University Press, Cambridge, U.K.
9. Ewing, W.H. 1986. *Edwards and Ewing's identification of Enterobacteriaceae*, 4th ed. Elsevier Science Publishing Co., Inc., New York.
10. Le Minor, L. 1972. *Le Diagnostic de Laboratoire des Bacilles a Gram Negatif Enterobacteries*. Tom. 1, 4th ed. Editions de La Tourelle, St. Mande-94, France.
11. MacFaddin, J.F. 2000. *Biochemical tests for identification of medical bacteria*, 3rd Ed. Lippincott, Williams & Wilkins, Baltimore.
12. Ferguson, W.W., and A.E. Hook. 1943. Urease activity of *Proteus* and *Salmonella* organisms. *J. Lab. Clin. Med.* 28:1715-1720.
13. Data on file at BD Diagnostics.
14. Simmons, J.S. 1926. A culture medium for differentiating organisms of typhoid-colon-aerogenes groups and for isolation of certain fungi. *J. Infect. Dis.* 39:209-214.
15. Moeller, V. 1955. Simplified tests for amino acid decarboxylases and for arginine dihydrolase system. *Acta Pathol. Microbiol. Scand.* 36:158-172.
16. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1999. *Manual of clinical microbiology*, 7th ed. American Society for Microbiology, Washington, D.C.


## Miscellaneous Gram-Negative Bacilli

Chart No. 1 (Motile by Peritrichous Flagella)


## Miscellaneous Gram-Negative Bacilli

Chart No. 2 (Motile by Polar Flagella)

## Miscellaneous Gram-Negative Bacilli

Chart No. 3 (Nonmotile)

- References:
1. Gilardi, G.L., *Identification of Glucose-Nonfermenting Gram-Negative Rods*, 1/90
  2. Manual of Clinical Microbiology, 5th Edition, American Society for Microbiology, Washington, D.C., 1991

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
## Miscelánea de bacterias gram-negativas

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1. Gilardi, G.L., *Identification of Glucose-Nonfermenting Gram-Negative Rods*, 1/90
2. Manual of Clinical Microbiology, 5th Edition, American Society for Microbiology, Washington, D.C., 1991

Table / Tableau / Tabelle / Tabella / Quadro / Tabla 1

Taxa In BBL Crystal E/NF ID System / Taxonomie dans le système BBL Crystal E/NF ID / Im BBL Crystal E/NF-ID-System gespeicherte Taxa / Unità tassonomiche nel sistema BBL Crystal E/NF ID / Grupos taxonómicos no Sistema de ID de E/NF BBL Crystal / Grupos taxonómicos del sistema BBL Crystal E/NF ID

<i>Acinetobacter baumannii</i>	<i>Flavobacterium breve</i>	<i>Salmonella arizone</i>
<i>Acinetobacter lwoffii</i>	<i>Flavobacterium gleum</i>	<i>Salmonella choleraesuis</i>
<i>Aeromonas caviae</i>	<i>Flavobacterium indologenes</i>	<i>Salmonella paratyphi A</i>
<i>Aeromonas hydrophila</i>	<i>Flavobacterium meningosepticum</i>	<i>Salmonella species</i>
<i>Aeromonas sobria</i>	<i>Flavobacterium odoratum</i>	<i>Salmonella typhi</i>
<i>Aeromonas veronii</i>	<i>Hafnia alvei</i>	<i>Serratia ficaria</i>
<i>Agrobacterium tumefaciens</i>	<i>Klebsiella ornithinolytica</i>	<i>Serratia fonticola</i>
<i>Burkholderia cepacia</i>	<i>Klebsiella oxytoca</i>	<i>Serratia liquefaciens</i>
<i>Burkholderia pseudomallei</i>	<i>Klebsiella ozaenae</i>	<i>Serratia marcescens</i>
<i>Cedecea davisae</i>	<i>Klebsiella pneumoniae</i>	<i>Serratia odorifera 1</i>
<i>Cedecea lapagei</i>	<i>Klebsiella rhinoscleromatis</i>	<i>Serratia odorifera 2</i>
<i>Cedecea neteri</i>	<i>Kluyvera ascorbata</i>	<i>Serratia plymuthica</i>
<i>Chromobacterium violaceum</i>	<i>Kluyvera cryocrescens</i>	<i>Serratia rubidaea</i>
<i>Chryseomonas luteola</i>	<i>Leclercia adecarboxylata</i>	<i>Shewanella putrefaciens</i>
<i>Citrobacter amalonaticus</i>	<i>Moellerella wisconsensis</i>	<i>Shigella dysenteriae</i>
<i>Citrobacter freundii</i>	<i>Morganella morganii</i>	<i>Shigella species (S. boydii, S. flexneri)</i>
<i>Citrobacter koseri</i>	<i>Pantoea agglomerans</i>	<i>Shigella sonnei</i>
<i>Edwardsiella hoshinae</i>	<i>Pasteurella aerogenes</i>	<i>Sphingobacterium multivorum</i>
<i>Edwardsiella tarda</i>	<i>Pasteurella haemolytica</i>	<i>Stenotrophomonas maltophilia</i>
<i>Enterobacter aerogenes</i>	<i>Pasteurella multocida</i>	<i>Tatumella pytseos</i>
<i>Enterobacter asburiae</i>	<i>Plesiomonas shigelloides</i>	<i>Vibrio alginolyticus</i>
<i>Enterobacter cloacae</i>	<i>Proteus mirabilis</i>	<i>Vibrio cholerae</i>
<i>Enterobacter gergoviae</i>	<i>Proteus penneri</i>	<i>Vibrio damsela</i>
<i>Enterobacter sakazakii</i>	<i>Proteus vulgaris</i>	<i>Vibrio fluvialis</i>
<i>Enterobacter taylorae</i>	<i>Providencia alcalifaciens</i>	<i>Vibrio hollisae</i>
<i>Escherichia coli</i>	<i>Providencia rettgeri</i>	<i>Vibrio metschnikovii</i>
<i>Escherichia coli</i> serogroup / séogroupe / Serogruppe / sierogruppo / serogruppo O111	<i>Providencia rustigianii</i>	<i>Vibrio mimicus</i>
<i>Escherichia coli</i> serogroup / séogroupe / Serogruppe / sierogruppo / serogruppo O157	<i>Providencia stuartii</i>	<i>Vibrio parahaemolyticus</i>
<i>Escherichia coli</i> AD	<i>Pseudomonas aeruginosa</i>	<i>Vibrio vulnificus</i>
<i>Escherichia fergusonii</i>	<i>Pseudomonas diminuta</i>	<i>Weeksella virosa/zoohelcum</i>
<i>Escherichia hermannii</i>	<i>Pseudomonas fluorescens</i>	<i>Yersinia enterocolitica</i> group / groupe / Gruppe / gruppo / grupo
<i>Escherichia vulneris</i>	<i>Pseudomonas gladioli</i>	( <i>Y. enterocolitica</i> , <i>Y. frederiksenii</i> , <i>Y. intermedia</i> , <i>Y. kristensenii</i> )
<i>Ewingella americana</i>	<i>Pseudomonas paucimobilis</i>	<i>Yersinia pseudotuberculosis</i>
<i>Flavimonas oryzzihabitans</i>	<i>Pseudomonas putida</i>	<i>Yokenella regensburgei</i>
	<i>Pseudomonas stutzeri</i>	Miscellaneous Gram-Negative Bacilli <sup>1</sup>
	<i>Pseudomonas vesicularis</i>	
	<i>Rahnella aquatilis</i>	

<sup>1</sup> "Miscellaneous Gram-Negative Bacilli" refers to a group of oxidase positive species that are relatively inactive and indistinguishable from each other in **BBL Crystal** Enteric/Nonfermenter ID System. Refer to Tables 1 and 2 provided in this insert for further identification when the first choice identification is "Miscellaneous Gram-Negative Bacilli."

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Benex Limited  
Rineanna House  
Shannon Free Zone  
Shannon, County Clare, Ireland