
Revisions


Rev from	Rev to	ECO #
0603	0604	2871-04

Notes:


- BD Cat. Number 240658, 240661, 240679, 240827, 240826, 240680
- Blank (Sheet) Size : Length: 11" Width: 34"
 Number of Pages: 16 Number of Sheets: 1
 Page Size: Length 11" Width 4.25" Final Folded Size: 1.375" x 4.25"
- Style (see illustrations below): # 4


#1


#2


#3


#4

#5

#6

#7

- See Specification Control Number N/A for Material Information
- Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# 2755, Blue
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 8810061		Category and Description Package Insert, BBL Coagulase Plasmas	Sheet: 1 of 17 Scale: 1:1	A

BD BBL™ Coagulase Plasmas
8810061
2004/06

English: pages 1 - 4 **Italiano: pagine 9 - 12**
Français : pages 4 - 6 **Español: páginas 12 - 15**
Deutsch: Seiten 7 - 9

See symbol glossary at end of insert. / Viz popis symbolů na konci příbalového letáku. / Se symbolglossaret i slutningen af indlægsedlen. / Zie lijst met symbolen aan het einde van de bijsluiter. / Vaadake sümbole seletust infolehe lõpus. / Katsopakkauselosteen lopussa olevaa kuvamerkkien sanastoa. / Voir le glossaire des symboles à la fin de la notice. / Siehe Symbol-Erklärungen am Ende der Packungsbeilage. / Δείτε το γλωσσάριο των συμβόλων στο τέλος του ένθετου. / A jelmagyarázat a használati utasítás végén található. / Vedere il glossario dei simboli alla fine del foglio illustrativo. / Žr. informacinio lapelio pabaigoje pateikiamą simbolių glosarijų. / Se i symbolforklaringen på slutten av produktvedlegget. / Zobacz objaśnienie symboli na końcu ulotki. / Consulte o glossário de símbolos no fim do folheto informativo. / Pozri slovník symbolov na konci letáka. / Consulte el glosario de símbolos al final del prospecto. / Se symbolförteckningen vid slutet av bipacksedeln.

Contact your local BD representative for instructions. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar.

INTENDED USE

BBL™ Coagulase Plasma, Rabbit and **BBL** Coagulase Plasma, Rabbit with EDTA are used to qualitatively determine the pathogenicity of staphylococci using the direct tube method.

SUMMARY AND EXPLANATION

Identification of staphylococci is based on microscopic examination, colonial morphology and cultural and biochemical characteristics. Staphylococci associated with acute infection (*Staphylococcus aureus* in humans and *S. intermedius* and *S. hyicus* in animals) can clot plasma. The most widely used and generally accepted criterion for identification of these pathogenic organisms is based on the presence of the enzyme coagulase.¹ The ability of *Staphylococcus* to produce coagulase was first reported by Loeb² in 1903.

Coagulase binds plasma fibrinogen, causing the organisms to agglutinate or plasma to clot. Two different forms of coagulase can be produced, free and bound. Free coagulase is an extracellular enzyme produced when the organism is cultured in broth. Bound coagulase, also known as clumping factor, remains attached to the cell wall of the organism. The tube test can detect the presence of both bound and free coagulase. Isolates that do not produce clumping factor must be tested for the ability to produce extracellular coagulase (free coagulase).

BBL Coagulase Plasma, Rabbit and **BBL** Coagulase Plasma Rabbit with EDTA are recommended for performing the direct tube test. The inoculum used for testing must be pure because a contaminant may produce false results after prolonged incubation. For the coagulase test, **BBL** Coagulase Plasma, Rabbit with EDTA is superior to citrated plasma because citrate-utilizing organisms such as *Pseudomonas* species, *Serratia marcescens*, *Enterococcus faecalis* and strains of *Streptococcus* will clot citrated plasma in 18 h.³

PRINCIPLES OF THE PROCEDURE

S. aureus produces two types of coagulase, free and bound. Free coagulase is an extracellular enzyme produced when the organism is cultured in broth. Bound coagulase, also known as the clumping factor, remains attached to the cell wall of the organism.

In the direct tube test, free coagulase liberated from the cell acts on prothrombin in the coagulase plasma to give a thrombin-like product. This product then acts on fibrinogen to form a fibrin clot.⁴

The tube test is performed by mixing an overnight broth culture or colonies from a non-inhibitory agar plate into a tube of rehydrated coagulase plasma. The tube is incubated at 37°C. The formation of a clot in the plasma indicates coagulase production.

REAGENTS

Coagulase Plasma, Rabbit is lyophilized rabbit plasma with 0.85% sodium citrate and 0.85% sodium chloride, approximately.

Coagulase Plasma, Rabbit with EDTA is lyophilized rabbit plasma with 0.15% EDTA (ethylenediaminetetraacetic acid) and 0.85% sodium chloride, approximately.

Warnings and Precautions

For *in vitro* Diagnostic Use.

The Packaging of This Product Contains Dry Natural Rubber.

Observe aseptic technique and established precautions against microbiological hazards throughout all procedures. After use, specimens, containers, slides, tubes and other contaminated material must be sterilized by autoclaving.

Directions for use should be followed carefully.

Storage

Store unopened lyophilized **BBL** Coagulase Plasma, Rabbit and **BBL** Coagulase Plasma, Rabbit with EDTA at 2-8°C.

Store reconstituted plasma at 2-8°C for up to 14 days, or aliquot and freeze promptly at -20°C for up to 30 days. Do not thaw and refreeze.

Expiration date applies to product in its intact container when stored as directed. Do not use if the product is caked, discolored or shows other signs of deterioration. Examine reconstituted reagents for evidence of contamination, evaporation or other signs of deterioration, such as cloudiness or partial clotting.

SPECIMEN COLLECTION AND PREPARATION

Collect specimens or samples in sterile containers or with sterile swabs and transport immediately to the laboratory according to recommended guidelines.^{1,4-9}

Process each specimen using procedures appropriate for that sample.^{1,4-9}

Select well-isolated colonies. The test described below requires the use of a pure test culture.

Suspicious growth, such as black colonies on Vogel and Johnson Agar or Tellurite Glycine Agar, or golden, hemolytic colonies from **Trypticase**[™] Soy Blood agar plates should be selected for testing.

Using a bacteriological loop, transfer a well-isolated colony from a pure culture into a tube of sterile Brain Heart Infusion Broth or **Trypticase** Soy Broth. Incubate for 18-24 h or until a dense growth is observed. Alternatively, 2-4 colonies (1 loopful) taken directly from a non-inhibitory agar plate such as **Trypticase** Soy Agar may be used as an inoculum instead of a broth culture.

PROCEDURE

Materials Provided: BBL Coagulase Plasma, Rabbit, BBL Coagulase Plasma, Rabbit with EDTA.

Materials Required But Not Provided: Bacteriological inoculating loop, Pipettes, Sterile purified water, Culture tubes, small (10 x 75 mm), Water bath or incubator (35-37°C), **Trypticase** Soy Broth or Brain Heart Infusion (BHI) Broth.

Reagent Preparation

Rehydrate **BBL** Coagulase Plasma, Rabbit and **BBL** Coagulase Plasma, Rabbit with EDTA by adding sterile purified water to the vial as indicated below. Mix by gentle end-over-end rotation of the vial.

Product Size	Sterile Purified Water	Approximate Number of Tests
3 mL	3 mL	6
15 mL	15 mL	30
25 mL	25 mL	50

Test Procedure

- Using a sterile 1 mL pipette, add 0.5 mL of rehydrated **BBL** Coagulase Plasma, Rabbit or **BBL** Coagulase Plasma, Rabbit with EDTA to a 10 x 75 mm test tube supported in a rack.
- Using a sterile 1 mL serological pipette, add approximately 0.05 mL of the overnight broth culture of the test organism to the tube of plasma. Alternatively, using a sterile bacteriological loop, thoroughly emulsify 2-4 colonies (1 loopful) from a non-inhibitory agar plate in the tube of plasma.
- Mix gently.
- Incubate in a water bath or incubator at 35-37°C for up to 4 h.
- Examine the tubes periodically by gently tipping the tube. Avoid shaking or agitating the tube, which could cause breakdown of the clot and, consequently, doubtful or false negative test results. Any degree of clotting in a 3-4 h period is regarded as a positive result. Many weak enzyme-producing strains will coagulate the plasma only after 24 h of incubation.
- Record results.

User Quality Control


At the time of use, test both positive and negative control cultures to check performance of the coagulase plasma, techniques and methodology. The following cultures listed are the minimum that should be used for performance testing.

Organism	ATCC [™]	Reaction
<i>Staphylococcus aureus</i>	25923	Clot in tube
<i>Staphylococcus epidermidis</i>	12228	No clot in tube

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

Results

Any degree of clotting in **BBL** Coagulase Plasma, Rabbit or **BBL** Coagulase Plasma, Rabbit with EDTA is considered a positive test. The following chart can be used as a guide in interpreting the reactions.

Negative	No evidence of fibrin formation
1 + Positive	Small unorganized clots
2 + Positive	Small organized clot
3 + Positive	Large organized clot
4 + Positive	Entire content of tube coagulates and is not displaced when tube is inverted.

LIMITATIONS OF THE PROCEDURE

- Some species of organisms utilize citrate in their metabolism and will yield false-positive reactions for coagulase activity. Normally, this does not cause problems since the coagulase test is performed almost exclusively on staphylococci. However, it is possible that bacteria that utilize citrate may contaminate *Staphylococcus* cultures on which the coagulase test is being performed. These

- contaminated cultures may, upon prolonged incubation, give false-positive results due to citrate utilization.⁴
- When checking results of the coagulase test, observe tubes hourly during the first 4 h of incubation. Some strains of *S. aureus* produce staphylokinase, which may lyse clots. If the tubes are not read until 24 h of incubation, false-negative results may occur.¹
 - Do not use plasma if a heavy precipitate or clot has formed before inoculation.

PERFORMANCE CHARACTERISTICS^{10,11}

The performance of **Bacto**[™] Coagulase Plasma (now **BBL** Coagulase Plasma) was compared to four other tests for the identification of *Staphylococcus aureus* in a study by Ad Luijenkijk, van Belkum, Verbrugh and Kluytmans.¹⁰ The free-coagulase (tube) test was performed. In addition, the bound coagulase (agar) test and three commercial latex agglutination tests were used to identify isolates.

Of the 330 staphylococcal isolates tested, 300 were *S. aureus* and 30 were non-*S. aureus*. All of the tests produced negative results for the 30 non-*S. aureus* isolates for a specificity of 100%. The table below summarizes the sensitivity of each test system for the 300 *S. aureus* isolates.

Test	MSSA* (222 isolates)		MRSA** (78 isolates)		Total (300 isolates)	
	No. of false-negative results	Test Sensitivity (%)	No. of false-negative results	Test Sensitivity (%)	No. of false-negative results	Test Sensitivity (%)
Free Coagulase	0	100	6	92.3	6	98.0
Bound Coagulase	0	100	3	96.1	3	99.0
Latex Agglutination Test #1	0	100	0	100	0	100
Latex Agglutination Test #2	2	99.1	12	84.6	14	95.3
Latex Agglutination Test #3	0	100	0	100	0	100

*Methicillin-susceptible *S. aureus*.

**Methicillin-resistant *S. aureus*.

In a second study by McDonald and Chapin,¹¹ the performance of **BBL** Coagulase Plasma was compared in a 2-h tube coagulase test (TCT) to two commercial latex agglutination tests for identifying *S. aureus* directly from blood culture broths and pellets obtained from supernatants of **BACTEC**[™] bottles. One hundred twelve (112) clinical blood culture isolates and 68 negative blood culture bottles seeded with a variety of gram-positive organisms were evaluated.

The table below gives the results of the coagulase test and the latex agglutination tests for both the seeded and clinical specimens.

Cultures and organism(s)	No. aerobic/ No. anaerobic (total)	No. of positive results					
		Direct			Pellet		
		Latex Test #1	Latex Test #2	TCT	Latex Test #1	Latex Test #2	TCT
Seeded							
<i>Staphylococcus aureus</i>	11/8 (19)	0	0	19	2	2	19
Coagulase-negative staphylococci	9/9 (18)	0	0	0	0	0	0
<i>Streptococcus pneumoniae</i>	5/5 (10)	0	0	0	0	0	0
<i>Enterococcus</i> spp.	6/5 (11)	5	5	0	5	5	0
<i>Streptococcus agalactiae</i>	3/5 (8)	0	0	0	0	0	0
<i>Streptococcus pyogenes</i>	1/1 (2)	0	0	0	0	0	0
Clinical							
Coagulase-negative staphylococci	70/3 (73)	0	0	0	0	0	0
<i>Staphylococcus aureus</i>	35/4 (39)	5	4	31	5	4	30

Among the 68 seeded blood culture bottles, the 2-h tube coagulase test using **BBL** Coagulase Plasma correctly identified 19 of 19 blood cultures seeded with *S. aureus*. There were no false positives with the 2-h tube coagulase test. Among the 112 clinical specimens tested, the 2-h tube coagulase test correctly identified 31 out of 39 *S. aureus* isolates directly in blood culture broth and 30 out of 39 *S. aureus* isolates in pelleted supernatants, for sensitivities of 79.5 and 76.9%, respectively. The specificity with both seeded and clinical isolates was 100% for the tube coagulase test.

AVAILABILITY

Cat. No. Description

- 240658 **BBL**[™] Coagulase Plasma, Rabbit, 10 X 3.0 mL
 240661 **BBL**[™] Coagulase Plasma, Rabbit, 10 X 15.0 mL
 240679 **BBL**[™] Coagulase Plasma, Rabbit, 10 X 25.0 mL
 240827 **BBL**[™] Coagulase Plasma, Rabbit with EDTA, 10 X 3.0 mL
 240826 **BBL**[™] Coagulase Plasma, Rabbit with EDTA, 10 X 15.0 mL
 240680 **BBL**[™] Coagulase Plasma, Rabbit with EDTA, 10 X 25.0 mL

REFERENCES

- Kloos, W. E., and T. L. Bannerman. 1999. *Staphylococcus* and *Micrococcus*, p. 264-282. In P.R. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover and R.H. Tenover, Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.
- Loeb, L. 1903. The influence of certain bacteria on the coagulation of the blood. *J. Med. Res.* 10:407-419.
- Bayliss, B.G. and E.R. Hall. 1965. Plasma coagulation by organisms other than *Staphylococcus aureus*. *J. Bacteriol.* 89:101-104.
- Pezzlo, M. (ed.). 1994. Aerobic bacteriology, p. 1.0.0.-1.20.47. In H. D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
- Baron, E.J., L.R. Peterson and S.M. Finegold. 1994. Bailey & Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis, MO.

6. Association of Official Analytical Chemists. 2000. Official methods of analysis of AOAC International, 17th ed. AOAC International, Arlington, VA.
7. Association of Official Analytical Chemists. 2001. FDA Bacteriological analytical manual online. <<http://www.cfsan.fda.gov/~ebam/bam-mm.html>>.
8. Downes, F.P. and K. Ito (ed.). 2001. Compendium of methods for the microbiological examination of foods, 4th ed. American Public Health Association, Washington, D.C.
9. Flowers, R.S., W. Andrews, C.W. Donnelly and E. Koenig. 1993. Pathogens in milk and milk products, p. 103-212. In R.T. Marshall (ed.), Standard methods for the examination of dairy products, 16th ed. American Public Health Association, Washington, D.C.
10. Luijendijk, A., A. van Belkum, H. Verbrugh and J. Kluytmans. 1996. Comparison of five tests for identification of *Staphylococcus aureus* from clinical samples. J. Clin. Microbiol. 34:2267-2269.
11. McDonald, C.L. and K. Chapin. 1995. Rapid Identification of *Staphylococcus aureus* from blood culture bottles by a classic 2-hour tube coagulase test. J. Clin. Microbiol. 33:50-52.

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DISPONIBILIDAD

N° de cat. Descripción

- 240658 **BBL** Coagulase Plasma, Rabbit, 10 X 3,0 mL
 240661 **BBL** Coagulase Plasma, Rabbit, 10 X 15,0 mL
 240679 **BBL** Coagulase Plasma, Rabbit, 10 X 25,0 mL
 240827 **BBL** Coagulase Plasma, Rabbit with EDTA, 10 X 3,0 mL
 240826 **BBL** Coagulase Plasma, Rabbit with EDTA, 10 X 15,0 mL
 240680 **BBL** Coagulase Plasma, Rabbit with EDTA, 10 X 25,0 mL

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ATCC is a trademark of the American Type Culture Collection.

Bacto is a trademark of Difco Laboratories, Inc., a subsidiary of Becton, Dickinson and Company.

BD, BD Logo, BACTEC, BBL and Trypticase are trademarks of Becton, Dickinson and Company. © 2004 BD.

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BENEX Limited