
## Revisions

Rev from	Rev to	ECO #
0304	1206	4258-06

**Notes:**

1. BD Cat. Number 244003
2. Blank (Sheet) Size: Length: Horizontal: N/A Vertical: N/A  
 Number of Pages: N/A Number of Sheets: N/A  
 Page Size: Length Horizontal N/A Vertical N/A Final Folded Size: N/A
3. Style (see illustrations below): N/A

4. See Specification Control Number 244003 for Material Information
5. Ink Colors: Printed two sides  Yes  No  
 No. of Colors: N/A
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	<b>Becton, Dickinson and Company</b> 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 88102025		Category and Description Package Insert, GonoGen II	Sheet: 1 of 9 <hr/> Scale: N/A	A

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## INTENDED USE

GonoGen™ II is a monoclonal antibody-based colorimetric test intended for the confirmatory identification of *Neisseria gonorrhoeae* from culture.

## SUMMARY AND EXPLANATION

Infection with *N. gonorrhoeae* requires different treatment than infections with other *Neisseria* species. Infections with nonpathogenic *Neisseria*, such as *N. flava*, *N. sicca* or *N. subflava*, usually require no treatment, whereas infections with pathogens such as *N. meningitidis* may require different antibiotic therapy than that for *N. gonorrhoeae*. In order to differentiate *N. gonorrhoeae* from other *Neisseria* species, it is necessary to perform a confirmatory test. The sugar utilization test is the standard reference method, but it requires prolonged incubation and pure, viable isolates before unequivocal results can be obtained. The GonoGen II assay for the confirmation of *N. gonorrhoeae* is more rapid than sugar utilization tests and does not require pure or viable isolates.

## PRINCIPLES OF THE PROCEDURE

The GonoGen II test for *N. gonorrhoeae* is comprised of a specific anti-gonococcal reagent. The specific antigonococcal reagent is composed of a pool of murine monoclonal antibodies (IgG) that have been prepared against a purified outer membrane protein, Protein I, of *N. gonorrhoeae*. Protein I is a major protein molecule that is exposed on the surface of the gonococcus and its epitopes are largely responsible for serotype specific reactions of the gonococcus.<sup>1-3</sup> By including monoclonal antibodies to the various serotypes of *N. gonorrhoeae*, maximum specificity and sensitivity are achieved. These antibodies are absorbed to suspended metal sol particles, which give the reagent its raspberry red color.

When a culture of *N. gonorrhoeae* is suspended in the solubilizing buffer, the outer membrane is stripped from the organism, releasing Protein I-containing complexes into solution, enabling these complexes to be bound by the antibody-sol particles. When the solution is then passed through the special matrix test device, the Protein I-sol particle complexes bind to the matrix, resulting in a color change. Sol particles which have not bound Protein I will yield a negative test (white to pale pink spot).

## REAGENTS

- Reagent 1** GonoGen II solubilizing buffer. Organisms from a suspected isolate of *N. gonorrhoeae* are suspended in this buffer prior to addition of GonoGen II reagent. Contains 0.09% sodium azide (preservative).
- Reagent 2** GonoGen II murine monoclonal antibodies to the Protein I antigens of *N. gonorrhoeae* that have been absorbed to metal sol particles. Contains 0.05% sodium azide (preservative).
- Control +** GonoGen II heat-killed *N. gonorrhoeae*. When mixed with the GonoGen II reagent and applied to a well on the reaction device, a positive reaction is visible as a pink to red dot. Contains 0.05% sodium azide and 0.01% gentamicin sulfate (preservatives).

**Control –** GonoGen II heat-killed *Neisseria* species other than *N. gonorrhoeae*. When mixed with the GonoGen II reagent and applied to a well on the test tray, no reaction will be visible. Contains 0.05% sodium azide and 0.01% gentamicin sulfate (preservatives).

**Test Tray:** Consists of wells with a special matrix and absorbent material. When the GonoGen II sample reactant is added, a positive (pink to red dot) or negative (white to pale pink dot) result is observed on the matrix.

## Warnings and Precautions:

For *in vitro* Diagnostic Use.

This product contains dry natural rubber.

**Reagents:** Do not use beyond the expiration date. Upon removal from the refrigerator, allow reagents to warm to room temperature before use. Do NOT mix reagents from different kit lot numbers. Prior to use reagents should be vigorously shaken or vortexed for 10 s prior to use.

**Warning:** Reagents contain sodium azide which may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide buildup.

**Controls:** Do not use the kit if positive and negative controls do not yield appropriate results.

To assure proper drop delivery, the dropper rod assembly, transfer pipet and reagent bottles should be held vertically, dispensing free-falling drops.

**Storage of Reagents:** Upon receipt, refrigerate reagents at 2 to 8°C. DO NOT FREEZE. Reagents should be recapped and returned to refrigeration when not in use, taking care not to mix color-coded caps.

Observe established precautions against micro-biological hazards when performing procedures, and during disposal of reagents/tests.

## SAMPLE PREPARATION

Clinical samples should be processed as quickly as possible. Specimens should be grown on selective, enriched media such as GC-Lect™, modified Thayer-Martin or Martin-Lewis media to ensure growth of isolates. Inoculated media should be incubated at 35 to 37°C in a humid 3 to 10% CO<sub>2</sub> atmosphere for 18 to 48 h. Suspected organisms are inspected for typical colonial morphology, Gram stain appearance and oxidase reaction (BBL™ Oxidase Reagent Droppers, see "Availability").

## SPECIMEN COLLECTION

Do not use calcium alginate swabs to transfer suspected organisms from the culture media to the buffer.

Colonies that have grown on selective or enriched plated media which are oxidase positive and appear as gram-negative diplococci can be considered to be presumptively identified as *Neisseria* species. These are then tested to confirm them as *N. gonorrhoeae* with the GonoGen II test.

If there is sufficient presumptive growth, the primary culture may be used to perform the GonoGen II test. Whether subcultures or primary subcultures are used for testing, viable cultures are not needed for testing. Any cultures incubated within 18 to 48 h can be tested with equal confidence.

**PROCEDURE**

Materials Provided:	#244003 (24 Determin.)	#244004 (40 Determin.)
<b>Reagent 1</b> GonoGen II Buffer	27 mL	27 mL
<b>Reagent 2</b> GonoGen II Reagent	1 mL	2 mL
<b>Control +</b> GonoGen II Positive Control	1 mL	2 mL
<b>Control -</b> GonoGen II Negative Control	1 mL	2 mL
Test Tray and accessories.	3	5

**Materials Required But Not Provided:** Test tubes (12 x 75 mm), test tube rack, cotton/polyester swab or loop, McFarland No. 1 turbidity standard.

Also required are the necessary equipment and labware used for preparation, storage and handling of clinical specimens.

**Test Procedure:** Review "Precautions", "Sample Preparation" and "Specimen Collection" prior to performing procedures. The testing area, reagents and test components should be at room temperature when used.

- Label a test tube (12 x 75 mm) for each specimen.
  - Using the dropper rod assembly provided, dispense 500 µL (demarcation line) of **Reagent 1** into each tube.
  - Using a cotton/polyester swab, make a suspension of test colonies (approximately 30 colonies) to match a McFarland No. 1 turbidity standard (barely visible turbidity).
- NOTE:** When using an inoculating loop to remove bacteria, dispense 300 µL of **Reagent 1** into a test tube, suspend bacteria to a McFarland No. 1 turbidity reference, vortex and perform the test as described, starting with step 6.
- Press the swab against the inside of the tube to express as much liquid as possible.
  - Discard the swab in disinfectant or appropriate biohazard container.
  - Vigorously shake or vortex **Reagent 2**, and add 1 drop into each of the tubes to be tested.
  - Mix well.
  - Allow tubes to set for at least 5 min.
  - With a transfer pipet, add 2 drops of each test suspension to be tested into a well of the test tray, using a separate well for each test.
  - Using a clean transfer pipet, add 1 drop of **Reagent 1** to each completed test well.

- Interpret results:  
Pink to red dot in well of test tray = *N. gonorrhoeae*.  
White to pale pink dot in well of test tray – NOT *N. gonorrhoeae*.

**NOTE:** A color reaction more intense than the negative control should be interpreted as positive.

If color reaction is questionable, reincubate tube at room temperature 3 min and repeat test.

**CAUTION:** If specimen suspension is made too turbid, a faint background color will occur. This should not be interpreted as a positive reaction.

- Properly dispose of all materials used.

**NOTE:** If all eight wells of the test tray are not used during a given test period, the unused wells can be used at a later time. Reacted test trays may be saved as a permanent record.

**User Quality Control:** Controls should be tested each day the GonoGen II test is performed for patients, to ensure the system is functioning properly. The controls may be performed along with the test specimen.

The **Control + (Positive)** and **Control - (Negative)** should be performed each day a test is used.

- Label a small test tube (12 x 75 mm) with **Control +**.
- Label a small test tube (12 x 75 mm) with **Control -**.
- Dispense 500 µL of **Reagent 1** into each of these tubes.
- Add 1 drop of well-mixed **Control +** into the tube marked **Control +**.
- Add 1 drop of well-mixed **Control -** into the tube marked **Control -**.
- Mix well by shaking vigorously.
- Add 1 drop of **Reagent 2** into the **Control +** tube, and 1 drop into the **Control -** tube.
- Mix and wait at least 5 min.
- Add 2 drops of the **Control +** into a well in the test tray.
- Add 2 drops of the **Control -** into a separate well in the test tray.
- Using a clean transfer pipet, add 1 drop of **Reagent 1** to each completed test well.

**Reading Results:**

Pink to red spot in the test tray well = **Positive**.

White to pale pink in test tray well = **Negative**.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI (formerly NCCLS) guidance and CLIA regulations for appropriate Quality Control practices.

**LIMITATIONS OF THE PROCEDURE**

No single diagnostic test result should be considered conclusive in diagnosing disease. The overall clinical and laboratory findings should be taken into consideration before a physician renders a diagnosis. Depending upon exposed antigenic sites and antigenic composition, some gonococci may not be identifiable with GonoGen II reagent, and others may vary in color intensity. In the rare case of an extremely weak or non-specific reaction with GonoGen II, confirmation by other methods, such as carbohydrate utilization, may be necessary.

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Becton, Dickinson and Company. ©2006 BD.
GonoGen is a trademark of New Horizons Diagnostics Corp.
BBL, BD, BD Logo and GC-Lect are trademarks of