
Revisions


Rev from	Rev to	ECO #
2009/11	2010/10	5571-10

Notes:


- BD Cat. Number 245127
- Blank (Sheet) Size : Length: 7 7/8" Width: 8 1/2"
 Number of Pages: 32 Number of Sheets: 8
 Page Size: Length 7 7/8" Width 4 1/4" Final Folded Size: 7 7/8" x 4 1/2"
- Style (see illustrations below): # 5


#1


#2


#3


#4

#5

#6

#7

- See Specification Control Number N/A for Material Information
- Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# Black
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 88-2041-1JAA		Category and Description Package Insert, BACTEC MGIT 960 SIRE EMB 7.5	Sheet: 1 of 33 Scale: N/A	A

BD BACTEC™ MGIT™ 960 SIRE Kit For the Antimycobacterial Susceptibility Testing of *Mycobacterium tuberculosis*

English: pages 1 – 6 Italiano: pagine 17 – 22
Français : pages 7 – 12 Español: páginas 22 – 27
Deutsch: Seiten 12 – 17

 88-2041-1JAA
2010/10

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / ΕΠΙΚΟΙΝΩΝΗΣΤΕ ΜΕ ΤΟΝ ΤΟΠΙΚΟ ΑΝΤΙΠΡΟΣΩΠΟ ΤΗΣ ΒD ΓΙΑ ΟΔΗΓΙΕΣ. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcję użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / **Свържете се с местния представител на BD за инструкции.** / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / **Обратите се свом локалном представнику компаније BD за упутства.** / Для получения инструкций свяжитесь с местным представителем компании BD. / **Өзіндік жергілікті БД өкіліне жүгініп нұсқау алыңыз.** / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

The **BACTEC™ MGIT™ 960 SIRE Kit** is used as a rapid qualitative procedure for susceptibility testing of *Mycobacterium tuberculosis*, from culture, to streptomycin, isoniazid, rifampin and ethambutol using the **BACTEC™ MGIT™ 960** and **BACTEC™ MGIT™ 320 Systems**.

SUMMARY AND EXPLANATION

Antimycobacterial susceptibility testing is necessary for the proper treatment of patients with tuberculosis. The treatment of tuberculosis is commonly through a multiple drug regimen which includes the antimycobacterial drugs streptomycin, isoniazid, rifampin and/or ethambutol. It is important that the antimycobacterial drugs prescribed show appropriate activity against *Mycobacterium tuberculosis*, i.e., susceptibility of the isolate to the drug.

Multidrug resistant *Mycobacterium tuberculosis* (MDR-TB) has recently become a serious public health problem.¹ Resistance to any of the four primary drugs, streptomycin (STR), isoniazid (INH), rifampin (RIF), and ethambutol (EMB), makes the disease more difficult and expensive to treat. The rapid detection of these strains is critical to the effective treatment of the patient.

Two methods have been widely used for antimycobacterial susceptibility testing. The first method, known as the Method of Proportion,² uses Middlebrook and Cohn 7H10 Agar. It compares colony counts on drug-containing and drug-free media. Resistance to a drug is detected when 1% or more of the bacterial population is resistant to the drug concentration under test. Results are generally available after 21 days of incubation. The second method, known as the **BACTEC™ 460TB** radiometric susceptibility method,³ generally takes from 4 to 12 days. It is based on the production of radioactive ¹⁴C-labeled carbon dioxide by the growing mycobacteria, manifested by a growth index increase in the system.

Historically, the Method of Proportion (MOP) procedure has included susceptibility testing of *Mycobacterium tuberculosis* using two concentrations of antimicrobials. The National Committee for Clinical Laboratory Standards (NCCLS) continues to recommend that the MOP test procedure include two concentrations of the primary drugs for testing except rifampin. The recommended low concentrations for the MOP procedure are generally considered to be the critical concentrations for these drugs. The critical concentration is defined as the drug concentration which inhibits the wild type sub-population while allowing sufficient growth of the mutant resistant sub-population in order to determine resistance at the critical proportion of 1%. The high drug concentration is not considered to be the critical concentration. However, resistance at the high concentration is evidence that resistance is generally distributed among the population of the test *M. tuberculosis* strain. Some clinicians use the susceptibility test results at the high concentrations to profile the degree of resistance of the test strain.

The **BACTEC™ MGIT™ 960 SIRE** test provides the susceptibility result in approximately the same time frame as the **BACTEC™ 460TB** system. In addition, this method is non-radiometric and allows appropriate antibiotic susceptibility results to be reported earlier, in most cases, than with the MOP procedure.

The **BACTEC MGIT** system has been developed to allow susceptibility testing at the critical concentration for streptomycin, isoniazid, rifampin and ethambutol, and at a higher concentration for streptomycin, isoniazid, and ethambutol. These concentrations correlate with the two concentrations used in the MOP procedure. A susceptible result at the critical concentration can be reported and no other tests are necessary. However, any strain determined to be resistant with the **BACTEC MGIT 960 SIRE** kit to streptomycin, isoniazid, and ethambutol at the critical concentration may be, at a minimum, tested at the high concentration. In this case, a final result of resistant at the critical concentration may be reported, with notification that an additional test at the high concentration is being performed.

PRINCIPLES OF THE PROCEDURE

The **BBL™ MGIT™ 7 mL Mycobacteria Growth Indicator Tube** is a tube containing a modified Middlebrook 7H9 Broth which supports the growth and detection of mycobacteria (see **BBL MGIT 7 mL** package insert). The **MGIT** tube contains a fluorescent compound embedded in silicone on the bottom of a 16 x 100 mm round-bottom tube. The fluorescent compound is sensitive to the presence of oxygen dissolved in the broth. The initial concentration of dissolved oxygen quenches the emission from the compound, and little fluorescence can be detected. Later, actively respiring microorganisms consume the oxygen which allows the compound to fluoresce.

The **BACTEC MGIT 960 SIRE Kit** is a 4 – 13 day qualitative test. The test is based on growth of the *Mycobacterium tuberculosis* isolate in a drug-containing tube compared to a drug-free tube (Growth Control). The **BACTEC MGIT** instrument monitors tubes for increased fluorescence. Analysis of fluorescence in the drug-containing tube compared to the fluorescence of the Growth Control tube is used by the instrument to determine susceptibility results.

The **BACTEC MGIT** instrument automatically interprets these results using predefined algorithms (which compare growth in the drug containing tube to that in the Growth Control tube) and the susceptible or resistant result is reported by the instrument.

REAGENTS

BACTEC MGIT 960 SIRE Kit contains one each lyophilized vials of streptomycin, isoniazid, rifampin and ethambutol, and eight vials of SIRE Supplement.

Approximate Formula* Per Vial Lyophilized drug: Streptomycin (STR) 332 µg

Approximate Formula* Per Vial Lyophilized drug: Isoniazid (INH) 33.2 µg

Approximate Formula* Per Vial Lyophilized drug: Rifampin (RIF) 332 µg

Approximate Formula* Per Vial Lyophilized drug: Ethambutol (EMB) 1660 µg

BACTEC MGIT STR 4.0 Kit contains one vial lyophilized streptomycin and two vials of SIRE Supplement.

Approximate Formula* Per Vial Lyophilized drug: Streptomycin 664 µg

BACTEC MGIT INH 0.4 Kit contains one vial lyophilized isoniazid and two vials of SIRE Supplement.

Approximate Formula* Per Vial Lyophilized drug: Isoniazid 66.4 µg

BACTEC MGIT EMB 7.5 Kit contains one vial lyophilized ethambutol and two vials of SIRE Supplement.

Approximate Formula* Per Vial Lyophilized drug: Ethambutol 1245 µg

BACTEC MGIT 960 SIRE Supplement contains 20 mL Middlebrook OADC enrichment

Approximate Formula* Per L Purified Water

Bovine albumin.....50.0 g Catalase.....0.03 g

Dextrose20.0 g Oleic Acid.....0.6 g

*Adjusted and/or supplemented as required to meet performance criteria.

Storage and reconstitution of reagents: **BACTEC MGIT 960 SIRE Drug vials** – on receipt, store the lyophilized drug vials at 2 – 8°C. Once reconstituted, the antibiotic solutions may be frozen and stored at -20°C or colder up to six months, not to exceed the original expiration date. Once thawed, use immediately. Discard unused portions.

BACTEC MGIT SIRE Supplement – On receipt, store in dark at 2 – 8°C. Avoid freezing or overheating. Open and use prior to the expiration date. Minimize exposure to light.

Directions For Use:

Reconstitute each **BACTEC MGIT 960 SIRE Kit Streptomycin** lyophilized drug vial with **4 mL** of sterile distilled/deionized water to make a stock solution of 83 µg/mL.

Reconstitute each **BACTEC MGIT 960 SIRE Kit Isoniazid** lyophilized drug vial with **4 mL** of sterile distilled/deionized water to make a stock solution of 8.3 µg/mL.

Reconstitute each **BACTEC MGIT 960 SIRE** Kit Rifampin lyophilized drug vial with 4 mL of sterile distilled/deionized water to make a stock solution of 83 µg/mL.

Reconstitute each **BACTEC MGIT 960 SIRE** Kit Ethambutol lyophilized drug vial with 4 mL of sterile distilled/deionized water to make a stock solution of 415 µg/mL.

NOTE: The following are reconstituted with a different volume. Failure to use the appropriate volume of sterile distilled water for reconstitution of the higher drug concentrations will invalidate those test results.

Reconstitute each **BACTEC™ MGIT™ 960 STR 4.0** Kit streptomycin lyophilized drug vial with 2 mL of sterile distilled/deionized water to make a stock solution of 332 µg/mL.

Reconstitute each **BACTEC™ MGIT™ 960 INH 0.4** Kit isoniazid lyophilized drug vial with 2 mL of sterile distilled/deionized water to make a stock solution of 33.2 µg/mL.

Reconstitute each **BACTEC™ MGIT™ 960 EMB 7.5** Kit ethambutol lyophilized drug vial with 2 mL of sterile distilled/deionized water to make a stock solution of 622.5 µg/mL.

WARNINGS:

For *in vitro* Diagnostic Use.

POTENTIALLY INFECTIOUS TEST SPECIMEN: Observe “Universal Precautions”⁴ and institutional guidelines when handling and disposing of infectious materials.

Working with *M. tuberculosis* growth in culture requires Biosafety Level (BSL) 3 practices, containment equipment and facilities.

Read and follow directions contained in all appropriate package inserts including the **BBL™ MGIT™ 7 mL Mycobacteria Growth Indicator Tube**.

Prior to use, the user should examine the tubes and vials for evidence of contamination or damage. Discard any tubes or vials if they appear unsuitable or if **MGIT** tubes exhibit fluorescence prior to use.

In the event of tube breakage: 1) Close the instrument drawers; 2) Turn off the instrument; 3) Vacate the area immediately; 4) Consult your facility/CDC guidelines. An inoculated leaking or broken tube may produce an aerosol of mycobacteria; appropriate handling should be observed.

Autoclave all inoculated **MGIT** tubes prior to disposal.

SPECIMEN PREPARATION

All preparations detailed below are from cultures of *Mycobacterium tuberculosis*. The laboratory should confirm, by appropriate identification techniques, that the isolate to be tested is a pure culture.

Preparation of the Isolate from Solid Media:

1. Add 4 mL of **BBL™ Middlebrook 7H9 Broth** (or **BBL MGIT** broth) to a 16.5 x 128 mm sterile tube with cap containing 8 – 10 glass beads.
2. Scrape with a sterile loop as many colonies as possible from growth no more than 14 days old, trying not to remove any solid medium. Suspend the colonies in the Middlebrook 7H9 Broth. The suspension should exceed in turbidity a 1.0 McFarland standard.
3. Vortex the suspension for 2 – 3 min to break up the larger clumps.
4. Let the suspension sit for 20 min without disturbing.
5. Transfer the supernatant fluid to another 16.5 x 128 mm sterile tube with cap (avoid transferring any of the sediment) and let sit for another 15 min.
6. Transfer the supernatant fluid (it should be smooth, free of any clumps) to a third 16.5 x 128 mm sterile tube.
7. Adjust the suspension to a 0.5 McFarland standard by visual comparison to a 0.5 McFarland turbidity standard.
8. Dilute 1 mL of the adjusted suspension in 4 mL of sterile saline (1:5 dilution).

Preparation from a Positive BACTEC MGIT Tube:

1. For the preparation of the test inoculum, a positive 7 mL **MGIT** tube should be used the day after it first becomes positive on the **BACTEC MGIT** instrument (Day 1), up to and including the fifth day (Day 5) after instrument positivity. A tube which has been positive longer than five days should be subcultured to a fresh 7 mL **MGIT** tube containing **BACTEC MGIT** Growth Supplement and tested on the **BACTEC MGIT** instrument until positive, and used from one to five days following positivity.
2. If the tube is a Day 1 or Day 2 positive, proceed to “Inoculation Procedure for Susceptibility Test.”

- If the tube is a Day 3, Day 4, or Day 5 positive, then dilute 1 mL of the positive broth in 4 mL of sterile saline (1:5 dilution). Use the diluted suspension for the inoculation procedures. Proceed to "Inoculation Procedure for Susceptibility Test."

PROCEDURE

Materials Provided: BACTEC MGIT 960 SIRE Kit containing one vial each lyophilized drug and eight vials of SIRE Supplement (approximately 40 tests per drug per kit). BACTEC MGIT 960 STR 4.0 Kit containing one vial lyophilized drug and two vials of SIRE Supplement (approximately 20 tests per kit), BACTEC MGIT 960 INH 0.4 Kit containing one vial lyophilized drug and two vials of SIRE Supplement (approximately 20 tests per kit), and BACTEC MGIT 960 EMB 7.5 Kit containing one vial lyophilized drug and two vials of SIRE Supplement (approximately 20 tests per kit).

Materials Required But Not Provided: BBL MGIT 7 mL Mycobacteria Growth Indicator Tubes, ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

Inoculation Procedure for BACTEC MGIT 960 SIRE Kit Susceptibility Test:

- Label five 7 mL MGIT tubes for each test isolate. Label one as GC (Growth Control), one as STR, one as INH, one as RIF, and the last as EMB. Place the tubes in the correct sequence in the appropriate size AST set carrier (see BACTEC MGIT Instrument User's Manual).
- Aseptically add 0.8 mL of BACTEC MGIT SIRE Supplement to each tube. NOTE: It is important to use the correct supplement.
- Aseptically pipet, using a micropipet, 100 µL of 83 µg/mL MGIT STR solution to the appropriately labeled MGIT tube. Aseptically pipet 100 µL of 8.3 µg/mL MGIT INH solution to the appropriately labeled MGIT tube. Aseptically pipet 100 µL of 83 µg/mL MGIT RIF solution to the appropriately labeled MGIT tube. Aseptically pipet 100 µL of 415 µg/mL MGIT EMB solution to the appropriately labeled MGIT tube. No drug solutions should be added to the MGIT GC tube.

Drug	Concentration of Drug after Reconstitution**	Volume Added to MGIT Tubes for Test	Final Concentration in MGIT Tubes
MGIT STR	83 µg/mL	100 µL	1.0* µg/mL
MGIT INH	8.3 µg/mL	100 µL	0.1* µg/mL
MGIT RIF	83 µg/mL	100 µL	1.0* µg/mL
MGIT EMB	415 µg/mL	100 µL	5.0* µg/mL

*Equivalent to CDC⁴ recommended critical drug concentrations.

** These drugs must be reconstituted using 4 mL sterile/deionized water to achieve concentrations indicated.

- Growth Control tube preparation and inoculation:** Aseptically pipet 0.1 mL of the organism suspension (see "Specimen Preparation") into 10 mL of sterile saline to prepare the 1:100 Growth Control suspension. Mix the Growth Control suspension thoroughly. Inoculate 0.5 mL of the 1:100 Growth Control suspension into the MGIT tube labeled "GC."
- Drug-containing tube inoculation:** Aseptically pipet 0.5 mL of the organism suspension (see "Specimen Preparation") into each of the FOUR remaining drug tubes (STR, INH, RIF, EMB).
- Tightly recap the tubes and mix well.
- Enter the AST set into the BACTEC MGIT 960 SIRE Kit using the AST set entry feature (refer to the BACTEC MGIT Instrument User's Manual). Ensure that the order of the tubes in the AST Set Carrier conforms to the Set Carrier definitions selected when performing the AST set entry feature.
- Streak 0.1 mL of the organism suspension to a **Trypticase™** Soy Agar with 5% Sheep Blood (TSA II) plate. Enclose in a plastic bag. Incubate at 35 – 37°C.
- Check the blood agar plate at 48 h for bacterial contamination. If the blood agar plate shows no growth, then allow AST testing to proceed. If the blood agar plate shows growth, discard the AST set (refer to the BACTEC MGIT Instrument User's Manual) and repeat testing with pure culture.

Inoculation Procedure for BACTEC MGIT 960 STR 4.0, INH 0.4 and EMB 7.5 Kits Susceptibility Test:

If resistance occurs at the critical concentration, an optional profile test may be performed which, at minimum, tests the high concentration of the drug to which it was originally resistant:

Isolate Source: The isolate used for this testing must have been prepared as described in "Specimen Preparation." A seed tube may be prepared from the drug free Growth Control tube

from the previously tested AST set of the isolate, from which 0.5 mL may be inoculated to a fresh **MGIT** 7 mL tube containing **BACTEC MGIT** Growth Supplement. Once the seed tube is instrument positive, proceed as described in “Specimen Preparation: Preparation from a Positive **MGIT** tube.”

1. Label enough **MGIT** 7 mL tubes for the test isolate to have a **MGIT** GC (Growth Control) and a **MGIT** drug tube for each antimicrobial tested.
2. Aseptically add 0.8 mL of **BACTEC MGIT** SIRE Supplement to each tube. NOTE: it is important to use the correct supplement.
3. Aseptically pipet, using a micropipet, 100 µL of the drug solution to the appropriately labeled **MGIT** tube. No antibiotics should be added to the **MGIT** GC tube.

Drug	Concentration of Drug after Reconstitution**	Volume Added to MGIT Tubes for Test	Final Concentration in MGIT Tubes
MGIT STR 4.0	332 µg/mL	100 µL	4.0* µg/mL
MGIT INH 0.4	33.2 µg/mL	100 µL	0.4* µg/mL
MGIT EMB 7.5	622.5 µg/mL	100 µL	7.5* µg/mL

* Equivalent to MOP high drug concentrations.⁴

** These drugs must be reconstituted using 2 mL sterile/deionized water to achieve concentrations indicated.

4. **Growth Control tube preparation and inoculation:** Aseptically pipet 0.1 mL of the organism suspension (see “Specimen Preparation”) into 10 mL of sterile saline to prepare the 1:100 Growth Control suspension. Mix the Growth Control suspension thoroughly. Inoculate 0.5 mL of the 1:100 Growth Control suspension into the **MGIT** tube labeled “GC.”
5. **Drug-containing tube inoculation:** Aseptically pipet 0.5 mL of the organism suspension (see “Specimen Preparation”) into each of the drug tubes.
6. Tightly recap the tubes and mix well.
7. Enter the AST set into the **BACTEC MGIT** instrument using the AST set entry feature (refer to the **BACTEC MGIT** Instrument User’s Manual). Ensure that the order of the tubes in the AST Set Carrier conforms to the Set Carrier definitions selected when performing the AST set entry feature.
8. Streak 0.1 mL of the organism suspension to a **Trypticase** Soy Agar with 5% Sheep Blood (TSA II) plate. Enclose in a plastic bag. Incubate at 35 – 37°C.
9. Check the blood agar plate at 48 h for bacterial contamination. If the blood agar plate shows no growth, then allow AST testing to proceed. If the blood agar plate shows growth, discard the AST set (refer to the **BACTEC MGIT** Instrument User’s Manual) and repeat testing with pure culture.

NOTE: The susceptibility test may be configured in a variety of formats. For example, a five tube carrier set containing only the critical concentrations can be configured into the system. A variety of other tube carrier sets can be configured depending on the optional profile tests being run (refer to **BACTEC MGIT** Instrument User’s manual).

User Quality Control: Upon receipt of a new shipment or lot number of **BACTEC MGIT** 960 SIRE Kit vials, it is suggested that the control organism shown below be tested (see “Inoculation Procedure for Susceptibility Test”). Observation of the proper results, as shown below, within 4 – 13 days indicates that the **BACTEC MGIT** 960 SIRE Kit is ready for use in testing patient isolates. If the proper results are not observed, repeat the test. If, after repeating the test, the proper results are still not observed, do not use the product until you have contacted Technical Services at (800) 638-8663 (United States only).

The same control organism should be run as batch QC once each week when susceptibility testing is performed. If the batch QC fails, do not report patient results for the drug (s) that failed for that testing period. Repeat the QC for the drug(s) and patient isolates affected by the initial QC failure. If the repeat QC does not perform as expected, do not report patient results. Do not use the product until you have contacted Technical Services at (800) 638-8663 (United States only).

Strains	GC	MGIT STR	MGIT INH	MGIT RIF	MGIT EMB
<i>M. tuberculosis</i> ATCC™ 27294	Positive	Susceptible	Susceptible	Susceptible	Susceptible

Strains	GC	MGIT STR 4.0	MGIT INH 0.4	MGIT EMB 7.5
<i>M. tuberculosis</i> ATCC™ 27294	Positive	Susceptible	Susceptible	Susceptible

REPORTING OF RESULTS

The **BACTEC MGIT** instrument will monitor susceptibility test sets until a susceptible or resistant determination is made. Once the set testing is completed, the results are reported by the **BACTEC MGIT** instrument (refer to the **BACTEC MGIT** Instrument User's Manual).

LIMITATIONS OF THE PROCEDURE

Suspensions made from solid media must be allowed to settle for the prescribed times prior to standardization. Inoculum preparations made from solid media without the use of a 0.5 McFarland turbidity standard and the appropriate dilutions may give inaccurate results.

Use only pure cultures of *M. tuberculosis*. Cultures which are contaminated or which contain multiple strains of mycobacteria may give erroneous results.

Failure to reconstitute the drugs with the appropriate volume of sterile deionized water may give inaccurate results.

Failure to use a 1:100 dilution of the isolate for the inoculation of the Growth Control tube may give inaccurate results.

Failure to load the tubes of the AST set into the AST Set Carrier in the proper sequence may give inaccurate results.

Failure to use the SIRE Supplement in the AST set may give inaccurate results. Do not add **BACTEC MGIT** Growth Supplement to the AST set.

EXPECTED VALUES

See addendum

PERFORMANCE CHARACTERISTICS

See addendum

AVAILABILITY

Cat. No.	Description
245123	BACTEC™ MGIT™ 960 SIRE Kit, carton of 4 lyophilized drug vials, and 8 SIRE Supplements.
245125	BACTEC™ MGIT™ 960 STR 4.0 Kit, carton of 1 lyophilized drug vial and 2 SIRE Supplements.
245126	BACTEC™ MGIT™ 960 INH 0.4 Kit, carton of 1 lyophilized drug vial and 2 SIRE Supplements.
245127	BACTEC™ MGIT™ 960 EMB 7.5 Kit, carton of 1 lyophilized drug vial and 2 SIRE Supplements.

REFERENCES

1. Barenfanger, J. 1993. Making your lab safe against multi-drug resistant *Mycobacterium tuberculosis*. Clin. Microbiol. Newsl. 15: 76-80.
2. Clinical and Laboratory Standards Institute. 2003. Approved Standard: M24-A. Susceptibility testing of mycobacteria, nocardiae, and other aerobic actinomycetes. CLSI, Wayne, Pa.
3. BD Diagnostic Systems. **BACTEC™** 460TB System Product and Procedure Manual.
4. Kent, P.T., and G.P. Kubica. 1985. Public health mycobacteriology: a guide for the level III laboratory. USDHHS. Centers for Disease Control, Atlanta.
5. Data on file at BD Diagnostic Systems.

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ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo, BACTEC, BBL, MGIT and Trypticase are trademarks of Becton, Dickinson and Company. © 2010 BD.
Benex Limited