Revisions

SO 0191-5

Rev from	Rev to	ECO #
2010/01	2011/07	5864-11

Notes:

1. BD Cat. Number 221602, 221606, 221607, 221608, 221609, 221625
2. Blank (Sheet) Size : Length: 9" Width: 28.875"
 Number of Pages: 14 Number of Sheets: 1
 Page Size: Length 9" Width 4.125" Final Folded Size: 2.25" x 4.125"
3. Style (see illustrations below): # 4

4. See Specification Control Number 8830161 for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# 347 Green
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 8830161JAA		Category and Description Package Insert, Port-A-Cul Specimen Collection and Transport Product	Sheet: 1 of 15 Scale: 1:1	A

BD BBL™ Port-A-Cul™ Specimen Collection and Transport Products

English:	pages	1 – 3	Italiano:	pagine	7 – 9	8830161JAA
Français :	pages	3 – 5	Español:	páginas	9 – 11	2011/07
Deutsch:	Seiten	5 – 7				

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзініздің жергілікті БД өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

BBL™ Port-A-Cul™ tubes, jars and vials contain a reduced transport medium and are intended to maintain the viability of anaerobic, facultative and aerobic microorganisms during transport from the patient to the laboratory. Sterile packages are for collection of specimens in clean areas; e.g., surgical suites.

SUMMARY AND EXPLANATION

Anaerobic microorganisms require special transport methods to prevent even brief exposure to oxygen.^{1,2} Maintaining aerobic, anaerobic and facultative organisms in an anaerobic environment provides an ideal method of transport.³⁻⁶

BBL Port-A-Cul Tubes are for swab specimens. Swab specimens are inserted into a reduced solidified holding medium, in contrast to other methods in which swabs are placed into fluid media where organisms may be diluted or into an empty tube where organisms may dry out. **Port-A-Cul** Transport Jars are for tissue and biopsy specimens. The wide-mouth jar allows for easier insertion of the specimen into the reduced solidified holding medium. **Port-A-Cul** Vials are for fluid specimens. Fluid specimens are injected directly onto the solid agar surface. The **Port-A-Cul** Fluid Collection Set provides a sterile-packed syringe and needle.

PRINCIPLES OF THE PROCEDURE

Swab and tissue specimens are protected from air and remain moist in a column of solid medium under an anaerobic atmosphere. Agar inhibits oxygen diffusion. Reducing agents combine with free oxygen. Salts and buffers provide a protective environment. Resazurin indicates the presence of oxygen (pink to blue).

Fluid specimens are protected from air in a **Port-A-Cul** vial because of the anaerobic atmosphere and exchange of reducing agents from the solid medium to the specimen.

REAGENTS

Port-A-Cul medium - balanced formula of reducing agents and resazurin in a buffered isotonic agar base.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

This product contains dry natural rubber.

 **BBL Port-A-Cul** is for single use only; reuse may cause a risk of infection and/or inaccurate results.

Exercise care in opening container with tight cap.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"⁷⁻¹⁰ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids. Autoclave contaminated materials after use and discard in appropriate biohazard waste container.

Storage Instructions: On receipt, store at 20 – 25°C. Media stored as labeled until just prior to use may be inoculated up to the expiration date and transported for the recommended transportation times.

Product Deterioration: Do not use if medium shows evidence of contamination, drying, cracking or other signs of deterioration (flocculation and precipitation are not indications of deteriorated medium) or if integrity of Sterile Pack packaging has been compromised. Do not use if a pink or blue band at the top of the medium is deeper than 3 mm, as this is indicative of excessive oxidation.

SPECIMEN COLLECTION AND TRANSPORT

Follow accepted specimen collection practices.¹¹ Fluid specimens are preferred. Do not culture body sites normally contaminated with indigenous flora. Transport to the laboratory at ambient temperatures (20 – 25°C) within 72 h.

PROCEDURE

Material Provided: BBL Port-A-Cul products (see "Availability").

Materials Required But Not Provided: Syringes and needles for Cat. Nos. 221608/221609; swabs for Cat. No. 221606; sterile forceps and materials required for organism isolation, identification and other microbiological procedures.

Test Procedure:

SWAB SPECIMENS – Obtain specimen, loosen screw cap, insert swab into medium within approximately 5 mm from bottom, break shaft evenly with tube lip and quickly replace and tighten cap. (H₂S may be smelled upon opening tube; this does not adversely affect medium.) Transport to laboratory promptly (within 72 h).

Remove swab by grasping applicator shaft with sterile forceps; streak onto an appropriate growth medium.

FLUID SPECIMENS – Remove green flip cap exposing rubber stopper. Swab rubber stopper with disinfectant. Expel air from syringe and needle. Obtain specimen, push needle through stopper and **slowly** inject on surface of agar. Transport to laboratory promptly (within 72 h).

Prior to aspirating specimen for culture, swab stopper with disinfectant.

TISSUE SPECIMENS – Obtain specimen, loosen screw cap, insert specimen into medium within approximately 5 mm from bottom and quickly replace and tighten cap. Transport to laboratory promptly (within 72 h).

To remove specimen, loosen cap and remove specimen with sterile forceps.

STERILE PACK – Peel envelope apart outside clean area without touching contents of envelope. Wearing sterile gloves, remove sterile tube and swabs, jar, vial or vial and syringe. Take into clean area. See directions for use on envelope for swab, fluid, or tissue transport.

Quality Control:

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI (formerly NCCLS) guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Little or no growth of microorganisms should occur within the specimen while it is being transported in the **Port-A-Cul** medium. However, organisms in logarithmic growth phase or in the nutritious environment of the specimen itself may continue to grow until endogenous or exogenous nutrients are depleted.

LIMITATIONS OF THE PROCEDURE

For best results, specimens should be transported to the laboratory as quickly as possible, but no longer than 72 h. Overgrowth may occur with polymicrobial infections. Avoid temperature extremes during transport. Organisms in small numbers (< 100 CFU/mL) may not survive longer than 24 h.

Viability of *Neisseria gonorrhoeae* decreases rapidly at temperatures below 35°C and may not survive more than 8 h even in high concentration.

Transport media, staining reagents, immersion oil, glass slides and specimens themselves sometimes contain dead organisms visible upon Gram staining.

PERFORMANCE CHARACTERISTICS

The survival of *Streptococcus pyogenes* and *Haemophilus influenzae* was found to be better in **Port-A-Cul** than in an aerobic transport device.⁶ **Port-A-Cul** was evaluated with stock cultures of a variety of clinically significant anaerobes and found to maintain viability for at least 72 h at 25°C.⁴ In another study, using simulated specimens stored up to 48 h at ambient temperature, recovery of ten specimens of anaerobic bacteria was greater in **Port-A-Cul** than in other controlled conditions.⁵

AVAILABILITY

Cat. No.	Description
221606	BBL™ Port-A-Cul™ Tube - pkg. of 10 size D tubes (11 mL/tube). ☹☹
221607	BBL™ Port-A-Cul™ Tube and Swabs Sterile Pack - pkg. of 10 envelopes, each containing one size D tube (11 mL/tube) and two swabs. ☹☹
221608	BBL™ Port-A-Cul™ Vial - pkg. of 10 vials (5 mL/vial). ☹☹
221609	BBL™ Port-A-Cul™ Vial Sterile Pack - pkg. of 10 envelopes, each containing one vial (5 mL/vial). ☹☹
221602	BBL™ Port-A-Cul™ Transport Jar Sterile Pack - pkg. of 10 envelopes, each containing one transport jar (20 mL/jar). ☹☹
221625	BBL™ Port-A-Cul™ Fluid Collection Set (Sterile) - pkg. of 10 envelopes, each containing one vial (5 mL/vial) and one BD 5 cc Luer-Lok™ syringe with 22 G needle.

REFERENCES

1. Loesche, W.J. 1969. Oxygen sensitivity of various anaerobic bacteria. *Appl. Microbiol.* 18:723-727.
2. Chow, A.W., R.D. Leake, T. Yamauchi, B.F. Anthony, and L.B. Guze. 1974. The significance of anaerobes in neonatal bacteremia: analysis of 23 cases and review of the literature. *Pediatrics* 54:736-745.
3. Finegold, S.M., V.L. Sutter, H.R. Attebery, and J.E. Rosenblatt. 1974. Isolation of anaerobic bacteria, p. 365-375. In E.H. Lennette, E.H. Spaulding, and J.R. Truant (ed.), *Manual of clinical microbiology*, 2nd ed. American Society for Microbiology, Washington, D.C.
4. Abramson, I.J., G.L. Evans, and P. Bathurst. 1977. Laboratory evaluation of Port-A-Cul™ with stock cultures of anaerobes, aerobes, and facultative organisms, abstr. C 197, p. 68. *Abstr. 77th Annu. Meet. Am. Soc. Microbiol.* 1977.
5. Mena. E., F.S. Thompson, A.Y. Armfield, V.R. Dowell, Jr., and D.J. Reinhardt. 1978. Evaluation of Port-A-Cul transport system for protection of anaerobic bacteria. *J. Clin. Microbiol.* 8:28-35.
6. Hamilton, R.L., G.L. Evans, and T.E. O'Neill. 1979. Survival of facultatively anaerobic bacteria in Port-A-Cul™ tubes and an aerobic transport device, abstr. C 173, p. 338. *Abstr. 79th Annu. Meet. Am. Soc. Microbiol.* 1979.
7. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
8. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
9. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
10. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262, 17/10/2000, p. 0021-0045.*
11. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenover (ed.). 2003. *Manual of clinical microbiology*, 8th ed. American Society for Microbiology, Washington, D.C.

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Benex Limited