BD BBL™ Sensi-Disc™ Antimicrobial Susceptibility Test Discs

8840621
2004/09

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See symbol glossary at end of insert.

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD representant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar.

Refer to product label for CE mark availability. / Informace o dostupnosti značky CE najdete na etiketě na výrobku / Produktet er tilgjengelig, hvis CE-mærkingen finnes på etiketten / Toote saadavuse teadaasaamiseks kontrollige CE-märgi olemasolu etiketil / Pour vérifier le marquage CE du produit, se référer à l'étiquette / Siehe Produktetikett für die Verfügbarkeit CE markierter Produkte / Το προϊόν είναι διαθέσιμο εάν η ετικέτα φέρει τη σήμανση CE / Lásd a termék címkéjén levő CE jel meglétét / CE žymės ieškote ant produkto etiketės / O dostopnosti produktu informuje znak CE na etykietce / Prodotto disponibile se l'etichetta reca il marchio CE / Verificar a disponibilidade do produto pela marca CE no rótulo / Pre kontrolu prítomnosti CE značky nahľadnite do etikety výrobku / Ver disponibilidade del marcado CE en la etiqueta del producto / Produkten finns tillgänglig om dess etikett är CE-märkt

INTENDED USE – These discs are used for semi-quantitative *in vitro* susceptibility testing by the agar disc diffusion test procedure of common, rapidly growing and certain fastidious bacterial pathogens. These include the *Enterobacteriaceae*, *Staphylococcus* spp., *Pseudomonas* spp., *Acinetobacter* spp., *Enterococcus* spp., *Vibrio cholerae* and, by modified procedures, *Haemophilus influenzae*, *Neisseria gonorrhoeae*, *Streptococcus pneumoniae* and other streptococci. NOTE: Special procedures are required for testing pneumococci, enterococci and methicillin/oxacillin-resistant staphylococci, for performing β -lactamase tests and for screening and confirmatory tests for ESBLs; see the "RESULTS" section.

For zone diameter interpretive criteria adopted in France, refer to the instructions in the French language section of this insert.

SUMMARY AND EXPLANATION – Agar diffusion methods employing dried filter paper discs impregnated with specific concentrations of antimicrobial agents were developed in the 1940s. In order to eliminate or minimize variability in this testing, Bauer et al. developed a standardized procedure in which Mueller Hinton Agar was selected as the test medium.^{1,2}

Various regulatory agencies and standards-writing organizations subsequently published standardized reference procedures based on the Bauer-Kirby method. Among the earliest and most widely accepted of these standardized procedures were those published by the U.S. Food and Drug Administration (FDA)³ and the World Health Organization (WHO).^{4,5} The procedure was adopted as a consensus standard by the National Committee for Clinical Laboratory Standards (NCCLS) and is periodically updated.^{6,7} The latest NCCLS documents should be consulted for current recommendations.

PRINCIPLES OF THE PROCEDURE – Discs containing a wide variety of antimicrobial agents are applied to the surface of Mueller Hinton Agar plates (or Haemophilus Test Medium Agar for *H. influenzae*, GC II Agar with IsoVital[™] Enrichment for *N. gonorrhoeae* or Mueller Hinton Agar with 5% Sheep Blood for *S. pneumoniae*, β -hemolytic and viridans group streptococci) that have been inoculated with pure cultures of clinical isolates. Following incubation, the plates are examined and the zones of inhibition surrounding the discs are measured and compared with established zone size ranges for individual antimicrobial agents in order to determine the agent(s) most suitable for use in antimicrobial therapy.

REAGENTS – Sensi-Disc[™] brand discs are 6-mm discs prepared by impregnating high quality absorbent paper with accurately determined amounts of antibiotic or other chemotherapeutic agents. Discs are clearly marked on both sides with letters and numbers designating the agent and the drug content. (See chart giving concentrations of reactive ingredients.) The drug content of discs is assayed by the methods established by the FDA or by methods similar or comparable to those published in the United States Federal Register.

Sensi-Disc agents are furnished in cartridges containing 50 discs each. The last disc in each cartridge is marked "X" and contains the drug as coded. Cartridges are for use in BBL[™] Sensi-Disc[™] Dispensers; these include a Single Disc Dispenser, an 8-Place Dispenser for 100 mm-style Petri dishes, 6- and 8-Place Self-Tamping Dispensers for 100 mm-style dishes and a Self-Tamping 12-Place Dispenser for 150 mm-style plates.

Precautions: For *in vitro* Diagnostic Use.

Follow directions for use; disc performance depends not only on disc potency, but on use of proper inoculum and control cultures, functional pretested plates, proper storage temperature and other factors.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. Sterilize cultures, containers and other contaminated materials after use.

Storage Instructions:

- On receipt, store discs at -20 – +8°C. If the laboratory refrigerator is frequently opened and closed, and a suitable temperature is not maintained, place there a supply sufficient only for use within a week. Some discs (e.g., β -lactams) should preferably be kept frozen at -20°C.
- Allow containers to come to room temperature before opening. Return unused discs to the refrigerator when application of discs has been completed.
- Use the oldest discs first.
- Discard expired discs. Also, cartridges from which discs have been frequently removed during a week and discs left out overnight in the laboratory should be discarded, or else the discs should be tested for acceptable performance prior to continued use.
- If the discs form incorrect zones with the recommended control organisms, the entire procedure should be checked; faulty zone size may be due to the disc, the inoculation, the preparation or depth (about 4 mm) of medium, or other factors.

The expiration date applies only to discs in intact containers, stored as directed.

SPECIMENS – Specimens should not ordinarily be employed in this test. See Directions, which include preparation of inoculum. If possible, cultures should be derived from specimens obtained from patients prior to the initiation of antimicrobial therapy.

PROCEDURE

Material Provided: Sensi-Disc[™] susceptibility test discs as labeled.

Materials Required But Not Provided: Ancillary culture media, reagents, quality control organisms and laboratory equipment required to perform disc diffusion susceptibility testing by the standardized procedure. Prepare a 0.5 McFarland turbidity standard by adding 0.5 mL of 0.048 M BaCl₂ [1.175% (wt/vol) BaCl₂·2H₂O] to 99.5 mL of 0.18 M [0.36N] H₂SO₄ [1% (vol/vol)]. Verify by using a spectrophotometer with a 1-cm light path and matched cuvette; absorbance at 625 nm should be 0.08 – 0.10.

Directions, Including User Controls:⁵

- Preparation of inoculum with test and control cultures.
 - Perform a Gram stain. Use only pure cultures.
 - Select three to five similar colonies and transfer with inoculation needle or loop into 4 – 5 mL of a suitable broth such as **Trypticase[™]** Soy Broth (or Mueller Hinton Broth for fastidious organisms).
 - Incubate the broth cultures at 35°C for 2 – 6 h, if necessary, to develop a turbidity equivalent to the 0.5 McFarland turbidity standard (approximately 1 to 2 x 10⁸ CFU/mL). Alternatively, make a direct broth or saline suspension of colonies selected from an agar plate incubated overnight (a nonselective medium such as blood agar, or chocolate agar for *H. influenzae* and *N. gonorrhoeae*, should be used). The direct colony suspension method is preferred for *Staphylococcus* spp., *S. pneumoniae* and other streptococci, *Haemophilus* spp. and *N. gonorrhoeae*.⁶
 - Dilute, if required, to obtain turbidity equivalent to the 0.5 McFarland turbidity standard. For diluent, use sterile broth or saline. Alternatively, standardize the inoculum photometrically; to facilitate inoculum adjustment of rapidly growing organisms, the **Prompt[™]** Inoculation System (volumetric inoculum preparation device) may be used.⁸
- Overnight broth cultures should not be used as inoculum.
- Inoculation.
 - Within 15 min, dip a sterile cotton swab into the properly adjusted inoculum and rotate it firmly several times against the upper inside wall of the tube to express excess fluid.
 - Streak the entire agar surface of a Mueller Hinton Agar (or other appropriate agar) plate three times, turning the plate 60° between streakings to obtain even inoculation.
 - The lid may be left ajar for 3 – 5 min, but no more than 15 min, to allow for any surface moisture to be absorbed before applying the drug-impregnated discs.
- Select appropriate discs (such as recommended in reference 7, Tables 1 and 1A of M100-513 [M2]).
- Apply the discs by means of a BBL[™] dispenser, using aseptic precautions. Deposit discs so that the centers are at least 24 mm apart. It is preferable to deposit penicillin and cephalosporin discs so that they are no less than 10 mm from the edge of the Petri dish, and their centers are at least 30 mm apart. Avoid placing such discs adjacent to one another. With *H. influenzae*, *N. gonorrhoeae* and *S. pneumoniae*, use no more than nine discs per 150 mm plate or four discs per 100 mm plate. If discs have been placed on the agar with other than the Self-Tamping Dispensers, press them down with a sterile needle or forceps to make contact with the surface.
- Within 15 min, place the plates agar side up in a 35°C incubator. *Haemophilus* spp., *N. gonorrhoeae*, *S. pneumoniae* and other streptococci should be incubated in an atmosphere enriched with 5% CO₂.
- Examine the plates after 16 – 18 h of incubation (20 – 24 h for *N. gonorrhoeae*, *S. pneumoniae* and other streptococci). A full 24 h of incubation is recommended for *Staphylococcus* spp. to detect methicillin/nafcillin/oxacillin-resistant staphylococci and *Enterococcus* spp. for vancomycin resistance. The diameters of the zones of complete inhibition are measured, as determined by gross visual inspection. Zones are measured to the nearest whole millimeter. For further details in measuring zones of inhibition, consult the reference.⁶ If only isolated colonies grow, the inoculum is too light and the test should be repeated. Zones around discs containing different drugs are not comparable for the purpose of comparing activity of drugs. See the Zone Diameter Interpretive Chart, which gives expected values from testing common aerobes. Zone measurement may be simplified by using a BBL[™] Sensi-Disc[™] Zone Interpretation Set.
- Control tests using prescribed cultures should be included each day susceptibility testing is performed or weekly if satisfactory performance can be documented according to the NCCLS standard.⁶ Typical zone sizes of *E. coli* ATCC[™] 25922, *S. aureus* ATCC 25923, *P. aeruginosa* ATCC 27853, *H. influenzae* ATCC 49247, *H. influenzae* ATCC 49766, *N. gonorrhoeae* ATCC 49226, *S. pneumoniae* ATCC 49619, *E. coli* ATCC 35218 (β -lactamase-producing strain), *E. faecalis* ATCC 29212 (for quality control testing of gentamicin 120 μ g and streptomycin 300 μ g discs) and *Klebsiella pneumoniae* ATCC 700603 (for screening and confirmatory tests for ESBLs) are given in the chart (or footnotes) and indicate the correct performance of the entire procedure. *E. faecalis* ATCC 29212 (or 33186) is also recommended for evaluating new lots of Mueller Hinton Agar for low thymine and thymidine content (see footnote t). *H. influenzae* ATCC 10211 is recommended as a useful additional quality control strain to verify the growth promotion properties of Haemophilus Test Medium Agar.⁷

RESULTS^{6,7} – NOTE: Recommended interpretive criteria are based on usual dosage regimens and routes of administration in the U.S.

Compare recorded zone diameters with those in the chart; results with a specific organism may be reported as Resistant, Intermediate or Susceptible. For some organism/antimicrobial combinations, the absence of resistant strains precludes defining any results categories other than "Susceptible." For strains yielding results suggestive of a "nonsusceptible" category, organism identification and antimicrobial susceptibility test results should be confirmed. Subsequently, the isolates should be saved and submitted to a reference laboratory that will confirm results using a NCCLS reference dilution method.

For organisms excluded from the Zone Diameter Interpretive Chart (e.g., *Campylobacter*, *Corynebacterium* spp., *Bacillus* spp.), studies are not yet adequate to develop reproducible definitive standards for interpretation of results. If necessary, a dilution method usually will be the most appropriate testing method, which may require submitting the organism to a reference laboratory.⁷

For fecal isolates of *Salmonella* and *Shigella* spp., only ampicillin, a quinolone and trimethoprim/sulfamethoxazole should be tested and reported routinely. In addition, chloramphenicol and a third generation cephalosporin should be tested and reported for extraintestinal isolates of *Salmonella* spp. For *Salmonella* and *Shigella* spp., aminoglycosides and first and second generation cephalosporins may appear active *in vitro* but are not effective clinically and should not be reported as susceptible.⁷

Enterobacter, *Citrobacter*, and *Serratia* may develop resistance during prolonged therapy with third generation cephalosporins. Therefore, isolates that are initially susceptible may become resistant within 3 to 4 days after initiation of therapy. Testing of repeat isolates may be warranted.

Non-*Enterobacteriaceae* other than *P. aeruginosa* and *Acinetobacter* spp. should be tested by the dilution method (see M7-A6⁹).

P. aeruginosa may develop resistance during prolonged therapy with all antibiotics. Isolates that are initially susceptible may become resistant within 3 to 4 days after initiation of therapy and testing of repeat isolates may be warranted.

The susceptibility of *Pseudomonas aeruginosa* isolated from patients with cystic fibrosis can be reliably determined by the disc method, but may require extended incubation up to 24 h before reporting as susceptible.

Enterococci may be resistant to penicillin and ampicillin because of the production of low-affinity, penicillin-binding proteins (PBPs), or the production of β -lactamase. The disc diffusion test can accurately detect isolates with altered PBPs, but it will not reliably detect β -lactamase producing strains. The latter strains are best detected by using a direct β -lactamase test;⁶ e.g., with **Cefinase[™]** nitrocefin discs or chromogenic cephalosporin discs.

For *Enterococcus* spp., cephalosporins, aminoglycosides (except for high level resistance screening), clindamycin and trimethoprim/sulfamethoxazole may appear active *in vitro* but are not effective clinically and isolates should not be reported as susceptible.

Only results of testing with ampicillin, one of the third-generation cephalosporins, chloramphenicol and meropenem should be reported routinely with cerebrospinal fluid isolates of *H. influenzae*.

Amoxicillin/clavulanic acid, azithromycin, clarithromycin, ceftazidime, cefeprozil, loracarbef, cefdinir, cefixime, cefpodoxime and cefuroxime axetil are oral agents that may be used as empiric therapy for respiratory tract infections due to *Haemophilus* spp. The results of susceptibility tests with these antimicrobial agents are often not useful for management of individual patients. However, susceptibility testing of *Haemophilus* spp. with these compounds may be appropriate for surveillance or epidemiologic studies.

Extended-spectrum β -lactamases (ESBLs) are enzymes produced by gram-negative bacilli that arise by mutation in genes for common plasmid-mediated β -lactamases. Strains of *Klebsiella* spp. and *E. coli* that produce ESBLs may be clinically resistant to therapy with penicillins, cephalosporins, or aztreonam, despite apparent *in vitro* susceptibility to some of these agents. Some of these strains will show zones of inhibition below the normal susceptible population but above the standard breakpoints for certain extended-spectrum cephalosporins or aztreonam; such strains should be screened for potential ESBL production by using the ESBL screening breakpoints before reporting results for penicillins, extended-spectrum cephalosporins or aztreonam. Other strains may test intermediate or resistant by standard breakpoints to one or more of these agents. In all strains with ESBLs the zone diameters for one or more of the extended-spectrum cephalosporins or aztreonam should increase in the presence of clavulanic acid as determined in phenotypic confirmatory testing. For all confirmed ESBL-producing strains, the test interpretation should be reported as resistant for all penicillins, cephalosporins, and aztreonam. See footnote t for ESBL screening and confirmatory tests. The decision to perform ESBL screening tests on all urine isolates should be made on an institutional basis, considering prevalence, therapy and infection control issues.⁷

For recognition of methicillin-resistant staphylococci, the oxacillin disc test is more likely to detect resistance than the use of methicillin or nafcillin discs. Therefore, the 1 μ g oxacillin disc should be used to test for methicillin/oxacillin resistance. Any zone surrounding the oxacillin disc should be inspected carefully using transmitted light for small colonies or a light "film" of growth within the zone of inhibition after a full 24 h incubation. Methicillin-resistant staphylococci are often resistant to multiple classes of antimicrobial agents including aminoglycosides, macrolides, clindamycin, phenolics, quinolones, sulfonamides and tetracycline. The observation of multiple resistance should be a clue to the possibility of methicillin-resistance. However, strains of methicillin-resistant *S. aureus* that do not exhibit resistance to other classes of antimicrobial agents have been isolated from both inpatient and outpatient populations. If the disc diffusion test result is in doubt with a possible methicillin-resistant *Staphylococcus* spp., perform additional confirmatory tests as outlined in NCCLS document M7.⁹ Methicillin/oxacillin-resistant *S. aureus* (MRSA) and coagulase negative staphylococci (MRS) should be reported as resistant (or not reported at all) to all penicillins, cepheps, carbapenems and other β -lactams, such as amoxicillin/clavulanic acid, ampicillin/sulbactam, ticarcillin/clavulanic acid, piperacillin/tazobactam and imipenem, regardless of the *in vitro* test results with those agents. This is because most cases of documented infections due to methicillin-resistant staphylococci have responded poorly to β -lactam therapy and convincing clinical data have yet to be presented that document clinical efficacy for those agents.⁶ Isolates of staphylococci that are shown to carry the *mecA* gene, or the produce PBP 2a, the *mecA* gene product, should be reported as oxacillin resistant.

Interpretive criteria for coagulase-negative staphylococci correlate with the presence or absence of the gene encoding methicillin resistance (*mecA*) for *S. epidermidis*. These interpretive criteria may overcall resistance for other coagulase-negative staphylococci, e.g., *S. lugdunensis* or *S. saprophyticus*. For serious infections with coagulase-negative staphylococci other than *S. epidermidis*, testing for *mecA* or the protein expressed by *mecA*, the penicillin binding protein 2a (PBP 2a, "also known as" PBP 2') may be appropriate for strains having zone diameters in the intermediate or resistant range. Isolates that are not shown to carry *mecA* or do not produce PBP 2a should be reported as oxacillin-susceptible.

It has been reported that disc susceptibility testing is not an accurate method for the determination of methicillin (oxacillin) susceptibility for coagulase-negative staphylococci (i.e., *S. saprophyticus*).¹⁰ Routine testing of urine isolates of *S. saprophyticus* is not advised, because infections respond to concentrations achieved in urine of antimicrobial agents commonly used to treat acute, uncomplicated urinary tract infections (e.g., nitrofurantoin, trimethoprim/sulfamethoxazole, or a fluoroquinolone).

SIZE: 5.375" W x 19" L
COLOR: Standard Black

A rapid β -lactamase test (e.g., using **Cefinase** discs) may yield clinically relevant information earlier than results of a disc diffusion test with *Haemophilus* spp., *N. gonorrhoeae* and *Moraxella catarrhalis*; it is the only reliable test for detecting β -lactamase-producing *Enterococcus* spp. A positive β -lactamase test predicts resistance to penicillin, ampicillin and amoxicillin among *Haemophilus* spp., *N. gonorrhoeae* and *M. catarrhalis* and resistance to penicillin, including acylamino-, carboxy- and ureido-penicillins among staphylococci and enterococci. A negative β -lactamase test does not rule out resistance due to other mechanisms. Do not test members of the *Enterobacteriaceae*, *Pseudomonas* spp. and other aerobic gram-negative bacilli because the results may not be predictive of susceptibility to the β -lactams most often used for therapy. Accurate detection of β -lactamase in staphylococci may require induction of the enzyme and incubation of a nitrocefin-based test for up to 1 h. Induction can be easily accomplished by testing the growth from the zone margin surrounding an oxacillin disc test. Care must be exercised to ensure accurate results, including testing of known positive and negative control strains at the time clinical isolates are examined.⁶

Susceptibility testing of penicillins and other β -lactams approved by the U.S. Food and Drug Administration for treatment of Group A and B streptococci is not necessary for clinical purposes and need not be done routinely, since as with vancomycin, resistant strains have not been recognized. However, some strains of *S. agalactiae* may give penicillin-intermediate results. Interpretive criteria are provided for pharmaceutical development, epidemiology or monitoring for emerging resistance. Any strain found to be intermediate or resistant should be referred to a reference laboratory for confirmation.

Disc diffusion tests with ampicillin, penicillin, and rifampin for *Neisseria meningitidis* are unreliable. Minimal inhibitory concentration (MIC) tests should be used for these organisms.⁷

LIMITATIONS OF THE PROCEDURE

- The test as herein described applies primarily to rapidly growing aerobic pathogens. Fastidious bacteria, other than *H. influenzae*, *N. gonorrhoeae*, *S. pneumoniae* and other streptococci, should be tested by a dilution method.⁹ Testing of anaerobes requires special procedures.¹¹
- The classifications of Resistant, Intermediate and Susceptible vary only by one millimeter, which is within normal laboratory error. Some cultures may give a borderline zone that varies from day to day or from laboratory to laboratory; such cultures are relatively uncommon.
- For detecting pneumococcal and enterococcal resistance, strictly adhere to NCCLS recommended methods.⁶
- Antimicrobial agents other than those listed in the Chart may be in current use. Susceptibility tests employing these agents should be interpreted on the basis of presence or absence of a definite zone of inhibition and should be considered as only qualitative until such time as interpretive zones have been established. All zone diameters should be recorded.
- ESBL confirmatory testing is only valid when the four discs (cefotaxime, cefotaxime/clavulanic acid, ceftazidime, ceftazidime/clavulanic acid) are used simultaneously. Individual usage of these discs is not recommended by NCCLS.^{6,7}

REFERENCES

- Bauer, A.W., W.M.M. Kirby, J.C. Sherris, and M. Turck. 1966. Antibiotic susceptibility testing by a standardized single disk method. *Am. J. Clin. Pathol.* 45:493-496.
- Ryan, K.J., F.D. Schoenknecht, and W.M.M. Kirby. 1970. Disc sensitivity testing. *Hospital Practice* 5:91-100.
- Federal Register. 1972. Rules and regulations. Antibiotic susceptibility discs. *Fed. Regist.* 37:20525-20529. Erratum, 38:2756, 1973.
- Ericsson, H.M., and J.C. Sherris. 1971. Antibiotic sensitivity testing. Report of an international collaborative study. *Acta Pathol. Microbiol. Scand. Sec. B. Suppl.* 217:1-90.
- World Health Organization Expert Committee on Biological Standardization. 1977. Technical report series 610. W.H.O., Geneva.
- National Committee for Clinical Laboratory Standards. 2003. Approved standard M2-A8. Performance standards for antimicrobial disk susceptibility tests, 8th ed. NCCLS, Wayne, Pa.
- National Committee for Clinical Laboratory Standards. 2003. M100-S13 (M2). Disk Diffusion Supplemental Tables, NCCLS, Wayne Pa.
- Baker, C.N., C. Thornsberry, and R.W. Hawkinson. 1983. Inoculum standardization in antimicrobial susceptibility testing: evaluation of overnight agar cultures and the rapid inoculum standardization system. *J. Clin. Microbiol.* 17:450-457.
- National Committee for Clinical Laboratory Standards. 2003. Approved standard M7-A6. Methods for dilution antimicrobial susceptibility tests for bacteria that grow aerobically, 6th ed. NCCLS, Wayne, Pa.
- York, M.K., L. Gibbs, F. Chehab, and G.F. Brooks. 1996. Comparison of PCR detection of *mecA* with standard susceptibility testing methods to determine methicillin resistance in coagulase-negative staphylococci. *J. Clin. Microbiol.* 34: 249-253.
- National Committee for Clinical Laboratory Standards. 2001. Approved standard M11-A5. Methods for antimicrobial susceptibility testing of anaerobic bacteria, 5th ed. NCCLS, Wayne, Pa.
- Bushby, S.R.M. 1973. Trimethoprim-sulfamethoxazole: in vitro microbiological aspects, p. 10-30. *In* M. Finland and E.H. Kass (ed.), Trimethoprim-sulfamethoxazole: microbiological, pharmacological, and clinical considerations. University of Chicago Press, Chicago.

Zone Diameter Interpretive Chart †												
Antimicrobial Agent	Code	Disc Potency	Zone Diameter Interpretive Standards (mm)			Control Zone Diameter Limits (mm)						
			Resistant	Intermediate ^a	Susceptible ^b	<i>E. coli</i> ATCC 25922	<i>S. aureus</i> ATCC 25923	<i>P. aeruginosa</i> ATCC 27853	<i>H. influenzae</i> ATCC 49247 ^c	<i>H. influenzae</i> ATCC 49766 ^c	<i>N. gonorrhoeae</i> ATCC 49226 ^d	<i>S. pneumoniae</i> ATCC 49619 ^e
Amdinocillin ^f <i>Enterobacteriaceae</i>	AMD-10	10 µg	≤15	—	≥16	23 – 29	—	—	—	—	—	—
Amikacin <i>Enterobacteriaceae</i> , <i>P. aeruginosa</i> , <i>Acinetobacter</i> and staphylococci	AN-30	30 µg	≤14	15 – 16	≥17	19 – 26	20 – 26	18 – 26	—	—	—	—
Amoxicillin/Clavulanic Acid ^{g,h,i} <i>Enterobacteriaceae</i> <i>Staphylococcus</i> spp. ^j <i>Haemophilus</i> spp. ^{c,k}	AmC-30	20/10 µg	≤13	14 – 17	≥18	18 – 24 ⁱⁱ	28 – 36	—	—	—	—	—
Ampicillin ^{h,i} <i>Enterobacteriaceae</i> ⁱⁱ and <i>V. cholerae</i> ^m Staphylococci spp. ^{i,iii} <i>Enterococcus</i> spp. ^{n,o,ii} <i>Listeria monocytogenes</i> ^f <i>Haemophilus</i> spp. ^{c,k,p} Streptococci (non- <i>S. pneumoniae</i> , β -hemolytic only) ^{e,i,aaa,ccc}	AM-10	10 µg	≤13	14 – 16	≥17	16 – 22	27 – 35	—	—	15 – 23 ^c	—	30 – 36 ^e
Ampicillin/Sulbactam ^{g,h,i} <i>Enterobacteriaceae</i> , <i>P. aeruginosa</i> , <i>Acinetobacter</i> ^q and staphylococci ^j <i>Haemophilus</i> spp. ^{c,k}	SAM-20	10/10 µg	≤11	12 – 14 ⁱⁱ	≥15 ⁱⁱ	19 – 24 ⁱⁱ	29 – 37	—	—	—	—	—
Azithromycin <i>Staphylococcus</i> spp. ^r <i>Haemophilus</i> spp. ^c <i>S. pneumoniae</i> and other streptococci ^{e,r,s}	AZM-15	15 µg	≤13	14 – 17	≥18	—	21 – 26	—	—	14 – 22 ^c	—	19 – 25 ^e
Azlocillin <i>P. aeruginosa</i>	AZ-75	75 µg	≤17	—	≥18	—	—	24 – 30	—	—	—	—
Aztreonam <i>Enterobacteriaceae</i> , ^t <i>P. aeruginosa</i> & <i>Acinetobacter</i> <i>Haemophilus</i> spp. ^c	ATM-30	30 µg	≤15	16 – 21	≥22	28 – 36	—	23 – 29	—	30 – 38 ^c	—	—
Bacitracin ^f <i>Enterobacteriaceae</i> and <i>Acinetobacter</i> <i>P. aeruginosa</i>	B-10	10 U	≤ 8	9 – 12	≥13	—	12 – 22	—	—	—	—	—
Carbencillin <i>Enterobacteriaceae</i> and <i>Acinetobacter</i> <i>P. aeruginosa</i>	CB-100	100 µg	≤19 ⁱⁱ	20 – 22 ⁱⁱ	≥23	23 – 29	—	18 – 24	—	—	—	—
Cefaclor ^{h,i} <i>Enterobacteriaceae</i> ^u and staphylococci ^j <i>Haemophilus</i> spp. ^{c,k}	CEC-30	30 µg	≤14	15 – 17	≥18	23 – 27	27 – 31	—	—	25 – 31 ^c	—	—
Cefamandole <i>Enterobacteriaceae</i> and staphylococci ^j	MA-30	30 µg	≤14	15 – 17	≥18	26 – 32	26 – 34	—	—	—	—	—
Cefazolin <i>Enterobacteriaceae</i> ^u and staphylococci ^j	CZ-30	30 µg	≤14	15 – 17	≥18	21 – 27 ⁱⁱ	29 – 35	—	—	—	—	—
Cefdinir ^h <i>Enterobacteriaceae</i> ^{kk} and methicillin-susceptible staphylococci ^j <i>Haemophilus</i> spp. ^c	CDR-5	5 µg	≤16	17 – 19	≥20	24 – 28	25 – 32	—	—	24 – 31 ^c	—	—
Cefepime ^{h,i} <i>Enterobacteriaceae</i> , <i>P. aeruginosa</i> , <i>Acinetobacter</i> and staphylococci ^j <i>Haemophilus</i> spp. ^c <i>N. gonorrhoeae</i> ^d Viridans Streptococci (non- <i>S. pneumoniae</i>) ^{e,ccc} Streptococci (non- <i>S. pneumoniae</i> , β -hemolytic only) ^{aaa,ccc}	FEP-30	30 µg	≤14	15 – 17	≥18	31 – 37 ⁱⁱⁱ	23 – 29	24 – 30	—	25 – 31 ^c	—	37 – 46 ^{d,ii} 28 – 35 ^e
Cefixime ^h <i>Enterobacteriaceae</i> ^v <i>Haemophilus</i> spp. ^c <i>N. gonorrhoeae</i> ^d	CFM-5	5 µg	≤15	16 – 18	≥19	23 – 27	—	—	—	25 – 33 ^c	—	37 – 45 ^d
Cefmetazole <i>Enterobacteriaceae</i> and staphylococci ^j <i>N. gonorrhoeae</i> ^d	CMZ-30	30 µg	≤12	13 – 15	≥16	26 – 32	25 – 34	—	—	—	—	31 – 36 ^d
Cefonicid <i>Enterobacteriaceae</i> and staphylococci ^j <i>Haemophilus</i> spp. ^{c,k}	CID-30	30 µg	≤14	15 – 17	≥18	25 – 29	22 – 28	—	—	30 – 38 ^c	—	—
Cefoperazone <i>Enterobacteriaceae</i> , <i>P. aeruginosa</i> , <i>Acinetobacter</i> and staphylococci ^j	CFP-75	75 µg	≤15	16 – 20	≥21	28 – 34	24 – 33	23 – 29	—	—	—	—
Cefotaxime ^h <i>Enterobacteriaceae</i> , ^{tx} <i>P. aeruginosa</i> , <i>Acinetobacter</i> and staphylococci ^j <i>Haemophilus</i> spp. ^c <i>N. gonorrhoeae</i> ^d Viridans Streptococci (non- <i>S. pneumoniae</i>) ^{e,i,ccc} Streptococci (non- <i>S. pneumoniae</i> , β -hemolytic only) ^{aaa,ccc}	CTX-30	30 µg	≤14	15 – 22 ⁱⁱ	≥23 ⁱⁱ	29 – 35	25 – 31	18 – 22	—	31 – 39 ^c	—	38 – 48 ^d 31 – 39 ^e
Cefotaxime/Clavulanic Acid ^t	CTX/CLA	30/10 µg	—	—	—	—	—	—	—	—	—	—
Cefotetan <i>Enterobacteriaceae</i> and staphylococci ^j <i>N. gonorrhoeae</i> ^d	CTT-30	30 µg	≤12	13 – 15	≥16	28 – 34	17 – 23	—	—	—	—	30 – 36 ^d
Cefoxitin <i>Enterobacteriaceae</i> and staphylococci ^j <i>N. gonorrhoeae</i> ^d	FOX-30	30 µg	≤14	15 – 17	≥18	23 – 29	23 – 29	—	—	—	—	33 – 41 ^d
Cefpodoxime ^{h,i} <i>Enterobacteriaceae</i> ^{t,u,v} and staphylococci ^j <i>Haemophilus</i> spp. ^c <i>N. gonorrhoeae</i> ^d	CPD-10	10 µg	≤17	18 – 20	≥21	23 – 28	19 – 25	—	—	25 – 31 ^c	—	35 – 43 ^d
Cefprozil ^{h,i} <i>Enterobacteriaceae</i> ^{u,v} and staphylococci ^j <i>Haemophilus</i> spp. ^{c,k}	CPR-30	30 µg	≤14	15 – 17	≥18	21 – 27	27 – 33	—	—	20 – 27 ^c	—	—
Ceftazidime <i>Enterobacteriaceae</i> , ^t <i>P. aeruginosa</i> , <i>Acinetobacter</i> and staphylococci ^j <i>Haemophilus</i> spp. ^c <i>N. gonorrhoeae</i> ^d	CAZ-30	30 µg	≤14	15 – 17	≥18	25 – 32	16 – 20	22 – 29	—	27 – 35 ^c	—	35 – 43 ^d
Ceftazidime/Clavulanic Acid ^t	CAZ/CLA	30/10 µg	—	—	—	—	—	—	—	—	—	—
Ceftibuten ^{h,i} <i>Enterobacteriaceae</i> ^{z,ii} <i>Haemophilus</i> spp. ^c	CTB-30	30 µg	≤17	18 – 20	≥21	27 – 35	—	—	—	29 – 36 ^{c,ii}	—	—
Ceftizoxime ^{h,i} <i>Enterobacteriaceae</i> , ^{tx} <i>P. aeruginosa</i> , ⁱⁱ <i>Acinetobacter</i> and staphylococci ^j <i>Haemophilus</i> spp. ^c <i>N. gonorrhoeae</i> ^d	ZOX-30	30 µg	≤14	15 – 19	≥20	30 – 36	27 – 35	12 – 17	—	29 – 39 ^c	—	42 – 51 ^d
Ceftriaxone ^h <i>Enterobacteriaceae</i> , ^{tx} <i>P. aeruginosa</i> , <i>Acinetobacter</i> and staphylococci ^j <i>Haemophilus</i> spp. ^c <i>N. gonorrhoeae</i> ^d Viridans Streptococci (non- <i>S. pneumoniae</i>) ^{e,i,ccc} Streptococci (non- <i>S. pneumoniae</i> , β -hemolytic only) ^{aaa,ccc}	CRO-30	30 µg	≤13	14 – 20 ⁱⁱ	≥21 ⁱⁱ	29 – 35	22 – 28	17 – 23	—	31 – 39 ^c	—	39 – 51 ^d 30 – 35 ^e
Cefuroxime (sodium) ^{h,i} CXM-30 <i>Enterobacteriaceae</i> ^u and staphylococci ^j (parenteral) <i>Haemophilus</i> spp. ^{c,k} <i>N. gonorrhoeae</i> ^d	CXM-30	30 µg	≤14	15 – 17	≥18	20 – 26	27 – 35	—	—	28 – 36 ^c	—	33 – 41 ^d
Cephalothin <i>Enterobacteriaceae</i> ^u and staphylococci ^j	CF-30	30 µg	≤14	15 – 17	≥18	15 – 21	29 – 37	—	—	—	—	—
Chloramphenicol <i>Enterobacteriaceae</i> , ^r <i>P. aeruginosa</i> , ^r <i>Acinetobacter</i> , ^r staphylococci, ^r enterococci ^{ty} and <i>V. cholerae</i> ^{bb} <i>Haemophilus</i> spp. ^c <i>S. pneumoniae</i> ^{e,r} Streptococci (non- <i>S. pneumoniae</i>) ^{e,r}	C-30	30 µg	≤12	13 – 17	≥18	21 – 27	19 – 26	—	—	31 – 40 ^c	—	23 – 27 ^e
Cinoxacin <i>Enterobacteriaceae</i>	CIN-100	100 µg	≤14	15 – 18	≥19	26 – 32	—	—	—	—	—	—

Zone Diameter Interpretive Chart †												
Antimicrobial Agent	Code	Disc Potency	Zone Diameter Interpretive Standards (mm)			Control Zone Diameter Limits (mm)						
			Resistant	Intermediate ^a	Susceptible ^b	<i>E. coli</i> ATCC 25922	<i>S. aureus</i> ATCC 25923	<i>P. aeruginosa</i> ATCC 27853	<i>H. influenzae</i> ATCC 49247 ^c	<i>H. influenzae</i> ATCC 49766 ^c	<i>N. gonorrhoeae</i> ATCC 49226 ^d	<i>S. pneumoniae</i> ATCC 49619 ^e
Ciprofloxacin <i>Enterobacteriaceae</i> , ^{ddd} <i>P. aeruginosa</i> , <i>Acinetobacter</i> , staphylococci ^{aa} and enterococci <i>Haemophilus</i> spp. ^c <i>N. gonorrhoeae</i> ^d	CIP-5	5 µg	≤15	16 – 20	≥21	30 – 40	22 – 30	25 – 33	34 – 42 ^c	—	48 – 58 ^d	—
Clarithromycin <i>Staphylococcus</i> spp. ^r <i>Haemophilus</i> spp. ^c <i>S. pneumoniae</i> and other streptococci ^{e,f,s}	CLR-15	15 µg	≤13	14 – 17	≥18	—	26 – 32	—	11 – 17 ^c	—	—	25 – 31 ^e
Clindamycin ^r <i>Staphylococcus</i> spp. ⁱⁱ <i>S. pneumoniae</i> and other streptococci ^e	CC-2	2 µg	≤14	15 – 20	≥21	—	24 – 30	—	—	—	—	19 – 25 ^e
Colistin ^{f,dd}	CL-10	10 µg	≤ 8	9 – 10	≥11	11 – 15	—	—	—	—	—	—
Doxycycline ^{ee} <i>Enterobacteriaceae</i> , <i>P. aeruginosa</i> , <i>Acinetobacter</i> , staphylococci and enterococci ^{yy}	D-30	30 µg	≤12	13 – 15	≥16	18 – 24	23 – 29	—	—	—	—	—
Enoxacin <i>Enterobacteriaceae</i> ^{ddd} and staphylococci ^{aa,ff} <i>N. gonorrhoeae</i> ^{d,cc}	ENX-10	10 µg	≤14	15 – 17	≥18	28 – 36	22 – 28	22 – 28	—	—	43 – 51 ^d	—
Ertapenem ^{h,i} <i>Enterobacteriaceae</i> and <i>Staphylococcus</i> spp. ^j <i>Haemophilus</i> spp. ^c	ETP-10	10 µg	≤15	16 – 18	≥19	29 – 36	24 – 31	13 – 21	20 – 28	27 – 33	—	28 – 35
Erythromycin <i>Staphylococcus</i> spp. ^r and enterococci ^{r,yy} <i>S. pneumoniae</i> and other streptococci ^{e,f,s}	E-15	15 µg	≤13	14 – 22 ⁱⁱ	≥23 ⁱⁱ	—	22 – 30	—	—	—	—	25 – 30 ^e
Fosfomicin ^z <i>E. coli</i> and <i>E. faecalis</i> only	FOS-200	200 µg	≤12	13 – 15	≥16	22 – 30 ⁱⁱ	25 – 33	—	—	—	—	—
Gatifloxacin <i>Enterobacteriaceae</i> ^{ddd} and <i>Staphylococcus</i> spp. <i>P. aeruginosa</i> , <i>Acinetobacter</i> spp. and enterococci ² <i>H. influenzae</i> ^c and <i>H. parainfluenzae</i> ^c <i>N. gonorrhoeae</i> ^d <i>S. pneumoniae</i> and other streptococci (non- <i>S. pneumoniae</i> , β-hemolytic only) ^e	GAT-5	5 µg	≤14	15 – 17	≥18	30 – 37	27 – 33	20 – 28 ^{mm}	33 – 41 ^c	—	45 – 56 ^d	24 – 31 ^e
Gentamicin Testing enterococci for high level resistance ^{n,9,99}	GM-120	120 µg	6	7 – 9 ^{hh}	≥10	—	—	—	—	—	—	—
<i>Enterobacteriaceae</i> , <i>P. aeruginosa</i> , <i>Acinetobacter</i> and staphylococci	GM-10	10 µg	≤12	13 – 14	≥15	19 – 26	19 – 27	16 – 21	—	—	—	—
Grepafloxacin <i>Enterobacteriaceae</i> ^{ddd} and staphylococci ^{aa} <i>Haemophilus</i> spp. ^c <i>N. gonorrhoeae</i> ^d <i>S. pneumoniae</i> and other streptococci ^e	GRX-5	5 µg	≤14	15 – 17	≥18	28 – 36	26 – 31	20 – 27 ⁱⁱ	32 – 39 ^c	—	44 – 52 ^d	21 – 28 ^e
Imipenem ^{h,i} <i>Enterobacteriaceae</i> , <i>P. aeruginosa</i> , <i>Acinetobacter</i> and staphylococci ^j <i>Haemophilus</i> spp. ^c	IPM-10	10 µg	≤13	14 – 15	≥16	26 – 32	—	20 – 28	—	—	—	—
Kanamycin <i>Enterobacteriaceae</i> and staphylococci	K-30	30 µg	≤13	14 – 17	≥18	17 – 25	19 – 26	—	—	—	—	—
Levofloxacin <i>Enterobacteriaceae</i> ^{ddd} , <i>P. aeruginosa</i> , <i>Acinetobacter</i> , staphylococci ^{aa} and enterococci <i>Haemophilus</i> spp. ^c <i>S. pneumoniae</i> and other streptococci (non- <i>S. pneumoniae</i> , β-hemolytic only) ^e	L VX-5	5 µg	≤13	14 – 16	≥17	29 – 37	25 – 30	19 – 26	32 – 40 ^c	—	—	20 – 25 ^e
Linezolid <i>Staphylococcus</i> spp. <i>Enterococcus</i> spp. <i>S. pneumoniae</i> and other streptococci ^e	LZD-30	30 µg	—	—	≥21	—	25 – 32 ⁱⁱ	—	—	—	—	25 – 34 ^{e,ii}
Lomefloxacin <i>Enterobacteriaceae</i> ^{ddd} , <i>P. aeruginosa</i> , <i>Acinetobacter</i> and staphylococci ^{aa} <i>Haemophilus</i> spp. ^c <i>N. gonorrhoeae</i> ^{d,cc}	LOM-10	10 µg	≤18	19 – 21	≥22	27 – 33	23 – 29	22 – 28	33 – 41 ^c	—	45 – 54 ^d	—
Loracarbef ^{h,i} <i>Enterobacteriaceae</i> ^{u,kk} and staphylococci ^j <i>Haemophilus</i> spp. ^{c,k}	LOR-30	30 µg	≤14	15 – 17	≥18	23 – 29	23 – 31	—	—	26 – 32 ^c	—	—
Meropenem ^{h,i} <i>Enterobacteriaceae</i> , <i>P. aeruginosa</i> , <i>Acinetobacter</i> and staphylococci ^j <i>Haemophilus</i> spp. ^c	MEM-10	10 µg	≤13	14 – 15	≥16	28 – 34	29 – 37 ⁱⁱ	27 – 33	—	—	—	—
Methicillin ^{ll}	—	—	—	—	≥20	—	—	—	20 – 28 ^c	—	—	—
Mezlocillin ⁱⁱ <i>Enterobacteriaceae</i> and <i>Acinetobacter</i> <i>P. aeruginosa</i>	MZ-75	75 µg	≤17	18 – 20	≥21	23 – 29	—	19 – 25	—	—	—	—
Minocycline ^{ee} <i>Enterobacteriaceae</i> , <i>P. aeruginosa</i> , <i>Acinetobacter</i> , staphylococci and enterococci ^{yy}	MI-30	30 µg	≤14	15 – 18	≥19	19 – 25	25 – 30	—	—	—	—	—
Moxalactam <i>Enterobacteriaceae</i> , <i>P. aeruginosa</i> , <i>Acinetobacter</i> and staphylococci ^j	MOX-30	30 µg	≤14	15 – 22	≥23	28 – 35	18 – 24	17 – 25	—	—	—	—
Moxifloxacin <i>Enterobacteriaceae</i> and <i>Staphylococcus</i> spp. <i>H. influenzae</i> ^c and <i>H. parainfluenzae</i> ^c <i>S. pneumoniae</i> ^e	MXF-5	5 µg	≤15	16 – 18	≥19	28 – 35	28 – 35	—	31 – 39 ^c	—	—	25 – 31 ^e
Nafcillin <i>Staphylococcus aureus</i> ^{j,nn}	NF-1	1 µg	≤10	11 – 12	≥13	—	16 – 22	—	—	—	—	—
Nalidixic Acid <i>Enterobacteriaceae</i> ^z	NA-30	30 µg	≤13	14 – 18	≥19	22 – 28	—	—	—	—	—	—
Neomycin ^f	N-30	30 µg	≤12	13 – 16	≥17	17 – 23	18 – 26	—	—	—	—	—
Netilmicin <i>Enterobacteriaceae</i> , <i>P. aeruginosa</i> , <i>Acinetobacter</i> and staphylococci	NET-30	30 µg	≤12	13 – 14	≥15	22 – 30	22 – 31	17 – 23	—	—	—	—
Nitrofurantoin <i>Enterobacteriaceae</i> , staphylococci and enterococci	F/M-300	300 µg	≤14	15 – 16	≥17	20 – 25	18 – 22	—	—	—	—	—
Norfloxacin ⁱⁱ <i>Enterobacteriaceae</i> ^{ddd} , <i>P. aeruginosa</i> , <i>Acinetobacter</i> , staphylococci ^{aa} and enterococci	NOR-10	10 µg	≤12	13 – 16	≥17	28 – 35	17 – 28	22 – 29	—	—	—	—
Novobiocin ^f (Mueller Hinton agar with sheep blood for veterinary use)	NB-30	30 µg	≤17	18 – 21	≥22	—	22 – 31	—	—	—	—	—
Ofloxacin <i>Enterobacteriaceae</i> ^{ddd} , <i>P. aeruginosa</i> , <i>Acinetobacter</i> and staphylococci ^{aa} <i>Haemophilus</i> spp. ^c <i>N. gonorrhoeae</i> ^d <i>S. pneumoniae</i> ^e and other streptococci (non- <i>S. pneumoniae</i> , β-hemolytic only) ^e	OFX-5	5 µg	≤12	13 – 15	≥16	29 – 33	24 – 28	17 – 21	31 – 40 ^c	—	43 – 51 ^d	16 – 21 ^e
Oxacillin <i>Staphylococcus aureus</i> ^{j,nn,oo} Staphylococci, coagulase-negative ^{j,nn,oo} <i>S. pneumoniae</i> (for penicillin G susceptibility) ^{e,h}	OX-1	1 µg	≤10	11 – 12	≥13	—	18 – 24	—	—	—	—	≤12 ^{e,bbb}
Oxolinic Acid ^f <i>Staphylococcus</i> spp. ^{j,pp} <i>Enterococcus</i> spp. ^{n,o} <i>L. monocytogenes</i> ^f <i>N. gonorrhoeae</i> ^{d,qq,ii} Streptococci (non- <i>S. pneumoniae</i> , β-hemolytic only) ^{e,i,r,t,aaa,ccc}	OA-2	2 µg	≤10	—	≥11	20 – 24	10 – 13	—	—	—	—	—
Penicillin ^h <i>Staphylococcus</i> spp. ^{j,pp} <i>Enterococcus</i> spp. ^{n,o} <i>L. monocytogenes</i> ^f <i>N. gonorrhoeae</i> ^{d,qq,ii} Streptococci (non- <i>S. pneumoniae</i> , β-hemolytic only) ^{e,i,r,t,aaa,ccc}	P-10	10 U	≤28	—	≥29	—	26 – 37	—	—	—	—	26 – 34 ^d
Piperacillin <i>Enterobacteriaceae</i> and <i>Acinetobacter</i> <i>P. aeruginosa</i>	PIP-100	100 µg	≤17	18 – 20	≥21	24 – 30 ⁹	—	25 – 33	—	—	—	—
Piperacillin/ Tazobactam ⁹ <i>Enterobacteriaceae</i> and <i>Acinetobacter</i> ⁱⁱ <i>Staphylococcus</i> spp. ^{j,ii} and <i>P. aeruginosa</i> ⁱⁱ	TZP-110	100/10 µg	≤17	18 – 20	≥21	24 – 30 ⁹	27 – 36	25 – 33 ⁱⁱ	—	—	—	—
Polymyxin B ^{f,dd}	PB-300	300 U	≤ 8	9 – 11	≥12	12 – 16	—	—	—	—	—	—
Quinupristin/Dalfopristin SYN-15 4.5/10.5 µg <i>Staphylococcus</i> spp., <i>Enterococcus faecium</i> and <i>S. pyogenes</i> ^e only	SYN-15	4.5/10.5 µg	≤15	16 – 18	≥19	—	21 – 28 ⁱⁱ	—	—	—	—	19 – 24 ^{e,ii,cc}
Rifampin <i>Staphylococcus</i> spp. ^{ss} and <i>Enterococcus</i> spp. ^{yy} <i>Haemophilus</i> spp. ^c <i>S. pneumoniae</i> ^{e,ss}	RA-5	5 µg	≤16	17 – 19	≥20	8 – 10	26 – 34	—	22 – 30 ^c	—	—	25 – 30 ^e
Sparfloxacin <i>Staphylococcus</i> spp. ^{aa} <i>S. pneumoniae</i> ^e	SPX-5	5 µg	≤15	16 – 18	≥19	30 – 38	27 – 33	21 – 29 ⁱⁱ	—	—	—	21 – 27 ^e
Spectinomycin <i>N. gonorrhoeae</i> ^d	SPT-100	100 µg	≤14	15 – 17 ^w	≥18	—	—	—	—	—	—	23 – 29 ^d
Streptomycin Testing enterococci for high level resistance ^{n,9,99}	S-300	300 µg	6	7 – 9 ^{hh}	≥10	—	—	—	—	—	—	—
<i>Enterobacteriaceae</i>	S-10	10 µg	≤11	12 – 14	≥15	12 – 20	14 – 22	—	—	—	—	—
Sulfisoxazole ^{tt} <i>Enterobacteriaceae</i> , <i>P. aeruginosa</i> , <i>Acinetobacter</i> , staphylococci and <i>V. cholerae</i> ^{tn}	G-.25	250 µg	≤12	13 – 16	≥17	15 – 23	23 – 34	—	—	—	—	—
Tetracycline ^{ee} <i>Enterobacteriaceae</i> , <i>P. aeruginosa</i> , <i>Acinetobacter</i> , staphylococci, enterococci ^{yy} and <i>V. cholerae</i> ^{tn} <i>Haemophilus</i> spp. ^c <i>N. gonorrhoeae</i> ^{d,uu} <i>S. pneumoniae</i> and other streptococci ^e	Te-30	30 µg	≤14	15 – 18	≥19	18 – 25	24 – 30	—	14 – 22 ^c	—	30 – 42 ^d	27 – 31 ^e
Ticarcillin <i>Enterobacteriaceae</i> and <i>Acinetobacter</i> <i>P. aeruginosa</i>	TIC-75	75 µg	≤14	15 – 19	≥20	24 – 30	—	21 – 27 ⁱⁱ	—	—	—	—
Ticarcillin/ Clavulanic Acid ⁹ <i>Enterobacteriaceae</i> and <i>Acinetobacter</i> <i>P. aeruginosa</i> <i>Staphylococcus</i> spp. ^j	TIM-85	75/10 µg	≤14	15 – 19	≥20	24 – 30 ^{9,ii}	29 – 37	20 – 28	—	—	—	—
Tobramycin <i>Enterobacteriaceae</i> , <i>P. aeruginosa</i> , <i>Acinetobacter</i> and staphylococci	NN-10	10 µg	≤12	13 – 14	≥15	18 – 26	19 – 29	19 – 25	—	—	—	—
Trimethoprim <i>Enterobacteriaceae</i> and staphylococci	TMP-5	5 µg	≤10	11 – 15	≥16	21 – 28	19 – 26	—	—	—	—	—
Trimethoprim/ Sulfamethoxazole ^{tt} <i>Enterobacteriaceae</i> , <i>P. aeruginosa</i> , <i>Acinetobacter</i> , staphylococci and <i>V. cholerae</i> ^{tn} <i>Haemophilus</i> spp. ^c <i>S. pneumoniae</i> ^e	SXT	1.25 µg 23.75 µg	≤10	11 – 15	≥16	23 – 29 ⁱⁱ	24 – 32	—	24 – 32 ^c	—	—	20 – 28 ^e
Trovafloxacin <i>Haemophilus</i> spp. ^c <i>N. gonorrhoeae</i> ^d <i>S. pneumoniae</i> and other streptococci ^e	TVA-10	10 µg	≤13 ^f	14 – 16 ^f	≥17 ^f	29 – 36	29 – 35	21 – 27	32 – 39 ^c	—	42 – 55 ^d	25 – 32 ^e

aaa Strains of β -hemolytic streptococci with ampicillin, cefepime, cefotaxime, ceftriaxone or penicillin zone diameters of less than 24 mm have not been observed; submit such strains to a reference laboratory.

bbb Deterioration in oxacillin disc content is best assessed with *S. aureus* ATCC 25293, with an acceptable zone diameter of 18 - 24 mm.

ccc For ampicillin, cefepime, cefotaxime, ceftriaxone and penicillin, **Streptococci, β -hemolytic only** includes the large-colony-forming pyogenic strains of streptococci with Group A (*S. pyogenes*), C or G antigens and strains with Group B (*S. agalactiae*) antigen. For cefepime, cefotaxime and ceftriaxone, **Viridans Streptococci** includes small-colony-forming β -hemolytic strains with Group A, C, F or G antigens (*S. anginosus*, previously termed *S. milleri*) as well as *S. mitis*, *S. oralis*, *S. sanguis*, *S. salivarius*, *S. intermedius*, *S. constellatus*, *S. mutans* and *S. bovis*.

ddd Fluoroquinolone-susceptible strains of *Salmonella* that test resistant to nalidixic acid may be associated with clinical failure or delayed response in fluoroquinolone-treated patients with extraintestinal salmonellosis. Testing of extraintestinal *Salmonella* isolates for nalidixic acid resistance may be considered.

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-  In Vitro Diagnostic Medical Device
-  Temperature limitation
-  Catalog number
-  Consult Instructions for Use
-  Authorized Representative in the European Community
-  Batch Code (Lot)
-  Manufacturer
-  Use by
YYYY-MM-DD / YYYY-MM (MM = end of month)

 Becton, Dickinson and Company
7 Loveton Circle
Sparks, Maryland 21152 USA
800-638-8663

 BENEX Limited
Bay K 1a/d, Shannon Industrial Estate
Shannon, County Clare, Ireland
Tel: 353-61-47-29-20
Fax: 353-61-47-25-46

ATCC is a trademark of the American Type Culture Collection. Prompt is a trademark of 3M. BD, BD Logo, BBL, Cefinase, IsoVitalEX, Sensi-Disc and Trypicase are trademarks of Becton, Dickinson and Company © 2004 BD.