Rev From	Rev To	ECO #	Date	Appr.
0499	0703	1998-03		

Notes


- BD Cat. No. 221550, 221553
- Blank (Sheet) Size : Length: 13.75" Width: 16.125"
 Number of Pages: 6 Number of Sheets: 1
 Page Size: Length 13.75" Width 5.375" Final Folded Size: 1.71875" x 5.375"
- Style (see illustrations below): 4


#1

#2

#3


#4


#5

#6

#7

- See Specification Control No. NA for Material Information
- Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS # Black
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

VS Controlled by BD Caribe, Ltd.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 250 Schilling Circle Cockeysville, MD. 21030-0243 USA
Proofer	Date		
Checked By	Date		
Part Number:	8840671	Category and Description	Sheet: 1 of 7
		Package Insert Taxo Nystatin	Scale:
			A

BBL™ Taxo™ Nystatin Discs for Inhibition of Yeasts

8840671
2003/07

English: pages 1 Italiano: pagine 3 - 4
 Français: pages 1 - 2 Español: páginas 4
 Deutsch: Seiten 2 - 3

See symbol glossary at end of insert. / Se symbolglossaret i slutningen af indlægssedlen. / Voir le glossaire des symboles à la fin de la notice. / Siehe Symbol-Erklärungen am Ende der Packungsbeilage. / Vedere il glossario dei simboli alla fine del foglio illustrativo. / Consulte o glossário de símbolos no fim do folheto informativo. / Consulte el glosario de símbolos al final del prospecto. / Se symbolförteckningen vid slutet av bipacksedeln.

Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Contacte o seu representante local da BD para obter instruções. / Kontakta lokal Becton Dickinson-representant för anvisningar.

INTENDED USE

Taxo Nystatin Discs are for use in a qualitative procedure for the differentiation and isolation of bacteria from old specimens and those likely to contain a mixed or large population of bacteria and yeasts.¹ They are not for susceptibility testing.

SUMMARY AND EXPLANATION

Nystatin (fungicidin) was isolated from *Streptomyces noursei* by Hazen and Brown and reported in 1951. They found that a crude preparation protected mice from lethal doses of *Candida albicans*. The drug was later refined and used for clinical fungal infections.²

Nystatin is an antimicrobial drug which affects cell membranes and is useful for the suppression of yeasts, especially *Candida*, and some bacteria.

The drug is impregnated into paper discs and applied to agar surfaces. The drug diffuses into the medium inhibiting susceptible organisms.

PRINCIPLES OF THE PROCEDURE

The discs have been found to be particularly valuable in detection of *Neisseria* in specimens. For *Neisseria*, the **Taxo** Nystatin Discs may be placed on the area of initial (heaviest) inoculation of **BBL** Thayer-Martin media or Transgrow Medium to supplement the nystatin already present in these media, which is insufficient to inhibit some strains of yeasts.³ Nystatin Discs are also recommended for use in or on other selective and nonselective media for examination of miscellaneous specimens.

Plates are observed for bacterial growth within zones where nystatin discs have inhibited the growth of yeasts or other undesirable microorganisms.

REAGENTS

Taxo Nystatin Discs are made from high quality absorbent paper impregnated with about 100 units of nystatin* per disc. *Mycostatin™: Bristol-Myers-Squibb.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Directions for use should be read and followed carefully.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures.

Storage Instructions: On receipt, store at -20 to +8°C. After use, store cartridge to protect product integrity at 2 to 8°C. Apply control tests at time of use.

Use the oldest discs first and discard expired discs. Also discard cartridges from which discs have been frequently removed during a week or so. Discard containers left out overnight in the laboratory, or else test the discs for performance.

The expiration date applies to produce in intact container stored as directed. Do not open until ready to use.

PROCEDURE

Material Provided: **Taxo** Nystatin Discs

Materials Required But Not Provided:

1. Suitable agar medium, such as **Trypticase™** Soy Agar with 5% Sheep Blood, Modified Thayer-Martin (MTM) Agar or Transgrow Medium.
2. Sterile forceps or Single Disc Dispenser.

Test Procedure:

The recommended method is to streak or swab the specimen on an agar medium suitable for isolation of the desired bacteria. Place one, or preferably two, **Taxo** Nystatin Discs on the area of initial (heaviest) inoculation. When two discs are used, place them approximately 15 mm apart.

User Quality Control: Control organisms producing known positive and negative results should be included at time of test. Cultures of *Candida albicans* and *Neisseria gonorrhoeae* are recommended.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Observe for bacterial growth within zones where the nystatin discs have inhibited growth of yeasts or other undesirable microorganisms.

PERFORMANCE CHARACTERISTICS

Prior to release, all lots of **Taxo** Nystatin Discs are tested to specific product specifications. Samples of the lot are assayed for nystatin content using *Candida (pseudotropicalis) kefyr* (ATCC™ 8553) as the test organism. The assay method is based on the protocol for antimicrobial agents in the United States *Code of Federal Regulations*, Title 21, Part 460, Subpart A (1998). The acceptable average potency of the **Taxo** Nystatin Discs is 80% to 140% of the theoretical disc potency.

AVAILABILITY

Cat. No.	Description
231550	BBL Taxo Nystatin Discs, Cartridge of 50 discs
231553	BBL Taxo Nystatin Discs, Pkg. of 10 cartridges

REFERENCES

1. Vera, H.D. 1971. Quality control in diagnostic microbiology. Health Lab. Sci. 8:175-189.
2. Hazen, E.L., and R. Brown. 1951. Fungicidin, an antibiotic produced by a soil actinomycete. Proc. Soc. Exp. Biol. Med. 76:93.
3. Thayer, J.D., and J.E. Martin, Jr. 1966. Improved medium selective for cultivation of *N. gonorrhoeae* and *N. meningitidis*. Public Health Rep. 81:559-562.

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© 2003 BD.
ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo, BBL, Taxo and Trypticase are trademarks of Becton, Dickinson and Company.
BENEX Limited