
Revisions

Rev from	Rev to	ECO #
1003	0804	3022-04

Notes:

1. BD Cat. Number 235201, 240764, 240765
2. Blank (Sheet) Size : Length: 14" Width: 35"
 Number of Pages: 14 Number of Sheets: 1
 Page Size: Length 14" Width 5" Final Folded Size: 3.5" x 2.5"
3. Style (see illustrations below): #4

4. See Specification Control Number L000135 for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS #2755 Blue
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	<p style="font-size: small;">COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION</p> <p style="font-weight: bold; font-size: large;">BD</p> <p>Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA</p>	<p style="font-weight: bold; font-size: x-large;">A</p>	
Proofer	Date			
Checked By	Date			
Part Number: L000135		Category and Description Package Insert VDRL Antigen	Sheet: 1 of 15	
			Scale: 1 : 1	

BD VDRL Antigen with Buffered Saline VDRL Test Control Serum Set

English: pages 1 – 3 Italiano: pagina 9 – 11
Français: page 3 – 6 Español: página 11 – 14
Deutsch: Seite 6 – 9

L000135
2004/08

See symbol glossary at end of insert. / Viz popis symbolů na konci příbalového letáku. / Se symbolglossaret i slutningen af indlægssedlen. / Zie lijst met symbolen aan het einde van de bijsluiter. / Vaadake sümboleite seletust infolehe lõpus. / Katso pakkausteosten lopussa olevaa kuvamerkkien sanastoa. / Voir le glossaire des symboles à la fin de la notice. / Siehe Symbol-Erklärungen am Ende der Packungsbeilage. / Δείτε το γλωσσάριο των συμβόλων στο τέλος του ένθετου. / A jelmagyarazát a használati utasítás végén található. / Vedere il glossario dei simboli alla fine del foglio illustrativo. / Zr. informacinio lapelio pabaigoje pateikiama simbolių glosarijų. / Se i symbolforklaringen på slutten av produktvedlegget. / Zobacz objaśnienie symboli na końcu ulotki. / Consulte o glossário de símbolos no fim do folheto informativo. / Pozri slovník symbolov na konci letáku. / Consulte el glosario de símbolos al final del prospecto. / Se symbolförteckningen vid slutet av bipacksedeln.

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukities vietos BD įgaliotojo atstovo. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcje ziskate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar.

INTENDED USE

VDRL Antigen with Buffered Saline is recommended for use in the Venereal Disease Research Laboratory (VDRL)¹ test for the detection of reagin, an antibody-like substance, by the qualitative and quantitative slide flocculation tests.

VDRL Test Control Serum Set is recommended for use in the quality control testing of VDRL Antigen by the slide flocculation test.

SUMMARY AND EXPLANATION

Treponema pallidum is the causative agent of syphilis. Syphilis is a chronic infection with many clinical manifestations which occur in distinct stages. Specific laboratory tests are recommended for the detection of each stage of the disease.²⁻⁴

The VDRL Antigen is a nontreponemal antigen composed of cardiolipin cholesterol and lecithin. The nontreponemal tests measure anti-lipid antibodies, which are formed by the host in response to lipids released from damaged host cells early in infection with *T. pallidum*, and lipid-like material from the treponemal cell surface.⁵ During syphilis infection, an antibody-like substance called reagin can be detected in the patient's serum. In syphilis infection of the central nervous system, reagin can be detected in the cerebrospinal fluid (CSF).

Reactive nontreponemal tests confirm the diagnosis in the presence of early or late lesion syphilis. They offer a clue in latent subclinical syphilis, and are effective tools for detecting cases in epidemiologic investigations. Nontreponemal tests are superior to the treponemal test for following the response to therapy.²

Nontreponemal antigen tests are not entirely specific for syphilis, nor do they have satisfactory sensitivity in all stages of syphilis. Whenever the results of a nontreponemal antigen test disagree with the clinical impression, a treponemal antigen test such as the FTA-ABS^{2,3} should be performed. Nontreponemal tests such as the VDRL are used to screen patient serum, while treponemal tests such as the FTA-ABS are used for confirmation. The likelihood of obtaining a reactive VDRL test result in various stages of untreated syphilis has been reported as follows.³

Stages of Untreated Syphilis	% Reactive VDRL Test
Primary	78
Secondary	100
Latent	96
Late	71

PRINCIPLES OF THE PROCEDURE

In the VDRL Test procedures, the patient's serum is heat-inactivated and mixed with a buffered saline suspension of VDRL Antigen containing cardiolipin, lecithin and cholesterol. The combination of reagin and VDRL Antigen forms microscopic clumping called flocculation. With modification, the serum procedure can be used for testing CSF.¹

REAGENTS

VDRL Antigen is 0.03% cardiolipin and 0.9% cholesterol dissolved in absolute alcohol with sufficient lecithin (approximately 0.18 - 0.20%) to produce standard reactivity. It is prepared with modifications according to the directions given by Harris, Rosenberg and Riedel.⁶ Cardiolipin and lecithin are prepared according to directions given by Pangborn.^{7,8,9}

VDRL Buffered Saline is a 1% sodium chloride solution, pH 6.0 ± 0.1, with 0.5% formaldehyde as a preservative.

Nontreponemal Antigen Reactive Serum is a lyophilized human serum with 0.02% thimerosal as a preservative which is standardized to provide a reactive reading when tested according to the USR or VDRL test procedure.

VDRL Weakly Reactive Serum is a lyophilized human serum with 0.02% thimerosal as a preservative which is standardized to provide a weakly reactive reading when tested according to the VDRL test procedure.

Nontreponemal Antigen Nonreactive Serum is a lyophilized human serum with 0.02% thimerosal as a preservative which is standardized to provide a nonreactive reading when tested according to the USR or VDRL test procedure.

Precautions:

- For *in vitro* Diagnostic Use.
- WARNING: POTENTIAL BIOHAZARDOUS REAGENTS.** Each donor unit used in preparation of VDRL Test Control Serum Set was tested by an FDA licensed method for the presence of the antibody to human immunodeficiency virus (HIV) as well as for hepatitis B surface antigen (HbsAg) and found to be negative (were not repeatedly reactive). Because no test method can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent, these reagents should be handled at the Biosafety Level 2 as recommended for any potentially infectious human serum or blood specimen in the Centers for Disease Control/National Institutes of Health Manual *Biosafety in Microbiological and Biomedical Laboratories*, 1999.
- Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"¹⁰⁻¹³ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.
- VDRL Antigen
HIGHLY FLAMMABLE. IRRITATING TO EYES, RESPIRATORY SYSTEM AND SKIN (US) POSSIBLE RISK OF IRREVERSIBLE EFFECTS (US), POSSIBLE RISK OF HARM TO THE UNBORN CHILD (US) TARGET ORGANS; Blood, Intestines, Liver, Muscle, Nerves.
Avoid contact with skin and eyes. Do not breathe mist. Wear suitable protective clothing. Keep container tightly closed. Keep away from sources of ignition. No smoking.
- VDRL Test Control Serum Set
This Product Contains Dry Natural Rubber

Storage Instructions:

Store VDRL Antigen at room temperature (15 – 30°C) in the dark.

Store VDRL Buffered Saline at 15 – 30°C. After bottle is opened, store at 2 – 8°C.

Store lyophilized control sera in the VDRL Test Control Serum Set at 2 – 8°C.

Upon re-hydration of the control sera, store at 2 – 8°C if contents will be consumed within one day. Alternatively, divide the control sera into aliquots sufficient for one day of testing and store at ≤ -20°C.

Do not use if reagents show signs of contamination, evaporation, precipitation or other signs of deterioration.

SPECIMEN COLLECTION AND PREPARATION

- Collect 5 – 8 mL of blood by aseptic venipuncture into clean, dry tube without an anticoagulant.
- Allow blood to clot at room temperature, and centrifuge to obtain serum.
- Store serum specimens at room temperature for up to 4 h; after 4 h, store at 2 – 8°C. Serum specimens may be refrigerated for up to 5 days, then frozen at < -20°C. Avoid repeat freezing – thawing of specimens.
- Serum specimens must be clear, free of hemolysis and show no visible evidence of bacterial contamination (turbidity, hemolysis or particulate matter). Refer to appropriate references for more information on collection of specimens.^{1,3,14}
- Before testing, heat the test sera at 56°C for 30 min. Specimens that are not tested within 4 h must be reheated for 10 min at 56°C.
- Specimens must be at 23 – 29°C when tested.

PROCEDURE

Materials Provided: VDRL Antigen with Buffered Saline, VDRL Test Control Serum Set.

Materials Required But Not Provided:

0.9% Saline

Nondisposable syringe, 1cc

Nondisposable calibrated needles without bevel:

Serum test: 18 gauge

CSF test: 21 or 22 gauge

Bottles, 30 mL, round, narrow-mouthed, 35 mm in diameter with glass stoppers and a flat inner bottom surface

NOTE: Properly washed 25 mL glass-stoppered Erlenmeyer flasks may be substituted if 30 mL bottles are not available.

Micropipettor, 50 µL

Pipettes, serologic, graduated to tip:

1.0 mL graduated in 1/100 mL,

5.0 mL, graduated in 1/10 mL,

10.0 mL graduated in 1/10 mL;

Slides:

Serum test: 2 x 3 inches with paraffin or ceramic rings approximately 14 mm in diameter and high enough to prevent spillage during rotation.

CSF test: Kline concavity slides, 3 x 2 inches x 3 mm thick, 12 concavities measuring 16 mm in diameter and 1.75 mm in depth.

Slide holder for 2 x 4 inch slides

Mechanical rotator adjustable to 180 ± 2 rpm, circumscribing a circle 19 mm in diameter on a horizontal plane

Water bath, 56°C

Light microscope with 10x ocular and 10x objective

Sterile distilled or deionized water

Absolute alcohol

Acetone

Timer

Reagent Preparation

- To rehydrate control sera in the VDRL Test Control Serum Set, add 3 mL sterile purified water and rotate gently to completely dissolve the contents.
- VDRL Antigen and VDRL Buffered Saline are ready to use in preparing VDRL Antigen suspension.

Preparation of Specific Glassware

Wash syringes with needles and suspension bottles by hand in the following manner:

- Rinse with tap water.
- Soak and wash thoroughly in a glassware detergent solution.
- Rinse with tap water 6 – 8 times.
- Rinse with unused distilled or deionized water.
- Rinse with absolute alcohol.
- Rinse with acetone.
- Air dry until acetone odor is completely eliminated.
- Remove needles from syringes for storage.

Ceramic-ringed slides:

- Rinse with tap water.
- Wash with a glassware detergent solution.
- Rinse with tap water 3 – 4 times.
- Rinse with unused distilled or deionized water.
- Wipe dry with lint-free cloth. If, after cleaning, the slide does not allow serum to spread evenly within the inner surface of the circle, treat the slide as follows.

- Scrub the slides with a nonabrasive cleanser.
- Rinse, dry and polish with clean lint-free cloth.

Avoid prolonged soaking of ceramic-ringed slides in detergent solution because the ceramic rings will become brittle and flake off.

Test Procedure

Prepare the Antigen Suspension

Check the pH of VDRL Buffered Saline before preparing VDRL antigen suspension. Discard if outside the range of pH 6.0 ± 0.1. Allow VDRL Antigen and VDRL Buffered Saline to reach 23 – 29°C before preparing VDRL antigen suspension.

Use only suspension bottle with flat inner-bottom surfaces that allow the initial VDRL buffered saline to evenly cover the inner-bottom surface of the bottle. If the VDRL buffered saline beads or does not spread evenly to cover the bottom of the bottle, rewash the bottle as described above. (See NOTE above regarding use of 25 mL glass stoppered flasks).

For reproducible results, the VDRL Antigen suspension must be checked daily for proper reactivity with VDRL Test Control Serum Set. Only those VDRL suspensions producing the established reactivity pattern of the control serum should be used.

- Prepare a fresh VDRL antigen suspension each testing day. The temperature of the buffered saline, antigen and equipment should be between 23 – 29°C at the time the antigen suspension is prepared.
- Pipette 0.4 mL of VDRL Buffered Saline to the bottom of a round, 30 mL glass-stoppered bottle with a flat inner-bottom surface. Gently tilt bottle so that VDRL buffered saline will cover the entire inner-bottom surface of the bottle.
- Add 0.5 mL of VDRL Antigen (from the lower half of a 1.0 mL pipette graduated to the tip) directly into the saline while continuously but gently rotating the bottle on a flat surface. Add antigen drop by drop at a rate allowing approximately 6 s for 0.5 mL of antigen. Keep the pipette tip in the upper third of the bottle. Do not splash saline onto the pipette. The proper speed of rotation is obtained when the center of the bottle circumscribes a 5 cm diameter circle approximately 3 times per s.
- Expel the last drop of antigen from the pipette without touching the pipette to the saline and continue rotation of the bottle for 10 s.
- Add 4.1 mL of buffered saline from a 5 mL pipette. Do not drop the saline directly onto the antigen; allow it to flow down the side of the bottle.
- Cap the bottle and shake it from bottom to top and back approximately 30 times in 10 s. The antigen suspension is ready for use and may be used during that day (8 h).
- Mix the VDRL Antigen suspension by gently swirling it each time it is used. Do not mix the suspension by forcing it back and forth through the syringe and needle, since this may cause breakdown of particles and loss of reactivity.

Test Accuracy of Antigen Suspension Needle for Serum Test

- The accuracy of the test depends on the amount of antigen suspension used. Check the calibration needle periodically to ensure delivery of the correct volume of VDRL antigen suspension.
- For the qualitative and quantitative tests on serum, dispense antigen suspension from a syringe fitted with an 18 gauge needle without bevel, which will deliver 60 ± 2 drops for antigen suspension per mL when held vertically.
- Place the needle on a 1 mL syringe. Fill the syringe with VDRL Antigen suspension. Holding the syringe in a vertical position, count the number of drops delivered in 0.5 mL. The needle is correctly calibrated if 30 ± 1 drops are delivered in 0.5 mL.
- Replace the needle if it does not meet this specification. Repeat calibration on the new needle.

VDRL Qualitative Slide Test on Serum

- Slide flocculation tests for syphilis are affected by the room temperature. For reliable and reproducible test results, the VDRL Antigen suspensions, controls and test specimens must be at room temperature, 23 – 29°C, when tests are performed.
- Pipette 50 µL of serum into one ring of paraffin or ceramic-ringed slide using a safety pipetting device. Do not use a glass slide with concavities, wells or glass rings. Spread the serum with a circular motion of the pipette tip so that the serum covers the entire inner surface of the paraffin or ceramic ring.
- Gently resuspend the VDRL Antigen suspension.
- Holding the VDRL Antigen suspension dispensing needle and syringe in a vertical position, dispense several drops to clear the needle of air. Then add exactly 1 free-falling drop (17 µL) of antigen suspension to each circle containing serum. Do not allow the needle to touch the serum.
- Place the slide on the mechanical rotator. Rotate the slide for 4 min at 180 ± 2 rpm. When performing the test in a dry climate, cover the slides with a moist, humidifying cover during rotation to prevent excessive evaporation.
- Immediately after rotating the slide, remove it from the rotator and read the test results microscopically using 100x magnification.

VDRL Quantitative Test on Serum

- To quantitate serum samples to endpoint titer, prepare serum dilutions on the slide at 1:1, 1:2, 1:4, and 1:8.
- Dispense 50 µL of 0.9% saline in circles numbered 2 – 4. Do not spread the saline.
- Dispense 50 µL of serum in circle 1 and 50 µL of serum in circle 2.
- Mix the saline and the serum in circle 2 by drawing the mixture up and down in the pipette approximately 8 times. Avoid forming bubbles.
- Transfer 50 µL from circle 2 (1:2) to circle 3 (1:4) and mix.
- Transfer 50 µL from circle 3 (1:4) to circle 4 (1:8), mix and then discard the last 50 µL.
- Gently resuspend the antigen suspension.
- Holding the antigen suspension dispensing needle and syringe in a vertical position dispense several drops to clear the needle of air. Then add exactly 1 free-falling drop (17 µL) of antigen suspension to each circle.
- Place the slide on the mechanical rotator. Rotate the slide for 4 min at 180 ± 2 rpm. When performing the test in a dry climate, place the slide under a moist, humidifying cover during rotation to prevent excessive evaporation.
- Immediately after rotation, read the test microscopically using 100x magnification.
- If the highest dilution tested (1:8) is reactive:
 - Prepare a 1:8 dilution of the test specimen by adding 0.1 mL of serum to 0.7 mL of 0.9% saline. Mix thoroughly.
 - Place 50 µL of 0.9% saline into the second, third and fourth paraffin rings in a row on the slide. Prepare additional serial dilutions for strongly reactive specimens.
 - Add 50 µL of the 1:8 dilution of the test specimen to paraffin rings 1 and 2.
 - Prepare serial twofold dilutions beginning with ring 2 and complete the test described above.

User Quality Control of VDRL Antigen Suspension

- Prepare a fresh antigen suspension each testing day. Once prepared, it should be used within 8 h.
- Store the prepared antigen suspension at 23 – 29°C.
- Test antigen suspension reactivity with control sera (Reactive, Weakly Reactive and Nonreactive). Test serum dilutions within 1 h after heat inactivation.
- Use the antigen suspension only if it produces the expected reactivity with the control sera (Reactive, Weakly Reactive and Nonreactive) comparable to results obtained with the reference antigen.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Qualitative

- Read and record results as follows:
Medium to large clumps-Reactive (R)
Small clumps-Weakly Reactive (WR)
NO clumping or very slight roughness - Nonreactive (NR)
- Verify that the control sera results are as expected (Reactive, Weakly Reactive, Nonreactive). If reactions are not expected, the test is invalid and results cannot be reported.
- Perform a quantitative test to endpoint on all serum specimens that produce Reactive, Weakly Reactive or "rough" Nonreactive results in the qualitative slide test.

Quantitative

Report the titer as the highest dilution that produces a Reactive (not Weakly Reactive) result. For example:

Undiluted (1:1)	Serum Dilutions					Report
	1:2	1:4	1:8	1:16	1:32	
R	W	N	N	N	N	Reactive, undiluted
R	R	W	N	N	N	Reactive, 1:2 dilution
R	R	R	W	N	N	Reactive, 1:4 dilution
W	W	R	R	W	N	Reactive, 1:8 dilution
N (rough)	W	R	R	R	N	Reactive, 1:16 dilution
W	N	N	N	N	N	Weakly Reactive, undiluted

If reactive results are obtained through dilution 1:32, prepare further twofold serial dilutions in 0.9% saline (1:64, 1:128 and 1:256) and retest using the quantitative test procedure.

Interpretation

- The results of the serum VDRL Test must be confirmed by a treponemal test.
- The diagnosis of syphilis depends upon the results of the VDRL, treponemal confirmatory test, clinical signs and symptoms and risk factors.
- A reactive VDRL Test may indicate past or present infection with a pathogenic treponeme. However, it may be a false positive reaction. A false positive is determined if the confirmatory treponemal test is negative.
- A nonreactive VDRL Test with clinical evidence of syphilis may indicate early primary syphilis, a prozone reaction in secondary syphilis, or late syphilis.
- A nonreactive VDRL Test with no clinical evidence of syphilis indicates no current infection or an effectively treated infection.
- A quantitative VDRL Test detects changes in reagent titer. Therefore, a serum specimen showing a fourfold increase in titer on a repeat specimen may indicate an infection, a reinfection or a treatment failure. Likewise, a fourfold decrease during treatment indicates adequate syphilis therapy.

VDRL Test on Cerebrospinal Fluid

Refer to the appropriate reference for the procedure when testing cerebrospinal fluids with the VDRL Test.¹

LIMITATIONS OF THE PROCEDURE

A prozone reaction may occur in which reactivity with undiluted serum is inhibited. The prozone phenomenon often gives Weakly Reactive or "rough" Nonreactive results in the qualitative test. Therefore, all specimens with these results must be tested in the quantitative test.

Biological false positive reactions can occur with nontreponemal tests in persons who abuse drugs, have diseases such as lupus erythematosus, mononucleosis, malaria, leprosy or viral pneumonia, or have recently been immunized.¹

Cross reactions may occur with other treponemal infections, such as yaws, pinta, bejel, or siti.

In manufacturing, VDRL Antigen with Buffered Saline is tested only with serum. For a modification of the serum test products and procedures for CSF testing, consult the appropriate reference.¹ The user is responsible for modifying the products and procedures and for the required quality control standards.

If the temperature of the testing area, specimens or reagents is less than 23°C, test reactivity is decreased; if the temperature is greater than 29°C, test reactivity is increased.¹

Test results are unpredictable when testing hemolyzed, contaminated or extremely turbid serum specimens. Plasma is not an acceptable specimen.

For correct test results, strictly adhere to correct speed and length of time for rotating the specimens and antigen.

PERFORMANCE CHARACTERISTICS^{15,16}

The performance of the VDRL slide test was compared to two rapid plasma reagin (RPR) card tests and a fluorescent treponemal antibody absorption (FTA-ABS) test for the diagnosis of syphilis in a study by Perryman, Larsen, Hambie, Pettit, Mullally and Whittington.¹⁵ Fresh serum samples obtained from 505 individuals were tested using these four test procedures.

Of the 505 fresh serum samples, 57 sera were reactive in at least one of the nontreponemal tests. The three sera giving borderline FTA-ABS test results were not entered into the tabulation. See Table 1 for distribution of reactive sera among the three nontreponemal tests and the FTA-ABS test.¹⁵

Table 1

Nontreponemal test pattern of reactivity ^a			No. of sera with pattern	No. of sera with following result by FTA-ABS test ^a		
RPR Card Test		VDRL Slide		R		
RPR Test #1	Standard RPR Test #2			R	B	N
R	R	R	43	38	2	3
N	R	R	1	1		
N	N	R	4	2		2
N	R	N	1	1		
R	R	N	5	5		
R	N	N	1	1		
R	N	R	2	1	1	

^aR, Reactive; B, borderline; N, nonreactive

Table 2 shows the diagnoses associated with reactive nontreponemal tests.¹⁵

Table 2

Stage of Syphilis	No. of sera with indicated reaction by following test ^a									
	RPR card test				VDRL slide test			FTA-ABS Test		
	Test #1		Standard Test #2							
	R	N	R	N	R	W	N	R	B	N
Primary										
Untreated	5	1	6		4	2		6		
Treated	7	1	7	1	5	1	2	8		
Secondary										
Untreated	6		6		6			6		
Treated	5		5		5			5		
Latent										
Untreated	9		8	1	8		1	8	1	
Treated	15	1	15	1	10	3	3	14		2
Past history of syphilis										
	1	1	1	1		2		2		
Other sexually transmitted diseases										
	3	2	2	3	3	2			2	3

^aR, Reactive; N, nonreactive; W, weakly reactive (considered reactive in tabulation); B, borderline

The sensitivity of all three nontreponemal tests, based on the 21 untreated cases of syphilis was 95.2% (20 out of 21 results in agreement). With the 31 sera from individuals with treated syphilis, the sensitivity of the VDRL slide test was 83.9% (26 out of 31 results in agreement). The overall sensitivity (treated and untreated) was 88.5% for the VDRL slide test, as compared to 92.3% for both RPR card tests.

The specificities of the three nontreponemal tests based on the 453 sera drawn from individuals presumed not to have syphilis were 99.3% for RPR Test #1, 99.6% for the standard RPR Test #2 and 98.9% for the VDRL slide test.¹⁵

In another study conducted by Pope, Hunter and Feeley¹⁶ at the Centers for Disease Control in Atlanta, Georgia, 297 sera were tested by an enzyme-linked immunosorbent assay (ELISA), VDRL slide test, FTA-ABS test and the microhemagglutination assay for *T. pallidum* antibodies (MHA-TP). Table 3 below lists the sensitivities and specificities of each method published in this study.¹⁶

Table 3

Test	Sensitivity	Specificity
ELISA	89.3%	98.5%
VDRL	93.3%	92.7%
FTA-ABS	100 %	97.8%
MHA-TP	76 %	98.2%

AVAILABILITY

Cat. No.	Description
240764	VDRL Antigen with Buffered Saline; 1 x 5 mL Ampules with 1 x 60 mL Saline
240765	VDRL Antigen with Buffered Saline; 10 x 0.5 mL Ampules with 1 x 60 mL Saline
235201	VDRL Test Control Serum Set, 1 set of 3 vials: Nontreponemal Antigen Reactive Serum, one 3 mL vial VDRL Weakly Reactive Serum, one 3 mL vial Nontreponemal Antigen Nonreactive Serum, one 3 mL vial

REFERENCES

- Larsen, S.A., V. Pope, R. Johnson, and E.J. Kennedy, Jr. (ed.). 1998. A manual of tests for syphilis, 9th ed. American Public Health Association, Washington, D.C.
- Creighton, E.T. 1990. Darkfield microscopy for the detection and identification of *Treponema pallidum*, p. 49-61. In S. A. Larsen, E. F. Hunter, and S. J. Kraus (ed.), Manual of tests for syphilis, 8th ed. American Public Health Association, Washington, D.C.
- Janda, W.M. (ed.). 1994. Immunology, p. 9.7.1-9.7.20. In H. D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 2. American Society for Microbiology, Washington, D.C.
- Norris, S.J., V. Pope, R.E. Johnson and S.A. Larsen. 2003. *Treponema* and other human host-associated spirochetes, p. 955-971. In P.R. Murray, E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenover (ed.), Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
- Matthews, H.M., T.K. Yang, and H.M. Jenkin. 1979. Unique lipid composition of *Treponema pallidum* (Nichols virulent strain). Infect. Immun. 24:713-719.
- Harris, A., A.A. Rosenberg, and L.M. Riedel. 1946. A microflocculation test for syphilis using cardioliipin antigen. J. Ven. Dis. Infor 27:169-174.
- Pangborn, M.C. 1941. A new serologically active phospholipid from beef heart. Proc. Soc. Exp. Biol. Med. 48:484-486.
- Pangborn, M.C. 1944. Acid cardioliipin and an improved method for the preparation of cardioliipin from beef heart. J. Biol. Chem. 153:343-348.
- Pangborn, M.C. 1945. A simplified preparation of cardioliipin, with a note on purification of lecithin for serologic use. J. Biol. Chem. 161:71-82.
- National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, PA.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80
- U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC) 4th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
- Thomson, R.B. and J.M. Miller. 2003. Specimen collection, transport and processing: bacteriology, p. 286-330. In Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller and R.H. Tenover (ed.), Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
- Perryman, M.W., S.A. Larsen, E.A. Hambie, D.E. Pettit, R.L. Mullally and W. Whittington. 1982. Evaluation of a new rapid plasma reagin card test as a screening test for syphilis. J. Clin. Microbiol. 16:286-290.
- Pope, V., E.F. Hunter and J.C. Feeley. 1982. Evaluation of the microenzyme-linked immunosorbent assay with *Treponema pallidum* antigen. J. Clin. Microbiol. 15:630-634.

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BENEX Limited