
**Revisions**


Rev from	Rev to	ECO #
0905	2010/06	5370-10

**Notes:**


- BD Cat. Number 231749
- Blank (Sheet) Size : Length: 12"      Width: 19.5"  
 Number of Pages: 12      Number of Sheets: 1  
 Page Size: Length 12"      Width 3.25"      Final Folded Size: N/A
- Style (see illustrations below): # 4


#1


#2


#3


#4

#5

#6

#7

- See Specification Control Number N/A for Material Information
- Ink Colors: Printed two sides  Yes       No  
 No. of Colors: 1      PMS# Black
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

VS Controlled by BD Caribe, LTD.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: L000147		Category and Description Package Insert, DrySlide Nitrocefim	Sheet: 1 of 13 Scale: N/A	<b>A</b>

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## INTENDED USE

**BBL™ DrySlide™ Nitrocefin** is used for detecting  $\beta$ -lactamase production by bacteria.

## SUMMARY AND EXPLANATION

The chromogenic cephalosporin beta-lactamase test is a qualitative procedure for detecting the production of beta-lactamase by bacteria. The reaction is based on the production of a colored compound when the substrate, nitrocefin, is exposed to a  $\beta$ -lactamase-producing culture.<sup>1</sup>

The enzyme, beta-lactamase, originally described by Abraham and Chain,<sup>2</sup> is produced by various organisms and is a mechanism of their resistance to penicillins and cephalosporins.<sup>3</sup> Test methods used to detect  $\beta$ -lactamase include iodometric, acidimetric and chromogenic cephalosporin procedures.<sup>4-7</sup>

Some staphylococci may require induction (exposure to a  $\beta$ -lactam agent) to increase production of  $\beta$ -lactamase to detectable levels. According to standard references, any negative  $\beta$ -lactamase test result from a non-induced staphylococcal isolate, despite the test method, should be confirmed by induction of the isolate and repeat testing.<sup>8</sup>

**BBL DrySlide Nitrocefin** employs nitrocefin, a cephalosporin compound first described by Glaxo Research (Middlesex, England), in the chromogenic cephalosporin test methodology.

## PRINCIPLES OF THE PROCEDURE

If beta-lactamase is produced by a culture, it hydrolyzes the  $\beta$ -lactam ring of nitrocefin, producing cephalosporanic acid.<sup>9</sup> A distinctive color change is associated with this reaction. The pale yellow nitrocefin turns pink after hydrolysis.<sup>1,4-7</sup>

Both aerobic and anaerobic  $\beta$ -lactamase-producing bacteria produce this color change.<sup>10,11</sup> Organisms that are  $\beta$ -lactamase negative do not alter the pale yellow color of nitrocefin within the time limits of the test.

## REAGENTS

**BBL DrySlide Nitrocefin** is a disposable slide having four filter paper reaction areas containing nitrocefin.

### Warnings and Precautions:

For *in vitro* Diagnostic Use.

Handle **BBL DrySlide Nitrocefin** only by the slide mount; do not touch the reaction area.

Follow proper established laboratory procedures in handling and disposing of infectious materials.

**Storage Instructions:** Store **BBL DrySlide Nitrocefin** at 15 – 30°C.

The expiration date applies to the product in its intact container when stored as directed.

**Product Deterioration:** Do not use a product if it fails to meet specifications for identity and performance.

## SPECIMEN COLLECTION AND PREPARATION

Collect specimens or samples in sterile containers or with sterile swabs and transport immediately to the laboratory according to recommended guidelines.<sup>5-7</sup>

Process each specimen using procedures appropriate for that sample.<sup>5-7</sup>

Obtain test organism from a primary isolation culture streaked on an appropriate solid medium (chocolate agar or tryptic soy blood agar). Select well-isolated colonies. Mixed cultures require subculturing for purification. Insufficient growth indicates the need for further incubation.

Pick isolated colonies or a sweep of confluent growth from the pure culture.

## PROCEDURE

**Material Provided:** **BBL DrySlide Nitrocefin.**

**Materials Required But Not Provided:** Inoculation loop (platinum or plastic) or wooden applicator stick, Pasteur pipette and bulb, distilled or deionized water and quality control organisms.

### Test Procedure:

1. Open the **BBL DrySlide Nitrocef**in pouch and remove the slide. Do not place used or partially used **BBL DrySlide Nitrocef**in (with organisms) back into the pouch.
2. Using a Pasteur pipette, dispense a small drop of distilled or deionized water onto a reaction area of **BBL DrySlide Nitrocef**in. Avoid adding excess water.
3. Using an inoculating loop, pick isolated colonies or a sweep of confluent growth from the culture to be tested.
4. Generously smear the specimen onto the moistened reaction area of **BBL DrySlide Nitrocef**in.
5. Examine the reaction area for appearance of a color change from yellow to pink in 5 to 60 min.

### User Quality Control:

**Identity Specifications** – Each pouch contains 1 slide mount and 1 desiccant container. Each slide mount has 4 yellow reaction areas.

**Performance Response** – Perform  $\beta$ -lactamase test using **BBL DrySlide Nitrocef**in following above procedure with the test organisms listed below.

Organism	ATCC™	$\beta$ -Lactamase Reaction	Color Response
<i>Haemophilus influenzae</i>	35056	+	Pink
<i>Staphylococcus aureus</i>	25923	–	Yellow

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

### RESULTS

Beta-lactamase-positive organisms change the color of the reaction area from yellow to pink. A positive result will develop within 5 min for most bacterial strains. However, positive reactions may take up to 60 min to develop for some staphylococci.

Beta-lactamase-negative organisms do not change the color of the reaction area.

### LIMITATIONS OF THE PROCEDURE

**BBL DrySlide Nitrocef**in is intended as a rapid test procedure and should not replace conventional susceptibility testing. For additional information on susceptibility testing, consult appropriate references.<sup>5-7,10,11</sup>

Testing for resistance to penicillin G and ampicillin using the chromogenic cephalosporin method has been predictive of resistance for *Prevotella*, *Porphyromonas*, *Bacteroides* and other anaerobes.<sup>9,10</sup> Under selected situations, the chromogenic cephalosporin test for beta-lactamase does provide clinically relevant information quickly. However, resistance to  $\beta$ -lactam drugs is not always mediated by  $\beta$ -lactamase production as with *Bilophila wadsworthia*, *Bacteroides gracilis* and some strains of *B. distasonis* and *B. fragilis*.<sup>10</sup>

Chromogenic cephalosporin-based tests provide reliable results with *Neisseria gonorrhoeae*, *Haemophilus* spp., *Branhamella catarrhalis*, enterococci and staphylococci. More specifically, resistance to penicillin, ampicillin and amoxicillin is detected for *N. gonorrhoeae*, *Haemophilus* spp. and *B. catarrhalis*. Resistance to penicillin, ampicillin, amoxicillin, carbenicillin, mezlocillin, piperacillin and ticarcillin is detected for staphylococci and enterococci.<sup>11-15</sup>

$\beta$ -lactamase negative,  $\beta$ -lactam antibiotic resistant organisms have been reported. An organism's resistance due to these other mechanisms will not be detected with **BBL DrySlide Nitrocef**in.

### PERFORMANCE CHARACTERISTICS

Prior to release, all lots of **BBL DrySlide Nitrocef**in slides are tested for performance characteristics. Representative samples of the lot are inoculated in the reaction area of the slide (after it is wetted with a drop of purified water) with a loopful of 18 to 24 cultures of *Haemophilus influenzae* (ATCC 35056), *Staphylococcus aureus* (VSL 7969) and *S. aureus* (ATCC 25923). After inoculation, the slides are examined from 5 to 60 min for a color change from yellow to pink. *H. influenzae* is positive (pink) within 5 min. *S. aureus* (VSL 7969) is positive (pink) within 60 min. *S. aureus* (ATCC 25923) is negative (yellow) after 60 min.

### AVAILABILITY

#### Cat. No. Description

231749 **BBL™ Dryslide™ Nitrocef**in, 25 slides.

### REFERENCES

1. O'Callaghan, C.H., A. Morris, S.M. Kirby, and A.H. Shingler. 1972. Novel method for detection of  $\beta$ -lactamase by using a chromogenic cephalosporin substrate. *Antimicrob. Agents Chemother.* 1 (4): 283-288.
2. Abraham, E.P. and E. Chain. 1940. An enzyme from bacteria able to destroy penicillin. *Nature* 146: 837.
3. Bene, V.E., and W.E. Farrar, Jr. 1973. Cephalosporinase activity in *Bacteroides fragilis*. *Antimicrob. Agents Chemother.* 3 (3): 369-372.
4. Montgomery, K., L. Raymundo, Jr., and W.L. Drew. 1979. Chromogenic cephalosporin spot test to detect beta-lactamase in clinically significant bacteria. *J. Clin. Microbiol.* 9 2: 205-207.
5. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.). 1995. *Manual of clinical microbiology*, 6th ed. American Society for Microbiology, Washington, D.C.
6. Pezzlo, M. (ed.). 1994. Aerobic bacteriology, p. 1.0.1.-1.20.47. In H.D. Isenberg (ed.), *Clinical microbiology procedures handbook*, vol. 1. American Society for Microbiology, Washington, D.C.
7. Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. *Bailey & Scott's diagnostic microbiology*, 9th ed. Mosby-Year Book, Inc., St. Louis, MO.

8. Leitch, C., and S. Boonlayangoor. 1994.  $\beta$ -lactamase tests, p. 5.3.1 In H.D. Isenberg (ed.), *Clinical microbiology procedures handbook*, vol. 2. American Society for Microbiology, Washington, D.C.
9. Lee, D.T., and J.E. Rosenblatt. 1983. A comparison of four methods of detecting  $\beta$ -lactamase in anaerobic bacteria. *Diagn. Microbiol. Infect. Dis.* 1: 173-175.
10. National Committee for Clinical Laboratory Standards. 1997. Approved standard M11-A4. Methods for antimicrobial susceptibility testing of anaerobic bacteria, 4th ed. National Committee for Clinical Laboratory Standards, Wayne, Pa.
11. National Committee for Clinical Laboratory Standards. 1997. Approved standard M2-A6. Performance standards for antimicrobial disk susceptibility tests, 6th ed. National Committee for Clinical Laboratory Standards, Wayne, Pa.
12. Murray, B.E., B. Mederski-Samoraj, S.K. Foster, J.L. Brunton, and P. Harford. 1986. In vitro studies of plasmid-mediated penicillinase from *Streptococcus faecalis* suggest a staphylococcal origin. *J. Clin. Invest.* 77: 289-293.
13. Murray, B.E., D.A. Church, A. Wanger, K. Zscheck, M.E. Levison, M.J. Ingerman, E. Abrutyn, and B. Mederski-Samoraj. 1986. Comparison of two  $\beta$ -lactamase-producing strains of *Streptococcus faecalis*. *Antimicrob. Agents Chemother.* 30: 861-864.
14. Bush, L.M., J. Calmon, C.L. Cherney, M. Wendler, P. Pitakis, J. Poupard, M.E. Levison, and C. Johnson. 1989. High-level penicillin resistance among isolates of enterococci, *Ann. Intern. Med.* 110: 515-520.
15. Patterson, J.E., B.L. Masecar, and M.J. Zervos. 1988. Characterization and comparison of two penicillinase-producing strains of *Streptococcus (Enterococcus) faecalis*. *Antimicrob. Agents and Chemother.* 32: 122-124.

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Company. © 2010 BD.
ATCC is a trademark of the American Type Culture Collection.
DrySlide is a trademark of Difco Laboratories, a subsidiary of
BD, BD Logo and BBL are trademarks of Becton, Dickinson and
Benex Limited