Revisions

Rev from	Rev to	ECO #
0703	2010/06	5368-10

Notes:

1. BD Cat. Number 231725
2. Blank (Sheet) Size: Length: 13.75" Width: 16.125"
 Number of Pages: 6 Number of Sheets: 1
 Page Size: Length 13.75" Width 5.375" Final Folded Size: 1.25" x 5.375"
3. Style (see illustrations below): # 4

4. See Specification Control Number N/A for Material Information
5. Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# Black
6. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

VS Controlled by BD Caribe, Ltd.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: L000151		Category and Description Package Insert, Taxo ALA	Sheet: 1 of 7	A
			Scale: N/A	

BD BBL™ Taxo™ Differentiation Discs ALA

English: pages 1 Italiano: pagina 4
Français : page 2 Español: página 5
Deutsch: Seite 3

 L000151
2010/06

Pokyny vám poskytnú miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіндік жөргілікті BD өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

BBL™ Taxo™ Differentiation Discs ALA (delta-aminolevulinic acid) are used for rapidly detecting porphyrin as a means of speciating *Haemophilus* species.

SUMMARY AND EXPLANATION

The differentiation of *Haemophilus influenzae* from *H. parainfluenzae* and other members of this genus using delta-aminolevulinic acid was first described by White and Granick.¹ Enzymes that convert ALA to porphyrins in the biosynthesis of heme (X factor) are not produced by *H. influenzae*, *H. aegyptius* and *H. canis*. These enzymes are produced by *H. parainfluenzae*, *H. parahaemolyticus*, *H. gallinarum*, *H. parasuis*, *H. parahaemolyticus* and *H. aphrophilus*.² The porphyrin test can be used to determine an organism's X-factor (hemin) requirement while avoiding the problem of X-factor carryover from primary culture media.^{3,4}

PRINCIPLES OF THE PROCEDURE

Heme-producing strains of *Haemophilus* possess the enzyme porphobilinogen synthase. This enzyme converts delta-aminolevulinic acid (ALA) into porphyrins. Porphyrins are detected by a reddish-orange fluorescence under UV light (360 nm).

REAGENTS

BBL Taxo Differentiation Discs ALA are paper discs that contain buffered delta-aminolevulinic acid.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Follow established laboratory procedures in handling and disposing of infectious materials.

Storage Instructions: Store BBL Taxo Differentiation Discs ALA at 2 – 8°C in the dark. Allow about 10 – 15 min for the container to reach room temperature before opening. Promptly replace the package at 2 – 8°C after use.

The expiration date applies to the product in its intact container when stored as directed.

Product Deterioration: Do not use a product if it fails to meet specifications for identity and performance.

SPECIMEN COLLECTION AND PREPARATION

Collect specimens or samples in sterile containers or with sterile swabs and transport immediately to the laboratory in accordance with recommended guidelines.⁵⁻⁷

Process each specimen, using procedures appropriate for that sample.⁵⁻⁷

PROCEDURE

Material Provided: BBL Taxo Differentiation Discs ALA.

Materials Required But Not Provided: Glass slide 1" x 3", Petri dishes, inoculating loop, Bunsen burner or incinerator, forceps, absorbent cotton, incubator 35 ± 2°C, ultraviolet light (360 nm), filter paper and quality control organisms.

Test Procedure:

1. Obtain a pure culture of the organism to be tested. Select well-isolated colonies.
2. Place 3 BBL Taxo ALA Discs equally separated from each other on a 1" x 3" glass slide or on the inside of a Petri dish. If discs are on a slide, place the slide in a Petri dish.
3. Add .04 mL sterile distilled water to each disc.
4. Rub a loopful of the test organism into one moistened disc, holding it in place with flamed and cooled forceps. Repeat for the positive control and negative control cultures, using the other moistened discs.
5. Place a piece of filter paper or cotton saturated with water in the Petri dish as far away from the discs as possible.
6. Incubate the Petri dish in the dark at 35 ± 2°C.
7. Examine the discs at hourly intervals for 5 h. Remove the top of the Petri dish and expose the discs to UV (360 nm) light in a darkened room.

User Quality Control:

Identity Specifications – White 1/4" paper discs. The discs are imprinted with ALA.

Cultural Response – Inoculate ALA discs with a heavy loopful of test organisms following label instructions. Incubate at 35 ± 2°C for 5 h. Read reactions by exposing the discs to long wave UV light (360 nm) in a dark room. A positive reaction is indicated by an orange-red fluorescence.

Organism	ATCC™	Porphyrin Production
<i>Haemophilus influenzae</i>	35056	–
<i>Haemophilus parainfluenzae</i>	7901	+

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

A positive reaction is an orange-red fluorescence, which indicates the presence of porphyrin. A negative reaction will show no fluorescence.

LIMITATIONS OF THE PROCEDURE

Use this procedure for differentiating *Haemophilus* species only.

Best results are obtained when a heavy inoculum of an 18 – 24 h culture is used.

Delta-aminolevulinic acid (ALA) is light sensitive. Protect the discs from light.

PERFORMANCE CHARACTERISTICS

Prior to release, all lots of BBL Taxo Differentiation Discs ALA are tested for performance characteristics. Representative samples of the lot are each placed on a glass slide in a Petri dish and processed water is added to the disc. A loopful of *Haemophilus influenzae* (ATCC 35056), *H. parahaemolyticus* (ATCC 10014) or *H. parainfluenzae* (ATCC 7901) is rubbed onto each moistened disc. A piece of cotton saturated with water is placed inside the Petri dish as far away from the discs as possible. The Petri dishes are covered and incubated in the dark at 35 ± 2°C. The discs are examined at hourly intervals for up to 5 h by removing the Petri dish lid and exposing the discs to UV (360 nm) light in a darkened room. *H. parahaemolyticus* and *H. parainfluenzae* give a positive reaction (orange-red fluorescence) after 2 h incubation. *H. influenzae* gives a negative reaction (no fluorescence) after 5 h incubation.

AVAILABILITY

Cat. No. Description

231725 BBL™ Taxo™ Differentiation Discs ALA, 50 discs.

REFERENCES

1. White, D.C., and S. Granick. 1963. Hemin biosynthesis in *Haemophilus*. J. Bacteriol. 85: 842-850.
2. Biberstein, E.L., P.D. Mini, and M.G. Gills. 1963. Action of *Haemophilus* cultures on delta-aminolevulinic acid. J. Bacteriol. 86: 835-842.
3. Kilian, M. 1974. A rapid method for the differentiation of *Haemophilus* strains, the porphyrin test. Acta. Path. Microbiol. Scand. Sect. B. 82: 835-842.
4. Kilian, M. 1976. A taxonomic study of the genus *Haemophilus*, with the proposal of a new species. J. Gen. Microbiol. 93: 9-62.
5. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.). 1995. Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
6. Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. Bailey & Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis.
7. Pezzlo, M. (ed.). 1994. Aerobic bacteriology, p. 1.0.1-1.20.47. In H.D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.

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Benex Limited