
L-000159  
Revised: July 1999

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#### INTENDED USE

CO<sub>2</sub> Indicators are intended to be used in conjunction with carbon dioxide generating systems used for the growth of bacterial cultures. The CO<sub>2</sub> Indicator signify that the carbon dioxide concentration in the atmosphere is 3% or greater in jar or pouch incubation systems.

#### SUMMARY AND EXPLANATION

The stimulatory effect of carbon dioxide on the growth of aerobic and microaerophilic bacteria has been studied extensively.<sup>1-7</sup> An aerobic environment enriched with up to 10% carbon dioxide is essential for optimal growth of many *Haemophilus* spp. and *Neisseria* spp., as well as *Mycobacterium* spp. and *Streptococcus pneumoniae*.<sup>7-9</sup> A microaerophilic environment enriched with approximately 5 – 10% carbon dioxide is essential for the growth of *Helicobacter* spp. and *Campylobacter* spp.<sup>8,9</sup>

The disposable CO<sub>2</sub> Indicator provide a means for monitoring the production and maintenance of the enriched carbon dioxide atmosphere required for enhanced growth of capnophilic and microaerophilic bacteria. These indicators can be used in the following Becton Dickinson Microbiology Systems products:

**GasPak 100™** system or **GasPak 150™** system operated with **GasPak™** Disposable Carbon Dioxide Generator Envelopes or **CampyPak™** or **CampyPak Plus™** Envelopes or **GasPak CO<sub>2</sub> Pouch™** or **Campy Pouch™**; or **Bio-Bag™** Type Cj or **Bio-Bag™** Type C.

#### PRINCIPLES OF THE PROCEDURE

In an atmosphere enriched with 3% or greater carbon dioxide, a phenol red indicator solution in the activated device will respond to the drop of pH caused by absorbed carbon dioxide and change from red to yellow in color within 30 min. This indicates proper functioning of the carbon dioxide generating system.

At carbon dioxide concentrations of less than 3%, the indicator solution will remain pink to orange in color, indicating improper activation of the carbon dioxide generating system or a leak in the incubation system used.

The indicator solution reacts reversibly upon re-exposure to a non-CO<sub>2</sub> enriched atmosphere and changes from yellow to a pink or red color within 30 min.

Maintenance of the yellow color, during the culture incubation period, indicates that the incubation system integrity has been maintained.

#### REAGENTS

CO<sub>2</sub> Indicators contain 0.5 mL of buffered, aqueous 0.002% Phenol Red solution.

**Precautions:** For Laboratory Use.

Care should be taken when crushing ampules so that over manipulation does not lead to injury.

**Storage Instructions:** Store at controlled room temperature 15 – 30°C (59 – 86°F). Use indicator once and then discard.

**Product Deterioration:** This reagent is hermetically sealed in an ampule which affords protection of the solution from chemical instability until the expiration date. Do not use indicators if they show evidence of being punctured or crushed. Do not use indicators that are yellow before activation. Do not use after the expiration date.

#### PROCEDURE

**Material Provided:** Single use, CO<sub>2</sub> Indicators .

**Materials Required But Not Provided:** Ancillary culture media and laboratory equipment as required for the procedure.

#### Directions for Use:

1. Remove CO<sub>2</sub> Indicator from box; ensure that the label is on the clear tube creating a non-rolling device.
2. Grasp the middle of the clear tube with thumb and forefinger and squeeze gently to break ampule inside the tube. **Caution: Break ampule close to its center one time only. Do not manipulate indicator any further, as the plastic may puncture and injury may occur.**
3. Shake well to ensure glass dispersion and proper activation.
4. Place CO<sub>2</sub> Indicators on top of the stack of plates within a jar or next to plates in a pouch system. *Do not place CO<sub>2</sub> Indicators vertically along side a stack of plates in a jar.*
5. Activate CO<sub>2</sub> generator system according to manufacturer's instructions.
6. Observe indicator for change from red to yellow color within 30 min as an indicator of the carbon dioxide concentration in the atmosphere is 3% or greater.

#### User Quality Control:

Examine the indicator for signs of deterioration (see "Product Deterioration").

#### RESULTS

Depending upon the number of inoculated plates contained in the culture system, incubation temperature and gas generating system used, the indicator will change from red to yellow within 30 min of sealing the system, indicating an atmosphere enriched with 3% or greater of carbon dioxide.

#### LIMITATIONS OF THE PROCEDURE

Use indicators once and discard.

The carbon dioxide concentration in candle jars may be less than 3%. The color change of the indicator solution will not be complete at less than 3%. Above 1% and below 3% the indicator solution will be melon to orange in color.

Activated CO<sub>2</sub> Indicators left in contact with paper or other absorbent material will drain off liquid over time and any resulting color on paper or within solution should be ignored.

**BECTON  
DICKINSON**

SIZE: 4.5" W x 11" L  
COLOR: Standard Black

Original by: Translingua

ENG – 6.5 pt on 7.5 pt leading  
MLL – 6 pt on 7 pt leading

Edited by: Translingua  
Date: 8/25/98  
Translate & Typeset MLL

Indicators placed vertically in jars will convert to yellow more slowly than correctly placed indicators and the color change may not be complete.

Do not use CO<sub>2</sub> Indicators in Gono-Pak or JEMBEC® transport systems.

CO<sub>2</sub> Indicators demonstrate the presence or absence of CO<sub>2</sub> only. No color change occurs in response to changes in O<sub>2</sub> levels.

Do not use CO<sub>2</sub> Indicators in CO<sub>2</sub> incubators.

#### AVAILABILITY

##### Cat. No. Description

271055 CO<sub>2</sub> Indicators, packaged 50 droppers/carton.

#### REFERENCES

1. Rockwell, B.E. 1923. The influence of carbon dioxide on the growth of bacteria. *J. Infect. Dis.* 32: 98-104.
2. Valley, G., and L.F. Rettger. 1926. Carbon dioxide requirements of bacteria. *J. Bacteriol.* 11: 78-79.
3. Valley, G., and L.F. Rettger. 1927. The influence of carbon dioxide on bacteria. *J. Bacteriol.* 14: 110-137.
4. Rockwell, B.E., and J.H. Highberger. 1927. The necessity of carbon dioxide for the growth of bacteria, yeast and molds. *J. Infect. Dis.* 40: 438-445.
5. Valley, G. 1928. The effect of carbon dioxide on bacteria. *Quart. Rev. Biol.* 3: 209-234.
6. Gladstone, G.P., P. Fildes, and B.M. Richardson. 1935. Carbon dioxide as an essential factor in the growth of bacteria. *Br. J. Exp. Pathol.* 16: 335-348.
7. Kempner, W., and C. Schlayer. 1942. Effect of CO<sub>2</sub> on the growth rate of the pneumococcus. *J. Bacteriol.* 43: 387-396.
8. Ballows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy (ed.), 1991. Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.
9. Baron, E.J., and S.M. Finegold. 1990. Bailey & Scott's diagnostic microbiology, 8th ed. The C.V. Mosby Company, St. Louis.

**TECHNICAL INFORMATION:** In the United States, telephone Technical Services, toll free (800) 638-8663.

©1999 Becton Dickinson and Company.  
Bio-Bag, CampyPak, CampyPak Plus, Campy Pouch,  
GasPak, GasPak 100, GasPak 150, and GasPak CO<sub>2</sub> Pouch  
are trademarks of Becton Dickinson and Company.  
JEMBEC is a trademark of Miles Scientific.

**BECTON  
DICKINSON**

**Becton Dickinson Microbiology Systems**  
Becton Dickinson and Company  
7 Loveton Circle  
Sparks, Maryland 21152 USA

**Becton Dickinson France S.A.**  
11 rue Aristide Bergès  
38800 Le Pont de Claix, France

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