
**Revisions**


Rev from	Rev to	ECO #
0306	2010/06	5347-10

**Notes:**


- BD Cat. Number 261192, 261193
- Blank (Sheet) Size: Length: 11"      Width: 18"  
 Number of Pages: 8      Number of Sheets: 1  
 Page Size: Length 11"      Width 4.5"      Final Folded Size: 4.5" x 1 7/8"
- Style (see illustrations below): # 4


#1


#2


#3


#4

#5

#6

#7

- See Specification Control Number D11826RM for Material Information
- Ink Colors: Printed two sides  Yes     No  
 No. of Colors: 1      PMS# 2755 Blue
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: L001126		Category and Description Package Insert, Voges-Proskauer Reagent Droppers	Sheet: 1 of 9 Scale: N/A	<b>A</b>

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## INTENDED USE

BBL™ Voges-Proskauer Reagent Droppers are intended for use in the Voges-Proskauer test, one of several qualitative tests used to distinguish between members of the *Enterobacteriaceae*.<sup>1</sup>

## SUMMARY AND EXPLANATION

The Voges-Proskauer test is one of the series known as the IMViC tests. The letters stand for Indole, Methyl Red, Voges-Proskauer, and Citrate. (The “i” is inserted for euphony.) They are used primarily to distinguish between the coliforms, but may be used for other organisms in the *Enterobacteriaceae* family. Specifically, the Voges-Proskauer test detects the presence of acetylmethylcarbinol. This compound is produced by certain microorganisms during growth in a buffered peptone-glucose broth, specifically MR-VP broth.<sup>1</sup>

## PRINCIPLES OF THE PROCEDURE

Upon addition of Voges-Proskauer reagents A and B to a culture tube of MR-VP broth, any acetylmethylcarbinol present will be oxidized to diacetylmethylcarbinol. The latter compound reacts with creatine (from the broth) to form a red-colored compound. The presence of this compound is detected visually and is considered a positive test.<sup>1</sup>

## REAGENTS

BBL Voges-Proskauer Reagent A Droppers contain 0.5 mL of 5% wt/vol alpha-naphthol in absolute alcohol.

BBL Voges-Proskauer Reagent B Droppers contain 0.5 mL of 40% wt/vol potassium hydroxide in distilled water.

## Warnings and Precautions:

For *in vitro* Diagnostic Use

Voges-Proskauer Reagent B is caustic. Avoid contact with the skin. Rinse thoroughly with water if spilled.

Follow established laboratory procedures in handling and disposing of infectious materials.

**Storage Instructions:** Store at room temperature 15 – 30°C (59 – 86°F).

**Product Deterioration:** Reagents are hermetically sealed in an ampule, which affords protection of the solution from chemical instability until expiration date. The reagent droppers need no refrigeration. Opened reagent droppers should be disposed of within 24 h.

## PROCEDURE

**Material Provided:** Voges-Proskauer Reagent Droppers.

**Materials Required But Not Provided:** Ancillary culture media, quality control organisms and laboratory equipment as required for these procedures.

### Test Procedure

1. Hold dropper upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to break ampule inside the dropper. **Caution: Break ampule close to its center *one time only*. Do not manipulate the dropper any further as the plastic may puncture and injury may occur.** Tap bottom on tabletop a few times.
2. Invert for drop-by-drop dispensing according to the test procedure performed.

#### A. Standard Method (Barritt's Modification)<sup>1</sup>

Inoculate buffered peptone-glucose broth (MR-VP medium) with a pure culture of the organism in question. Incubate the culture at 37°C for no less than 48 h.

Empty contents (15 drops) from reagent A dropper and 5 drops from reagent B dropper into 1 mL of broth culture. Shake well after the addition of each reagent to aerate sample.

#### B. Other Methods

Voges-Proskauer Reagent Droppers can be used to achieve highly accurate tests with the commercially prepared bacterial identification systems. Follow the manufacturer's instructions for use.

## User Quality Control

Using either of the methods described above (Test Procedure), an authentic culture of *Enterobacter cloacae* or *Serratia marcescens* will yield a positive result, while *Escherichia coli* or *Salmonella enteritidis* will yield a negative result.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

## RESULTS

A positive reaction is indicated by the development of a distinct red color which occurs within 5 min. The development of a copper color in some tests should be disregarded.

## LIMITATIONS OF THE PROCEDURE

One should allow at least 15 min for color to develop before considering the Voges-Proskauer test negative.

## PERFORMANCE CHARACTERISTICS

Barritt<sup>2</sup> evaluated two Voges-Proskauer (VP) formulations: one containing  $\alpha$ -naphthol and another containing creatine. Three hundred sixty eight (368) strains were tested by both formulations. Two hundred fifty nine (259) strains were negative by both methods and were identified as *E. coli*. Thirty one (31) strains were positive by both methods and were identified as *E. aerogenes*. Seven (7) strains gave a negative reaction with the VP reagent containing creatine and a positive reaction with the VP reagent containing  $\alpha$ -naphthol. Reactions of these strains on EMB Agar along with methyl red and citrate tests identified all of these strains as *E. aerogenes*, indicating that the formulation with  $\alpha$ -naphthol was more sensitive. The remaining 71 strains known as the "intermediate group," showed the greatest discrepancy between the VP reagent containing creatine and the VP reagent containing  $\alpha$ -naphthol. All 71 of these strains were methyl red and citrate positive. Twenty-one (21) of the 71 strains (29%) gave a positive result with the VP reagent containing creatine while 54 strains (76%) gave a positive result with the VP reagent containing  $\alpha$ -naphthol. This result enhanced the observation that the VP reagent containing  $\alpha$ -naphthol was more sensitive than the VP reagent containing creatine and increased the correlation between acetylmethylcarbinol production and citrate utilization, which further helped to distinguish *E. coli* from *E. aerogenes*.

## AVAILABILITY

### Cat. No. Description

261192 **BBL™** Voges-Proskauer Reagent A Droppers, packaged 50 droppers/carton.  
261193 **BBL™** Voges-Proskauer Reagent B Droppers, packaged 50 droppers/carton.

## REFERENCES

1. Lennette, E.H., A. Balows, W.J. Hausler, Jr., and H.J. Shadomy (ed.) 1985. Manual of clinical microbiology, 4th ed. American Society for Microbiology, Washington, D.C.
2. Barritt, M.M. 1936. The intensification of the Voges-Proskauer reaction by the addition of  $\alpha$ -naphthol. J. Pathol. Bacteriol. 42:441-454.

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Benex Limited