



QUALITY CONTROL PROCEDURES

I INTRODUCTION

MacConkey II Agar is a selective and differential medium for the detection of coliform organisms and enteric pathogens.

II PERFORMANCE TEST PROCEDURE

- Inoculate representative samples with dilutions of the cultures listed below.
 - Streak the plates for isolation. Use an 18- to 24-h broth culture of *Enterococcus faecalis* diluted to 10⁴–10⁵ CFU/plate. For the remaining organisms, use an 18- to 24-h broth culture diluted to yield 10³– 10⁴ CFU/plate.
 - Incubate plates at 35 ± 2°C in an aerobic atmosphere.
 - Include **Trypticase™** Soy Agar with 5% Sheep Blood (TSA II) plates as nonselective controls for all organisms.
- Examine plates after 18–24 h for growth, colony size, pigmentation and selectivity.
- Expected Results

CLSI Organisms	ATCC™	Recovery	Colony Color
* <i>Escherichia coli</i>	25922	Growth	Pink
* <i>Proteus mirabilis</i>	12453	Growth Inhibition of swarming (partial)	Colorless
* <i>Salmonella choleraesuis</i> subsp. <i>choleraesuis</i> serotype Typhimurium	14028	Growth	Colorless
* <i>Enterococcus faecalis</i>	29212	Inhibition (partial)	May be pink
Additional Organisms			
<i>Pseudomonas aeruginosa</i>	10145	Growth	Pink to green
<i>Shigella dysenteriae</i>	9361	Growth	Colorless to pink

*Recommended organism strain for User Quality Control.

III ADDITIONAL QUALITY CONTROL

- Examine plates as described under "Product Deterioration."
- Visually examine representative plates to assure that any existing physical defects will not interfere with use.
- Determine the pH potentiometrically at room temperature for adherence to the specification of 7.1 ± 0.2.
- Note the firmness of plates during the inoculation procedure.
- Incubate uninoculated representative plates aerobically at 35 ± 2°C for 72 h and examine for microbial contamination.

PRODUCT INFORMATION

IV INTENDED USE

MacConkey II Agar is a selective and differential medium for the detection of coliform organisms and enteric pathogens.

V SUMMARY AND EXPLANATION

At the present time, many culture media are available to the laboratorian for the isolation, cultivation and identification of enteric bacteria. One of the earliest of these was developed by MacConkey and first described as a brief published note.¹ The landmark paper on MacConkey Agar was published in 1905 and contained detailed descriptions of the medium and the bacterial growth patterns obtained.² This formulation was devised in the knowledge that bile salts are precipitated by acids and certain enteric microorganisms ferment lactose whereas others do not possess this ability.

Since the publication of the early papers, the MacConkey Agar formula has been modified many times. A compilation of culture media published in 1930 lists ten modifications which were published up to that time.³ More recent modifications include use of additives (e.g., kanamycin) and the deletion of certain ingredients (e.g., crystal violet, and neutral red⁴).

MacConkey Agar is recommended for use with clinical specimens likely to contain mixed microbial flora, such as urine, respiratory and wound, because it allows a preliminary grouping of enteric and other gram-negative bacteria.^{5,6} MacConkey Agar is also used in the BAM (*Bacteriological Analytical Manual*) of the Food and Drug Administration (FDA) procedure for isolating *E. coli* from foods.⁷

The BBL MacConkey II Agar formulation was made available in 1983. It was specially designed to improve the inhibition of swarming *Proteus* species, to achieve more definitive differentiation of lactose fermenters and nonfermenters, and for the promotion of superior growth of enteric pathogens.

VI PRINCIPLES OF THE PROCEDURE

MacConkey II Agar is a selective and differential medium. It is only slightly selective since the concentration of bile salts, which inhibits gram-positive microorganisms, is low in comparison with other enteric plating media. Crystal violet also is included in the medium to inhibit the growth of gram-positive bacteria, especially enterococci and staphylococci.

Differentiation of enteric microorganisms is achieved by the combination of lactose and the neutral red indicator. Colorless or pink to red colonies are produced depending upon the ability of the isolate to ferment the carbohydrate.

VII REAGENTS

MacConkey II Agar

Approximate Formula* Per Liter Purified Water

Pancreatic Digest of Gelatin	17.0 g	Sodium Chloride	5.0 g
Pancreatic Digest of Casein	1.5 g	Neutral Red	0.03 g
Peptic Digest of Animal Tissue	1.5 g	Crystal Violet	0.001 g
Lactose	10.0 g	Agar	13.5 g
Bile Salts.....	1.5 g		

*Adjusted and/or supplemented as required to meet performance criteria.

Warnings and Precautions: For *in vitro* Diagnostic Use.

If excessive moisture is observed, invert the bottom over an off-set lid and allow to air dry in order to prevent formation of a seal between the top and bottom of the plate during incubation.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"⁸⁻¹¹ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids. After use, prepared plates, specimen containers and other contaminated materials must be sterilized by autoclaving before discarding.

Storage Instructions: On receipt, store plates in the dark at 2–8°C. Avoid freezing and overheating. Do not open until ready to use. Minimize exposure to light. Prepared plates stored in their original sleeve wrapping at 2–8°C until just prior to use may be inoculated up to the expiration date and incubated for recommended incubation times. Allow the medium to warm to room temperature before inoculation.

Product Deterioration: Do not use plates if they show evidence of microbial contamination, discoloration, drying, cracking or other signs of deterioration.

VIII SPECIMEN COLLECTION AND HANDLING

Specimens suitable for culture may be handled using various techniques. For detailed information, consult appropriate texts.^{12,13} Specimens should be obtained before antimicrobial therapy has been administered. Provision must be made for prompt delivery to the laboratory.

IX PROCEDURE

Material Provided: MacConkey II Agar

Materials Required But Not Provided: Ancillary culture media, reagents, quality control organisms and laboratory equipment as required.

Test Procedure: Observe aseptic techniques.

The agar surface should be smooth and moist, but without excessive moisture.

Streak the specimen as soon as possible after it is received in the laboratory. The streak plate is used primarily to isolate pure cultures from specimens containing mixed flora. A nonselective medium should also be streaked to increase the chance of recovery when the population of gram-negative organisms is low and to provide an indication of other organisms present in the specimen.

Alternatively, if material is being cultured directly from a swab, roll the swab over a small area of the surface at the edge; then streak from this inoculated area.

Incubate plates, protected from light, at 35 ± 2°C (do not use CO₂-enriched atmosphere with MacConkey II Agar) or other appropriate temperature for 18–24 h.

User Quality Control: See "Quality Control Procedures."

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI (formerly NCCLS) guidance and CLIA regulations for appropriate Quality Control practices.

X RESULTS

After incubation most plates will show an area of confluent growth. Because the streaking procedure is, in effect, a "dilution" technique, diminishing numbers of microorganisms are deposited on the streaked areas. Consequently, one or more of these areas should exhibit isolated colonies of the organisms contained in the specimen. Better isolation is obtained due to the inhibitory action of the medium.

Typical colonial morphology on MacConkey II Agar is as follows:

E. coliPink to rose-red (may be surrounded by a zone of precipitated bile)

Enterobacter/Klebsiella.....Mucoid, pink

Proteus.....Colorless, swarming in areas of isolated colonies is inhibited

SalmonellaColorless

ShigellaColorless

Pseudomonas.....Irregular, colorless to pink

Gram-positive bacteriaNo growth to slight growth

XI LIMITATIONS OF THE PROCEDURE

It has been reported that some *Enterobacteriaceae* and *Pseudomonas aeruginosa* are inhibited on MacConkey Agar when incubated in a CO₂-enriched atmosphere.¹³

Not all strains of *E. coli* ferment lactose.

For identification, organisms must be in pure culture. Morphological, biochemical, and/or serological tests should be performed for final identification. Consult appropriate texts for detailed information and recommended procedures.^{5,12,14-17}

A single medium is rarely adequate for detecting all organisms of potential significance in a specimen. Cultures of specimens grown on selective media should, therefore, be compared with specimens cultured on nonselective media to obtain additional information and help ensure recovery of potential pathogens.

XII PERFORMANCE CHARACTERISTICS

Prior to release, all lots of MacConkey II Agar are tested for performance characteristics. Representative samples of the lot are streak-inoculated with the following cultures: *Escherichia coli* (ATCC 25922), *Proteus mirabilis* (ATCC 12453), *Pseudomonas aeruginosa* (ATCC 10145), *Salmonella typhimurium* (ATCC 14028), *Shigella dysenteriae* (ATCC 9361) and *Enterococcus faecalis* (ATCC 29212). The inoculum for *E. faecalis* is diluted to yield 10⁴–10⁵ colony-forming units (CFU) per plate; the inocula for all other organisms is diluted to yield 10³–10⁴ CFU/plate. After inoculation, the plates are incubated at 35 ± 2°C in an aerobic atmosphere. After 18–24 h incubation, colonies of *E. coli* are rose-red and may be surrounded by precipitated bile; *P. mirabilis* exhibits fair to heavy growth of colorless colonies and swarming of the colonies is inhibited; *P. aeruginosa* shows areas of confluent growth which may exhibit green to yellow-green pigmentation while individual colonies show pink to green pigmentation; *S. typhimurium* gives fair to heavy growth of colorless colonies; *S. dysenteriae* shows growth of colorless to pink colonies; *E. faecalis* is completely to partially inhibited (fair growth) and the colonies may be pink in color.

XIII AVAILABILITY

Cat. No.	Description
221172	BBL™ MacConkey II Agar, Pkg. of 20 plates
221884	BBL™ MacConkey II Agar, Pkg. of 30 Space-saver plates
221270	BBL™ MacConkey II Agar, Ctn. of 100 plates
292539	BBL™ MacConkey II Agar, Ctn. of 150 Space-saver plates

XIV REFERENCES

1. MacConkey, A.T. 1900. Note on a new medium for the growth and differentiation of the *Bacillus coli communis* and the *Bacillus typhi abdominalis*. The Lancet, Part II:20.
2. MacConkey, A. 1905. Lactose-fermenting bacteria in faeces. J. Hyg. 5:333-379.
3. Levine, M., and H.W. Schoenlein. 1930. A compilation of culture media for the cultivation of microorganisms. The Williams & Wilkins Company, Baltimore.
4. MacFaddin, J.F. 1985. Media for isolation-cultivation-identification-maintenance of medical bacteria, vol. I. Williams & Wilkins, Baltimore.
5. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2002. Bailey & Scott's diagnostic microbiology, 11th ed. Mosby, Inc., St. Louis.
6. Farmer, J.J., III. 2003. *Enterobacteriaceae*: introduction and identification, p. 636-671. In P.R. Murray, E.J. Baron, J.H. Jorgensen, M.A. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
7. Hitchins, A., P. Feng, W. Watkins, S. Rippey, and L. Chandler. 1998. *Escherichia coli* and the coliform bacteria, p. 4.01-4.29. In FDA bacteriological analytical manual. Association of Official Analytical Chemistry International, Gaithersburg, Md.
8. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, PA.
9. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
10. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
11. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
12. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Tenover, and R.H. Tenover (ed.). 2003. Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
13. Mazura-Reetz, G., T.R. Neblett, and J.M. Galperin. 1979. MacConkey agar: CO₂ vs. ambient incubation, abstr. C 179, p. 339. Abstr. 79th Annu. Meet. Am. Soc. Microbiol. 1979.
14. Holt, J.G., N.R. Krieg, P.H.A. Sneath, J.T. Staley, and S.T. Williams (ed.). 1994. Bergey's Manual™ of determinative bacteriology, 9th ed. Williams & Wilkins, Baltimore.
15. MacFaddin, J.F. 2000. Biochemical tests for identification of medical bacteria, 3rd ed. Lippincott Williams & Wilkins, Baltimore.
16. Koneman, E.W., S.D. Allen, W.M. Janda, P.C. Schreckenberger, and W.C. Winn, Jr. 1997. Color atlas and textbook of diagnostic microbiology, 5th ed. Lippincott-Raven, Philadelphia.
17. Isenberg, H.D. (ed.). 2004. Clinical microbiology procedures handbook, vol. 1, 2 and 3, 2nd ed. American Society for Microbiology, Washington, D.C.

Becton, Dickinson and Company
7 Loveton Circle
Sparks, Maryland 21152 USA
800-638-8663

ATCC is a trademark of the American Type Culture Collection.

BD, BD Logo, BBL and Trypticase are trademarks of Becton, Dickinson and Company. ©2006 BD.