Revisions

Rev from	Rev to	ECO #
0108	2010/02	5262-10

Notes:

- BD Cat. Number 442260
- Blank (Sheet) Size: Length: N/A Width: N/A
 Number of Pages: N/A Number of Sheets: N/A
 Page Size: Length N/A Width N/A Final Folded Size: N/A
- Style (see illustrations below): # N/A

- See Specification Control Number VS-13-7220-0PR for Material Information
- Ink Colors: Printed two sides Yes No
 No. of Colors: N/A PMS# N/A
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

VS Controlled by BD Caribe, Ltd.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: PP105JAA		Category and Description Package Insert, BACTEC Standard/10 Aerobic/F	Sheet: 1 of 25 <hr/> Scale: N/A	A

BD BACTEC™ Standard/10 Aerobic/F Culture Vials

English:	pages	1 – 3	Español:	páginas	11 – 13
Français :	pages	3 – 6	Dansk:	side	13 – 15
Deutsch:	Seiten	6 – 8	Português:	páginas	16 – 18
Italiano:	pagine	8 – 11	Svenska:	sidan	18 – 20

PP105JAA
2010/02

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőtől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / **Свържете се с местния представител на BD за инструкции.** / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіңіздің жергілікті BD өкіліне жүгініп нұсқау алыңыз. / Kontaktarj lokalnog predstavnika BD za upute.

INTENDED USE

BACTEC™ Standard/10 Aerobic/F culture vials (enriched Soybean-Casein Digest broth with CO₂) are for aerobic blood cultures. Principal use is with the BACTEC fluorescent series instruments for the qualitative culture and recovery of aerobic microorganisms (bacteria and yeast) from blood.

SUMMARY AND EXPLANATION

The sample to be tested is inoculated into one or more vials which are inserted into the BACTEC fluorescent series instrument for incubation and periodic reading. Each vial contains a chemical sensor which can detect increases in CO₂ produced by the growth of microorganisms. The sensor is monitored by the instrument every ten minutes for an increase in its fluorescence, which is proportional to the amount of CO₂ present. A positive reading indicates the presumptive presence of viable microorganisms in the vial. Detection is limited to microorganisms that will grow in a particular type of medium.

PRINCIPLES OF THE PROCEDURE

If microorganisms are present in the test sample inoculated into the BACTEC vial, CO₂ will be produced when the organisms metabolize the substrates present in the vial. Increases in the fluorescence of the vial sensor caused by the higher amount of CO₂ are monitored by the BACTEC fluorescent series instrument. Analysis of the rate and amount of CO₂ increase enables the BACTEC fluorescent series instrument to determine if the vial is positive; i.e., that the test sample contains viable organisms.

REAGENTS

The BACTEC Standard/10 Aerobic/F culture vials contain the following active ingredients prior to processing:

List of Ingredients

(WTR) Processed Water.....	40 mL
(SCB) Soybean-Casein Digest Broth	3.0% w/v
(YEX) Yeast Extract	0.3% w/v
(ATD) Animal Tissue Digest	0.01% w/v
(SCR) Sucrose	0.1% w/v
(HEM) Hemin	0.0005% w/v
(MEN) Menadione	0.00005% w/v
(PXH) Pyridoxal HCl (Vitamin B ₆)	0.001% w/v
(SBC) Sodium Bicarbonate	0.04% w/v
(SPS) Sodium Polyanetholsulfonate	0.035% w/v

All BACTEC media are dispensed with added CO₂.

Warnings and Precautions

For *in vitro* Diagnostic Use.

This Product Contains Dry Natural Rubber.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"¹⁻⁴ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

Prior to use, each vial should be examined for evidence of damage, contamination or deterioration. Vials displaying evidence of damage or contamination such as leakage, cloudiness, discoloration (darkening), bulging or depressed septum should not be used.

A contaminated vial could contain positive pressure. If a contaminated vial is used for direct draw, contaminated culture media could be refluxed into the patient's vein. Vial contamination may not be readily apparent. When using direct draw procedures, monitor the process closely to avoid refluxing materials into patient.

On rare occasions, the glass bottle neck may be cracked and the neck may break during removal of the flip-off cap or in handling. Also, on rare occasions, a vial may not be sealed sufficiently. In both cases the contents of the vial may leak or spill. If the vial has been inoculated, treat the leak or spill with caution, as pathogenic organisms/agents may be present. Before discarding, sterilize all inoculated vials by autoclaving.

Positive culture vials for subculturing or staining, etc.: Before sampling it is necessary to release gas which often builds up due to microbial metabolism. Sampling should be performed in a biological safety cabinet if possible, and appropriate protective clothing, including gloves and masks, should be worn. See Procedure section for more information on subculturing.

To minimize the potential of leakage during inoculation of specimen into culture vials, use syringes with permanently attached needles or Luer-Lok™ brand tips.

Storage Instructions

The BACTEC vials are ready for use as received and require no reconstitution or dilution. Store at 2–25°C, in a dry place **out of direct light**.

SPECIMEN COLLECTION

The specimen must be collected using sterile techniques to reduce the chance of contamination. The typical specimen volume is 8 – 10 mL. It is recommended that the specimen be inoculated into the BACTEC vials at bedside. Most commonly, a 10 cc or 20 cc syringe with a Luer-Lok brand tip is used to draw the sample. If appropriate, a Vacutainer™ brand Needle Holder and a Vacutainer™ brand Blood Collection Set, Vacutainer™ Safety-Lok™ Blood Collection Set or other tubing "butterfly" set may be used. If using a needle and tubing set (direct draw), carefully observe the direction of blood flow when starting sample collection. The vacuum in the vial will usually exceed 10 mL, so the user should monitor the volume collected by means of the 5 mL graduation marks on the vial label. When the desired 8 – 10 mL has been drawn, the flow should be stopped by crimping the tubing and

removing the tubing set from the BACTEC vial. Sample volumes as low as 3 mL can be used, however, recovery will not be as great as with larger volumes. The inoculated BACTEC vial should be transported as quickly as possible to the laboratory.

PROCEDURE

Remove the flip-off cap from BACTEC vial top and inspect the vial for cracks, contamination, excessive cloudiness, and bulging or indented septums. **DO NOT USE** if any defect is noted. Before inoculating, swab the septum with alcohol (iodine is **NOT** recommended). Aseptically inject or draw directly 8 – 10 mL of specimen per vial. If sample volumes of 3 – 4 mL are used, recovery will not be as great as with larger volumes (see Limitations of the Procedure). **Inoculated aerobic vials should be placed in the BACTEC fluorescent series instrument as soon as possible** for incubation and monitoring. If placement of an inoculated vial into the instrument has been delayed and visible growth is apparent, it should not be tested in the BACTEC fluorescent series instrument, but rather it should be subcultured, Gram-stained and treated as a presumptively positive bottle.

Vials entered into the instrument will be automatically tested every ten minutes for the duration of the testing protocol period. Positive vials will be determined by the BACTEC fluorescent series instrument and identified as such (see the appropriate BACTEC fluorescent series instrument User's Manual). The sensor inside the bottle will not appear visibly different in positive and negative vials, however the BACTEC fluorescent series instrument can determine a difference in fluorescence.

If at the end of the testing period a negative BACTEC Standard/10 Aerobic/F vial appears visually positive (i.e., chocolate blood, bulging septum, lysed and/or very darkened blood), it should be subcultured, Gram-stained and treated as a presumptive positive.

Positive vials should be subcultured and a Gram-stained slide prepared. In a great majority of cases, organisms will be seen and a preliminary report can be made to the physician. Subcultures to selective media and a preliminary direct antimicrobial susceptibility test may be prepared from fluid in the BACTEC vials.

Subculturing: Prior to subculturing, put the vial in an upright position, and place an alcohol wipe over the septum. To release pressure in the vial, insert a sterile needle with an appropriate filter or pledget through the alcohol wipe and septum. The needle should be removed after the pressure is released and before sampling the vial for subculture. The insertion and withdrawal of the needle should be done in a straight-line motion, avoiding any twisting motions.

For maximum yield of isolates, negative cultures may be checked by stain and/or subcultured at some point prior to discarding as negative.

QUALITY CONTROL

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI (formerly NCCLS) guidance and CLIA regulations for appropriate Quality Control practices.

DO NOT USE culture vials past their expiration date.

DO NOT USE culture vials that exhibit any cracks or defects; discard the vial in the appropriate manner.

Quality Control Certificates are provided with each carton of media. Quality Control Certificates list test organisms, including ATCC™ cultures specified in the CLSI Standard, *Quality Control for Commercially Prepared Microbiological Culture Media*.⁵

The range of time-to-detection in hours was ≤ 72 hours for each of the organisms listed on the Quality Control Certificate for this medium:

<i>Streptococcus pyogenes</i> ATCC 19615	<i>Escherichia coli</i> ATCC 25922	<i>Streptococcus pneumoniae</i> * ATCC 6305
<i>Pseudomonas aeruginosa</i> ATCC 27853	<i>Candida albicans</i> ATCC 18804	<i>Neisseria meningitidis</i> ATCC 13090
<i>Alcaligenes faecalis</i> ATCC 8750	<i>Haemophilus influenzae</i> ATCC 19418	<i>Staphylococcus aureus</i> ATCC 25923

*CLSI Strain

For information on Quality Control for the BACTEC fluorescent series instrument, refer to the appropriate BACTEC fluorescent series instrument User's Manual.

LIMITATIONS OF THE PROCEDURE

Contamination

Care must be taken to prevent contamination of the sample during collection and inoculation into the BACTEC vial. A contaminated sample will give a positive reading, but will not indicate a relevant clinical result. Such a determination must be made by the user based on such factors as type of organisms recovered, occurrence of the same organism in multiple cultures, patient history, etc.

Recovery of SPS Sensitive Organisms from Blood Samples

Because blood can neutralize the toxicity of SPS toward organisms sensitive to SPS (such as some *Neisseria* species), the presence of optimum volumes of blood (8 – 10 mL) is a benefit in the recovery of these organisms.

Some fastidious organisms, such as certain *Haemophilus* species, require growth factors, such as NAD, or factor V, which are provided by the blood specimen. If the blood specimen volume is 3.0 mL or less, an appropriate supplement may be required for recovery of these organisms. BACTEC™ FOSTM, Fastidious Organism Supplement, may be used as a nutritional supplement.

Nonviable Organisms

A Gram-stained smear from culture medium may contain small numbers of nonviable organisms derived from media constituents, staining reagents, immersion oil, glass slides, and specimens used for inoculation. In addition, the patient specimen may contain organisms that will not grow in the culture medium or in media used for subculture. Such specimens should be subcultured to special media as appropriate.⁶

Recovery of *Streptococcus pneumoniae*

In aerobic media, *S. pneumoniae* will typically be visually and instrument positive, but in some cases no organism will be seen on Gram stain or recovered on routine subculture. If an anaerobic vial was also inoculated, the organism can usually be recovered by performing an aerobic subculture of the anaerobic vial, since this organism has been reported to grow well under anaerobic conditions.⁷

General Considerations

Optimum recovery of isolates will be achieved by adding 8 – 10 mL of blood. Use of lower or higher volumes may adversely affect recovery and/or detection times. Blood may contain antimicrobials or other inhibitors which may slow or prevent the growth of microorganisms. False negative readings may result when certain organisms are present which do not produce enough CO₂ to be detected by the system or significant growth has occurred before placing the vial into the system. False positivity may occur when the white blood cell count is high.

EXPECTED RESULTS

Seeded culture studies were performed using inoculum levels targeted at 10 to 50 CFU per culture vial of a combination of ATCC and wild microbial strains. The following is a list of organisms which were detected as positive in BACTEC Standard/10 Aerobic/F medium within a five (5) day period.

<i>Acinetobacter anitratus</i>	<i>Cryptococcus neoformans</i>	<i>Proteus mirabilis</i>	<i>Staphylococcus epidermidis</i>
<i>Acinetobacter lwoffii</i>	<i>Enterobacter cloacae</i>	<i>Providencia stuartii</i>	<i>Streptococcus pneumoniae</i>
<i>Alcaligenes faecalis</i>	<i>Enterococcus faecalis</i>	<i>Pseudomonas aeruginosa</i>	<i>Streptococcus pyogenes</i> (grp. A)
<i>Candida (Torulopsis) glabrata</i>	<i>Escherichia coli</i>	<i>Salmonella typhimurium</i>	<i>Xanthomonas maltophilia</i>
<i>Candida albicans</i>	<i>Haemophilus influenzae</i>	<i>Serratia marcescens</i>	
<i>Corynebacterium J-K</i>	<i>Klebsiella pneumoniae</i>	<i>Staphylococcus aureus</i>	

PERFORMANCE CHARACTERISTICS

A clinical study was conducted comparing BACTEC Standard/10 Aerobic/F blood culture medium and BACTEC Standard Aerobic/F blood culture medium. A total of 1384 paired vials were evaluated. No false positive vials were observed. Three false negative BACTEC Standard/10 Aerobic/F blood culture medium vials were identified during this study. The false negative rate was calculated as 0.2%.

A total of 113 organisms were recovered. Table 1 shows the isolates recovered by media type. Of these 76 (67.3%) were judged to be clinically significant. Of 76 clinically significant isolates, 53 isolates (69.7%) were recovered in both media, 8 (10.5%) were recovered only in BACTEC Standard/10 Aerobic/F blood culture medium and 15 isolates (19.7%) were recovered only in BACTEC Standard Aerobic/F blood culture medium. The difference in recovery only in BACTEC Standard/10 Aerobic/F blood culture medium compared to recovery only in BACTEC Standard Aerobic/F blood culture medium was not statistically significant. The average time to detection for all groups was 20.4 hours in BACTEC Standard/10 Aerobic/F blood culture medium compared to 21.5 hours in BACTEC Standard Aerobic/F blood culture medium.

TABLE 1: Clinical Study Isolate Recovery – Media Type

Organism	Recovered in Standard/10 Aerobic/F ONLY	Recovered in Standard Aerobic/F ONLY	Recovered in BOTH
GRAM NEGATIVES	7	11	30
GRAM POSITIVES	0	4	22
YEAST	1	0	1

List of organisms recovered in BACTEC Standard/10 Aerobic/F media:

<i>Candida albicans</i>	<i>Escherichia coli</i>	<i>Pseudomonas aeruginosa</i>	<i>Streptococcus pneumoniae</i>
<i>Candida tropicalis</i>	<i>Haemophilus influenzae</i>	<i>Staphylococcus aureus</i>	<i>Streptococcus</i> species
<i>Citrobacter diversus</i>	<i>Klebsiella pneumoniae</i>	<i>Staphylococcus</i> coag. negative	
<i>Enterobacter aerogenes</i>	<i>Morganella morganii</i>	<i>Staphylococcus haemolyticus</i>	
<i>Enterococcus faecalis</i>	<i>Proteus mirabilis</i>	<i>Streptococcus constellatus</i>	

AVAILABILITY

Cat. No. Description

442260 BACTEC™ Standard/10 Aerobic/F Culture Vials, case of 50 vials.

REFERENCES:

1. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
3. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262, 17/10/2000, p. 0021-0045.*
5. Clinical and Laboratory Standards Institute. 2004. Approved Standard M22-A3. Quality control of commercially prepared microbiological culture media, 3rd ed. CLSI, Wayne, Pa.
6. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry and M.A. Pfaller (ed.). 2007. *Manual of clinical microbiology*, 9th ed. American Society for Microbiology, Washington, D.C.
7. Howden, R.J., *J. Clin. Path.* 1976, 29:50-53.

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Benex Limited