**Revisions**

Rev from	Rev to	ECO #
0108	2010/02	5262-10

**Notes:**

- BD Cat. Number 442003, 442288
- Blank (Sheet) Size : Length: N/A      Width: N/A  
 Number of Pages: N/A    Number of Sheets: N/A  
 Page Size: Length N/A    Width N/A      Final Folded Size: N/A
- Style (see illustrations below): # N/A

- See Specification Control Number VS-13-7270-0PR for Material Information
- Ink Colors: Printed two sides  Yes     No  
 No. of Colors: 1      PMS #199 Red
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

VS Controlled by BD Caribe, Ltd.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: PP162JAA		Category and Description Package Insert, BACTEC Myco/F Lytic Blood Culture	Sheet: 1 of 39 Scale: N/A	<b>A</b>

# **BD BACTEC™ Myco/F Lytic Culture Vials**

## Supplemented Middlebrook 7H9 and Brain Heart Infusion Broth

### For Use with BACTEC Fluorescent Series Instruments

English: pages 1 – 5    Italiano: pagine 14 – 18    Português : páginas 27 – 31  
 Français : pages 5 – 9    Español: páginas 18 – 22    Svenska: sidan 32 – 36  
 Deutsch : Seiten 10 – 14    Dansk: side 23 – 27

 PP162JAA  
2010/02

Pokyny vám poskytnú miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőtől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcije získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіндіздің жергілікті БД өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

#### INTENDED USE

BACTEC™ Myco/F Lytic culture medium when used with the BACTEC fluorescent series instruments is a nonselective culture medium to be used as an adjunct to aerobic blood culture media for the recovery of mycobacteria, yeast and fungi from blood. This medium may also be used for the culture of sterile body fluids when yeast or fungi are suspected.

#### SUMMARY AND EXPLANATION

Since the mid-1980s and expanding size of the immunocompromised patient population, the incidence of septicemia caused by opportunistic pathogens such as yeast, fungi and mycobacteria has risen. *Mycobacterium tuberculosis* (MTB) and mycobacteria other than tuberculosis (MOTT), especially *Mycobacterium avium* complex (MAC), have become resurgent. From 1985 to 1992, the number of MTB cases reported increased 18%. Between 1981 and 1987, AIDS case surveillances indicated that 5.5% of the patients with AIDS had disseminated nontuberculous mycobacterial infections, e.g., MAC. By 1990, the increased cases of disseminated nontuberculous mycobacterial infections had resulted in a cumulative incidence of 7.6%.<sup>1</sup> It has also been noted that the incidence of fungemia has steadily increased since the early 1980s. This has increased the need for the clinical laboratory to have effective diagnostic procedures for fungemia and mycobacteremia.

The Centers for Disease Control and Prevention (CDC) have recommended that every effort must be made for laboratories to use the most rapid methods available for diagnostic mycobacteria testing. These recommendations include the use of a liquid medium for mycobacterial culture.<sup>2-4</sup>

The BACTEC fluorescent series instruments are designed for the rapid detection of microorganisms in clinical specimens. BACTEC Myco/F Lytic Culture medium is a Middlebrook 7H9 and Brain Heart Infusion broth formulation for the recovery of mycobacteria from blood specimens and yeast and fungi from blood and sterile body fluids. Specific modifications were made to enhance the growth and recovery of mycobacteria, yeast and fungi. These modifications include ferric ammonium citrate to provide an iron source for specific strains of mycobacteria and fungi, the addition of saponin as a blood lysing agent and the addition of specific proteins and sugars to provide nutritional supplements. Each vial contains a sensor which can detect decreases in oxygen concentration in the vial resulting from microorganism metabolism and growth. The sensor is monitored by the BACTEC fluorescent series instrument for increasing fluorescence which is proportional to the decrease in oxygen. A positive determination indicates the presumptive presence of viable microorganisms in the vial.

#### PRINCIPLES OF THE PROCEDURE

The BACTEC Myco/F Lytic culture vial is designed for the rapid detection of mycobacteria in blood, and yeast and fungi in blood and sterile body fluids. Specimens are inoculated into the BACTEC Myco/F Lytic vial either with a syringe or direct draw with a needle and tubing. The vial is placed into the BACTEC fluorescent series instrument and is continuously agitated and incubated at 35°C for maximum recovery. The default testing protocol is 42 days. The recommended testing protocol for the following organisms are 7 days for yeast, 30 days for fungi and 42 days for mycobacteria. Each vial contains a sensor which can detect decreases in oxygen concentration in the vial resulting from microorganism metabolism and growth. The sensor is monitored by the BACTEC fluorescent series instrument every ten minutes. Analysis of the rate of oxygen decrease as measured by increasing fluorescence enables the BACTEC fluorescent series instrument to determine if the vial is instrument positive. A positive determination indicates the presumptive presence of viable microorganisms in the vial. Detection is limited to microorganisms that will grow in the medium at 35°C. The medium is not selective and will support the growth of other aerobic organisms including bacteria which may interfere, if present, with the recovery of slower growing mycobacteria, yeast and fungi. Culture vials which remain negative after the completion of protocol and which show no visible sign of positivity are removed from the instrument and sterilized prior to discarding.

#### REAGENTS

Each BACTEC Myco/F Lytic culture vial contains the following active ingredients prior to processing:

##### List of Ingredients

Processed Water .....	40 mL qs
7H9 Middlebrook Broth Base without phosphate salts .....	0.12% w/v
Brain Heart Infusion.....	0.5% w/v
Casein Hydrolysate.....	0.10% w/v
Supplement H.....	0.10% w/v
Inositol.....	0.05% w/v
Glycerol.....	0.10% w/v
Sodium Polyanetholsulfonate .....	0.025% w/v
Polysorbate 80.....	0.0025% w/v
Pyridoxal HCl .....	0.0001% w/v
Ferric Ammonium Citrate .....	0.006% w/v
Potassium Phosphate .....	0.024% w/v
Saponin .....	0.24% w/v
Antifoam.....	0.01% w/v

Composition may have been adjusted to meet specific performance requirements.

This BACTEC medium is dispensed with added CO<sub>2</sub> and O<sub>2</sub>.

BACTEC Myco/F Lytic medium requires no supplement addition. Each 40 mL vial of BACTEC Myco/F Lytic is ready for use when received. The appearance of the media upon receipt should be clear and light amber in color.

## Warnings and Precautions

For *in vitro* Diagnostic Use.

This Product Contains Dry Natural Rubber.

**Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"<sup>5-8</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.**

If recovery of mycobacteria is intended, CDC-NIH guidelines strongly recommend that the test instrument be placed in the mycobacteria laboratory where the additional safety issues that the recovery of mycobacteria present can be addressed.<sup>9</sup>

**BACTEC Myco/F Lytic vials will accept more than the recommended maximum of 5 mL of specimen volume. Monitoring of fill should be conducted.**

For activities involving the propagation and manipulation of *Mycobacterium tuberculosis* or *Mycobacterium bovis* grown in culture, Biosafety Level 3 practice, containment equipment and facilities are recommended.<sup>9</sup>

Prior to use, each vial should be examined for evidence of contamination such as cloudiness, bulging or depressed septum, or leakage. **DO NOT USE** any vial showing evidence of contamination, leakage or damage. Vial contamination may not be readily apparent. A contaminated vial could contain positive pressure. If a contaminated vial is used for direct draw, gas or contaminated culture media could be refluxed into the patient's vein. On rare occasions, the glass bottle neck may be cracked and the neck may break during removal of the flip-off cap or in handling. Also, on rare occasions, a vial may not be sealed sufficiently. In both cases the contents of the vials may leak or spill, especially if the vial is inverted.

To minimize the potential of leakage during inoculation by syringe of specimen into culture vials, use syringes with Luer-Lok™ brand tips. A one-handed inoculation technique and a suitable vial holder should be employed to prevent accidental needle stick injury.

Before discarding, sterilize all inoculated BACTEC Myco/F Lytic vials by autoclaving.

**Positive culture vials for subculturing or staining, etc.:** Before sampling it is necessary to release gas which often builds up due to microbial metabolism. Sampling and venting of vials must be performed in a biological safety cabinet, and appropriate protective clothing, including gloves and masks, should be worn. See PROCEDURE Section for more information on subculturing.

## Leaking or Broken Vials

**CAUTION: Because an inoculated leaking or broken vial may produce an aerosol of mycobacteria, including *M. tuberculosis*, or other bacteria, appropriate handling should be observed.**

If an inoculated vial is found to be leaking or is accidentally broken during collection or transport, use the established procedure in your facility for dealing with mycobacterial spills. As a minimum, "Standard Precautions" should be employed. Vials should be discarded in an appropriate manner.

If a vial is found to have leaked contents into the instrument proper, or if a vial is accidentally broken, turn off the instrument immediately. Vacate the affected area. Contact your facility's Safety or Infection Control Officer(s). Determine the necessity of turning off or modifying the settings of the air handling units serving the affected area. Do not return to the area until any potential aerosols have settled or have been removed by appropriate ventilation. BD Diagnostics should be notified by calling 1-800-638-8663 in the USA or the appropriate BD representative in your area. Guidelines for proper handling of accidental mycobacterial contamination due to breakage of culture tubes or broth suspensions have been issued by the CDC.<sup>9</sup>

## Storage Instructions

Store at 2 – 25°C in a dry location out of direct light.

DO NOT use after expiration date.

## SPECIMEN COLLECTION

**NOTE: It is recommended that this procedure be reviewed with the appropriate personnel prior to use of medium to ensure proper specimen collection techniques as described in this section.**

The specimen must be collected using sterile technique to reduce the chance of contamination. The range of blood volume which can be cultured is 1 mL to 5 mL, with optimum recovery obtained at 3 mL to 5 mL. It is recommended that the specimen be inoculated at bedside. Most commonly, a syringe with a Luer-Lok brand tip is used to draw the specimen. If appropriate, a Vacutainer™ brand Needle Holder and a Vacutainer brand Blood Collection Set, Vacutainer Safety-Lok™ Blood Collection Set or other tubing "butterfly" set may be used. If using a needle and tubing (direct draw), carefully observe the direction of the blood flow when starting sample collection. Prior to inoculation, the medium fill volume should be noted on the label with a pen or marker to indicate the starting point of specimen collection. The vacuum in the bottle will usually exceed 5 mL, so the user should monitor the volume collected by means of the 5 mL graduation marks on the vial label. When the desired 1 – 5 mL of blood has been drawn, the flow should be stopped by crimping the tubing and removing the needle from the BACTEC vial. The BACTEC vial should be transported as quickly as possible to the laboratory and placed in the BACTEC instrument. A yellow-top Vacutainer Brand Tube containing SPS may also be used to collect the blood sample from the patient. The tube should be transported to the laboratory as quickly as possible for transfer into the BACTEC culture vial.

## PROCEDURE

**Materials Provided:** BACTEC Myco/F Lytic Culture Vials.

**Materials Required But Not Provided:** Biological Safety Cabinet, autoclave, venting unit, mycobacterial disinfectant, 70% isopropyl alcohol, Quality Control Organisms (*Mycobacterium intracellulare*, ATCC™ 13950; *Candida glabrata*, ATCC 15545; *Cryptococcus neoformans*, ATCC 13690), microscope and materials for staining slides and subculturing vials.

Inoculation of BACTEC Myco/F Lytic Culture Vials

1. Remove the flip-off cap from the BACTEC vial top and inspect the vial for cracks, leaks, contamination, excessive cloudiness and bulging or indented septum. DO NOT USE if any defect is noted.
2. Label culture vial with specimen identification and mark medium fill graduation line on vial label.
3. Before inoculating, swab the septum with alcohol. Aseptically inject with a syringe or draw directly with the aid of the graduation lines on the vial label 1 – 5 mL of specimen per vial (see the section on Limitations of the Procedure). Inoculated vials should be placed into the BACTEC fluorescent series instruments as soon as possible for incubation and monitoring.
4. Vials entered into the instruments will be automatically tested for the duration of the testing protocol. The default testing protocol is 42 days. The recommended testing protocol for the following organisms are 7 days for yeast, 30 days for fungi and 42 days for mycobacteria. See the appropriate BACTEC User's Manual, to set protocol length. If at the end of the protocol, a negative BACTEC Myco/F Lytic vial appears visually positive (i.e., bulging septum), it should be subcultured, AFB and Gram stained and treated as a presumptive positive.
5. Positive vials will be identified by the BACTEC fluorescent series instrument. The sensor inside the vial may not appear visibly different in positive or negative vials; however, the BACTEC fluorescent series instrument can determine a difference in sensor fluorescence. **All positive vials should be handled in a Biological Safety Cabinet. Biosafety Level 3 practices, containment equipment and facilities are recommended.**

Positive vials should be subcultured and an appropriate smear prepared.

Processing an instrument-positive vial

- a) Remove the vial from the instrument and place in a biological safety cabinet.
- b) Invert vial to mix contents.
- c) Vent the vial to equilibrate vial pressure with atmosphere.
- d) Remove aliquot from vial ( approx. 0.1 mL) for stain preparations (AFB and Gram).
- e) Inspect smear and report preliminary results only after smear evaluation.

**Subculturing:** Subculturing should be performed in a biological safety cabinet, and appropriate clothing, including gloves and masks, should be worn. Prior to subculturing, place the vial in an upright position, and place an alcohol wipe over the septum. To release any positive pressure in the vial which could be caused by growth of possible contaminants, insert a sterile 25-gauge (or smaller) needle equipped with an appropriate filter or pledget through the alcohol wipe and septum. The needle should be removed after any pressure is released and before sampling the vial for subculture. The insertion and withdrawal of the needle should be done in a straight-line motion, avoiding any side-to-side motions which could permanently damage the septum. **Do not re-cap the needle. Discard needles and syringes in a puncture-resistant biohazard container.**

#### QUALITY CONTROL

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI (formerly NCCLS) guidance and CLIA regulations for appropriate Quality Control practices.

It is recommended that each new shipment or lot of BACTEC Myco/F Lytic media be tested with the ATCC control organisms identified in the chart below as a positive control, and an uninoculated vial as a negative control.

Organism	Range of Time-to-Detection (days)
<i>Mycobacterium intracellulare</i> , ATCC 13950	8 to 16
<i>Candida glabrata</i> , ATCC 15545	< 3
<i>Cryptococcus neoformans</i> , ATCC 13690	< 3

The positive vials should be inoculated using a 1:100 dilution of a McFarland #1 suspension grown on solid medium. Inoculate the vial with 0.1 mL of the diluted culture. The vials and an uninoculated control vial should be scanned into the instrument and tested. The inoculated vial should be detected as positive by the instrument within the test protocol. The negative control should remain negative. If expected results for Quality Control are not obtained, do not use the medium and contact BD Technical Services (in the U.S only: 1-800-638-8663) or your local BD Representative for further assistance.

For information on quality control for the BACTEC System, refer to the appropriate BACTEC User's Manual.

#### Reporting of RESULTS

An instrument positive vial must be confirmed by acid-fast smear or Gram stain. A positive result indicates the presumptive presence of viable microorganisms in the vial.

If AFB smear or Gram stain positive, subculture to solid media and report as: instrument-positive, AFB or Gram stain positive, ID pending.

If no microorganisms are present on the smears, subculture to solid media, re-enter the vial into the instrument as an ongoing negative vial within 3 h of removal and allow to complete test protocol. No reportable result.

Perform subcultures from the BACTEC Myco/F Lytic vial for identification and susceptibility testing.

#### LIMITATIONS OF THE PROCEDURE

BACTEC Myco/F Lytic vials are not selective and will support the growth of other aerobic organisms besides mycobacteria, yeast and fungi. Positive vials may contain one or more species of mycobacteria and/or other non-mycobacterial species. If present, fast growing organisms may mask the detection of slower growing mycobacteria, yeast and fungi. Subculture and additional procedures are required. The consistency of microscopic morphology in BACTEC Myco/F Lytic has not been established.

Inoculation of blood volumes of 1 – 5 mL are acceptable, but optimum recovery is obtained with 3 – 5 mL. During internal studies with less than 3 mL of blood, *M. intracellulare*, *M. malmoense*, *M. haemophilum* and *M. xenopi* exhibited detection delays and/or compromised recovery with BACTEC Myco/F Lytic. False positivity most likely will increase when the blood volume is above 5 mL.

Care must be taken to prevent contamination of the sample during collection and inoculation into the BACTEC vial. A contaminated vial will give a positive instrument reading, but will not indicate a relevant clinical result. Such a determination must be made by the user, based on such factors as stain results, type of organism recovered, occurrence of the same organism in multiple cultures, patient history, etc.

Mycobacteria may vary in acid-fastness depending on strain, age of culture and other variables.

Blood may contain antimicrobials or other inhibitors which may slow or prevent the growth of microorganisms.

BACTEC Myco/F Lytic vials are incubated at 35°C potentially precluding the recovery of mycobacteria requiring other incubation temperatures (e.g., *M. marinum*, *M. ulcerans*, *M. haemophilum*). Recovery of such organisms may require additional culture methods.

*Penicillium purpurescens* and *Blastomyces dermatitidis* were not detectable in the BACTEC Myco/F Lytic culture medium. *Hansenula anomala*, *Exophiala jeamselmei*, *Actinomyces bovis*, *Rhodotorula rubra*, and *Mucor ramosissimus* exhibited inconsistent results at low inoculum levels (<10 CFU/vial) with seeded culture studies. Recovery of such organisms may require additional culture methods.

The following isolates were detected as positive in the BACTEC 9240 instrument using BACTEC Myco/F Lytic medium during internal seeded studies and/or clinical trials:

<i>Mycobacterium terrae</i>	<i>Mycobacterium abscessus</i>	<i>Cryptococcus neoformans</i>
<i>Mycobacterium tuberculosis</i>	<i>Mycobacterium xenopi</i>	<i>Histoplasma capsulatum</i>
<i>Mycobacterium avium</i>	<i>Mycobacterium malmoense</i>	<i>Aspergillus flavus</i>
<i>Mycobacterium kansasii</i>	<i>Mycobacterium haemophilum</i>	<i>Aspergillus fumigatus</i>
<i>Mycobacterium fortuitum</i>	<i>Candida albicans</i>	<i>Noctardia asteroides</i>
<i>Mycobacterium intracellulare</i>	<i>Candida glabrata</i>	<i>Malassezia furfur</i>
<i>Mycobacterium gordonae</i>	<i>Candida parapsilosis</i>	<i>Trichophyton rubrum</i>
<i>Mycobacterium szulgai</i>	<i>Candida krusei</i>	
<i>Mycobacterium simiae</i>	<i>Candida tropicalis</i>	
<i>Mycobacterium celatum</i>	<i>Ajellomyces dermatitidis</i>	

## EXPECTED RESULTS

1,488 blood cultures obtained from patients suspected of mycobacterial, yeast or fungal infections were evaluated in the BACTEC Myco/F Lytic culture vial with the BACTEC 9240 Blood Culture System. There were 315 positive cultures with 243 clinically significant organisms recovered, of which 131 (53.9%) were mycobacteria, 35 (14.4%) were yeast or fungi and 77 (31.7%) were other bacteria. Of the 1,488 blood specimens tested in the clinical study, eleven BACTEC Myco/F Lytic culture vials (0.7%) were determined to be false positive (instrument-positive, smear and/or subculture-negative). Of the 315 instrument positive Myco/F Lytic vials, 11 (3.5%) were determined to be false positive. Of the 1,488 blood specimens tested in the clinical study, one (1) BACTEC Myco/F Lytic culture vial (0.07%) was determined to be false negative (instrument-negative, smear and/or subculture-positive). Of the 1,173 instrument negative BACTEC Myco/F Lytic culture vials, one (0.08%) was determined to be false negative. The contamination rate during this evaluation was 3.3%.

Frequency distribution of clinical trial specimens positive in the BACTEC Myco/F Lytic culture vials with the BACTEC 9000 Blood Culture System are illustrated in FIGURE 1 for times to detection (TTD) for mycobacteria and FIGURE 2 for TTD for yeast and fungi.

FIGURE 1



FIGURE 2

## PERFORMANCE CHARACTERISTICS

The BACTEC Myco/F Lytic medium was evaluated with the BACTEC 9240 instrument at two clinical sites considered large tertiary care teaching hospitals in geographically diverse areas. The site populations included patients infected with HIV, immunocompromised patients, transplant patients, and patients suspected of a mycobacterial infection. The BACTEC Myco/F Lytic medium was compared to the BACTEC 13A medium for the recovery and detection of mycobacteria from blood. A total of 1,100 blood culture specimens were tested during the evaluation. The total number of pathogenic mycobacteria isolates recovered in the study was 111 (See TABLE 1). Of these positives, ten (9%) were recovered in the BACTEC Myco/F Lytic medium only and three (3%) were recovered by BACTEC 13A medium only.

TABLE 1: SUMMARY OF MYCO/F LYTIC CULTURE MEDIUM ISOLATE RECOVERY DURING CLINICAL TRIALS

Organism	Total Isolates	Myco/F Lytic Medium Only	13A Medium Only	Both
<b>All Pathogenic Mycobacteria:</b>				
<i>Mycobacterium avium</i>	108	10	3	95
<i>Mycobacterium tuberculosis</i>	2	0	0	2
<i>Mycobacterium celatum</i>	1	0	0	1
<b>Total</b>	<b>111</b>	<b>10</b>	<b>3</b>	<b>98</b>

The BACTEC Myco/F Lytic medium was evaluated with the BACTEC 9240 instrument at four clinical sites considered large tertiary care teaching hospitals. The site populations included patients infected with HIV, immunocompromised patients, transplant patients, and patients suspected of a fungal infection. The BACTEC Myco/F Lytic medium was compared to the ISOLATOR™ System (Wampole Laboratories, Cranbrook, NJ) for the recovery and detection of yeast and fungi from blood. BACTEC Myco/F Lytic vials were inoculated with 1 – 5 mL of blood and ISOLATOR tubes were inoculated with 3 – 10 mL of blood. The ISOLATOR sediment was plated to Chocolate Agar, Brain Heart Infusion Agar with 5% sheep blood, and Sabouraud Dextrose Agar. A total of 748 specimens were tested during the evaluation. The total number of pathogenic yeast and fungal isolates recovered in the study was 32 (See TABLE 2). Of these positives, seven (22%) were recovered in the BACTEC Myco/F Lytic medium only and six (19%) were recovered in the ISOLATOR system only.

**TABLE 2: SUMMARY OF MYCO/F LYTIC MEDIUM ISOLATE RECOVERY DURING CLINICAL TRIAL**

Organism	Total Isolates	Myc/F Lytic Medium Only	Isolator Only	Both
<b>All Pathogenic Fungi:</b>				
<i>Candida albicans</i>	10	3	3	4
<i>Candida glabrata</i>	5	0	1	4
<i>Candida krusei</i>	2	2	0	0
<i>Candida parapsilosis</i>	1	1	0	0
<i>Candida tropicalis</i>	1	1	0	0
<i>Cryptococcus neoformans</i>	1	0	0	1
<i>Fusarium</i> species	1	0	1	0
<i>Histoplasma capsulatum</i>	11	0	1	10
Total	32	7	6	19

**AVAILABILITY**

Cat. No.	Description
442003	BACTEC™ Myco/F Lytic Culture Vials, case of 25 vials
442288	BACTEC™ Myco/F Lytic Culture Vials, case of 50 vials

**REFERENCES**

- Horsburg Jr., C.R. 1991. *Mycobacterium avium* Complex Infection in the acquired immunodeficiency syndrome. *New England Journal of Medicine* 324:1332-1338.
- Tenover, F.C., et al, 1993. The resurgence of Tuberculosis: Is Your Laboratory Ready? *Journal of Clinical Microbiology* 31:767-770.
- Isolation and Identification of *Mycobacterium tuberculosis*: A Guide for the Level II Laboratory. U.S. Department of Health and Human Services, Atlanta, GA, 1981.
- Kent, P.T., and G.P. Kubica. 1985. *Public Health Mycobacteriology: A Guide for the Level III Laboratory*. U.S. Department of Health and Human Services, Centers for Disease Control and Prevention, Atlanta, GA.
- Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17: 53-80.
- U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262, 17/10/2000, p. 0021-0045.*
- Biosafety in Microbiological and Biomedical Laboratories. U.S. Department of Health and Human Services Public Health Service/Centers for Disease Control, Atlanta, GA, May 1988 pp. 63-64.

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Benex Limited