Rev From	Rev To	DADS #	Date	Appr.
1274	0703	1953-03		

Notes


- BD Cat. No. 223091, 223101, 225851
- Blank (Sheet) Size : Length: 9" Width: 32
 Number of Pages: 16 Number of Sheets: 1
 Page Size: Length 9" Width 4" Final Folded Size: 2.25" x 4"
- Style (see illustrations below): #4


#1

#2

#3


#4


#5

#6

#7

- See Specification Control No. L0001620 for Material Information
- Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS #2755 Blue
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 250 Schilling Circle Cockeysville, MD. 21030-0243 USA	
Proofer	Date			
Checked By	Date			
Part Number:	S1220	Category and Description Package Insert, Difco Bordetella Antiserum and Antigen	Sheet: 1 of 17 Scale:	A

BD Difco™ Bordetella Antisera
Difco™ Bordetella Antigens
S1220
2003/07

English: pages 1 - 3 **Italiano: pagine 10 - 12**
Français: pages 4 - 6 **Español: páginas 12 - 15**
Deutsch: Seiten 7 - 9

See symbol glossary at end of insert. / Se symbolglossaret i slutningen af indlægsedlen. / Voir le glossaire des symboles à la fin de la notice. / Siehe Symbol-Erklärungen am Ende der Packungsbeilage. / Vedere il glossario dei simboli alla fine del foglio illustrativo. / Consulte o glossário de símbolos no fim do folheto informativo. / Consulte el glosario de símbolos al final del prospecto. / Se symbolförteckningen vid slutet av bipacksedeln.

Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Contacte o seu representante local da BD para obter instruções. / Kontakta lokal Becton Dickinson-representant för anvisningar.

INTENDED USE

Difco™ Bordetella Pertussis Antiserum and **Difco** Bordetella Parapertussis Antiserum are recommended for use in slide agglutination tests for the identification of *Bordetella pertussis* and *Bordetella parapertussis*.

Difco Bordetella Pertussis Antigen is used to demonstrate a positive quality control test in the slide agglutination test.

SUMMARY AND EXPLANATION

B. pertussis and *B. parapertussis* are two species of *Bordetella* that have been isolated from humans. These organisms adhere to, multiply among and remain localized in the ciliated epithelial cells of the respiratory tract. *B. pertussis* is the major cause of whooping cough or pertussis. *B. parapertussis* is associated with a milder, less frequently occurring form of the disease.¹ Data from molecular genetic studies show that these two species may be insufficiently different to be classified as separate species. Toxin production remains the major distinction between *B. pertussis* and *B. parapertussis*.

Classic pertussis is described as a respiratory tract infection in which the patient develops a cough with a characteristic "whoop." The term "whooping cough" was abandoned when it was determined that not all patients with clinical pertussis developed this characteristic cough.

Bordetella species are tiny, gram-negative, strictly aerobic coccobacilli occurring singly or in pairs and may exhibit a bipolar appearance. While some species are motile, *B. pertussis* and *B. parapertussis* are nonmotile. They do not produce acid from carbohydrates. *B. pertussis* will not grow on common blood agar bases or chocolate agar, while *B. parapertussis* will grow on blood agar and sometimes on chocolate agar. Media for primary isolation must include starch, charcoal, ion-exchange resins or a high percentage of blood to inactivate inhibitory substances.² *B. pertussis* may be recovered from secretions collected from the posterior nasopharynx, bronchoalveolar lavage and transbronchial specimens.

PRINCIPLES OF THE PROCEDURE

Identification of *Bordetella* species includes the isolation of the microorganism, biochemical identification and serological confirmation. Serological confirmation involves the reaction in which the microorganism (antigen) reacts with its corresponding antibody. The *in vitro* reaction produces macroscopic clumping called agglutination. The desired homologous reaction should be rapid, tightly binding (high affinity) and not dissociate (high avidity).

Because a microorganism (antigen) may agglutinate with an antibody produced in response to some other species, heterologous reactions are possible. These are characterized as forming slowly and being weak, and such unexpected and perhaps unpredictable reactions may lead to some confusion in serological identification. A positive homologous agglutination reaction should support the morphological and biochemical identification of the microorganism.

REAGENTS

Difco Bordetella Pertussis Antiserum and **Difco** Bordetella Parapertussis Antiserum are lyophilized, polyclonal rabbit antisera containing approximately 0.04% thimerosal as a preservative.

Difco Bordetella Pertussis Antigen is a suspension of a selected strain of *B. pertussis* containing approximately 0.04% thimerosal as a preservative.

When rehydrated and used as described, each vial of **Difco** Bordetella Pertussis and **Difco** Bordetella Parapertussis antisera diluted 1:10 contains sufficient reagent for approximately 200 slide tests.

When used as described, each vial of **Difco** Bordetella Pertussis Antigen contains sufficient reagent for approximately 140 slide tests.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

The packaging of this product contains dry natural rubber.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. After use, specimens, containers, slides, tubes and other contaminated material must be sterilized by autoclaving. Directions for use should be followed carefully.

Storage: Store lyophilized and rehydrated **Difco** Bordetella Pertussis Antiserum and **Difco** Bordetella Parapertussis Antiserum at 2–8°C.

Store **Difco** Bordetella Pertussis Antigen at 2–8°C.

Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.

Product Deterioration: The expiration date applies to the product in its intact container when stored as directed. Do not use if the product is caked, discolored or shows other signs of deterioration.

Difco Bordetella Pertussis Antigen is stable through the expiration date on the label when stored as described. Exposure to temperatures below 2°C can result in autoagglutination. Antigens must be smooth uniform suspensions. Examine antigen vials for agglutination before use. Suspensions with agglutination are not usable and should be discarded.

SPECIMEN COLLECTION AND PREPARATION

Isolation of *Bordetella* requires the use of certain media, such as Bordet-Gengou Agar. For specific recommendations consult appropriate references.^{2,3} Determine that a pure culture of the microorganism has been obtained and that biochemical test reactions are consistent with the identification of the organism as *Bordetella* spp. After these criteria are met, serological identification can be performed.

The isolate for serological testing should be subcultured from selective media to a nonselective agar.

PROCEDURE

Materials Provided: **Difco** Bordetella Pertussis Antiserum, **Difco** Bordetella Parapertussis Antiserum and **Difco** Bordetella Pertussis Antigen.

Materials Required But Not Provided: Agglutination slides with 1 inch squares, applicator sticks, sterile 0.85% saline, purified water and inoculating loops.

Reagent Preparation

Equilibrate all materials to room temperature before performing the tests. Ensure that all glassware and pipettes are clean and free of residues such as detergents.

To rehydrate **Difco** Bordetella Pertussis Antiserum and **Difco** Bordetella Parapertussis Antiserum, add 1 mL distilled or deionized water and rotate gently to completely dissolve the contents. Dilute rehydrated antiserum 1:10 with sterile 0.85% saline. Dilute only enough for 1–2 days testing requirements.

Difco Bordetella Pertussis Antigen is ready for use.

Test Isolate for Autoagglutination

1. From the test culture on nonselective media, transfer a portion of a loopful of growth to a drop (approximately 35 µL) of sterile 0.85% saline on a clean slide and emulsify the organism.
2. Rotate the slide for one min, then observe for agglutination.
3. If agglutination (autoagglutination) occurs, the culture is rough and cannot be tested. Subculture to nonselective agar, incubate and test the organism again as described in steps 1 and 2.
4. If no agglutination occurs, proceed with testing the organism.

Test Procedure

1. On an agglutination slide dispense 2 separate drops (approximately 35 µL, each) of the **Difco** Bordetella Antiserum to be tested. The first square is to be used for the test isolate and the second for the positive control.
2. **Negative control** – Dispense 1 drop of sterile 0.85% saline.
3. **Test Isolate** – From a solid agar medium, transfer a portion of a loopful of an isolated colony to the first drop of antiserum and a second portion to the negative control (0.85% saline).
4. **Positive control** – Add 1 drop of **Difco** Bordetella Pertussis Antigen or a small amount of a known culture of *B. pertussis* or *B. parapertussis* to the second drop of antiserum.
5. Mix each test and control serum reaction area with an applicator stick and rock gently for 1 min. Read for agglutination.

User Quality Control

At the time of use, test both positive and negative antigen controls to check performance of the antisera, techniques and methodology.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

- Read and record results as follows:
 - 4+ 100% agglutination; background is clear to slightly hazy
 - 3+ 75% agglutination; background is slightly cloudy
 - 2+ 50% agglutination; background is moderately cloudy
 - 1+ 25% agglutination; background is cloudy
 - No agglutination.
- The positive control should show 3+ or greater agglutination.
- The negative control should show no agglutination.
- For the test isolates, a 3+ agglutination within 1 min is a positive result.

LIMITATIONS OF THE PROCEDURE

- Correct interpretation of serological reactions depends upon culture purity, morphological characteristics and biochemical reactions that are consistent with identification of the microorganism as *Bordetella* species.
- Serological methods alone cannot identify the isolate as a *Bordetella* species.
- Excessive heat from external sources (hot bacteriological loop, burner flame, light source, etc) may prevent making a smooth suspension of the microorganism or cause evaporation or precipitation of the test mixture. False positive reactions may occur.
- Rough culture isolates do occur and will agglutinate spontaneously, causing agglutination of the negative control (autoagglutination). Smooth colonies must be selected and tested in serological procedures.
- Shake the **Difco** Bordetella Pertussis Antigen vial well before use to obtain a smooth, uniform suspension. Occasionally a *Bordetella* suspension may settle out during storage.
- Difco** Bordetella Pertussis Antigen will display irreversible autoagglutination if at any time during shipment or storage it is subjected to freezing temperatures. Do not allow to freeze.

PERFORMANCE CHARACTERISTICS

Sensitivity of **Difco** Bordetella Antisera is determined by demonstrating appropriate reactivity, as defined in the "Results" section, against a battery of well-characterized homologous *B. pertussis* and *B. parapertussis* cultures with known antigenic reactivity (refer to the table below). Specificity is determined by demonstrating non-reactivity against non-related (heterologous) *B. pertussis* and *B. parapertussis* cultures.

Catalog No.	Product	Homologous Cultures Tested
223091	Bordetella Pertussis Antiserum	<i>B. pertussis</i> 10536 MDH <i>B. pertussis</i> ATCC™ 12743 <i>B. pertussis</i> Malkiel NIH
223101	Bordetella Parapertussis Antiserum	<i>B. parapertussis</i> 32472 MDH <i>B. parapertussis</i> 170326 MDH <i>B. parapertussis</i> 17903 <i>B. parapertussis</i> Wolinski

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ATCC is a trademark of the American Type Culture Collection.
Difco is a trademark of Difco Laboratories, Inc., a subsidiary of
BD and BD Logo are trademarks of Becton, Dickinson and Company. © 2003 BD.
BENEX Limited