
Revisions

SO 0046-2


Rev From	Rev To	ECO #	Date	Appr.
1279	0803			

Notes


- BD Cat. No. 224301, 224311, 224321
- Blank (Sheet) Size : Length: 9" Width: 36"
 Number of Pages: 18 Number of Sheets: 1
 Page Size: Length 9" Width 4" Final Folded Size: 2.25" X 4"
- Style (see illustrations below): #4


#1

#2


#3


#4


#5

#6

#7

- See Specification Control No. S1279 for Material Information
- Ink Colors: Printed two sides Yes No
 No. of Colors: 1 PMS# 2755 (blue)
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 250 Schilling Circle Cockeysville, MD. 21030-0243 USA
Proofer	Date		
Checked By	Date		
Part Number:	S1279	Category and Description Package Insert Difco Vibrio Cholerae Antisera	Sheet: 1 of 19 Scale: 1:1 <div style="font-size: 2em; font-weight: bold; text-align: center;">A</div>

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CE S1279
2003/08

See symbol glossary at end of insert.

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INTENDED USE

Difco™ Vibrio Cholerae antisera are recommended for use in slide agglutination tests for the identification and serotyping of *Vibrio cholerae*.

SUMMARY AND EXPLANATION

V. cholerae is the causative agent of cholera. The spread of the disease is primarily through contaminated water and the fecal-oral route. If not treated, patients with severe cholera may die within five hours from massive fluid and electrolyte loss.¹ Because *Vibrio* spp. are natural inhabitants of aquatic environments, food products, such as uncooked seafoods or seafoods incorrectly handled, can spread infection. Gastrointestinal symptoms of cholera include "rice water stools" caused by a potent enterotoxin. From humans, the primary specimen is feces. Seafood products are frequently tested as vehicles of human infection.

Gardner and Ventkatraman² defined the identification of *Vibrio cholerae* as the recognition of the somatic antigen O1 (O group 1) as evidenced by agglutination with O1 (polyvalent) antiserum.

V. cholerae of the O1 serogroup are divided into serotypes: Ogawa, Inaba and Hikojima.^{1,3} The antigenic factors are:

Serotype	O Antigen Factors
Ogawa	AB
Inaba	AC
Hikojima	ABC

V. cholerae are facultative anaerobes that by microscopic morphology are gram-negative curved or straight bacilli. The microorganisms either require sodium chloride or grow best in its presence. Various special media designed to select and cultivate this microorganism are used in clinical and food microbiology laboratories and for environmental testing.

PRINCIPLES OF THE PROCEDURE

Identification of *V. cholerae* includes the isolation of the microorganism, biochemical identification and serological confirmation.

Serological confirmation involves the reaction in which the microorganism (antigen) reacts with its corresponding antibody. This *in vitro* reaction produces macroscopic clumping called agglutination. The desired homologous reaction is rapid, binds strongly (high affinity) and does not dissociate (high avidity).

Because a microorganism (antigen) may agglutinate with an antibody produced in response to some other species, heterologous reactions are possible. These are characterized as weak in strength and slow in formation. Such unexpected and perhaps unpredictable reactions may lead to some confusion in serological identification. Therefore, a positive homologous agglutination reaction should support the morphological and biochemical identification of the microorganism. Homologous reactions occur rapidly and are strong. Heterologous reactions form slowly and are weak.

REAGENTS

Difco Vibrio Cholerae antisera are stable, lyophilized polyclonal rabbit *Vibrio cholerae* O1 antisera containing approximately 0.04% Thimerosal as a preservative. Difco Vibrio Cholerae Antiserum Ogawa and Difco Vibrio Cholerae Antiserum Inaba are monospecific absorbed antisera.

Difco Vibrio Cholerae Antiserum Poly and the absorbed Difco Vibrio Cholerae Antiserum Ogawa and Difco Vibrio Cholerae Antiserum Inaba contain the following antibodies:

Antiserum	O Antibodies
Difco Vibrio Cholerae Antiserum Poly	ABC
Difco Vibrio Cholerae Antiserum Ogawa	B
Difco Vibrio Cholerae Antiserum Inaba	C

When reconstituted and used as described, each vial of **Difco** Vibrio Cholerae antisera contains sufficient reagent for 20 slide tests.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

The Packaging of This Product Contains Dry Natural Rubber.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. After use, specimens, containers, slides, tubes and other contaminated material must be sterilized by autoclaving. Directions for use should be followed carefully.

Storage: Store lyophilized and reconstituted **Difco** Vibrio Cholerae antisera at 2-8°C.

Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.

Reconstituted **Difco** Vibrio Cholerae Antiserum that is cloudy or has a precipitate at any time during its period of use should be discarded.

Product Deterioration

Expiration date applies to product in its intact container when stored as directed. Do not use if the product is caked, discolored or shows other signs of deterioration.

SPECIMEN COLLECTION AND PREPARATION

Clinical Specimens:

V. cholerae can be recovered from clinical specimens on selective media, such as TCBS (thiosulfate-citrate-bile salts-sucrose) agar. For specific recommendations, consult appropriate references.⁴⁻⁶ Determine that a pure culture of the microorganism has been obtained and that biochemical test reactions are consistent with the identification of the organism as *V. cholerae*. After these criteria are met, serological identification can be performed.

Food Specimens:

V. cholerae can be recovered from various types of foods, particularly seafood, when samples are processed to prevent overgrowth of competing microorganisms and selective media are used to enhance the growth of the microorganism. Some isolation media have been specifically developed for the food industry. Consult appropriate references for recommended procedures when testing food samples.^{7,8} After following an established protocol, determine that a pure culture of the microorganism has been obtained and that biochemical test reactions are consistent with the identification of the organism as *V. cholerae*. After these criteria are met, the serologic identification can be performed.

The isolate for serological testing should be subcultured from selective media to a nonselective agar, such as Nutrient Agar.

PROCEDURE

Materials Provided: **Difco** Vibrio Cholerae Antiserum Poly, **Difco** Vibrio Cholerae Antiserum Ogawa and **Difco** Vibrio Cholerae Antiserum Inaba.

Materials Required But Not Provided: Agglutination slides, applicator sticks and sterile 0.85% saline.

Reconstitution/Reagent Preparation:

Equilibrate all materials to room temperature before performing the tests. Ensure that all glassware and pipettes are clean and free of residues, such as detergents.

To reconstitute **Difco** Vibrio Cholerae antiserum add 3 mL of sterile 0.85% saline and rotate gently to completely dissolve the contents.

Test isolate for autoagglutination:

1. From the test culture on nonselective media, transfer a loopful of growth to a drop (approximately 35 µL) of sterile 0.85% saline on a clean slide and emulsify the organism.
2. Rotate the slide for 1 min, then observe for agglutination.
3. If agglutination (autoagglutination) occurs, the culture is rough and cannot be tested. Subculture to nonselective agar, incubate and test the organism again as described in steps 1 and 2.
4. If no agglutination occurs, proceed with testing the organism.

Test Procedure

Use the **Difco** Vibrio Cholerae Antiserum Poly to screen possible *V. cholerae* isolates.

Continue testing with **Difco** Vibrio Cholerae Antiserum Ogawa and **Difco** Vibrio Cholerae Antiserum Inaba.

1. Dispense a drop of the **Difco** Vibrio Cholerae Antiserum to be tested on an agglutination slide.
2. Transfer a loopful of growth of the test organism to the drop of antiserum and mix thoroughly.
3. Rotate the slide for 1 min and read for agglutination.

User Quality Control

At the time of use, test both positive and negative antigen controls to check performance of the antisera, techniques and methodology.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

1. Read and record results as follows:
 - 4+ 100% agglutination; background is clear to slightly hazy
 - 3+ 75% agglutination; background is slightly cloudy
 - 2+ 50% agglutination; background is moderately cloudy
 - 1+ 25% agglutination; background is cloudy
 - No agglutination
2. The positive control should show 3+ or greater agglutination.
3. The negative control should show no agglutination.
4. For the test isolates, a 3+ or greater agglutination within 1 min is a positive result.

INTERPRETATION OF RESULTS

1. Positive agglutination using **Difco** Vibrio Cholerae Antiserum Poly with typical biochemical test results gives presumptive identification of *V. cholerae* O1.
2. Cultures with positive agglutination in **Difco** Vibrio Cholerae Antiserum Poly may be serotyped using the **Difco** Vibrio Cholerae Antiserum Ogawa and **Difco** Vibrio Cholerae Antiserum Inaba. Positive agglutination in both antisera is rare and, when it occurs, is usually interpreted as identifying *V. cholerae* serotype Hikojima.⁹ *V. cholerae* serotype Hikojima is a rare serotype and should be sent to a reference laboratory for further study.
3. Positive agglutination will be immediate and strong. The strongest and most rapid reaction should be used to identify the serotype, because *V. cholerae* O1 strains frequently cross-react slowly or weakly in monospecific antiserum for the other serotype.
4. Isolates that weakly or slowly agglutinate with **Difco** Vibrio Cholerae Antiserum Poly and not with **Difco** Vibrio Cholerae Antiserum Inaba or **Difco** Vibrio Cholerae Antiserum Ogawa are usually considered *V. cholerae* non-O1. The isolate may be sent to a reference laboratory for further study.

LIMITATIONS OF THE PROCEDURE

1. Correct interpretation of serological reactions depends upon culture purity, morphological characteristics and biochemical reactions that are consistent with identification of the microorganism as *V. cholerae*.
2. Serologic methods alone cannot identify the isolate as *V. cholerae*.
3. Excessive heat from external sources (hot bacteriological loop, burner flame, light source, etc.) may prevent a smooth suspension of the microorganism or cause evaporation or precipitation of the test mixture. False positive reactions may occur.
4. Rough culture isolates do occur and will agglutinate spontaneously causing agglutination of the negative control (autoagglutination). Smooth colonies must be selected and tested in serological procedures.
5. **Difco** Vibrio Cholerae antisera have been tested using undiluted cultures taken from agar media. These antisera have not been tested using antigen suspensions in saline or other diluents. If the user applies variations in the recommended steps each lot of antiserum must be tested with known control cultures to verify expected reactions under the modified procedure.
6. Reconstituted **Difco** Vibrio Cholerae Antiserum that is cloudy or has a precipitate at any time during its period of use should be discarded.

PERFORMANCE CHARACTERISTICS

Sensitivity of **Difco** Vibrio Cholerae antisera is determined by demonstrating appropriate reactivity, as defined in the "Results" section, against a battery of well-characterized homologous *Vibrio cholerae* cultures with known antigenic reactivity (refer to the table below). Specificity is determined by demonstrating non-reactivity against non-related (heterologous) *Vibrio cholerae* cultures.

Cat. No.	Product	Homologous Cultures Tested
224321	Vibrio Cholerae Antiserum Poly	Hikojima 18 Hikojima 18ss Ogawa Bu93 Ogawa 1321 CDC Ogawa 20-A-11 Ogawa 11 CDC Ogawa VRL 628 Inaba 6 Inaba 11 CDC Inaba 35 A3NIH Inaba 6 CDC
224311	Vibrio Cholerae Antiserum Ogawa	Ogawa Bu93 Ogawa 1321 CDC Ogawa 20-A-11 Ogawa 11 CDC Ogawa VRL 628
224301	Vibrio Cholerae Antiserum Inaba	Inaba 6 Inaba 11 CDC Inaba 35 A3NIH Inaba 6 CDC

AVAILABILITY

Cat. No.	Description
224321	Difco™ Vibrio Cholerae Antiserum Poly, 3 mL
224311	Difco™ Vibrio Cholerae Antiserum Ogawa, 3 mL
224301	Difco™ Vibrio Cholerae Antiserum Inaba, 3 mL

REFERENCES

- World Health Organization. 2001. Weekly Epidemiol. Record. 76: 117-124.
- Gardner, A.D., and K.V. Venkatraman. 1935. The antigen of the cholera group of vibrios. J. Hyg. 35:262-282.
- Nobechi, K. 1923. Contributions to the knowledge of *V. cholerae*. 3. Immunological studies upon the types of *V. cholerae*. Sci. Repr. Govt. Inst. Infect. Dis. 2:43-87.
- Farmer III, J.J., J.M. Janda and K. Birkhead. 2003. *Vibrio*, p. 706-718. In P.R. Murray, E.J. Barron, J.H. Jorgensen, M.A. Pfaller and R.H. Tenover (ed.), Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
- Pezzlo, M. 1992. Aerobic bacteriology, p.1.0.1-1.20.47. In H. D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
- Forbes, B.A., D.F. Sahm, A.S. Weissfeld and E.A. Trevino. 2002. Bailey & Scott's diagnostic microbiology, 11th ed. Mosby, Inc., St. Louis, Mo.
- Kaysner, C.A., and A. DePaola, Jr. 2001. *Vibrio*, p. 405-420. In F.P. Downes and K. Ito (ed.), Compendium of methods for the microbiological examination of foods, 4th ed. American Public Health Association, Washington, D.C.
- AOAC International. 2001. FDA bacteriological analytical manual online. <<http://www.cfsan.fda.gov/~ebam/bam-mm.html>>.
- Kelly, M.T., F.W. Hickman-Brenner and J.J. Farmer III. 1991. *Vibrio*, p. 384-395. In A. Ballows, W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg and H.J. Shadomy (ed.), Manual of clinical microbiology, 5th ed. American Society of Microbiology, Washington, D.C.

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EC REP BENEX Limited