**Revisions**

Rev from	Rev to	ECO #
0707	E	5099-09

**Notes:**

- BD Cat. Number 240950
- Blank (Sheet) Size: Length: N/A      Width: N/A  
 Number of Pages: N/A      Number of Sheets: N/A  
 Page Size: Length N/A      Width N/A      Final Folded Size: N/A
- Style (see illustrations below): # N/A

- See Specification Control Number 240950 for Material Information
- Ink Colors: Printed two sides  Yes  No  
 No. of Colors: N/A      PMS# N/A
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: X5424		Category and Description Package Insert, BBL Streptocard Enzyme Latex Test	Sheet: 1 of 31 Scale: N/A	<b>A</b>

# **BD BBL™ Streptocard™ Enzyme Latex Test**

**A latex agglutination test for the identification of streptococcal groups A, B, C, D, F and G.**

X5424E  
2009/04

English:	pages	1 – 6	Italiano:	pagine	16 – 21
Français :	pages	6 – 11	Español:	páginas	21 – 26
Deutsch:	Seiten	11 – 16	Flow Chart		30

Pokyny vám poskytnú miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselététől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinizi danışın. / Обратитесь к своему локальному представителю компании BD за указания. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіңіздің жергілікті БД өкіліне жүгініп нұсқау алыңыз.

## INTENDED USE

The **BBL™ Streptocard™ Enzyme Latex Test** is a latex test system for the qualitative identification of Lancefield streptococcal groups A, B, C, D, F and G.

## SUMMARY AND EXPLANATION

Lancefield<sup>1</sup> showed that the majority of pathogenic streptococci possess specific carbohydrate antigens, which permit the classification of streptococci into groups. These streptococcal group antigens are extracted from the streptococcal cell wall in a liquid form, and reacted with group specific antibodies.

## PRINCIPLES OF THE PROCEDURE

The **BBL Streptocard** Enzyme Latex Test latex particles are sensitized with group specific antibody and will agglutinate in the presence of homologous antigen. In the absence of such antigen, the latex particles will remain in a smooth suspension. The use of a patented enzymatic extraction in the **BBL Streptocard** Enzyme Latex Test procedure considerably shortens the time required for antigen extraction and improves the antigen yield.

## REAGENTS

Test Latex A	1 x 2.5 mL	
Test Latex B	1 x 2.5 mL	Test Latex A, B, C, D, F and G consist of blue latex particles sensitized with rabbit antibody to appropriate group specific antigen, suspended in buffer containing 0.1% sodium azide (preservative).
Test Latex C	1 x 2.5 mL	
Test Latex D	1 x 2.5 mL	
Test Latex F	1 x 2.5 mL	
Test Latex G	1 x 2.5 mL	
Extraction Enzyme	2 x 2.1 g	When reconstituted with 12 mL of distilled water, contains 0.01% thimerosal (Lyophilized) as a preservative. The Extraction Enzyme is also used as the Negative Control.
Control +	1 x 2.5 mL	Positive Control contains a mixture of extracted antigen from streptococcal groups A, B, C, D, F and G with 0.1% sodium azide as a preservative.
Reaction Cards	50	Disposable; 6 reaction circles per card.
Mixing Sticks	250	Disposable.

## Warnings and Precautions:

For *in vitro* Diagnostic Use.

**Caution:** This product contains natural rubber latex which may cause allergic reactions.

Do not use test components beyond the expiration date.

Observe established precautions against microbiological hazards throughout all procedures. Extraction Enzyme does not always render bacteria nonviable. After use, contaminated materials must be sterilized by autoclaving.

Do not allow reagents to become contaminated by allowing the dropper tip to touch the specimen on the reaction card. Ensure caps are securely fitted on reagent bottles after each use to prevent contamination and drying out of the reagents.

**WARNING:** Extraction Enzyme contains Mercury. Don't put in trash. Recycle or dispose as hazardous waste.

Reagents contain sodium azide. Very toxic by inhalation, in contact with skin, and if swallowed. Contact with acids liberates very toxic gas. After contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

**CARDS:** Care should be taken not to finger-mark the test areas, since this may result in an oily deposit and improper test results.

## Storage:

**LATEX REAGENTS:** Reagents are ready for use. Vials should be stored at 2 – 8°C. **DO NOT FREEZE.** Under these conditions the reagents will retain their activity until the date shown on the bottle labels. After use return the kit to the refrigerator storing bottles in an upright position.

**EXTRACTION ENZYME:** The lyophilized **BBL™** Extraction Enzyme should be stored at 2 – 8°C. Under these conditions it will retain its activity until the expiration date shown on the bottle labels. The reconstituted **BBL** Extraction Enzyme when stored at 2 – 8°C will retain its activity for three months. However, the reconstituted **BBL** Extraction Enzyme must be used within the expiration date on the vial.

**CONTROL +:** Ready for use; store at 2 – 8°C. Under these conditions activity will be retained until the date shown on the bottle label.

**Product Deterioration:** Do not use the kit if the Control + does not yield appropriate results. Refer to “User Quality Control.” Examine the Control + and **BBL** Extraction Enzyme for evidence of contamination, evaporation or other signs of deterioration, such as turbidity. Examine the Test Latex, and do not use if not a homogenous suspension.

**Preparation of Cultures:** Samples for identification should be grown on a blood agar plate 16 – 24 h at 35 ± 2°C. Note the hemolytic reaction of suspect colonies. It is also advisable to perform a Gram stain and catalase test to confirm the presence of gram-positive, catalase-negative cocci. For further details, please consult standard references.<sup>2</sup>

**Reconstitution:** Prior to use, reconstitute the lyophilized vial of **BBL** Extraction Enzyme with 12 mL of distilled water. Agitate vial until extraction enzyme is completely dissolved. Store at 2 – 8°C. Use within 3 months of reconstitution, but prior to the expiration date on the vial.

## PROCEDURES

Review “Warnings and Precautions,” “Product Deterioration,” “Preparation of Cultures” and “Test Method” prior to performing procedures. The test area, all reagents, test specimens and test components should be at room temperature (15 – 30°C) when used.

**Materials Provided:** All materials as listed under “Reagents” and accessories.

**Materials Required But Not Provided:** Pipette (to measure 0.4 mL), microbiological loop, Pasteur pipettes, glass or plastic tubes and water bath or heating block at 37 ± 2°C. Also required are the necessary equipment and

labware used for preparation, storage and handling of serologic specimens.

**Test Method:** If fewer than six tests are to be performed, the card may be cut with scissors and the unused portion saved for future use.

1. Label a test tube appropriately and dispense 0.4 mL of **BBL** Extraction Enzyme into the tube for each specimen to be tested.
2. Select 2 – 5 similar colonies with a microbiological loop and emulsify in the **BBL** Extraction Enzyme. If the culture is mixed, avoid obvious contamination. If the colonies are small, use more than five and ensure that at least a slightly turbid suspension is obtained.
3. Incubate the tube for a total of 10 min at  $37 \pm 2^{\circ}\text{C}$  in a water bath or heat block (see “Limitations”). After 5 min incubation remove each tube and mix by shaking for 2 – 3 sec, then continue the incubation. Remove and allow to cool to room temperature. The extract is now ready for use.
4. Ensure that the Test Latex has warmed to room temperature. Make sure the Test Latex suspensions are mixed by shaking and expel any Test Latex from the dropper for complete mixing.
5. Dispense 1 drop from each Test Latex to be tested onto a separate circle on the reaction card.
6. Using a Pasteur pipette add 1 drop of extract to each of the six test circles.
7. With the mixing sticks provided, spread the mixture over the entire area of the circle, using a separate stick for each.
8. Gently rock the card manually for up to 1 min and observe for agglutination under normal lighting conditions. Read macroscopically: do not use magnification to aid reading.
9. Dispose of the reaction card in an appropriate biohazard container.

**User Quality Control:** Quality control requirements must be performed in accordance with applicable local, and/or federal regulations or accreditation requirements and your laboratory’s standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices. Initially, upon receipt, the laboratory should check each shipment or lot of material prior to use to verify the performance of the product.

Streptococci of known group reactivity should be subjected to the complete test procedure. The performance of the test is assessed by the presence of agglutination in one latex suspension only, with the other five suspensions showing a negative (no agglutination) reaction for each reference strain tested. This will evaluate both the efficiency of the extraction procedure and the specificity of each reagent.

The positive control procedure and negative control procedure should be tested each day of use. Local, regional, or other laboratory regulations may apply which supersede package insert directions for frequency of testing the positive and/or negative control.

**POSITIVE CONTROL** – Shake each Test Latex and dispense 1 drop onto a separate circle on the card. Dispense 1 drop of Control + onto each of six circles. Spread each mixture over the entire area of the circle, using a separate stick for each Test Latex. Rock the card manually for 1 min. Each of the six Test Latex should demonstrate obvious agglutination.

**NEGATIVE CONTROL** – Shake each Test Latex and dispense 1 drop onto a separate circle on the card. Dispense 1 drop of reconstituted **BBL** Extraction Enzyme onto each of the six circles. Spread each mixture over the entire area of the circle, using a separate stick for each Test Latex. Rock the card manually for 1 min. No obvious agglutination should be evident for any Test Latex suspension.

The **BBL™ Streptocard™** Positive Control (Cat. No. 240965) contains a mixture of extracted antigen from streptococcal groups A, B, C, D, F and G with sodium azide as a preservative. This control may also be used to perform quality control when using the individual test latex reagents (see “Availability”).

## **INTERPRETATION OF TEST RESULTS**

**Positive Result:** A positive result is obtained if obvious agglutination of the blue latex particles occurs within 1 min with a single Test Latex. Any weaker reaction which occurs in the presence of a substantially stronger positive reaction should be ignored.

**Negative Result:** A negative result is obtained if no agglutination occurs within 1 min. Faint traces of granular material may be observed in negative reactions and should be ignored.

**Uninterpretable Results:** If more than one Test Latex strongly agglutinates, then the possibility exists that a mixed culture of streptococcal groups are present. Examine the plate and carefully select organisms of like morphology and retest. Subculture, if the suspected organism is overgrown or insufficient. If the reaction pattern is unaltered, reisolate the organism or perform additional biochemical tests.

Refer to Flow Charts 1 and 2 (pages 30 – 31), "Recommended Scheme for Grouping Streptococci." When carrying out a serological identification of streptococci, the following initial observations should be made: (i) note hemolysis, (ii) note cell morphology, (iii) assess colonial growth for purity and quality.

- a. Rule out *Streptococcus pneumoniae*. This streptococcus is  $\alpha$ -hemolytic, bile soluble and optochin susceptible. Other streptococci are not bile soluble and are optochin resistant.<sup>2</sup>
- b. Aerococci are non- $\beta$ -hemolytic, grow in 6.5% NaCl broth and give variable reactions in the bile-esculin test and pyrrolidonyl arylamidase (PYR) test. They can be differentiated from enterococci by their arrangements in tetrads or as single cells, whereas enterococci are arranged as diplococci or short chains.<sup>3</sup>
- c. Staphylococci and *Listeria monocytogenes* are  $\beta$ -hemolytic and can be distinguished from streptococci by their cellular morphology and positive catalase reaction.<sup>4,5</sup>

#### LIMITATIONS OF THE PROCEDURE

False negative results can occur if an inadequate amount of the culture is used for extraction. Some streptococci, notably group F, produce minute colonies. When this occurs, use more colonies to prepare extract.

Nearly all the  $\beta$ -hemolytic streptococci isolated from human infections possess specific carbohydrate antigens, which can be recognized by serological reactivity. Attempts to extend these procedures to non- $\beta$ -hemolytic streptococci have been unsuccessful except for groups B and D.<sup>2</sup>

*Streptococcus pneumoniae* share common antigenic determinant(s) with group C streptococci<sup>6</sup> and therefore may react with Test Latex C. *S. pneumoniae* colonies are typically  $\alpha$ -hemolytic.

Certain strains of *Streptococcus milleri* possess A, C, F or G antigens and may therefore react with one or more of these Test Latex. *S. milleri* typically form minute and usually non-hemolytic colonies on blood agar plates. Identification of *S. milleri* may be performed using a scheme such as that described by Lawrence et al.<sup>7</sup>

*Streptococcus porcinus*, which is usually associated with swine, may react with Test Latex B. Differentiation from group B *Streptococcus* may be based on *S. porcinus* typically giving more pronounced zones of hemolysis and a positive PYR test.<sup>8</sup>

A water bath or heating block should be used to ensure adequate heating during the extraction procedure. Heat transfer in an air incubator may be insufficient. However, an equilibrated beaker of water within a 37°C air incubator will suffice.

Group D streptococci may require further biochemical tests to distinguish between the enterococcal and non-enterococcal species. In such cases the bile-esculin test and 6.5% NaCl test can be used to differentiate. Refer to Flow Charts 1 and 2 (pages 30 – 31).

The group D Test Latex may not easily react with some *S. bovis* strains. These strains may require further testing for identification.

Strains have been found which appear to have both D and G antigen.<sup>9,10</sup>

## PERFORMANCE CHARACTERISTICS

A total of 555 clinical isolates and stock cultures of streptococci were tested with the **BBL Streptocard** Enzyme Latex Test and compared to a "commercially available latex assay." Of these 193 were group A, 146 group B, 37 group C, 31 group D, 3 group F, 70 group G, 75 of the isolates did not belong to any of these groups and gave negative reactions with all 6 Test Latex (Table 1).

Table 1

Group	No. of Strains	BBL™ Streptocard™ Enzyme Latex Test	Commercial Latex Test	Correlation
A	193	193	193	100%
B	146	146	146	100%
C	37	37	37	100%
D	31	31	30	96.8%
F	3	3	3	100%
G	70	70	70	100%
Non A,B,C,D,F,G	75	75	75	100%

In a further study of 22 known group D streptococci strains also possessing group G antigen, the **BBL Streptocard** Enzyme Latex Test correctly identified all 22 strains as group D.

In a study of the **BBL™ Streptocard™** Enzyme Latex Test with 55 non-streptococcal strains from 14 different genera, one strain of *Escherichia coli* gave a positive reaction with the group C Test Latex. All other reactions were negative.

## AVAILABILITY

### Cat. No. Description

240950	<b>BBL™ Streptocard™</b> Enzyme Latex Kit, 50 Tests.
240951	<b>BBL™ Streptocard™</b> Acid Latex Kit, 50 Tests.
240960	<b>BBL™ Streptocard™</b> Test Latex A, one 2.5 mL bottle.
240961	<b>BBL™ Streptocard™</b> Test Latex B, one 2.5 mL bottle.
240865	<b>BBL™ Streptocard™</b> Test Latex C, one 2.5 mL bottle.
240962	<b>BBL™ Streptocard™</b> Test Latex D, one 2.5 mL bottle.
240866	<b>BBL™ Streptocard™</b> Test Latex F, one 2.5 mL bottle.
240867	<b>BBL™ Streptocard™</b> Test Latex G, one 2.5 mL bottle.
240964	<b>BBL™ Streptocard™</b> Extraction Enzyme, one 2.1 g bottle.
240965	<b>BBL™ Streptocard™</b> Positive Control, one 2.5 mL bottle.
240966	<b>BBL™ Streptocard™</b> Test Cards, box of 50.

## REFERENCES

1. Lancefield R.C. 1938. Proc. Soc. Exp. Bio. Med. 38: 473-478.
2. Ruoff, K.L. 1995. *Streptococcus*, p. 299-307. In Murray, P.R., Baron, E.J., Tenover, F.C. and Tenover, R.H. (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
3. Ruoff, K.L. 1995. *Leuconostoc*, *Pediococcus*, *Stomatococcus*, and Miscellaneous Gram-Positive Cocci that grow Aerobically, p. 315-323. In Murray, P.R., Baron, E.J., Tenover, F.C. and Tenover, R.H. (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
4. Kloos, W.E. and Bannerman, T.L. 1995. *Staphylococcus* and *Micrococcus*, p. 282-298. In Murray, P.R., Baron, E.J., Tenover, F.C. and Tenover, R.H. (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
5. Swaminathan, B., Rocourt, J. and Bille, J. 1995. *Listeria*, p. 341-348. In Murray, P.R., Baron, E.J., Tenover, F.C. and Tenover, R.H. (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
6. Lee, P.C. and Wetherall, B.L. 1987. J. Clin. Microbiol. 25: 152-153.
7. Lawrence, J., Yaiko, D.M. and Hadley, W.K. 1985. J. Clin. Microbiol. 22: 772-777.
8. Thompson, T., and Facklam, R. 1997. J. Clin. Microbiol. 35: 1885-1886.
9. Birch, B.R., Keaney, M.G.L. and Ganguli, L.A. 1984. J. Clin. Pathol. 37: 1289-1292.
10. McMurray, M.D. 1984. Lancet II: 1353.

---

---

# **BD BBL Streptocard Enzyme Latex Test**

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BENEX Limited