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REFERENCES

1. Bopp, C.A., Brenner, F.W., Fields, P.I., Wells, J.G., and N.A. Strockbine. 2003. Escherichia, Shigella, and Salmonella. In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
2. Muller, L. 1923. Un nouveau milieu d'enrichissement pour la recherche du bacille typhique et des paratyphiques. C. R. Soc. Biol. 89:434. Paris.
3. Kauffmann, F. 1930. Ein kombiniertes Anreicherungsverfahren für Typhus- und Paratyphusbacillen. Zentralb. Bakteriol. Parasitenkde. Infektionskr. Hyg. Abt. I Orig. 113:148.
4. Kauffmann, F. 1935. Weitere Erfahrungen mit den kombinierten Anreicherungsverfahren für Salmonellabacillen. Z. Hyg. Infektionskr. 117:26.
5. Jones, F. T., R. C. Axtell, D. V. Rives, S. E. Scheideler, F. R. Tarver, Jr., R. L. Walker, and M. J. Wineland. 1991. A survey of Salmonella contamination in modern broiler production. J. Food Prot. 54:502-507.
6. Barnhart, H. M., D. W. Dressen, R. Bastien, and O. C. Pancorbo. 1991. Prevalence of Salmonella enteritidis and other serovars in ovaries of layer hens at time of slaughter. J. Food Prot. 54:488-492.
7. Eckner, K. F., W. A. Dustman, M. S. Curiale, R. S. Flowers, and B. J. Robison. 1994. Elevated-temperature, colorimetric, monoclonal, enzyme-linked immunosorbent assay for rapid screening of Salmonella in foods: collaborative study. J. Assoc. Off. Anal. Chem. 77:374-383.
8. Andrews, W. H., G. A. June, P. S. Sherrod, T. S. Hammack, and R. M. Amaguana. 1995. Salmonella. p 5.01-5.20. In Bacteriological analytical manual, 8th ed. AOAC International. Gaithersburg, MD.
9. Russell, S. F., J.-Y. D'Aoust, W. H. Andrews, and J. S. Bailey. 1992. Salmonella, p.371-422. In Vanderzant, C. and D. F. Splittstoesser (ed.). Compendium of methods for the microbiological examination of food, 3rd ed. American Public Health Association, Washington, D.C.
10. Flowers, R. S., W. Andrews, C. W. Donnelly, and E. Koenig. 1993. Pathogens in milk and milk products, p. 103-212. In R. T. Marshall (ed.) Standard methods for the examination of dairy products. 16th ed., American Public Health Association, Washington, D.C.
11. United States Pharmacopeial Convention. 1995. The United States Pharmacopeia, 23rd ed. The United States Pharmacopeial Convention. Rockville, MD.
12. Federal Register. 1991. Animal and plant health inspection service: chicken affected by Salmonella enteritidis, final rule. Fed. Regist. 56:3730-3743.
13. MacFaddin, J. D. 1985. Media for isolation-cultivation-identification- maintenance of medical bacteria, p. 751-754, Williams & Wilkins, Baltimore, MD.
14. Isenberg, H. D. (ed.). 1992. Clinical microbiology procedures handbook, vol. 1, American Society for Microbiology, Washington, D. C.
15. Kist, M. et al. 2000. Infektionen des Darmes. In: Mauch, H., R. Lüttiken, and S. Gattermann (ed.). MIQ – Qualitätsstandards in der mikrobiologisch-infektiologischen Diagnostik. Vol 9. Urban und Fischer, München, Jena.
16. Knox, R., P. H. Gell, and M. R. Pollack. 1942. Selective media for organisms of the Salmonella group. J. Pathol. Bacteriol. 54:469-483.
17. Hinsley, A.P., and B.C. Berks. 2002. Specificity of respiratory pathways involved in the reduction of sulfur compounds by Salmonella enterica. Microbiology 148: 3631-3638.
18. Hammack, T.S., et al. 1999. Relative effectiveness of selenite cystine broth, tetrathionate broth, and Rappaport-Vassiliadis medium for the recovery of Salmonella spp. from foods with a low microbial load. J. Food Prot. 62: 16-21.

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BD and BD logo are trademarks of Becton, Dickinson and Company.
Bacto is a trademark of Difco Laboratories, division of Becton, Dickinson and Company.
ATCC is a trademark of the American Type Culture Collection.
© 2003 Becton, Dickinson and Company