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BIBLIOGRAFIA

1. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, PA.
2. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17: 53-80.
3. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
4. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
5. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry and M.A. Pfaller. 2007. Manual of clinical microbiology. 9th ed. American Society for Microbiology, Washington, D.C.
6. Gleaves, C.A., R.L. Hodinka, S.L.G. Johnston, and E.M. Swierkosz. 1994. Cumitech 15A. Laboratory diagnosis of viral infections. ASM, Washington, DC.
7. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2002. Bailey and Scott's diagnostic microbiology. 11th ed. Mosby, St. Louis, MO.
8. Wardford, A., M. Chernesky, and E. M. Peterson. 1999. Cumitech 19A, Laboratory diagnosis of *Chlamydia trachomatis* infections. ASM, Washington, DC.
9. Miller, J. M. 1999. A guide to specimen management in clinical microbiology, 2nd ed. ASM, Washington, DC.
10. Isenberg, H. D., 2004. Clinical microbiology procedures handbook, 2nd ed. ASM, Washington, DC.
11. Isenberg, H.D., 1998. Essential procedures for clinical microbiology. Chapter 14.12, Page 787. Packaging and shipping infectious substances.
12. Maass, M. and K. Dalhoff. 1995. Transport and storage conditions for cultural recovery of *Chlamydia pneumoniae*. J. Clin. Microbiol. 33:1793-1796.
13. Johnson, F. 2005. Transport of viral specimens. Clin. Microbiol. Rev 3:120-121.
14. 42CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. Interstate Shipment of Etiologic Agents.
15. Clinical and Laboratory Standards Institute. 2003. Quality control of microbiological transport systems. Approved Standard M40-A, CLSI, Wayne, PA.
16. Clinical and Laboratory Standards Institute. 2006. Viral culture; Approved Guideline M41-A, CLSI, Wayne, PA.

Do not reuse / Nepoužívejte opakovaně / Må ikke genbruges / Niet opnieuw gebruiken / Mitte kasutada korduvalt / Ei saa kasutada uudele / Usage unique / Nicht wiederverwenden / Μην το ξαναχρησιμοποιείτε / Egyszer használatos / Non riutilizzare / Tik vienkartiniam naudojimui / Må ikke gjenbrukes / Nie stosować powtórnie / Não reutilizar / Nepouzívajte opakovane / No reusar / Får ej återanvändas / Не използвайте отново / A nu se reutiliza / Tekrar kullanmayın / Ne upotrebljavajte ponovo / Не использовать повторно / Пайдаланбаңыз / Ne koristiti ponovo

STERILE

Method of sterilization: irradiation / Způsob sterilizace: záření / Sterilisationsmåde: Bestråling / Sterilisatiewijze: bestraling / Steriliseerimismeetod: kiirgus / Steriloitimenetelmä: säteilytys / Méthode de stérilisation: irradiation / Sterilisationsmethode: Bestrahlung / Μέθοδος αποστείρωσης: ακτινοβολία / Sterilizálás módszere: besugárzás / Metodo di sterilizzazione: irradiazione / Sterilizavimo būdas: radiacija / Steriliseringsmetode: bestrålning / Metoda sterylizacji: napromienianie / Método de esterilização: irradiação / Metóda sterilizácie: ožiarenie / Método de esterilización: irradiación / Steriliseringsmetod: strålning / Метод на стерилизация: ирадиация / Metodă de sterilizare: iradiere / Sterilizasyon yöntemi: ırradyasyon / Metoda sterilizacije: označavanje / Метод стерилизации: облучение / Стерилизация әдісі – сәуле түсімі / Metoda sterilizacije: zračenje (iridacija)

Manufacturer / Výrobce / Producent / Fabrikant / Tootja / Valmistaja / Fabricant / Hersteller / Κατασκευαστής / Gyártó / Ditta produttrice / Gamintojas / Producent / Fabricante / Výrobca / Tillverkare / Производител / Producător / Üretici / Proizvođač / Производител / Аткарушы

Use by / Spotřebujte do / Användes för / Houdbaar tot / Kasutada enne / Viimeinkäyttöpäivä / A utiliser avant / Verwendbar bis / Ημερομηνία λήξης / Felhasználhatóság dátuma / Usare entro / Naudokite iki / Brukes før / Stosować do / Utilizar em / Použite do / Usar antes de / Använd före / Използвайте до / A se utiliza până la / Son kulanma tarihi / Upotrebiti do / Использовать до / дейін пайдалануға / Upotrijebiti do / YYYY-MM-DD / YYYY-MM (MM = end of month) / RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutning af måned) / JJJJ-MM-DD / JJJJ-MM (MM = einde maand) / ÁÁÁÁ-KK-PP / ÁÁÁÁ-KK (KK = kuu lõpp) / VVVV-KK-PP / VVVV-KK (kuukauden loppuun mennessä) / ÁÁÁÁ-MM-JJ / ÁÁÁÁ-MM (MM = fin du mois) / JJJJ-MM-TT / JJJJ-MM (MM = Monatsende) / EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα) / ÉÉÉÉ-HH-NN / ÉÉÉÉ-HH (HH = hónap utolsó napja) / ÁÁÁÁ-MM-GG / ÁÁÁÁ-MM (MM = fine mese) / MMMM-MM-DD / MMMM-MM (MM = mêsio pabaiga) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutten av måneden) / RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = fim do mês) / RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca) / áaaa-mm-dd / áaaa-mm (mm = fin del mes) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutet på månaden) / ГГГГ-MM-ДД / ГГГГ-MM (MM = края на месеца) / ÁÁÁÁ-LL-ZZ / ÁÁÁÁ-LL (LL = sfârșitul lunii) / YYYY-AA-GG / YYYY-AA (AA = ayın sonu) / GGGG-MM-DD / GGGG-MM (MM = kraj meseca) / ГГГГ-MM-ДД / ГГГГ-MM (MM = конец месяца) / ЖЖЖЖ-АА-КК / ЖЖЖЖ-АА (АА = айдың соңы) / GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)

REF

Catalog number / Katalogové číslo / Katalognummer / Catalogusnummer / Catalogoogi number / Tuotenumero / Numéro catalogue / Bestellnummer / Αριθμός καταλόγου / Katalógusszám / Numero di catalogo / Katalogo numeris / Numer katalogowy / Número do catálogo / Katalogové číslo / Número de catálogo / Каталоген номер / Număr de catalog / Katalog numarası / Kataloški broj / Номер по каталогу / Каталог нөмірі

EC REP

Authorized Representative in the European Community / Autorizovaný zástupce pro Evropskou unii / Autoriseret repræsentant i EU / Erkend vertegenwoordiger in de Europese Unie / Volitatud esindaja Euroopa Nõukogus / Valtuutettu edustaja Euroopan yhteisössä / Représentant agréé pour la C.E.E. / Autorisierte EG-Vertretung / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Hivatalos képviselő az Európai Unióban / Rappresentante autorizzato nella Comunità europea / Įgaliotasis atstovas Europos Bendrijoje / Autoriseret representant i EU / Autorizowane przedstawicielstwo w Unii Europejskiej / Reprezentante autorizado na União Europeia / Autorizovaný zástupca v Európskom spoločenstve / Representante autorizado en la Comunidad Europea / Auktoriserad representant i EU / Оторизован представител в EU / Repräsentant autorizat în Uniunea Europeană / Автура Преплугуу Yetkilil Temsilcisi / Ovlašteni predstavnik u Evropskoj zajednici / Уполномоченный представитель в Европейском сообществе / Европа қауымдастығындағы уәкілетті екіні / Autorizuirani predstavnik u EU

IVD

In Vitro Diagnostic Medical Device / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medisch hulpmiddel voor in vitro diagnose / In vitro diagnostiska meditsiniaparatuur / Lääkinnällinen in vitro -diagnostikkalite / Dispositif médical de diagnostic in vitro / Medizinisches In-vitro-Diagnostikum / In vitro διαγνωστική ιατρική συσκευή / In vitro diagnostikai orvosi eszköz / Dispositivo medico diagnostico in vitro. / In vitro diagnostikos prietais / In vitro diagnostisk medisinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Medicinska pomůcka na diagnostiku in vitro / Dispositivo médico de diagnóstico in vitro / Medicinsk anordning för in vitro-diagnostik / Mediciniskni uredj za diagnostika in vitro / Aparatură medicală de diagnosticare in vitro / In Vitro Djanostik Tibbi Cihaz / Medicinski uredaj za in vitro dijagnostiku / Медицинский прибор для диагностики in vitro / Жасанды жагдайда жүргізетін медициналық диагностика аспабы / Medicinska pomagala za In Vitro Dijagnostiku

LOT

Batch Code (Lot) / Kód (číslo) šarže / Batch code (Lot) / Chargennummer (lot) / Partii kood / Eräkoodi (LOT) / Code de lot (Lot) / Chargencode (Chargenbezeichnung) / Κωδικός παρτίδας (Παρτίδα) / Tétel száma (Lote) / Codice del lotto (Partita) / Partijos numeris (Lot) / Batch-code (Serie) / Kod partii (seria) / Código do lote (Lote) / Kód série (šarža) / Código de lote (Lote) / Satskod (parti) / Код (Партида) / Număr lot (Lotul) / Parti Kodu (Lot) / Код партии (лот) / Топтама коды / Lot (kod)

Contains sufficient for <n> tests / Dostatečné množství pro <n> testů / Indeholder tilstrækkeligt til <n> test / Voldoende voor <n> tests / Küllaldane <n> testide jaoks / Sisältöön riittävä <n> testejä varten / Contenu suffisant pour <n> tests / Ausreichend für <n> Tests / Περιέχει επαρκή ποσότητα <n> εξετάσεις / <n> teszthez elegendő / Contenuto sufficiente per <n> test / Pakankamas kiekis atlikti <n> testų / Innholder tilstrekkelig for <n> tester / Zawiera ilość wystarczającą do <n> testów / Contêmo suficiente para <n> testes / Obsah vystačí na <n> testov / Contenido suficiente para <n> pruebas / Räckertill <n> antal tester / Съдържанието е достатъчно за <n> теста / Conține suficient pentru <n> teste / <n> testleri için yeterli miktarda içerir / Sadržaj dovoljan za <n> testova / Достаточно для <n> тестов(а) / <n> тесттері үшін жеткілікті / Sadržaj za (n) testova

Consult Instructions for Use / Prostudujte pokyny k použití / Læs brugsanvisningen / Raadpleeg gebruiksaanwijzing / Lugeda kasutusjuhendit / Tarkista käyttöohjeista / Consulter la notice d'emploi / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Skaitykite naudojimo instrukcijas / Se i bruksanvisningen / Zobacz instrukcja użytkowania / Consulte as instruções de utilização / Pozri Pokyny na používanie / Consultar las instrucciones de uso / Se bruksanvisningen / Направление справка в инструкциите за употреба / Consultați instrucțiunile de utilizare / Kullanım Talimatları'na başvurun / Pogledajte upute za upotrebu / См. руководство по эксплуатации / Пайдалану нұсқаулығымен танысып алыңыз / Koristi upute za upotrebu

Do not use if package damaged / Nepoužívejte, je-li obal poškozený / Må ikke anvendes hvis emballager beskadiget / Niet gebruikindie verpakking beschadigd is / Mitte kasutada, kui pakend on kahjustatud / Ei saa kasutada, jos pakkaus on vahingoittunut / Ne pas utiliser si l'emballage est endommagé / Inhal beschädigter Packung nicht verwenden / Μην το χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιό / Ne használja, ha a csomagolás sérült / Non usare se la confezione è danneggiata / Jei pakuotė pažeista, nenaudoti / Må ikke brukes hvis pakker skadet / Nie używać, jeśli opakowanie jest uszkodzone / Não usar se a embalagem estiver danificada / Nepoužívajte, ak je obal poškozený / No usar si el paquete está dañado / Använd ej om förpackningen är skadad / Не използвайте, ако опаковката е повредена / A nu se folosi dacă pachetul este deteriorat / Ambalaj hasar görmüşse kullanmayın / Ne koristite ako je pakovanje oštećeno / Не использовать при повреждении упаковки / Ерер пакет бұзылган болса, пайдаланба / Ne koristiti ako je oštećeno pakiranje

Becton, Dickinson and Company
7 Loveton Circle
Sparks, Maryland 21152 USA
800-638-8663
www.bd.com/ds

Becton, Dickinson France S.A.S.
11 rue Aristide Bergès
38800 Le pont de Claix, France
Tel: +33 (0) 476 68 36 36

Made in Italy by Copan

Bacti Disks is a trademark of Remel, Inc.

BD and BD Logo are trademarks of Becton, Dickinson and Company © 2010 BD.