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1. Loesche, W.J. 1969. Oxygen sensitivity of various anaerobic bacteria. *Appl. Microbiol.* 18:723-727.
2. Chow, A.W., R.D. Leake, T. Yamauchi, B.F. Anthony, and L.B. Guze. 1974. The significance of anaerobes in neonatal bacteremia: analysis of 23 cases and review of the literature. *Pediatrics* 54:736-745.
3. Finegold, S.M., V.L. Sutter, H.R. Attebery, and J.E. Rosenblatt. 1974. Isolation of anaerobic bacteria, p. 365-375. *In* E.H. Lennette, E.H. Spaulding, and J.R. Truant (ed.), *Manual of clinical microbiology*, 2nd ed. American Society for Microbiology, Washington, D.C.
4. Abramson, I.J., G.L. Evans, and P. Bathurst. 1977. Laboratory evaluation of Port-A-Cul with stock cultures of anaerobes, aerobes, and facultative organisms, abstr. C 197, p. 68. *Abstr. 77th Annu. Meet. Am. Soc. Microbiol.* 1977.
5. Mena, E., F.S. Thompson, A.Y. Armfield, V.R. Dowell, Jr., and D.J. Reinhardt. 1978. Evaluation of Port-A-Cul transport system for protection of anaerobic bacteria. *J. Clin. Microbiol.* 8:28-35.
6. Hamilton, R.L., G.L. Evans, and T.E. O'Neill. 1979. Survival of facultatively anaerobic bacteria in Port-A-Cul tubes and an aerobic transport device, abstr. C 173, p. 338. *Abstr. 79th Annu. Meet. Am. Soc. Microbiol.* 1979.
7. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
8. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
9. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
10. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021-0045.

11. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenkay (ed.). 2003. Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.

	Manufacturer / Výrobce / Producent / Fabrikant / Tootja / Valmistaja / Fabricant / Hersteller / Κατασκευαστής / Gyártó / Ditta produttrice / Gamintojas / Producent / Fabricante / Výrobca / Tillverkare / Производител / Producător / Üretici / Proizvođač / Производител / Атқарушы
	Use by / Spotřebujte do / Användes för / Houdbaar tot / Kasutada enne / Viimeinkäyttöpäivä / A utiliser avant / Verwendbar bis / Ημερομηνία λήξης / Felhasználhatóság dátuma / Usare entro / Naudokite iki / Brukes før / Stosować do / Utilizar em / Použite do / Usar antes de / Använd före / Используйте до / A se utiliza până la / Son kulanma tarihi / Upotrebiti do / Исползовать до / дейін пайдалануға / Upotrijebiti do / YYYY-MM-DD / YYYY-MM (MM = end of month) / RRRR-MM-DD / RRRR-MM (MM = konec měsíce) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutning af måned) / JJJJ-MM-DD / JJJJ-MM (MM = einde maand) / AAAA-KK-PP / AAAA-KK (KK = kuu lõpp) / VVVV-KK-PP / VVVV-KK (kuukauden loppuun mennessä) / AAAA-MM-JJ / AAAA-MM (MM = fin du mois) / JJJJ-MM-TT / JJJJ-MM (MM = Monatsende) / EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα) / ÉÉÉÉ-HH-NN / ÉÉÉÉ-HH (HH = hónap utolsó napja) / AAAA-MM-GG / AAAA-MM (MM = fine mese) / MMMM-MM-DD / MMMM-MM (MM = mensesio pabaiga) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = sluttet av måneden) / RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca) / AAAA-MM-DD / AAAA-MM (MM = fin du mès) / RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca) / aaaa-mm-dd / aaaa-mm (mm = fin del mes) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutet på månaden) / ГГГГ-MM-ДД / ГГГГ-MM (MM = края на месеца) / AAAA-LL-ZZ / AAAA-LL (LL = sfârșitul lunii) / YYYY-AA-GG / YYYY-AA (AA = ayin sonu) / GGGG-MM-DD / GGGG-MM (MM = kraj meseca) / ГГГГ-MM-ДД / ГГГГ-MM (MM = конец месяца) / ЖЖЖЖ-АА-КК / ЖЖЖЖ-АА (АА = айдың соңы) / GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)
	Catalog number / Katalogové číslo / Katalognummer / Catalogusnummer / Kataloogi number / Tuotenumero / Numéro catalogue / Bestellnummer / Αριθμός καταλόγου / Katalógusszám / Numero di catalogo / Katalogo numeris / Numer katalogowy / Número do catálogo / Katalogové číslo / Número de catálogo / Каталоген номер / Număr de catalog / Katalog numerasi / Kataloški broj / Номер по каталогу / Каталог нөмірі
	Authorized Representative in the European Community / Autorizovaný zástupce pro Evropskou unii / Autoriseret repræsentant i EU / Erkend vertegenwoordiger in de Europese Unie / Volitadud esindaja Euroopa Nõukogus / Valtuutettu edustaja Euroopan yhteisössä / Représentant agréé pour la C.E.E. / Autorisierte EG-Vertretung / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Hivatalos képvislet az Európai Unióban / Rappresentante autorizzato nella Comunità europea / Įgaliojasis atstovas Europos Bendrijoje / Autoriseret representant i EU / Autoryzowane przedstawicielstwo w Unii Europejskiej / Representante autorizado na União Europeia / Autorizovaný zástupca v Európskom spoločenstve / Representante autorizado en la Comunidad Europea / Auktoriserad representant i EU / Оторизирован представитель в EU / Reprezentant autorizat în Uniunea Europeană / Avrupa Topluluğu Yetkili Temsilcisi / Ovlašteni predstavnik u Evropskoj zajednici / Уполномоченный представитель в Европейском сообществе / Европа қауымдастығындағы уәкілетті өкіл / Autorizuirani predstavnik u EU
	In Vitro Diagnostic Medical Device / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medisch hulpmiddel voor in vitro diagnose / In vitro diagnostika meditsiniaparatuur / Lääkinnällinen in vitro -diagnostiikkalaitte / Dispositif médical de diagnostic in vitro / Medizinisches In-vitro-Diagnostikum / In vitro διαγνωστική ιατρική συσκευή / In vitro diagnosztikai orvosi eszköz / Dispositivo medico diagnostico in vitro. / In vitro diagnostikos prietais / In vitro diagnostisk medisinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Medicinska pomôcka na diagnostikę in vitro / Dispositivo médico de diagnóstico in vitro / Medicinsk anordning för in vitro-diagnostik / Медицински уред за диагностика ин витро / Aparatură medicală de diagnosticare in vitro / In Vitro Diagnostik Tibbi Cihaz / Medicinski uređaj za in vitro dijagnostiku / Медицинский прибор для диагностики ин витро / Жасанды жағдайда жүргізілетін медициналық диагностика аспабы / Medicinska pomagala za In Vitro Dijagnostiku
	Temperature limitation / Teplotní omezení / Temperaturbegrænsning / Temperatuurlimiet / Temperatuuri piirang / Lämpötilarajoitus / Température limite / Zulässiger Temperaturbereich / Οπο έρροοκρσσις / Hőmérsékleti határ / Temperatura limite / Laikymo temperatūra / Temperaturbegrænsning / Ograniczenie temperatury / Limitação da temperatura / Ohraničenje teploty / Limitación de temperatura / Temperaturbegrænsning / Температурни ограничения / Limitare de temperatură / Sicaklık sinirlaması / Ograničenje temperature / Ограничение температуры / Температураны шектеу / Dozvoljena temperatura
	Batch Code (Lot) / Kód (číslo) řarže / Batch kode (Lot) / Chargennummer (lot) / Partii kood / Eräkoodi (LOT) / Code de lot (Lot) / Chargencode (Chargenbezeichnung) / Кодифікаційний номер / Кодифікаційний номер (Партида) / Tétel száma (Lot) / Codice del lotto (partita) / Partijos numeris (Lot) / Batch-kode (Serie) / Kod partii (serie) / Código do lote (Lote) / Kód řerie (řarža) / Código de lote (Lote) / Satskod (parti) / Код (Партида) / Număr lot (Lotul) / Parti Kodu (Lot) / Kod řerie / Код партии (not) / Топтама коды / Lot (kod)
	Consult Instructions for Use / Prostudujte pokyny k použití / Læs brugsanvisningen / Raadpleeg gebruiksaanwijzing / Lugeda kasutusjuhendit / Tarkista käyttöohjeista / Consulter la notice d'emploi / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Skaitykite naudojimo instrukcijas / Se i bruksanvisningen / Zobacz instrukcja użytkowania / Consulte as instruções de utilização / Pozri Pokyny na používanie / Consultar las instrucciones de uso / Se bruksanvisningen / Направте справка в инструкциите за употреба / Consultați instrucțiunile de utilizare / Kullanim Talimatları'na bagvurun / Pogledajte uputstvo za upotrebu / См. руководство по эксплуатации / Пайдалану нұсқаулығымен танысып алыңыз / Koristi upute za upotrebu
	Do not reuse / Nepoužívejte opakovaně / Må ikke genbruges / Niet opnieuw gebruiken / Mitte kasutada korduvalt / Ei saa käyttää uudelleen / Usage unique / Nicht wiederverwenden / Μην το ξαναχρησιμοποιείτε / Egyszer használatos / Non riutilizzare / Tik vienkartiniam naudojimui / Må ikke gjenbrukes / Nie stosować powtórnie / Não reutilizar / Nepoužívať opakovane / No reusar / Får ej återanvändas / Не използвайте отново / A nu se reutiliza / Tekrar kullanmayın / Ne upotrebljavajte ponovo / Не использовать повторно / Пайдаланбаңыз / Ne koristiti ponovo
	Method of sterilization: ethylene oxide / Způsob sterilizace: etylenoxid / Sterilisationsmåde: Etylenoxid / Sterilisatiemethode: ethyleenoxide / Steriliseerimismeetod: etüleenoksiid / Sterilointimenetelmä: etyleeniksiidi / Méthode de stérilisation : oxyde d'éthylène / Sterilisationsmethode: Etylenoxid / Μέθοδος αποστείρωσης: αιθυλενοξείδιο / Sterilizálás módszere: etilén-oxid / Metodo di sterilizzazione: ossido di etilene / Sterilizavimo būdas: etileno oksidas / Steriliseringsmetode: etylenoksid / Metoda sterylizacji: tielen etylu / Método de esterilização: óxido de etileno / Metóda sterilizácie: etilénoksid / Método de esterilización: óxido de etileno / Steriliseringsmetod: etylenoxid / Метод на стерилизация: етиленов оксид / Metodă de sterilizare: oxid de etilenă / Sterilizasyon yöntemi: etilen oksit / Metoda sterilizacije: etilen oksid / Метод стерилизации: этиленоксид / Стерилизация әдісі – этилен тотығы

 Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA
800-638-8663
www.bd.com/ds

 Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

BD, BD Logo, BBL, Luer-Lok and Port-A-Cul are trademarks of Becton, Dickinson and Company. © 2013 BD.