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XIV REFERENCE

1. Brewer, J.H. 1940. Clear liquid medium for the "aerobic" cultivation of anaerobes. *J. Am. Med. Assoc.* 115:598-600.
2. Vera, H.D. 1944. A comparative study of materials suitable for the cultivation of clostridia. *J. Bacteriol.* 47:59-70.
3. U.S. Pharmacopeial Convention, Inc. 2009. The U.S. pharmacopeia 32/The national formulary 27. U.S. Pharmacopeial Convention, Inc., Rockville, Md.
4. MacFaddin, J.F. 1985. Media for isolation-cultivation-identification-maintenance of medical bacteria, vol. I. Williams & Wilkins, Baltimore.
5. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
6. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
7. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
8. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
9. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry and M.A. Pfaller (ed.) 2007. Manual of clinical microbiology, 9th ed. American Society for Microbiology, Washington, D.C.
10. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey & Scott's diagnostic microbiology, 12th ed. Mosby, Inc., St. Louis.
11. Horwitz, W. (ed.). 2000. Official methods of analysis of AOAC International, 17th ed, vol.1. AOAC International, Gaithersburg, Md.
12. Reischelderfer and Mangels. 1992. In Isenberg (ed.), Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
13. Holt, J.G., N.R. Krieg, P.H.A. Sneath, J.T. Staley, and S.T. Williams (ed.). 1994. Bergey's Manual of determinative bacteriology, 9th ed. Williams & Wilkins, Baltimore.

Tehnički servis i podrška na BD Diagnostics: obratite se lokalnom predstavniku tvrtke BD ili posjetite www.bd.com/ds.

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo, and all other trademarks are property of Becton, Dickinson and Company. © 2015 BD