

Rev from	Rev to	JOB #
02	03	8548-16

NOTES:

1. BD Catalog Number: 261192, 261193
2. Blank (Sheet) Size: Length: 11" Width: 22.5"
3. Number of Pages: 10 Number of Sheets: 1
4. Page Size: Length: 11" Width: 4.5" Final Folded Size: 4.5" x 1.875"
5. Ink Colors: No. of Colors: 3 PMS#: 2755 Blue, 032 Red, Standard Black
6. Printed two sides: Yes No
7. Style (see illustrations below): # 4

8. Vendor Printed Online/In House Printed Web
9. See specification control no. D11826RM for material information.
10. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	REVISED BY <small>By Mary Schmidt at 3:29 pm, May 31, 2016</small>	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	PROOFING APPROVED BY <small>By Helaine Jones at 3:37 pm, May 31, 2016</small>		
Checked By	THIRD EYE BY <small>By Kim Stone at 8:19 am, Jun 21, 2016</small>		
Part Number:	L001126	Category and Description Package Insert, Voges-Proskauer Reagent Droppers	Sheet: 1 of 11 Scale: N/A

A

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteys lähipään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítás kérje a BD helyi képviseletétől. / Нұсқаулар үшін жерпілкіті BD екілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukiteis vietas BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o representante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získeť u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontaktta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasla geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

BD BBL™ Voges-Proskauer Reagent Droppers are intended for use in the Voges-Proskauer test, one of several qualitative tests used to distinguish between members of the *Enterobacteriaceae*.¹

SUMMARY AND EXPLANATION

The Voges-Proskauer test is one of the series known as the IMViC tests. The letters stand for Indole, Methyl Red, Voges-Proskauer, and Citrate. (The “i” is inserted for euphony.) They are used primarily to distinguish between the coliforms, but may be used for other organisms in the *Enterobacteriaceae* family. Specifically, the Voges-Proskauer test detects the presence of acetyl methyl carbinol. This compound is produced by certain microorganisms during growth in a buffered peptone-glucose broth, specifically MR-VP broth.¹

PRINCIPLES OF THE PROCEDURE

Upon addition of Voges-Proskauer reagents A and B to a culture tube of MR-VP broth, any acetyl methyl carbinol present will be oxidized to diacetyl methyl carbinol. The latter compound reacts with creatine (from the broth) to form a red-colored compound. The presence of this compound is detected visually and is considered a positive test.¹

REAGENTS

BD BBL™ Voges-Proskauer Reagent A Droppers contain 0.5 mL of 5% wt/vol alpha-naphthol in absolute alcohol.

BD BBL™ Voges-Proskauer Reagent B Droppers contain 0.5 mL of 40% wt/vol potassium hydroxide in distilled water.

Warnings and Precautions:

For *in vitro* Diagnostic Use

Follow established laboratory procedures in handling and disposing of infectious materials.

Cat. No. Description

261192 **BD BBL™** Voges-Proskauer A Reagent Droppers, 50.

Danger

H225 Highly flammable liquid and vapour. **H318** Causes serious eye damage. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P303+P361+P353** IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P370+P378** In case of fire: Use for extinction: CO₂, powder or water spray. **P403+P235** Store in a well-ventilated place. Keep cool. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Cat. No. Description

261193 **BD BBL™** Voges-Proskauer B Reagent Droppers, 50.

Danger

H302 Harmful if swallowed. **H314** Causes severe skin burns and eye damage. **P103** Read label before use. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P303+P361+P353** IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage Instructions: Store at room temperature 15 – 30 °C (59 – 86 °F).

Product Deterioration: Reagents are hermetically sealed in an ampule, which affords protection of the solution from chemical instability until expiration date. The reagent droppers need no refrigeration. Opened reagent droppers should be disposed of within 24 h.

PROCEDURE

Material Provided: Voges-Proskauer Reagent Droppers.

Materials Required But Not Provided: Ancillary culture media, quality control organisms and laboratory equipment as required for these procedures.

Test Procedure

1. Hold dropper upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to break ampule inside the dropper. **Caution: Break ampule close to its center one time only. Do not manipulate the dropper any further as the plastic may puncture and injury may occur.** Tap bottom on tabletop a few times.
2. Invert for drop-by-drop dispensing according to the test procedure performed.

A. Standard Method (Barritt's Modification)

Inoculate buffered peptone-glucose broth (MR-VP medium) with a pure culture of the organism in question. Incubate the culture at 37 °C for no less than 48 h.

Empty contents (15 drops) from reagent A dropper and 5 drops from reagent B dropper into 1 mL of broth culture. Shake well after the addition of each reagent to aerate sample.

B. Other Methods

Voges-Proskauer Reagent Droppers can be used to achieve highly accurate tests with the commercially prepared bacterial identification systems. Follow the manufacturer's instructions for use.

User Quality Control

Using either of the methods described above (Test Procedure), an authentic culture of *Enterobacter cloacae* or *Serratia marcescens* will yield a positive result, while *Escherichia coli* or *Salmonella enteritidis* will yield a negative result.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

A positive reaction is indicated by the development of a distinct red color which occurs within 5 min. The development of a copper color in some tests should be disregarded.

LIMITATIONS OF THE PROCEDURE

One should allow at least 15 min for color to develop before considering the Voges-Proskauer test negative.

PERFORMANCE CHARACTERISTICS

Barritt² evaluated two Voges-Proskauer (VP) formulations: one containing α-naphthol and another containing creatine. Three hundred sixty eight (368) strains were tested by both formulations. Two hundred fifty nine (259) strains were negative by both methods and were identified as *E.coli*. Thirty one (31) strains were positive by both methods and were identified as *E. aerogenes*. Seven (7) strains gave a negative reaction with the VP reagent containing creatine and a positive reaction with the VP reagent containing α-naphthol. Reactions of these strains on EMB Agar along with methyl red and citrate tests identified all of these strains as *E. aerogenes*, indicating that the formulation with α-naphthol was more sensitive. the remaining 71 strains known as the "intermediate group," showed the greatest discrepancy between the VP reagent containing creatine and the VP reagent containing α-naphthol. All 71 of these strains were methyl red and citrate positive. Twenty-one (21) of the 71 strains (29%) gave a positive result with the VP reagent containing creatine while 54 strains (76%) gave a positive result with the VP reagent containing α-naphthol. This result enhanced the observation that the VP reagent containing α-naphthol was more sensitive than the VP reagent containing creatine and increased the correlation between acetyl methyl carbinol production and citrate utilization, which further helped to distinguish *E. coli* from *E. aerogenes*.

AVAILABILITY

Cat. No. Description

261192 BD BBL™ Voges-Proskauer A Reagent Droppers, 50.

261193 BD BBL™ Voges-Proskauer B Reagent Droppers, 50.

REFERENCES

1. Lennette, E.H., A. Balows, W.J. Hausler, Jr., and H.J. Shadomy (ed.) 1985. Manual of clinical microbiology, 4th ed. American Society for Microbiology, Washington, D.C.
2. Barritt, M.M. 1936. The intensification of the Voges-Proskauer reaction by the addition of α-naphthol. J. Pathol. Bacteriol. 42:441-454.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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Benex Limited
Becton Dickinson Pty Ltd.