Part Number:	L001203	BALTSO0191 Version 13.0 Template 4 Inserts	
Category and Description:	Package Insert, India Ink Reagent Droppers	Rev from: 02	Rev to: 03
		Job Number:	788-18

Catalog Number: 261194

Blank (Sheet) Size: Length: 11" Width: 22.5"


Number of Pages: 10 Number of Sheets: 1

Page Size: Length: 11" Width: 4.5" Final Folded Size: 4 1/2" X 1 7/8"

Ink Colors: Number of Colors: 1 PMS #: 2755 Blue

Printed Two Sides: Yes: ☒ No: ☐

Style (see illustrations below): # 4

Vendor Printed: ☒ Online / In House Printed: ☐ Web Printed: ☐

See Specification control no. N/A for material information.

BD Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Company confidential. This document is the property of Becton, Dickinson and Company and is not to be used outside the company without written permission. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Revised By:	REVISED BY By Tori Pagani at 3:16 pm, Oct 25, 2018
Proofing Approved By:	PROOFING APPROVED BY By Nichole Graham at 7:20 pm, Oct 31, 2018
Third Eye By:	THIRD EYE BY By Kim Stone at 9:16 am, Nov 01, 2018

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INTENDED USE

BD BBL™ India Ink Reagent Droppers are used to enhance the microscopic detection of *Cryptococcus* spp. in wet preparations.

SUMMARY AND EXPLANATION

India ink facilitates the visualization of cryptococcal polysaccharide capsules. The capsular material of cryptococci displaces the colloidal carbon particles of the ink so that the capsule appears as a clear halo around the microorganisms against a black background.¹

BD BBL India Ink Reagent Droppers may be used to stain pus, exudate, tissue, sputum and sediment of centrifuged urine and cerebral spinal fluid (CSF) specimens.² The India ink smear is particularly useful for CSF specimens because *Cryptococcus neoformans* is the most common cause of fungal central nervous system (CNS) infections throughout the world.³ The appearance of encapsulated yeast forms in an India ink smear from CSF is diagnostic for cryptococcal meningitis.⁴

PRINCIPLES OF THE PROCEDURE

India ink is a negative stain resulting in a dark background and unstained organisms. Cryptococci are spherical (sometimes oval) yeast forms which range from 4 µm to 20 µm in diameter. The cryptococcal polysaccharide capsule is present on nearly all *Cryptococcus* spp. and is resistant to coloration by India ink.

REAGENTS

BD BBL India Ink Reagent Droppers contain 0.5 mL of an aqueous suspension of India ink with suspending agent and preservative.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Follow proper laboratory procedures in handling and disposing of infectious materials.

Storage Instructions: Store at controlled room temperature 15–30 °C (59–86 °F). Each dropper is good for one day's use after breaking ampule.

Product Deterioration: Reagents are hermetically sealed in an ampule which affords protection of the solution from chemical instability until expiration date. Do not use if extraneous material or organisms are present in stain control slide (see "User Quality Control"). Do not use after the expiration date.

PROCEDURE

Material Provided: BD BBL India Ink Reagent Droppers.

Materials Required But Not Provided: Ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

Test Procedure

Note: CSF, and other body fluids, should be centrifuged at 1,500 x *g* for 15 min prior to preparation of slide if greater than 1 mL of specimen is available.²

1. Hold reagent dropper upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to break the ampule inside the dropper. **Caution: Break ampule close to its center one time only. Do not manipulate dropper any further as the plastic may puncture and injury may occur.**
2. Tap bottom of dropper on tabletop a few times. Then invert for convenient drop-by-drop dispensing of reagent.
3. To a small loopful of sterile saline, water or broth on a precleaned microscope slide, add a small amount of growth from a young agar culture or specimen sediment, or add a drop of CSF directly to the slide.
4. Add one drop of India ink and mix.
5. Cover with a thin, large coverslip. Press gently to allow the fluid to spread as a thin film beneath the coverglass.
6. Scan slide under low power and examine under high power or oil immersion for encapsulated yeast forms.

User Quality Control

1. Examine the reagent for signs of deterioration (see "Product Deterioration").
2. Positive (a 48-h culture of *Cryptococcus neoformans* ATCC™ 36556), negative (*Candida albicans* ATCC 10231) and stain (1 drop India ink + 1 drop sterile saline) controls should be run simultaneously with the organism or specimen being tested.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

In positive smears, cryptococcal capsules appear as clear, refractile halos either surrounding the yeast cells or positioned between the well-defined cell wall and the surrounding black mass of India ink particles. Capsules may be broad or narrow. The yeast cells may be round, oval, or elongate; buds may be absent, single, or, rarely, multiple. The buds may be detached from the mother cell but enclosed in a common capsule.

LIMITATIONS OF THE PROCEDURE

- 1. Care must be taken not to misinterpret fat droplets or lymphocytes as yeast forms. Fat droplets will not have a well-defined cell wall. Lymphocytes have a shaggy border and an eccentric nucleus. The capsules of cryptococci are usually sharply outlined and the well-defined yeast cell wall is centrally located within the capsule.³
- 2. Organisms in AIDS patients may not have capsules.¹
- 3. India ink smears performed on CSF sediment are positive in only 60% of all cases of cryptococcal meningitis.³ All requests for India ink smears should be supplemented with culture or antigenic detection.³
- 4. *Rhodotorula* spp. and rare *Toryopsis* spp. will produce capsules.⁴
- 5. Positive India ink smears may be used as a presumptive diagnosis only. Definitive diagnosis must be accomplished by culture and/or antigenic detection.^{1,3}
- 6. If smear is too dark, dilute the India ink 1:2 with sterile, distilled water and make a new smear. Specimens treated with 20% KOH are not to be used for India ink preparations.

PERFORMANCE CHARACTERISTICS

Cryptococcus neoformans is a causative agent of meningitis. Several studies ⁵⁻¹⁰ point to the utility of the India Ink stain as one of the tools for diagnosing *C. neoformans*. these articles show that the India Ink stain is useful in identifying *C. neoformans* in cerebrospinal fluid (CSF) samples.

AVAILABILITY

Cat. No.	Description
261194	BD BBL™ India Ink Reagent Droppers, 50.

REFERENCES

1. Baron, E.J., and S.M. Finegold. 1990. Bailey & Scott's diagnostic microbiology, 8th ed. The C.V. Mosby Co., St. Louis.

2. Kwon-Chung, K.J., and J.E. Bennett. 1992. Medical mycology. Lea & Febiger, Philadelphia.

3. Dalton, H.P., and H.C. Nottebart, Jr. 1986. Interpretive medical microbiology. Churchill Livingstone, New York.

4. Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy (ed.). 1991. Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.

5. McGinnis, M.R.: Detection of fungi in cerebrospinal fluid. Am. J. Med.; July 28, 1983: 129-138.

6. Wang, W., and A.R. Carm: Clinical manifestations of AIDS with cryptococcal meningitis. China Med J (Engl); Aug, 2001: 114(8): 841-843.

7. Monaco, L.S., Nieto F. Silva, E. Warley, M.R. Cervelli: AIDS-associated meningeal cryptococcosis in the Hospital Diego Paroissien from 1996-1999. Rev. Argent. Microbiol. (Rev); Apr-June, 2001: 33(2): 118-121.

8. Koneman, E.W.: Should an India Ink test be used to detect cryptococcal meningitis in a low-risk population – primarily in a general acute care setting – or are better stat tests available? CAP Today; Dec, 2000: 14(12): 62.

9. Imwidthaya, P. and N. Pongvarin: Cryptococcosis in AIDS. Postgrad. Med. J.; Feb, 2000: 76(892): 85-88.

10. Menezes, E.A., M.N. Monteiro, M.R. Angelo, C.D. Santos, C.C. Freire and F.A. Cunha: *Cryptococcus neoformans* causing meningitis in AIDS patients. Rev. Soc. Bras. Med. Trop; Sep-Oct, 2002: 35(5): 537-539.

Technical Information: In the United States contact BD Technical Service and Support at 800.638.8663 or www.bd.com.

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ATCC is a trademark of the American Type Culture Collection.
© 2018 BD. BD, the BD Logo and BBL are trademarks of Becton, Dickinson and Company.
Benex Limited
Becton Dickinson Pty Ltd.