Part Number:	L001206	BALTS00191 Version 13.0 Template 4 Inserts	
Category and Description: Package Insert, Calcofluor White Reagent Droppers	Rev from: 03 Rev to: 04		
	Job Number: 788-18		

Catalog Number: 261195

Blank (Sheet) Size: Length: 11" Width: 27"


Number of Pages: 12 Number of Sheets: 1

Page Size: Length: 11" Width: 4.5" Final Folded Size: 4 1/2" X 1 7/8"

Ink Colors: Number of Colors: 1 PMS #: 2755 Blue

Printed Two Sides: Yes: ☒ No: ☐

Style (see illustrations below): # 4

Vendor Printed: ☒ Online / In House Printed: ☐ Web Printed: ☐

See Specification control no. N/A for material information.

BD Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Company confidential. This document is the property of Becton, Dickinson and Company and is not to be used outside the company without written permission. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Revised By:	REVISED BY By Maria Zacharias at 3:35 pm, Oct 15, 2018
Proofing Approved By:	PROOFING APPROVED BY By Natalie Morio at 10:46 am, Oct 16, 2018
Third Eye By:	THIRD EYE BY By Nichole Graham at 3:34 pm, Oct 17, 2018

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INTENDED USE

BD BBL™ Calcofluor White Reagent Droppers are used in the rapid fluorescent microscopic detection of fungi in direct smears.

SUMMARY AND EXPLANATION

Calcofluor white is a nonspecific fluorochrome stain that binds to fungi and, depending upon the filter system employed, fluoresces either an apple green or blue white color when exposed to light of the appropriate wavelength.¹ It may be used on fresh, frozen, fixed, paraffin-embedded, and clinical specimens.² It has been reported that calcofluor white can be used in the detection of parasitic *Pneumocystis jiroveci* (formerly *Pneumocystis carinii*)³ and other opportunistic fungal organisms in bronchoalveolar lavage (BAL) specimens and aspirates from immunosuppressed individuals.^{4,5}

PRINCIPLES OF THE PROCEDURE

Calcofluor white is a disodium salt of 4,4'-bis-[4-anilino-bis-diethylamino-5-triazin-2-ylamino]-2,2'-stilbene-disulfonic acid, a colorless dye that is used in the textile and paper industry as a whitening agent. It has the ability to bind to β 1-3, β 1-4 polysaccharides (i.e., cellulose and chitin), and exhibits fluorescence when exposed to long-wavelength ultraviolet and short-wavelength visible light. It has been used as a biological marker to stain the cell wall of plants, and is therefore valuable in delineating fungal elements.⁶

REAGENTS

BD BBL Calcofluor White Reagent Droppers contain 0.5 mL of a 0.05% solution of calcofluor white in distilled water.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"⁷⁻¹⁰ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids. Prior to discarding, specimen containers and other contaminated materials must be sterilized by autoclaving.

Storage Instructions: Store at room temperature 15–30 °C (59–85 °F). Protect from light.

Product Deterioration: This reagent is hermetically sealed in an ampule, which affords protection of the solution from chemical instability until the expiration date. The staining solution should be clear with a slight green tint. Do not use if a heavy white precipitate is evident. Each dropper is good for one day's use after ampule has been broken. Do not use after the expiration date.

PROCEDURE

Material Provided: BD BBL Calcofluor White Reagent Droppers. Each dropper contains sufficient reagent to stain one slide.

Materials Required But Not Provided: Ancillary culture media, reagents, 10% Potassium Hydroxide Reagent Dropper, Cat. No. 261191, quality control organisms, a fluorescent microscope and other laboratory equipment as required for this procedure.

Test Procedure:

GENERAL SPECIMENS (hair, nails, skin, other tissue, culture, etc.)

1. Slides used should be clean and free of oils and debris.
2. Prepare a smear of the specimen to be stained.
3. Add 1 to 2 drops of a 10% potassium hydroxide solution and gently mix to spread.
4. Hold Calcofluor White Reagent dropper upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to break ampule inside the dropper. **Caution: Break ampule close to its center one time only. Do not manipulate the dropper any further as the plastic may puncture and injury may occur.** Tap the bottom of the dropper on the tabletop a few times. Invert for convenient drop-by-drop dispensing of the reagent.
5. Add 1 to 2 drops of Calcofluor White Reagent. Wait for 1 to 2 min.
6. Mount specimen with a cover slip.
7. Using a fluorescent microscope, examine smears initially at 10X, then confirm findings at a higher magnification.
8. Tissue sections showing fungi can be rinsed in distilled water and restained with periodic acid-Schiff stain or by the Gomori methenamine silver (GMS) technique.⁶

ASPIRATE OR BRONCHOALVEOLAR LAVAGE (BAL) SPECIMENS

1. Slides used should be clean and free of oils and debris.
2. Specimens may be fresh or centrifuged and resuspended to the appropriate concentration in Phosphate Buffered Saline (PBS), pH 7.4.
3. Mark off a well-defined square on the slide with a wax marking pencil.
4. Apply the specimen to the area within the wax square on the slide; allow to air dry.
5. Cover the area within the square with 100% methanol; allow to air dry.
6. Add sufficient 10% potassium hydroxide solution to just cover the specimen. Do not saturate. Wait for 5 s.

7. Hold Calcofluor White Reagent dropper upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to break ampule inside the dropper. **Cautions: Break ampule close to its center *one time only*. Do not manipulate the dropper any further as the plastic may puncture and injury may occur.** Tap the bottom of the dropper on the tabletop a few times. Invert for convenient drop-by-drop dispensing of the reagent.
8. Add 1 to 2 drops of Calcofluor White Reagent to the area within the square. Wait for 2 min.
9. Drain excess Calcofluor White Reagent from the slide.
10. Using a fluorescent microscope view the specimen:
 - a. To view as a wet mount slide, add Calcofluor White Reagent, as needed, a cover slip and view with appropriate objective.
 - or -
 - b. To view under oil immersion, after air-drying, view directly.

User Quality Control:

1. Examine the staining solution for signs of deterioration (see "Product Deterioration").
2. Check the performance of the stain with fresh cultures of *Candida albicans* ATCC™ 10231 (positive) and *Escherichia coli* ATCC 25922 (negative).

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

When exposed to light of the suggested wavelength (300–412 nm), fungi will fluoresce a bright apple green or blue white color. Bacteria will fluoresce weakly or not at all.

LIMITATIONS OF THE PROCEDURE

Many clinical laboratories now use fluorescence microscopy for the examination of other stains. For best results with Calcofluor White Reagent, a mercury vapor rather than a quartz halogen lamp, and UV/BV excitation are recommended.¹ Microscopes fitted with selective filters for the excitation of fluorescein (KP490 interference filter) are not recommended since these filters prevent the proper wave length from striking the specimen.¹

Certain debris may fluoresce and should be distinguished from fungal elements on the basis of morphology.

BD BBL Calcofluor White Reagent Droppers are used for the presumptive identification of fungi. All positive smears should be confirmed by culture.

Although the cells themselves of *Cryptococcus* will stain with calcofluor white, it has been reported that the capsule will not; further testing for confirmation is recommended.

Bacteria may fluoresce, although to a lesser degree than fungi and may be confused with calcofluor white-stained artifacts. The presence of fungi should be determined on the basis of morphology.

Fluorescent staining intensity may vary depending on the organism being stained.

PERFORMANCE CHARACTERISTICS

BD conducted a study of the BD BBL Calcofluor White Reagent using the following organisms:

<i>Aspergillus</i> sp. (ATCC 36607)	<i>Candida albicans</i> (ATCC 10231)
<i>Candida guilliermondii</i> (ATCC 14242)	<i>Cladosporium</i> sp. (ATCC 13026)
<i>Microsporium canis</i> (ATCC 10214)	<i>Penicillium roquefortii</i> (ATCC 9295)
<i>Candida glabrata</i> (ATCC 2001)	

All slides were prepared in triplicate with the wet mount method as described in the "GENERAL SPECIMENS" section of the Test Procedure (including the use of 10% potassium hydroxide reagent). There was one high positive control slide (cellulose via cotton fibers) and two negative control slides (*Escherichia coli* ATCC 25922 and an uninoculated slide).

Slides were examined with a fluorescence microscope (Olympus, BH-2) equipped with a 100w mercury lamp, exciter filter (490 nm) and barrier filter (545 nm).

The organisms evaluated in this study exhibited appropriate fluorescent staining characteristics.

AVAILABILITY

Cat. No. Description

261195	BD BBL™ Calcofluor White Reagent Droppers, 50.
261191	BD BBL™ 10% Potassium Hydroxide Reagent Droppers, 50.

REFERENCES

1. Harrington, B.J., and G.J. Hageage, Jr. 1991. Calcofluor white: tips for improving its use. Clin. Microbiol. Newsl. 13:3-5.
2. Monheit, J.E., D.F. Cowan, and D.G. Moore. 1984. Rapid detection of fungi in tissues using calcofluor white and fluorescence microscopy. Arch. Pathol. Lab. Med. 108:616-618.
3. Stringer, J.R., C.B. Beard, R.F. Miller, and A.E. Wakefield. 2002. A new name (*Pneumocystis jiroveci*) for *Pneumocystis* from humans. Emerg. Infect. Dis. 8(9):891-6.
4. Baselski, V.S., M.K. Robinson, L.W. Pifer, and D.R. Woods. 1990. Rapid detection of *Pneumocystis carinii* in bronchoalveolar lavage samples by using Cellufluor staining. J. Clin. Microbiol. 28:393-394.
5. Kim, Y.K., S. Parulekar, P.K.W. Yu, R.J. Pisani, T.F. Smith, and J.P. Anhalt. 1990. Evaluation of calcofluor white stain for detection of *Pneumocystis carinii*. Diagn. Microbiol. Infect. Dis. 13:307-310.
6. Hageage, G.J., and B.J. Harrington. 1984. Use of calcofluor white in clinical mycology. Lab. Med. 15:109-112.
7. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
8. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
9. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
10. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EECP). Official Journal L262, 17/10/2000, p. 0021-0045.

Technical Information: In the United States contact BD Technical Service and Support at 800.638.8663 or www.bd.com.

BD BBL Calcofluor White Reagent Droppers

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ATCC is a trademark of the American Type Culture Collection.
© 2018 BD. BD, the BD Logo and BBL are trademarks of Becton, Dickinson and Company.
Benex Limited
Becton Dickinson Pty Ltd.