
Revisions

BALTSO0191 Version 8.0 Template 4
Inserts

Rev from	Rev to	JOB #
01	02	7043-14

NOTES:

- BD Catalog Number: 261196
- Blank (Sheet) Size: Length: 11" Width: 22.5"
- Number of Pages: 10 Number of Sheets: 1
- Page Size: Length: 11" Width: 4.5" Final Folded Size: 4 1/2" x 1 7/8"
- Ink Colors: No. of Colors: 3 PMS#: 2755 Blue; 032 Red; Standard Black
- Printed two sides: Yes ☒ No ☐
- Style (see illustrations below): # 4

- See specification control no. N/A for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	REVISED BY By Nancy Carlsen at 4:34 pm, Apr 30, 2015	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.		 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	<i>Pharm 05/01/15</i>	Category and Description Package Insert, PYR Reagent Droppers		Sheet: 1 of 11	A
Checked By	<i>ICStmo 05-05-15</i>			Scale: N/A	
Part Number: L001210					

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INTENDED USE

BD BBL™ PYR Reagent Droppers are used in the rapid presumptive identification of group A β -hemolytic streptococci and group D enterococci.^{1,2}

SUMMARY AND EXPLANATION

PYR Reagent Droppers are used in a rapid test for determining the ability of bacteria to produce pyrrolidonyl peptidase. The production of pyrrolidonyl peptidase (PYRase) has been widely used as an aid in differentiating certain streptococci species.¹⁻³ Following incubation of organisms with a PYR substrate, the reagent (cinnamaldehyde) is added. Color development indicates a positive reaction. The color change is both rapid and long lasting, making results easy to interpret.

PRINCIPLES OF THE PROCEDURE

Pyrrolidonase peptidase activity can be evaluated by the hydrolysis of pyroglutamyl- β -naphthylamide (L-pyrrolidone- β -naphthylamide [PYR]) into L-pyrrolidone and β -naphthylamine, which combines with a PYR reagent to produce a red color.^{4,5}

PYR Reagent Droppers can be used with any PYR substrate. The substrate may be impregnated onto a paper disc or incorporated into agar or broth. After incubation with the substrate, the reagent is added. The reaction can be read after 2 min.

REAGENTS

PYR Reagent Droppers contain 0.5 mL of 0.015% *p*-dimethylaminocinnamaldehyde in acetic acid.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Follow proper laboratory procedures in handling and disposing of infectious materials.

Danger

H360 May damage fertility or the unborn child. **P101** If medical advice is needed, have product container or label at hand. **P281** Use personal protective equipment as required. **P202** Do not handle until all safety precautions have been read and understood. **P308+P313** IF exposed or concerned: Get medical advice/attention. **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage Instructions: Store at room temperature 15 – 30 °C (59 – 86 °F). Protect from light. Each dropper is good for one day's use after breaking ampule.

Product Deterioration: This reagent is hermetically sealed in an ampule, which affords protection of the solution from chemical instability until the expiration date. It is normal for cinnamaldehyde to change from yellow (yellow-orange) to orange to red upon storage. This does not affect the performance of the reagent.

PROCEDURE

Material Provided: PYR Reagent Droppers.

Materials Required But Not Provided: Ancillary culture media, PYR substrate discs, broth or agar, quality control organisms and laboratory equipment as required for these procedures.

Test Procedures:

To use PYR Reagent Droppers, hold upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to break ampule inside the dropper. **Caution: Break ampule close to its center one time only. Do not manipulate the dropper any further as the plastic may puncture and injury may occur.** Tap bottom of dropper on the tabletop a few times. When ready for use, invert for convenient drop-by-drop dispensing of reagent.

PYR Broth: Dispense 1 or 2 drops (depending on sample size) of reagent into a PYR Broth which has been incubated at 35 °C for 24 h. Observe for a red color development within 1 – 2 min indicating a positive test. A negative test is indicated by a yellow color development.

PYR/Esculin Disc: Place the disc on the agar surface of a primary isolation plate of colonies to be tested. Isolated colonies from an 18 – 24 h culture should be selected for testing. Inoculate a small area of the disc surface with a sterile loop or wooden applicator stick. Replace the lid and incubate the plate aerobically at 35 – 37 °C for 15 min. After interpreting the esculin reaction, dispense 1 or 2 drops of reagent onto the disc. Observe for a red color development within 1 – 2 min indicating a positive test. No color development indicates a negative test.

User Quality Control:

Examine the reagent for signs of deterioration (see "Product Deterioration"). Positive (*Streptococcus pyogenes* ATCC™ 19615 and *Enterococcus faecalis* ATCC 29212) and negative (*Streptococcus agalactiae* ATCC 10386) controls should be performed along with the organisms to be tested.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures.

It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Positive: Pink to red. **Negative:** No color change or yellow-orange. Consult appropriate references for PYR reactions with specific organisms.^{4,5}

LIMITATIONS OF THE PROCEDURE

- 1. Allow at least 2 min for the color to develop before considering the PYR test negative.
- 2. This test should only be used on β-hemolytic streptococci and *Enterococcus* species. *Staphylococcus* species, most *Corynebacterium haemolyticum*, and some *Micrococcus* species are PYR positive, as well as some *Enterobacteriaceae* and other gram-negative bacilli.^{6,7}
- 3. Use only pure cultures. Mixed cultures may produce misleading results.

PERFORMANCE CHARACTERISTICS

Godsey, Schulman and Eriquez reported that 100% of group A streptococci and 100% of group D enterococci hydrolyzed PYR producing a positive test result.⁸ They also reported that there were no false positive results obtained when testing group D nonenterococci or streptococci belonging to groups B,C,F, or G.⁸ False positive PYR reactions have been reported by Bosley, Facklam and Grossman when aerococci were test.⁹ Oberhofer found that *S. haemolyticus*, *S. intermedius* and *S. warneri* gave strong positive PYR reactions while all other species of staphylococci were found to be negative.¹⁰

An internal study was conducted comparing **BD BBL** PYR Reagent and another commercially available PYR Reagent.¹¹ As summarized in the table, all reactions were identical for the 6 strains tested using both PYR reagents.

Organism	BD BBL Reagent	Commercial Reagent
<i>Streptococcus pyogenes</i> ATCC 10389	+	+
<i>Streptococcus agalactiae</i> ATCC 4768	-	-
<i>Enterococcus faecalis</i> ATCC 29212	+	+
<i>Streptococcus pyogenes</i> ATCC 19615	+	+
<i>Streptococcus bovis</i> ATCC 9809	-	-
<i>Streptococcus</i> sp. Group B ATCC 12386	-	-

AVAILABILITY

Cat. No.	Description
261196	BD BBL™ PYR Reagent Droppers, 50.

REFERENCES

1. Gordon, D.B., P.C. DeGirolami, S. Bolivar, G. Karafotias, and K. Eichelberger. 1988. A comparison of the identification of group A streptococci and enterococci by two rapid pyrrolidonyl aminopeptidase methods. *Am. J. Clin. Pathol.* 90:210-212.

2. Facklam, R.R., L.G. Thacker, B. Fox, and L. Eriquez. 1982. Presumptive identification of streptococci with a new test system. *J. Clin. Microbiol.* 15:987-990.

3. Ellner, P.D., D.A. Williams, M.E. Hosmer, and M.A. Cohenford. 1985. Preliminary evaluation of a rapid colorimetric method for the presumptive identification of group A streptococci and enterococci. *J. Clin. Microbiol.* 22:880-881.

4. Baron, E.J., and S.M. Finegold. 1990. *Bailey and Scott's diagnostic microbiology*, 8th ed. The C.V. Mosby Co., St Louis.

5. Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy (ed.) 1991. *Manual of clinical microbiology*, 5th ed. American Society for Microbiology, Washington, D.C.

6. Facklam, R.R. 1987. Specificity study of kits for detection of group A streptococci directly from throat swabs. *J. Clin. Microbiol.* 25:504-508.

7. Lawrence, J., D.M. Yajko, and W.K. Hadley. 1984. Comparison of Strep-A-Check (a 10 min chromogenic substrate test) with bacitracin for presumptive identification of *Streptococcus pyogenes* (group A strep). Program Abstr. 24th Intersci. Conf. Antimicrob. Agents Chemoth., abstr 19.

8. Godsey, J., R. Schulman, and L. Eriquez. 1981. The hydrolysis of L-pyrrolidonyl-β-naphthylamide as an aid in the rapid identification of *Streptococcus pyogenes*, *S. aviium*, and group D-enterococci. Abstr. C-84, p. 276. Abstr. 81st Annu. Meet. Am. Soc. Microbiol. 1981.

9. Bosley, G.S., R.R. Facklam, and D. Grossman. 1983. Rapid identification of enterococci. *J. Clin. Microbiol.* 18:1275-1277.

10. Oberhofer, T.R. 1986. Value of the pyrrolidonyl-β-naphthylamide hydrolysis test for identification of select gram-positive cocci. *Diagn. Microbiol. Infect. Dis.* 4:43-47.

11. Data on file, BD Diagnostic Systems.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo and BBL are trademarks of Becton, Dickinson and Company. © 2015 BD.
Benex Limited
Becton Dickinson Pty Ltd.