






Revisions

BALTSO0191 Version 8.0 Template 4
Inserts

Rev from	Rev to	JOB #
01	02	7043-14

NOTES:

- BD Catalog Number: 261204
- Blank (Sheet) Size: Length: 11" Width: 22.5"
- Number of Pages: 10 Number of Sheets: 1
- Page Size: Length: 11" Width: 4.5" Final Folded Size: 4.5" x 1 7/8"
- Ink Colors: No. of Colors: 3 PMS#: 2755 Blue; 032 Red; Standard Black
- Printed two sides: Yes No ☒ ☐
- Style (see illustrations below): # 4

- See specification control no. N/A for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	REVISED BY By Nancy Carlsen at 4:01 pm, Apr 30, 2015	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.		Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	<i>Sharon 05/04/15</i>			
Checked By	<i>ICSM 05-04-15</i>	Category and Description		Sheet: 1 of 11
Part Number:	L001217	Package Insert, Methylene Blue Loeffler Reagent Droppers		Scale: N/A
				A

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INTENDED USE

BD BBL™ Methylene Blue Loeffler Stain Droppers are used for presumptive identification of *Corynebacterium diphtheriae* and as an adjunct to the Gram stain.

SUMMARY AND EXPLANATION

Methylene Blue is recognized as a simple stain that can be used for a variety of purposes, such as determining bacterial morphology and aiding in the interpretation of the Quellung test.¹ It is the stain of choice for the presumptive identification of *Corynebacterium diphtheriae*. Upon staining, this species is characterized by a beaded or banded appearance of the cell, with deeply stained metachromatic granules contrasted with a lighter blue cytoplasm.^{2,3}

Methylene Blue is used as an adjunct to the Gram stain, especially in the staining of gram-negative bacteria such as *Haemophilus influenzae* and *Neisseria* species.^{4,5} The stain can reveal the morphology of fusiform bacteria and spirochetes (from oral infections such as Vincent's angina) that may be difficult to see with the Gram stain.⁵ The stain is also used to detect the presence of fecal white blood cells that may show the presence of invasive bowel disease.⁵

BD BBL Methylene Blue Loeffler Stain Droppers can be used as a counterstain with **BD BBL** Flagella Stain Droppers, which are recommended for staining the flagella of bacteria.

BD BBL Methylene Blue Loeffler Stain Droppers are supplied in a hermetically sealed ampule for improved stability of the reagent prior to use.

PRINCIPLES OF THE PROCEDURE

Methylene Blue is a cationic dye that imparts a blue color to the negatively charged polyphosphate portions of cells (such as ribonucleic acid).¹ Metachromasia, which is observed with *C. diphtheriae*, is due to a high concentration of polymerized polyphosphates. The deeply blue stained polyphosphate granules, surrounded by a lighter blue stained cytoplasm, are known as metachromatic granules or Babes-Ernst bodies.

REAGENTS

BD BBL Methylene Blue Loeffler Stain Droppers contain: 0.3 g Methylene Blue, 30 mL Alcohol, 95 %, 0.1 mL Potassium Hydroxide, 10 %, and 100 mL Distilled Water.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Harmful by inhalation and if swallowed.

Follow established laboratory procedures in handling and disposing of infectious materials.

Warning

H226 Flammable liquid and vapour. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P233** Keep container tightly closed. **P303+P361+P353** IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. **P370+P378** In case of fire: Use for extinction: CO₂, powder or water spray. **P403+P235** Store in a well-ventilated place. Keep cool. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage Instructions: Store Methylene Blue Loeffler Stain Droppers at 15 – 30 °C. Methylene Blue Loeffler Stain Droppers are ready to use.

The expiration date applies to the product in its intact container when stored as directed.

Product Deterioration: Do not use a product if it fails to meet specifications for identity and performance.

SPECIMEN COLLECTION AND PREPARATION

Refer to appropriate references for specific procedures on processing specimens.⁴⁻⁷

The microscopic morphology of *C. diphtheriae* is particularly distinctive when grown on Loeffler Blood Serum.⁸

PROCEDURE

Materials Provided: Methylene Blue Loeffler Stain Droppers.

Materials Required But Not Provided: Glass slides, Bunsen burner, inoculating loop and quality control organisms.

Test Procedure:

1. Prepare a smear of the specimen to be stained on a clean glass slide.
2. Allow the smear to air dry thoroughly.
3. Heat or methanol fix the smear.
4. Hold dropper upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to break the inner ampule. **Caution: Break ampule close to its center one time only. Do not manipulate the dropper any further as the plastic may puncture and injury may occur.** Tap bottom on tabletop a few times. Invert for convenient drop by drop dispensing of solution.
5. Flood the slide with stain and allow to stain for 1 – 3 min.
6. Rinse the slide with tap or distilled water and blot dry.
7. Examine the smear with an oil-immersion objective.

User Quality Control:

Identity Specifications – Blue, clear solution, inside a glass ampule contained in a plastic container.

Cultural Response – Test the performance of Methylene Blue Loeffler Stain Droppers using the procedure described above. Observe *C. diphtheriae* for metachromatic granules.

Organism	ATCC™	Stain Reaction
<i>Escherichia coli</i>	25922	blue rods
<i>Streptococcus pyogenes</i>	19615	blue cocci
<i>Corynebacterium diphtheriae</i> biotype <i>gravis</i>	8028	deep blue granules

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Bacterial cells stain medium blue, while the background stains lighter blue. Stained *C. diphtheriae* appear as pleomorphic beaded rods with swollen ends, producing a club shape, and containing deeply stained metachromatic granules.⁶

LIMITATIONS OF THE PROCEDURE

Over-staining may reduce the contrast between granules and cytoplasm or between bacteria and background.⁷

The Loeffler staining test may be used in the presumptive identification of *Corynebacterium diphtheriae*. Additional biochemical testing is recommended for complete identification because some strains of *Actinomyces*, *Propionibacterium*, and pleomorphic forms of streptococci may demonstrate the characteristic stained appearance of *C. diphtheriae*.⁴

PERFORMANCE CHARACTERISTICS

Prior to release, all lots of **BD BBL™** Methylene Blue Loeffler Stain Droppers are tested to verify specific product characteristics. One drop of a heavy saline suspension (approximately McFarland standard no. 3) of *Escherichia coli* (ATCC 25922) and *Streptococcus pyogenes* (ATCC 19615) is placed on a microscope slide, allowed to air dry and heat fixed. A smear of *Corynebacterium diphtheriae* biotype *gravis* (ATCC 8028) is placed on a clean microscopic slide, allowed to air dry and heat fixed. The slides are flooded with **BD BBL** Methylene Blue Loeffler Stain and allowed to stand for 2 min. They are then rinsed with water and blotted dry. The slides are examined under the oil-immersion lens of a light microscope. Bacterial cells of *E. coli* and *S. pyogenes* stain blue although the color intensity may vary. Cell morphology is readily distinguishable. Metachromatic granules in the cells of *C. diphtheriae* stain a deeper blue than the surrounding cytoplasm.

AVAILABILITY

Cat. No.	Description
261204	BD BBL™ Methylene Blue Loeffler Stain Droppers, 50.

REFERENCES

1. Lillie, R. D. (ed.). 1977. H. J. Conn's biological stains, 9th ed. Williams and Wilkins, Baltimore.
2. Clark, G. (ed.). 1981. Staining procedures, 4th ed. Williams and Wilkins, Baltimore.
3. Saelinger, C. B., and F. F. Bonventre. 1983. *Corynebacterium diphtheriae*. In N. R. Rose and A. L. Baren (ed.), Microbiology: basic principles and clinical applications. Macmillan, New York.
4. Koneman, E. W., S. D. Allen, V. R. Dowell, Jr., and H. M. Sommers. 1983. Color atlas and textbook of diagnostic microbiology, 2nd ed. J. B. Lippincott Company, Philadelphia.
5. Murray, P. R., E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.). 1995. Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
6. Isenberg, H. D. (ed.). 1992. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
7. Baron, E. J., and S. M. Finegold. 1992. Bailey & Scott's diagnostic microbiology, 8th ed. The C. V. Mosby Company, St. Louis.
8. Lennette, E. H., A. Balows, W. J. Hausler, Jr., and J. P. Truant (ed.). 1980. Manual of clinical microbiology, 3rd ed. American Society for Microbiology, Washington, D.C.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

BD BBL Methylene Blue Loeffler Stain Droppers

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DISPONIBILIDAD

N° de cat. Descripción

261204 **BD BBL** Methylene Blue Loeffler Stain Droppers, 50.

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

Australian Sponsor:

Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113
Australia

ATCC is a trademark of the American Type Culture Collection.

BD, BD Logo and BBL are trademarks of Becton, Dickinson and Company. © 2015 BD.