






Revisions

BALTSO0191 Version 8.0 Template 4
Inserts

Rev from	Rev to	JOB #
01	02	7043-14

NOTES:

- BD Catalog Number: 261203
- Blank (Sheet) Size: Length: 11" Width: 27"
- Number of Pages: 12 Number of Sheets: 1
- Page Size: Length: 11" Width: 4.5" Final Folded Size: 4 1/2" x 1 7/8"
- Ink Colors: No. of Colors: 3 PMS#: 2755 Blue; 032 Red; Standard Black
- Printed two sides: Yes ☒ No ☐
- Style (see illustrations below): # 4

- See specification control no. N/A for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	REVISED BY By Nancy Carlsen at 3:40 pm, Apr 30, 2015	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.		 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	<i>[Signature]</i> 05/01/15			
Checked By	<i>[Signature]</i> 05-04-15	Category and Description		Sheet: 1 of 13
Part Number:	L001237	Package Insert, Catalase Reagent Droppers		Scale: N/A
				A

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INTENDED USE

BD BBL™ Catalase Reagent Droppers are used for determining the presence of catalase produced by bacteria.

SUMMARY AND EXPLANATION

Catalase Reagent Droppers are used in a qualitative procedure for determining catalase activity by bacteria. The enzyme catalase is present in most aerobic and facultatively anaerobic bacteria, with the major exception of streptococci and enterococci. The enzyme catalase decomposes hydrogen peroxide into water and oxygen.

Catalase Reagent Droppers contain a hydrogen peroxide solution of approximately 3% (2.5% to 3.5%). This reagent has been frequently cited as a standard method for performing the catalase test.¹⁻⁸

PRINCIPLES OF THE PROCEDURE

Hydrogen peroxide, an end product of aerobic carbohydrate metabolism, is extremely toxic to bacteria. Bacteria containing the enzyme catalase break down hydrogen peroxide into oxygen and water. Organisms that possess catalase demonstrate the reaction in 3% hydrogen peroxide by the rapid appearance of gas bubbles.

REAGENTS

Catalase Reagent Droppers contain: (approximately) 3% Hydrogen Peroxide and 0.05% or less of stabilizer.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Follow proper established laboratory procedures in handling and disposing of infectious materials.

Warning

H302 Harmful if swallowed. **H315+H320** Causes skin and eye irritation. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P302+P352** IF ON SKIN: Wash with plenty of soap and water. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage Instructions: Store Catalase Reagent Droppers at 15 – 30 °C. Catalase Reagent is light sensitive. Protect from light.

Catalase Reagent Droppers are ready for use. The expiration date applies to the product in its intact container when stored as directed.

After unsealing the tube, use the reagent within one working day. Discard any remaining reagent.

Product Deterioration: Do not use a product if it fails to meet specifications for identity and performance.

SPECIMEN COLLECTION AND PREPARATION

Collect samples or specimens in sterile containers or with sterile swabs and transport immediately to the laboratory according to recommended guidelines.¹⁻⁸

Process each specimen using procedures appropriate for that sample.¹⁻⁸

PROCEDURE

Material Provided: Catalase Reagent Droppers.

Materials Required But Not Provided: Sterile containers or sterile swabs, platinum inoculating loop or needle or applicator stick, glass slides, plated or tubed media, Bunsen burner or incinerator and quality control organisms.

Test Procedure:

Use of Reagent Dispenser

1. Holding the dispenser upright and with the tip pointing away from the user, twist and pull the cap to remove it from the tube. Avoid squeezing the dispenser while removing the cap.
2. Invert and squeeze slightly to dispense the reagent on a per drop basis.

Slide Test Method

1. Obtain a pure culture of the organism to be tested.
2. Using an inoculating needle or applicator stick, pick a well-isolated colony and transfer to a glass slide.
3. Add 1 or 2 drops of the Catalase Reagent to the smear.
4. Examine immediately for the rapid production of gas bubbles.

Tube or Agar Plate Method

1. Add a few drops of Catalase Reagent to the surface of an 18 – 24 h agar plate or slant that does not contain blood.
2. Examine immediately and after 5 min for the evolution of bubbles.

Broth Method

1. Express the contents of a Catalase Reagent Dropper into a 24 – 48 h broth culture.
2. Examine immediately and after 5 min for the evolution of bubbles.

User Quality Control:

Identity Specifications – The solution is colorless, clear and free of precipitate or foreign matter.

Cultural Response – Test the performance of Catalase Reagent Droppers using one of the methods described above.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

Organism	ATCC™	Reaction
<i>Staphylococcus aureus</i>	25923	+
<i>Streptococcus pyogenes</i>	19615	-

RESULTS

Positive catalase reactions are denoted by the production of bubbles that occur after the addition of the reagent to the culture.

LIMITATIONS OF THE PROCEDURE

The Catalase Reagent Droppers may be used as an aid in the identification of microorganisms. Additional biochemical testing using pure cultures is recommended for complete identification.

Culture media that contain blood are unacceptable for this test due to the presence of catalase in erythrocytes. When picking colonies from blood media, avoid carryover when performing the slide test. If a blood medium is used, a control slide catalase test should be performed where a small loopful of the blood-containing agar is tested with the reagent on the same slide as the organism. If the catalase reaction from the colony is much stronger than that from the agar alone, the test can be considered positive.²

Dirty glassware can cause false-positive results. Use properly cleaned tubes and slides.

Inoculating needles or loops containing iron may produce false-positive results. Avoid contact of the catalase reagent with iron-containing needles and loops.

Bacteria (e.g., lactobacilli) grown on media with low levels or no glucose may yield confusing reactions from pseudocatalase, a non-iron enzyme. The pseudocatalase reaction can be prevented by the addition of 1% glucose to the medium.^{9,10}

Anaerobic cultures grown on a blood-free medium must be exposed to air for a minimum of 30 min before testing.⁹

PERFORMANCE CHARATERISTICS¹¹

An internal study verified that the Catalase Reagent performed as expected biologically. Catalase Reagent is formulated to a concentration of approximately 3% hydrogen peroxide. Six organisms were tested in triplicate with testing performed on glass microscope slides and on Chocolate II Agar plates.

The acceptance criteria was established prior to beginning the experiment:

Organism	Slide Reaction	Plate Reaction
<i>N. meningitidis</i> ATCC 13090	+	+
<i>E. faecalis</i> ATCC 33186	-	-
<i>P. aeruginosa</i> ATCC 27853	+	+
<i>S. aureus</i> ATCC 25923	+	+
<i>S. epidermidis</i> ATCC 12228	+	+
<i>S. pyogenes</i> ATCC 19615	-	-

The results show that the reaction was acceptable at a concentration of 3% H₂O₂ for each organism tested.

Organism	% H ₂ O ₂	Slide Reaction	Plate Reaction
<i>N. meningitidis</i> ATCC 13090	3.0	+	+
<i>E. faecalis</i> ATCC 33186	3.0	-	-
<i>P. aeruginosa</i> ATCC 27853	3.0	+	+
<i>S. aureus</i> ATCC 25923	3.0	+	+
<i>S. epidermidis</i> ATCC 12228	3.0	+	+
<i>S. pyogenes</i> ATCC 19615	3.0	-	-

Availability

Cat. No. Description

261203 BD BBL™ Catalase Reagent Droppers, 50.

References

1. MacFaddin, J.F. 1980. Biochemical tests for identification of medical bacteria, 2nd ed. Williams & Wilkins, Baltimore MD.
2. Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. Bailey & Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis.
3. Vanderzant, C., and D.F. Splittstoesser (ed.). 1992. Compendium of methods for the microbiological examination of foods, 3rd ed. American Public Health Association, Washington, D.C.
4. Marshall, R.T. (ed.). 1992. Standard methods for the examination of dairy products. American Public Health Association, Washington, D.C.
5. Isenberg, H.D. (ed.). 1994. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
6. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.). 1995. Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
7. Greenberg, A.E., L.S. Clesceri, and A.D. Eaton (ed.). 1995. Standard methods to the examination of water and wastewater, 19th ed. American Public Health Association, Washington, D.C.
8. Cuniff, P. (ed.). 1995. Official methods of analysis of AOAC International, 16th ed., March 1996 supplement. AOAC International, Gaithersburg, Md.

9. Gerhardt, P.R., G.E. Murray, W.A. Wood, and N.R. Krieg (ed.). 1994. Methods for general and molecular bacteriology. American Society for Microbiology, Washington, D.C.

10. Doelle, H.W. 1969. Bacterial metabolism. Academic Press, New York.

11. Data on file, BD Diagnostic Systems.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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Benex Limited
Becton Dickinson Pty Ltd.