Part Number:	8800801JAA	BALTSO0191 Version 13.0 Template 4 Inserts	
Category and Description:	Package Insert BD BBL Cefinase Paper Discs	Rev from: 02	Rev to: 03
		Job Number:	706-18

Catalog Number: 231650

Blank (Sheet) Size: Length: 15.5" Width: 24"


Number of Pages: 10 Number of Sheets: 1

Page Size: Length: 15.5" Width: 4.8125" Final Folded Size: 4.8125" X 1.25"

Ink Colors: Number of Colors: 1 PMS #: Standard Black

Printed Two Sides: Yes: ☒ No: ☐

Style (see illustrations below): # 4

Vendor Printed: ☒ Online / In House Printed: ☐ Web Printed: ☐

See Specification control no. VS-88-0080-1PR for material information.

BD Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Company confidential. This document is the property of Becton, Dickinson and Company and is not to be used outside the company without written permission. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Revised By: **REVISED BY**
By Shanthi Sagayamary Ramu at 4:21 pm, Sep 20, 201

Proofing Approved By: **PROOFING APPROVED BY**
By Natalie Morio at 1:44 pm, Sep 25, 2018

Third Eye By: **THIRD EYE BY**
By Sonia F.Thompson at 3:48 pm, Sep 25, 2018

 BD BBL™ Paper Discs for the
Detection of β-Lactamase Enzymes
Cefinase Discs

English: pages 1 – 3 Italiano: pagine 6 – 7
Français: pages 3 – 4 Español: páginas 7 – 9
Deutsch: Seiten 4 – 6

8800801JAA(03)
2018-09

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyts lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőjétől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o representante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

BD BBL Cefinase™ discs are intended for use in rapid testing of isolated colonies of *Neisseria gonorrhoeae*, *Staphylococcus* species, *Haemophilus influenzae*, enterococci and anaerobic bacteria for the production of β-lactamase.

SUMMARY AND EXPLANATION

The ability of certain bacteria to produce enzymes that inactivate β-lactam antibiotics, i.e., penicillins and cephalosporins, has long been recognized. Abraham and Chain in 1940 first recognized an enzymatic activity in extracts of *Escherichia coli* that inactivated penicillin.¹ Since then a large number of similar enzymes have been isolated from a number of bacterial species with somewhat different substrate specificities. Some selectively hydrolyze penicillin-class antimicrobics (i.e., penicillin G, ampicillin, carbenicillin) and have been described as penicillinases. Others selectively hydrolyze the cephalosporin-class antimicrobics (i.e., cephalothin, cephalexin, cephadrine) and have been described as cephalosporinases. Still other enzymes hydrolyze both cephalosporins and penicillins.²

A large number of β-lactamase-resistant penicillin and cephalosporin class antimicrobics have been developed by various pharmaceutical companies. One group includes the semisynthetic penicillins; methicillin, oxacillin, nafcillin and others, which are resistant to the penicillinase enzymes produced by staphylococci.³ A large number of cephalosporins have also been developed which have varying degrees of resistance to β-lactamases. These include the second-generation cephalosporins (cefoxitin, cefamandole and cefuroxime) and third-generation cephalosporins (cefotaxime, moxalactam, cefoperazone and others).⁴

Several clinical tests have been developed for the detection of β-lactamases. These tests provide rapid information predictive of the development of resistance. Interpretation of β-lactamase test results must consider: the sensitivity of the test for different classes of β-lactamase enzymes, the types of β-lactamases produced by different taxonomic groups of organisms and the substrate specificities of the different β-lactamases.

The most commonly used clinical procedures include the iodometric method, the acidometric method, and a variety of chromogenic substrates.⁵ The iodometric and acidometric tests are generally performed using penicillin as a substrate and, therefore, can only detect enzymes which hydrolyze penicillin. One of the chromogenic cephalosporins, PADAC (Calbiochem-Behring) has been found effective in detecting most of the known β-lactamases except for some of the penicillinases produced by staphylococci, and some β-lactamases produced by anaerobic bacteria.⁶ Another chromogenic cephalosporin, nitrocefin (Glaxo Research), has been found effective in detecting all known β-lactamases including the staphylococcal penicillinases.⁷⁻⁹

For many taxonomic groups of organisms, e.g. *Enterobacteriaceae*, the β-lactamase test is of little value because a diversity of β-lactamase enzymes with different substrate specificities may be produced within the group, or even within a single strain.¹⁰

In other bacteria, for example, penicillin-resistant *Neisseria gonorrhoeae*,¹¹ *Staphylococcus aureus*,^{12,13} *Moraxella catarrhalis*,¹⁴ and ampicillin-resistant *Haemophilus influenzae*,^{5,9,15} only one class of enzyme is produced by resistant strains. The β-lactamase test performed with these organisms enables a prediction of resistance to be made immediately after primary isolation, 18–24 h prior to the time that growth-dependent susceptibility results would be available.

While the prevalence of β-lactamase-producing enterococci appears to be small, a low inoculum may result in strains going undetected by susceptibility-testing procedures, and routine screening by the nitrocefin disc procedure is recommended.¹⁶

With anaerobic bacteria, the relationship between the production of β-lactamase and resistance to β-lactam antimicrobics is complicated and somewhat similar to *Enterobacteriaceae*. β-lactamases are most commonly found within the *Bacteroides* species,¹⁷ however, β-lactamase-producing strains of *Clostridium butyricum*, *C. perfringens* and *Fusobacterium* sp. have been reported.^{18,19} Among the *Bacteroides* group, a variety of enzymes may be produced with different substrate specificities. The β-lactamases frequently found in strains of *Prevotella melanogenica* and *P. oralis* are usually specific for penicillins (penicillinase),²⁰ whereas the β-lactamases frequently found in the *B. fragilis* group are cephalosporinases.^{21,22} A variety of cephalosporinases have been reported in the *B. fragilis* group and they include some very active enzymes which can hydrolyze some of the reportedly β-lactamase-resistant cephalosporins such as cefotaxime.^{23,24} Rare strains have been reported which hydrolyze at high rates all known β-lactams including cefoxitin.^{24,25}

Even though the β-lactamases produced by the *B. fragilis* group are most active against cephalosporins, most strains are found to be resistant to penicillin, carbenicillin and ampicillin in growth-dependent susceptibility tests.^{17,26} This finding suggests that the *B. fragilis* group may be intrinsically resistant to penicillins through factors such as permeability barriers,²² or that the β-lactamase is produced in quantities sufficient to overcome the relatively slow hydrolysis rate of the enzyme with penicillins. Evidence which tends to support a contributory role for β-lactamase in the resistance to penicillins is found in reports that the combination of clavulanic acid (a β-lactamase inhibitor) and penicillins is many times more active against *B. fragilis* than is the penicillin alone.²⁷

Whatever the cause or causes of penicillin resistance in *B. fragilis*, all strains should probably be considered resistant.²⁸ The other gram-negative anaerobic strains are probably susceptible to penicillin so long as they are β-lactamase negative.²⁸

PRINCIPLES OF THE PROCEDURE

The BD BBL Cefinase disc is impregnated with the chromogenic cephalosporin, nitrocefin. This compound exhibits a very rapid color change from yellow to red as the amide bond in the β-lactam ring is hydrolyzed by a β-lactamase. When a bacterium produces this enzyme in significant quantities, the yellow-colored disc turns red in the area where the isolate is smeared.

Although other penicillins and cephalosporins may be used as substrates for specific enzymes, nitrocefin has the wide spectrum of susceptibility and sensitivity of the commercially available β-lactams. It is not known to react with other microbial enzymes.²⁹

Each disc is used to test one bacterial strain for the presence of β-lactamase.

REAGENTS

BD BBL Cefinase discs impregnated with nitrocefin.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

These discs are not for use in susceptibility testing.

Observe aseptic techniques and established precautions against microbiological hazards throughout all procedures. After use, prepared plates and other contaminated materials must be sterilized by autoclaving before discarding.

Nitrocefin induces mutations in certain strains of bacteria (Ames test) and may be sensitizing. Ingestion, inhalation or contact with the skin or eyes should be avoided.

Storage Instructions: Upon receipt, store unopened package at -20 to +8 °C. After use, the BD BBL Cefinase cartridge should be stored in any glass, air tight container containing desiccant and stored at -20 to +8 °C. Discard remaining BD BBL Cefinase discs 60 days after opening blister packaging. The expiration date on the cartridge applies only to discs intact in unopened blister packaging.

Indications of Deterioration: Do not use the cartridge if the discs appear orange or red in color.

SPECIMEN COLLECTION AND HANDLING

This procedure is not to be used directly with clinical specimens or other sources containing mixed microbial flora. The bacteria to be tested must first be isolated as separate colonies by streaking the specimen onto appropriate culture media plates.

PROCEDURE

Material Provided: BD BBL Cefinase discs, 50 discs per cartridge.

Materials Required But Not Provided: Ancillary reagents, quality control organisms and laboratory equipment as required for the procedure.

Test Procedure:

- Using a single disc dispenser, dispense the required number of discs from the cartridge into an empty Petri dish or onto a microscope slide.
- Moisten each disc with one drop of purified water.
- With a sterilized loop or applicator stick remove several well-isolated similar colonies and smear onto a disc surface.
- Observe disc for color change.
- Alternate procedure: Using forceps moisten disc with one drop of purified water and then wipe across colony.

User Quality Control: Control reference cultures should be run with each group of unknowns. The following organisms are recommended for use as test strains.

Test Strain	Expected Results
<i>Staphylococcus aureus</i> ATCC® 29213	Positive
<i>Haemophilus influenzae</i> ATCC 10211	Negative

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS AND INTERPRETATION

A positive reaction will show a yellow to red color change on the area where the culture was applied. Note: color change does not usually develop over the entire disc. A negative result will show no color change on the disc. For most bacterial strains a positive result will develop within 5 min. However, positive reactions for some staphylococci may take up to 1 h to develop.

Organism	Result	Approx. Reaction Time	Interpretation
<i>Staphylococcus aureus</i>	Positive	1 h	Resistant to penicillin, ampicillin, carbenicillin and ticarcillin. Probably susceptible to cephalothin, methicillin, oxacillin, nafcillin and other penicillinase-resistant penicillins.*
<i>Haemophilus influenzae</i>	Positive	1 min	Resistant to ampicillin. Susceptibile to cephalosporins.*
<i>Neisseria gonorrhoeae</i> and <i>Moraxella catarrhalis</i>	Positive	1 min	Resistant to penicillin.
<i>Enterococcus faecalis</i>	Positive	5 min	Resistant to penicillin and ampicillin
Anaerobic bacteria	Positive	30 min	Probable identification is <i>Bacteroides</i> species. Probably resistant to penicillin and may be resistant to cephalosporins including cefotaxime and rarely cefoxitin.

* Susceptibility should be confirmed by growth-dependent susceptibility tests. Negative results imply but do not guarantee susceptibility.

LIMITATIONS OF THE PROCEDURE

The efficacy of this test in predicting the β-lactam resistance of microorganisms other than *Neisseria gonorrhoeae*, *Haemophilus influenzae*, *Moraxella catarrhalis*, staphylococci, enterococci and certain anaerobic bacteria is unproven.

Resistance to β-lactam antibiotics has been on rare occasions reported in some of the above organisms without the production of β-lactamases.^{30,31} In these cases, resistance mechanisms such as permeability barriers have been postulated. Therefore, the β-lactamase test should be used as a rapid supplement and not a replacement for conventional susceptibility testing.

For some strains of staphylococci,¹³ particularly *S. epidermidis*, an inducible β-lactamase has been described that might result in a false-negative β-lactamase reaction with a strain which is resistant to penicillin or ampicillin.

SPECIFIC PERFORMANCE CHARACTERISTICS

In a comparative study of four methods for detecting β-lactamase activity in anaerobic bacteria, the following percentages of agreement with a “standard” employing nitrocefin-saturated filter paper were obtained: Cefinase, 100%; pyridine-2-azo-p-dimethylaniline cephalosporin, 96%; a penicillinase disc using bromcresol purple pH indicator, 72%; the slide iodometric technique, 78%.³²

AVAILABILITY

Cat. No.	Description
231650	BD BBL™ Cefinase™, 50

REFERENCES

- Abraham, E.P., and E. Chain. 1940. An enzyme from bacteria capable of destroying penicillin. *Nature* 146:837.
- McCarthy, L.R. 1980. β-lactamases. *Clin. Microbiol. Newsl.* 2 (2): 1-3. G.K. Hall and Co., Boston.
- Richmond, M.H. 1979. β-lactam antibiotics and β-lactamases: two sides of a continuing story. *Rev. Inf. Dis.* 1:30-36.
- Bush, K., and R.B. Sykes. 1982. Interaction of new β-lactams with β-lactamases and β-lactamases-producing gram-negative rods, p.47-63. *In* H.C. Neu (ed.), *New β-lactam antibiotics: review from chemistry to clinical efficacy of new cephalosporins*. College of Physicians of Philadelphia, Philadelphia.
- Thornsberry, C., T.L. Gavan, and E.H.Gerlach. 1977. Cumitech 6, New developments in antimicrobial agent susceptibility testing. Coordinating ed., J.C. Sherris. American Society for Microbiology, Washington, D.C.
- Jorgensen, J.H., S.A. Crawford, and G.A. Alexander. 1982. Pyridine-2-azo-p-dimethylaniline chromophore, a new chromogenic cephalosporin for rapid beta-lactamase testing. *Antimicrob. Agents Chemother.* 22:162-164.
- Montgomery, K., L. Raymundo, Jr., and W.L. Drew. 1979. Chromogenic cephalosporin spot test to detect beta-lactamase in clinically significant bacteria. *J. Clin. Microbiol.* 9:205-207.
- O'Callaghan, C.H., A. Morris, S.M. Kirby, and S.H. Shingler. 1972. Novel method for detection of β-lactamase by using a chromogenic cephalosporin substrate. *Antimicrob. Agents and Chemother.* 1:283-288
- Skinner, A., and R. Wise. 1977. A comparison of three rapid methods of β-lactamase activity in *Haemophilus influenzae*. *J. Clin. Pathol.* 30:1030-1032.
- Sykes, R.B., and M. Mathew. 1976. The β-lactamases of gram-netative bacteria and their role in resistance to β-lactam antibiotics. *J. Antimicrob. Chemother.* 2:115-157
- Ashford, W.A., R.G. Golash, and V.G. Hemming. 1976. Penicillinase-producing *Neisseria gonorrhoeae*. *Lancet ii*:657-658.
- Adam, A.P., A.L. Barry, and E. Benner. 1970. A simple rapid test to differentiate penicillin-susceptible from penicillin-resistant *Staphylococcus aureus*. *J. Infect. Dis.* 122:544-546
- Kirby, W.M.M. 1944. Extraction of a highly potent penicillin inactivator from penicillin resistant staphylococci. *Science* 99:452-453
- Malmvall, B.E., J.E. Brorsson, and J. Johnsson. 1977. *In vitro* sensitivity to penicillin V and β-lactamase production of *Branhamella catarrhalis*. *J. Antimicrob. Chemother.* 3:374-375.
- Khan, W., S. Ross, W. Rodriquez, G. Controni, and A.K. Saz. 1974. *Haemophilus influenzae* type b resistant to ampicillin. *J. Am. Med. Assoc.* 299:298-301.
- Neumann, M.A., D.F. Sahm, C. Thornsberry, and E.A. Kerr. 1991. Cumitech 6A, New developments in antimicrobial agent susceptibility testing: a practical guide. Coordinating ed., J.E. McGowan, Jr. American Society for Microbiology, Washington, D.C.
- Olsson, B., K. Dornbush, and C.E. Nord. 1977. Susceptibility testing of β-lactam antibiotics and production of β-lactamase in *Bacteroides fragilis*. *Med. Microbiol. Immunol.* 163:183-194.
- Hart, C.A., K. Barr, T. Makin, P. Brown, and R.W.I. Cooke. 1982. Characteristics of a β-lactamase produced by *Clostridium butyricum*. *J. Antimicrob. Chemother.* 10:31-35.
- Marrie, T.J., E.V. Haldane, C.A. Swantee, and E.A. Kerr. 1981. Susceptibility of anaerobic bacteria to nine antimicrobial agents and demonstration of decreased susceptibility of *Clostridium perfringens* to penicillin. *Antimicrob. Agents and Chemother.* 19:51-55.
- Salyers, A.A., J. Wong and T.D. Wilkins. 1977. β-lactamase activity in strains of *Bacteroides melaninogenicus* and *Bacteroides oralis*. *Antimicrob. Agents Chemother.* 11:142-146.
- Del Bene, V.E., and W.E. Farrar, Jr. 1973. Cephalosporinase activity in *Bacteroides fragilis*. *Antimicrob. Agents Chemother.* 3:369-372
- Timewell, R., E. Taylor, and I. Phillips. 1981. The β-lactamases of *Bacteroides* species. *J. Antimicrob. Chemother.* 7:137-146.
- Pechere, J.C., R. Guay, J. Dubois, and R. Letarte. 1980. Hydrolysis of cefotaxime by a β-lactamase from *Bacteroides fragilis*. *Antimicrob. Agents Chemother.* 17:1001-1003.
- Yotsuji, A., S. Minami, M. Inoue, and S. Mitsuhashi. 1983. Properties of novel β-lactamase produced by *Bacteroides fragilis*. *Antimicrob. Agents Chemother.* 24:925-929.
- Cuchural, G.J., F.P. Tally, N.V. Jacobus, P.K. Marsh, and J. W. Mayhew. 1983. Cefoxitin inactivation by *Bacteroides fragilis*. *Antimicrob. Agents Chemother.* 24:936-940
- Olsson, B., K. Dornbush, and C.E. Nord. 1979. Factors contributing to β-lactam antibiotics in *Bacteroides fragilis*. *Antimicrob. Agents Chemother.* 15:263-268.
- Lamonthe, F., F. Auger, and J.M. Lacroix. 1984. Effect of clavulanic acid on the activities of ten β-lactam agents against members of the *Bacteroides fragilis* group. *Antimicrob. Agents Chemother.* 25:662-665.
- Gabay, E.L., V.L. Sutter, and S.M. Finegold. 1981. Rapid β-lactamase testing in *Bacteroides*. *J. Antimicrob. Chemother.* 8:413-416.

29. Bush, K., and R.B. Sykes. 1984. β -lactamase (penicillinase, cephalosporinase), p. 280-285, 406, 407. *In* H.U. Bergmeyer (ed.) Methods of enzymatic analysis, 3rd ed, vol. IV. Verlag. Chemie, Deerfield Beach, Fla.

30. Sabath, L.D., F.F. Barrett, C. Wilcox, D.A. Gerstein, and M. Finland. 1969. Methicillin resistance of *Staphylococcus aureus* and *Staphylococcus epidermidis*, p. 302-306. *In* G.L. Hobby (ed.), Antimicrob. Agents Chemother. 1968. American Society for Microbiology, Washington, D.C.

31. Markowitz, S.M. 1980. Isolation of an ampicillin-resistant, non β -lactamase producing strain of *Haemophilus influenzae*. Antimicrob. Agents Chemother. 17:302-306.

32. Lee, D.T., and J.E. Rosenblatt. 1983. A comparison of four methods for detecting beta-lactamase activity in anaerobic bacteria, abstr. C302, p. 362. Abstr. Annu. Meet. Am. Soc. Microbiol. 1983.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

BD BBL Paper Discs for the Detection of β -Lactamase Enzymes Cefinase Discs

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Change History

Revision	Date	Change Summary
(03)	2018-09	Added Rx Only and IVD symbols.

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R_x Only This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."

	Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	Australian Sponsor: Becton Dickinson Pty Ltd. 4 Research Park Drive Macquarie University Research Park North Ryde, NSW 2113 Australia
	Benex Limited Pottery Road, Dun Laoghaire Co. Dublin, Ireland	

Nitrocefin is a product of Glaxo Research; distributed exclusively by BD Diagnostics.
ATCC is a trademark of American Type Culture Collection.
© 2018 BD. BD and the BD Logo are trademarks of Becton, Dickinson and Company.