Part Number:	L001148	BALTSO0191 Version 13.0 Template 4 Inserts	
Category and Description: Package Insert, BD BBL 10% Potassium Hydroxide Reagent Droppers	Rev from: 02	Rev to: 03	
	Job Number:		760-18

Catalog Number: 261191

Blank (Sheet) Size: Length: 11" Width: 22 1/2"


Number of Pages: 10 Number of Sheets: 1

Page Size: Length: 11" Width: 4 1/2" Final Folded Size: 4 1/2" X 1 7/8"

Ink Colors: Number of Colors: 3 PMS #: 2755 Blue; 032 Red; Standard Black

Printed Two Sides: Yes: ☒ No: ☐

Style (see illustrations below): # 4

Vendor Printed: ☒ Online / In House Printed: ☐ Web Printed: ☐

See Specification control no. NA for material information.

BD Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Company confidential. This document is the property of Becton, Dickinson and Company and is not to be used outside the company without written permission. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Revised By: **REVISED BY**
By Maria Zacharias at 5:12 pm, Sep 26, 2018

Proofing Approved By: **PROOFING APPROVED BY**
By Natalie Morio at 10:08 am, Sep 26, 2018

Third Eye By: **THIRD EYE BY**
By Nichole Graham at 12:43 pm, Sep 26, 2018

BD BBL™ 10% Potassium Hydroxide Reagent Droppers

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INTENDED USE

BD BBL™ 10% Potassium Hydroxide Reagent Droppers are intended for use in the examination of direct smears for fungal elements.

SUMMARY AND EXPLANATION

BD BBL 10% potassium hydroxide (KOH) facilitates the clearing of specimens for enhanced microscopic observation of fungi without altering the fungal elements present. It is most useful with skin, hair, nail and sputum specimens.¹

Each dropper is good for one day's use after the ampule has been broken.

PRINCIPLES OF THE PROCEDURE

When the smear is warmed, the KOH, aided by the solvent dimethyl sulfoxide, dissolves proteinaceous components, leaving intact the polysaccharide-containing fungal cell walls.^{1,2}

REAGENTS

BD BBL 10% Potassium Hydroxide Reagent Droppers contain 0.5 mL of 10% KOH with 1% dimethyl sulfoxide.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Danger

H314 Causes severe skin burns and eye damage. **P260** Do not breathe dust/fume/gas/mist/vapours/spray. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P303+P361+P353** IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage Instructions: Store at room temperature 15–30 °C (59–86 °F). Protect from light.

Product Deterioration: This reagent is hermetically sealed in an ampule, which affords protection of the solution from chemical instability until the expiration date. The reagent should be colorless but may contain clear glass-like particles floating in the solution. However, the reagent should not be used if KOH crystals are precipitated onto the wall of the ampule or are floating in the reagent.

PROCEDURE

Material Provided: BD BBL 10% Potassium Hydroxide Reagent Droppers.

Materials Required But Not Provided: Ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

Test Procedure

1. Hold reagent dropper upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to break ampule inside the dropper. **Caution: Break ampule close to its center *one time only*. Do not manipulate the dropper any further as the plastic may puncture and injury may occur.**
2. Tap bottom of dropper on tabletop a few times. Then invert for convenient drop-by-drop dispensing of reagent.
3. Dispense 1–3 drops of 10% KOH onto a clean microscope slide.
4. Add a small amount of specimen to be examined to the KOH.
5. Place coverslip over sample and tamp down lightly.
6. Warm slide gently to facilitate clearing of specimen.
7. Let smear sit for 5 min.
8. Examine under high dry or oil immersion objective with a brightfield or phase-contrast microscope.
9. If clearing is not complete, let smear sit an additional 5–10 min.

User Quality Control

1. Examine the reagent for signs of deterioration (see "Product Deterioration").
2. Smears of hair, skin or nails treated with 10% KOH should reveal partial cell disintegration. Water-treated smears will be unchanged.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CSLI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Fungal elements should become visible in thick, viscid or opaque specimens following treatment with 10% KOH.

LIMITATIONS OF THE PROCEDURE

Only experienced personnel should interpret direct examinations for fungal elements as background artifacts and animal cells may be misleading. It may be useful to compare preserved quality control specimens containing known fungi with the clinical specimen under examination.³

Clearing of thicker specimens may require an extended time.

Overheating of the smear may cause crystallization of the KOH.

PERFORMANCE CHARACTERISTICS

The Potassium Hydroxide (KOH) preparation (prep) has been documented as a reliable means of detecting fungi in clinical specimens. A study conducted during 1977-1982, involving 17,782 specimens, excluding those for dermatophyte cultures, were screened with the KOH prep and 7.6% were positive for yeast and or hyphae or pseudohyphae.⁴ Of the 11,938 specimens for dermatophyte culture, 17.1% had positive KOH preps.⁴ The percents of KOH positivity and culture positivity are summarized in the table below by specimen type.⁴

Specimen	KOH Positive (%)	Culture Positive (%)
Respiratory	55.9	98.7
Tissue	41.9	67
Urine	33.4	99.2
Vaginal	75	99.5
Sterile Body Fluids	20.3	98.4
Dermatophyte	42.7	90.5

AVAILABILITY

Cat. No. Description

261191 BD BBL™ 10% Potassium Hydroxide Reagent Droppers, 50.

REFERENCES

1. Coonrod, J.D., L.J. Kunz, and M.J. Ferraro (ed.). 1983. The direct detection of microorganisms in clinical samples. Academic Press, Inc., New York.
2. Baron, E.J., and S.M. Finegold. 1990. Bailey and Scott's diagnostic microbiology, 8th ed. The C.V. Mosby Co., St Louis.
3. Larsh, H.W., and N.L. Goodman. 1985. In E.H. Lennette, A. Balows, W.J. Hausler, Jr., and H.J. Shadomy (ed.), Manual of clinical microbiology, 4th ed. American Society for Microbiology, Washington D.C.
4. Quinn, P.R. and G.D. Roberts. 1984. Usefulness of the potassium hydroxide preparation for the detection of fungi in clinical specimens. Abstract F60, Abstracts of the Ann. Meet. of the Amer. Soc. Microbiol., St. Louis, 1984, p 302.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

 BD BBL 10% Potassium Hydroxide Reagent Droppers

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Benex Limited
Becton Dickinson Pty Ltd.