Part Number:	L001137	BALTSO0191 Version 13.0 Template 4 Inserts	
Category and Description: Package Insert, BD BBL Desoxycholate Reagent Droppers	Rev from: 02 Rev to: 03		
	Job Number: 760-18		

Catalog Number: 261183

Blank (Sheet) Size: Length: 11" Width: 22.5"


Number of Pages: 10 Number of Sheets: 1

Page Size: Length: 11" Width: 4.5" Final Folded Size: 4 1/2" X 1 7/8"

Ink Colors: Number of Colors: 1 PMS #: 2755 Blue

Printed Two Sides: Yes: ☒ No: ☐

Style (see illustrations below): # 4

Vendor Printed: ☒ Online / In House Printed: ☐ Web Printed: ☐

See Specification control no. NA for material information.

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Company confidential. This document is the property of Becton, Dickinson and Company and is not to be used outside the company without written permission. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Revised By: **REVISED BY**
By Maria Zacharias at 2:37 pm, Sep 27, 2018

Proofing Approved By: **PROOFING APPROVED BY**
By Natalie Morio at 10:00 am, Sep 27, 2018

Third Eye By: **THIRD EYE BY**
By Nichole Graham at 4:51 pm, Sep 27, 2018

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INTENDED USE

BD BBL™ Desoxycholate Reagent Droppers are used for the presumptive differentiation of pneumococci from other gram-positive cocci by the bile solubility test.

SUMMARY AND EXPLANATION

Sodium desoxycholate is a salt of purified bile acid and can be used in lower concentrations than other salts of bile acids. It is soluble in water and lyses pneumococcal cells treated with a 10% solution. The bile solubility test may be performed by adding the sodium desoxycholate solution to a cell suspension or directly to colonies on a blood agar plate.¹⁻³

Each dropper is good for one day's use after ampule has been broken.

PRINCIPLES OF THE PROCEDURE

The bile solubility test is based on the observation that cells of *Streptococcus pneumoniae* lyse when treated with a 10% solution of sodium desoxycholate, while other streptococci and gram-positive cocci are not bile soluble. Lysis occurs because bile soluble organisms contain an autolytic amidase that, when activated by bile salts, cleaves the bond between alanine and muramic acid in the cell wall.⁴

REAGENTS

BD BBL Desoxycholate Reagent Droppers contain 0.5 mL of 10% aqueous solution of sodium desoxycholate.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

This product contains dry natural rubber.

Follow proper laboratory procedures in handling and disposing of infectious materials. Avoid contact with the skin. Rinse thoroughly with water if spilled.

Storage Instructions: Store at room temperature 15–30 °C (59–86 °F). Protect from light.

Product Deterioration: This reagent is hermetically sealed in an ampule, which affords protection of the solution from chemical instability until the expiration date. Some flocculation of the reagent may occur below 23 °C. This does not indicate deterioration of the reagent; however, to eliminate flocculation ensure that the reagent is at a minimum of 23 °C prior to use.

PROCEDURE

Material Provided: BD BBL Desoxycholate Reagent Droppers.

Materials Required But Not Provided: Ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

Test Procedure and Results:

Hold reagent dropper upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to break ampule inside the dropper. **Caution: Break ampule close to its center one time only. Do not manipulate dropper any further as the plastic may puncture and injury may occur.** Tap bottom on tabletop a few times. When ready to use, invert for convenient drop-by-drop dispensing of reagent.

A. Direct Colony Method

1. Place 1 drop of reagent directly on a well-isolated 18–24 h colony.
2. Incubate plate at room temperature, agar side down, for 15 min.
3. Examine the colony for appearance. Colony flattening or deterioration is typical of a positive result.

B. Test Tube Method

1. Add 1.0 mL of sterile 0.85% saline to a 13 x 100 mm test tube.
2. Prepare a heavy suspension (approximately a 1.0–2.0 McFarland standard) of test organism. Vortex for 10 sec to obtain a uniform suspension. Adjust suspension to neutral pH, if necessary.
3. Divide the suspension into two 13 x 100 mm test tubes (0.5 mL/tube). Label one tube Test and the other Control.
4. Add 10 drops (0.5 mL) of sterile 0.85% saline to the tube labeled Control.
5. Add 10 drops of Desoxycholate Reagent to the tube labeled Test.
6. Gently mix each tube.
7. Incubate at room temperature for 15 min.
8. Compare Test result to Control. If Test appears clear or less turbid than Control, the result is positive. If Test appears equal to Control, the result is negative.

User Quality Control:

1. Positive (*Streptococcus pneumoniae* ATCC® 6303) and negative (*Enterococcus faecalis* ATCC 29212) controls should be run simultaneously with the organism to be tested.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

LIMITATIONS OF THE PROCEDURE

Use 18- to 24-h cultures to ensure accurate results. False negative results can occur when using older cultures.^{4,5}

False negative results may occur with the Direct Colony Method if the reagent contains flocculant material. Therefore, ensure that the reagent is at a minimum of 23 °C prior to use.

Failure to adjust the pH to neutrality prior to the addition of Desoxycholate Reagent in the Test Tube Method may result in false negative results due to precipitation of the reagent.^{1,2}

If only partial clearing occurs in the Test Tube Method, a second confirmatory test should be performed using a BD Pneumoslide™ test kit or BD Taxo™ P discs.

PERFORMANCE CHARACTERISTICS⁶

Downie, Stent and White tested 25 strains (19 smooth, 6 rough) of pneumococci with multiple dilutions of 10% Deoxycholate (Desoxycholate) acid using the tube method and reported the following results. Results were similar for all strains tested with the exception of one rough strain, which produced a positive result at a dilution of 1:20.

	Final Dilutions						
	1:100	1:200	1:400	1:800	1:1600	1:3200	1:6400
10% Desoxycholate acid	+	+	+	+	+	+	-
n = 24							

Availability

Cat. No. Description

261183 BD BBL™ Desoxycholate Reagent Droppers, 50.

REFERENCES

1. MacFaddin, J.F. 1980. Biochemical tests for the identification of medical bacteria, 2nd ed. Williams & Wilkins, Baltimore.
2. Lennett, E.H., A. Balows, W.J. Hausler, Jr., and H.J. Shadomy (ed.). 1985. Manual of clinical microbiology, 4th ed. American Society for Microbiology, Washington, D.C.
3. Baron, E.J., and S.M. Finegold. 1990. Bailey & Scott's diagnostic microbiology, 8th ed. The C.V. Mosby Company, St. Louis.
4. Joklik, W.K., H.P. Willett, D.B. Amos and C.M. Wilfert (ed.). 1988. Zinsser Microbiology, 19th ed. Appleton & Lange, Norwalk, Conn.
5. Washington, J.A. (ed.). 1981. Laboratory procedures in clinical microbiology. Springer-Verlag, New York.
6. Downie, A.W., L., Stent, and S.M., White, 1913. The bile solubility of pneumococcus, with special reference to the chemical structure of various bile salts. Br. J. Exp. Pathol., 12, 1.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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Benex Limited
Becton Dickinson Pty Ltd.