Part Number:	L001145	BALTSO0191 Version 13.0 Template 4 Inserts	
Category and Description:	Package Insert, BD BBL Ferric Chloride Reagent Droppers	Rev from: 02	Rev to: 03
		Job Number:	760-18

Catalog Number: 261190

Blank (Sheet) Size: Length: 11" Width: 18"


Number of Pages: 8 Number of Sheets: 1

Page Size: Length: 11" Width: 4 1/2" Final Folded Size: 4 1/2" X 1 7/8"

Ink Colors: Number of Colors: 3 PMS #: 2755 Blue; 032 Red; Standard Black

Printed Two Sides: Yes: ☒ No: ☐

Style (see illustrations below): # 4

Vendor Printed: ☒ Online / In House Printed: ☐ Web Printed: ☐

See Specification control no. NA for material information.

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Company confidential. This document is the property of Becton, Dickinson and Company and is not to be used outside the company without written permission. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Revised By: **REVISED BY**
By Maria Zacharias at 10:37 am, Sep 27, 2018

Proofing Approved By: **PROOFING APPROVED BY**
By Natalie Morio at 9:55 am, Sep 27, 2018

Third Eye By: **THIRD EYE BY**
By Nichole Graham at 4:37 pm, Sep 27, 2018

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INTENDED USE

BD BBL Ferric Chloride Reagent Droppers are intended for use in the differentiation of microorganisms capable of phenylalanine deamination.¹

SUMMARY AND EXPLANATION

BD BBL Ferric Chloride Reagent Droppers are used to determine if a specific microorganism is capable of producing phenylpyruvic acid from phenylalanine deamination.¹

Most members of the genera *Proteus* and *Providencia* are positive whereas most *Escherichia* and *Citrobacter* species are negative.²

PRINCIPLES OF THE PROCEDURE

Upon the addition of ferric chloride reagent to a culture of the microorganism, any phenylpyruvic acid present will react with the ferric salt to give a green color. The reaction is detected visually and is considered a positive test.¹

REAGENTS

BD BBL Ferric Chloride Reagent Droppers contain 0.5 mL of BD BBL 10% ferric chloride in aqueous solution.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Danger

H302 Harmful if swallowed. **H314** Causes severe skin burns and eye damage. **H318** Causes serious eye damage. **P260** Do not breathe dust/fume/gas/mist/vapours/spray. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P303 + P361 + P353** IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. **P305 + P351 + P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P405** Store locked up. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage Instructions: Store at room temperature 15–30 °C (59–86 °F).

Product Deterioration:

Reagents are hermetically sealed in an ampule, which affords protection of the solution from chemical instability until expiration date. The reagent droppers need no refrigeration. If the color of the reagent turns cloudy or if a definite precipitate forms, the reagent's effectiveness is not reduced.

PROCEDURE

Material Provided: BD BBL Ferric Chloride Reagent Droppers.

Materials Required But Not Provided: Ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

Test Procedure:

1. Hold upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to break ampule inside the dropper. **Caution: Break ampule close to its center one time only. Do not manipulate the dropper any further as the plastic may puncture and injury may occur.**
2. Tap bottom of ampule on tabletop. Then invert for convenient drop-by-drop dispensing of reagent directly onto a culture colony or into a test tube, depending upon the test performed and the recommended test procedure.
 - A. Phenylalanine Agar Method
Inoculate a slant of phenylalanine agar with the microorganisms in question and incubate for 24 h at 37 °C. Add 4 or 5 drops of ferric chloride test reagent to the growth on the slant. The presence of phenylpyruvic acid is indicated by the development of a green color in the syneresis fluid, diffusing into the agar slant within 5 min.
 - B. API – Follow directions supplied with API kit.

User Quality Control

Using the methods described above (Test Procedure), an authentic culture of *Proteus vulgaris* will yield a positive test within 5 min, while *Serratia marcescens* will yield a negative result.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

The test is positive if a visible green color develops upon addition of a solution of BD BBL 10% ferric chloride.

LIMITATIONS OF THE PROCEDURE

A heavy growth of organisms on the slant is essential. The reaction must be read promptly.

PERFORMANCE CHARACTERISTICS

Giammanco, Pignato and Agodi reported tests for the detection of phenylalanine deaminase and tryptophan deaminase for the differentiation of *Proteus* and *Providencia* several years ago.³ Out of 285 strains of different enterobacterial species, positive results were shown only by 75 *Proteus* and *Providencia* strains. The phenylalanine deaminase test is now used in the identification of *Enterobacteriaceae*. The deamination of phenylalanine is revealed by the formation of green color after addition of a ferric chloride solution to positive cultures, while deamination of tryptophan is revealed by a red color.

AVAILABILITY

Cat. No. Description

261190 BD BBL™ Ferric Chloride Reagent Droppers, 50.

REFERENCES

1. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 1998. Bailey & Scott's diagnostic microbiology, 10th ed. Mosby, Inc., St. Louis.
2. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.). 1999. Manual of clinical microbiology, 7th Edition, American Society for Microbiology, Washington, D.C.
3. Giammanco, G., S. Pignato, and A. Agodi. 1985. A simple chromogenic test for rapid screening of *Proteus* and *Providencia* bacteria. Microbiologica. 8:395-397.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

BD BBL Ferric Chloride Reagent Droppers

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Benex Limited
Becton Dickinson Pty Ltd.