Part Number:	L001122	BALTSO0191 Version 13.0 Template 4 Inserts	
Category and Description: Package Insert, BD BBL Acridine Orange Reagent Droppers	Rev from: 02 Rev to: 03		
	Job Number: 760-18		

Catalog Number: 261182

Blank (Sheet) Size: Length: 11" Width: 18"


Number of Pages: 8 Number of Sheets: 1

Page Size: Length: 11" Width: 4.5" Final Folded Size: 4.5" X 1.875"

Ink Colors: Number of Colors: 1 PMS #: 2755 Blue

Printed Two Sides: Yes: ☒ No: ☐

Style (see illustrations below): # 4

Vendor Printed: ☒ Online / In House Printed: ☐ Web Printed: ☐

See Specification control no. NA for material information.

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Company confidential. This document is the property of Becton, Dickinson and Company and is not to be used outside the company without written permission. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Revised By: REVISED BY
By Maria Zacharias at 10:20 am, Sep 27, 2018

Proofing Approved By: PROOFING APPROVED BY
By Natalie Morio at 9:53 am, Sep 27, 2018

Third Eye By: THIRD EYE BY
By Nichole Graham at 4:32 pm, Sep 27, 2018

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INTENDED USE

BD BBL™ Acridine Orange Reagent Droppers are used for fluorescent microscopic detection of microorganisms in direct smears.

SUMMARY AND EXPLANATION

BD BBL Acridine orange, a fluorochrome stain buffered at a low pH, differentiates between bright orange-stained bacteria and pale green to yellow-stained human cells and tissue.¹

Each dropper contains sufficient reagent to stain one slide.

PRINCIPLES OF THE PROCEDURE

BD BBL Acridine orange binds to nucleic acids of cells and bacteria. When viewed under UV light, single-stranded DNA and RNA fluoresce orange, whereas double-stranded DNA appears green. At low pH (3.5–4.0), bacteria and fungi stain bright orange. Cellular material stains pale green to yellow.² Nuclei of activated leukocytes may stain orange, yellow or red, depending upon the degree of increased RNA production. Erythrocytes either have no color or appear pale green.

REAGENTS

BD BBL Acridine Orange Reagent Droppers contain 0.5 mL of a 0.01% solution of [3,6-bis (dimethyl-amino) acridine hydrochloride] in a 0.5 M acetate buffer.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Follow proper laboratory procedures in handling and disposing of infectious materials. Avoid contact with skin. Rinse thoroughly with water if spilled.

Storage Instructions: Store at room temperature 15–30 °C (59–85 °F). Protect from light.

Product Deterioration: This reagent is hermetically sealed in an ampule, which affords protection of the solution from chemical instability until the expiration date. The staining solution should be orange, clear and without evidence of precipitation.

PROCEDURE

Material Provided: BD BBL Acridine Orange Reagent Droppers.

Materials Required But Not Provided: Ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

Test Procedure:

1. Prepare a smear of the specimen to be stained on a glass slide and air dry.
2. Fix the smear in methanol (50–100%) for 2 min, drain off excess methanol and allow smear to air dry.
3. Hold reagent dropper upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to break ampule inside the dropper.
Caution: Break ampule close to its center one time only. Do not manipulate the dropper any further as the plastic may puncture and injury may occur.
4. Tap bottom of dropper on tabletop a few times. Then invert for convenient drop-by-drop dispensing of reagent.
5. Flood the slide with acridine orange stain for 2 min.
6. Rinse thoroughly with tap water and allow to air dry.
7. Using a fluorescence microscope, examine smears initially at 100X to 400X, then confirm findings at 1000X using an oil immersion objective.

User Quality Control:

1. Examine the staining solution for signs of deterioration (see “Product Deterioration”).
2. Check the performance of the stain with fresh cultures of *Escherichia coli* and *Haemophilus influenzae*.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

Results

Microorganisms usually stain a brilliant orange on a black, light green or yellow background.

LIMITATIONS OF THE PROCEDURE

1. Acridine orange provides presumptive information on microorganisms which may be present in the specimen. Since microorganisms seen in smears may arise from external sources (i.e., collection devices, water used, slides, etc.), all positive smears should be confirmed by culture.
2. A minimum of approximately 10⁴ CFU/mL are required for detection by this method.³
3. Certain types of debris may fluoresce in acridine orange smears and can be distinguished on the basis of morphology when viewed at higher magnification.
4. Acridine orange staining does not distinguish between gram-positive and gram-negative organisms. Gram staining may be performed directly over the acridine orange-stained smear.

PERFORMANCE CHARACTERISTICS 4

Meseguer et al. conducted a study where a total of 1,592 blood cultures without macroscopic signs of bacterial growth in the first 12–24 hours of incubation were processed using both acridine orange stain and blind subculture. One hundred and twenty-one (7.6%) blood cultures were positive by either method; of these, 105 (86.8%) were detected by both acridine orange stain and subculture, 5 (4.1%) by subculture alone and 11 (9.1%) by acridine orange stain alone. The difference between the 116 blood cultures positive by acridine orange and the 110 blood cultures positive by subculture was not statistically significant ($p>0.1$).

Availability

Cat. No.	Description
261182	BD BBL™ Acridine Orange Reagent droppers, 50.

REFERENCES

1. Kronvall, G., and E. Myhre. 1977. Differential staining of bacteria in clinical specimens using acridine orange buffered at low pH. *Acta. Pathol. Microbiol. Scand. Sect. B* 85:249-254.
2. Lauer, B.A., L.B. Reller, and S. Mirrett. 1981. Comparison of acridine orange and Gram stains for detection of microorganisms in cerebrospinal fluid and other clinical specimens. *J. Clin. Microbiol.* 14:201-205.
3. McCarthy, L.R., and J.E. Senne. 1980. Evaluation of acridine orange stain for detection of microorganisms in blood cultures. *J. Clin. Microbiol.* 11:281-285.
4. Meseguer, M., L. de Rafael, M. Baquero, M. Martínez Ferrer, and M. Martínez López-Brea. 1984. Acridine orange stain in the early detection of bacteria in blood cultures. *Eur. J. Clin. Microbiol.* 3:113-115.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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Benex Limited
Becton Dickinson Pty Ltd.