
Revisions

BALTSO0191 Version 8.0 Template 4
Inserts

Rev from	Rev to	JOB #
01	02	7043-14

NOTES:

- BD Catalog Number: 261187
- Blank (Sheet) Size: Length: 11" Width: 27"
- Number of Pages: 12 Number of Sheets: 1
- Page Size: Length: 11" Width: 4.5" Final Folded Size: 4 1/2" x 1 7/8"
- Ink Colors: No. of Colors: 3 PMS#: 2755 Blue; 032 Red; Standard Black
- Printed two sides: Yes ☒ No ☐
- Style (see illustrations below): # 4

- See specification control no. N/A for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	REVISED BY By Nancy Carlsen at 1:32 pm, Apr 30, 2015	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.		 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	<i>Heckman</i> 04/01/15	Category and Description Package Insert, DMACA Indole Reagent Droppers		
Checked By	<i>KStone</i> 05-04-15			Scale: N/A
Part Number:	L001128			

BD BBL™ DMACA Indole Reagent Droppers

English: pages 1 – 2 Italiano: pagine 5 – 6
Français : pages 2 – 4 Português: páginas 7 – 8
Deutsch: Seiten 4 – 5 Español: páginas 8 – 9

L001128(02)
2015-02

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteys lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőjétől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

For the detection of indole production as an aid in the identification of aerobic, anaerobic, or facultatively anaerobic organisms.¹

SUMMARY AND EXPLANATION

BD BBL™ DMACA Indole Reagent Droppers are used in a rapid spot test for determining the ability of bacteria to produce indole by the deamination of tryptophan. The production of indole has been widely used as an aid in differentiating certain genera and species within the same genus. p-Dimethylaminocinnamaldehyde (DMACA) has been shown to be sensitive to the production of indole.² In addition, the color change with DMACA is both rapid and long lasting, making results easier to interpret than with more traditional reagents.³

Each dropper is good for one day's use after the ampule has been broken.

PRINCIPLES OF THE PROCEDURE

The reaction indicates the presence of the enzyme tryptophanase that reacts with tryptophan to produce indole. The indole produced reacts with DMACA to form a blue-green color.

REAGENTS

DMACA Indole Reagent Droppers contain 0.5 mL of aqueous 1% p-dimethylaminocinnamaldehyde in 10% hydrochloric acid.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Follow proper laboratory procedures in handling and disposing of infectious materials.

Danger

H314 Causes severe skin burns and eye damage. **H335** May cause respiratory irritation. **P103** Read label before use. **P260** Do not breathe dust/fume/gas/mist/vapours/spray. **P280** Wear protective gloves/protective clothing/eye protection/face protection **P303+P361+P353** IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. **P403+P233** Store in a well-ventilated place. Keep container tightly closed. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage Instructions: Store at room temperatures 15 – 30 °C (59 – 86 °F). Protect from light.

Product Deterioration: This reagent is hermetically sealed in an ampule, which affords protection of the solution from chemical instability until expiration date. Change in color of the reagent from light yellow to brown indicates improper storage, which may cause weaker reactions.

PROCEDURE

Material Provided: DMACA Indole Reagent Droppers.

Materials Required But Not Provided: Ancillary culture media, filter paper, reagents, quality control organisms and laboratory equipment as required for this procedure.

Test Procedure: Rapid Spot Method.

To use DMACA Indole Reagent Droppers, hold upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to break ampule inside the dropper.

Caution: Break ampule close to its center one time only. Do not manipulate dropper any further as the plastic may puncture and injury may occur. Tap bottom of the dropper on the tabletop a few times. When ready for use, invert for convenient drop-by-drop dispensing of reagent.

1. Moisten a piece of Whatman No. 1 filter paper (or equivalent) with a few drops of indole reagent.
2. Remove a well-isolated colony from an 18 – 24 h culture on a Blood Agar Plate with a sterile loop or wooden stick and smear on moistened filter paper.⁴
3. Observe for a blue to blue-green color within 2 min.

User Quality Control:

Examine the reagent for signs of deterioration (See "Product Deterioration").

Positive (*Escherichia coli* ATCC™ 25922) and negative (*Klebsiella pneumoniae* ATCC 33495) controls should be performed along with the organisms to be tested.

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CSLI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Positive: Blue to blue/green. Negative: No color to pink. Consult appropriate references for indole reactions with specific organisms.^{5,6}

LIMITATIONS OF THE PROCEDURE

Colonies to be tested should be isolated on media containing sufficient tryptophan, and no glucose⁷ (e.g., **BD Trypticase™** Soy Agar).

Allow at least 2 min for the absence of a reaction to develop before considering the indole reagent test negative.

Colonies should not be tested from media containing dyes (i.e., MacConkey or EMB agars). Dye carryover may give misleading results.

Colonies from mixed cultures should not be tested. Indole positive colonies have been shown to cause nearby (within 5 mm) indole negative colonies to appear weakly positive.⁸

False negative reactions may occur with some organisms, especially *Clostridium* spp. which form indole but break it down more rapidly.⁷

For *Flavobacterium* spp., a 48 h culture should be tested if 24 h culture was negative.⁴

PERFORMANCE CHARACTERISTICS

An internal study was performed comparing **BD BBL™** DMACA Indole Reagent and another commercially available DMACA Indole Reagent.⁹ As noted in the table, all reactions were identical for the 11 organisms tested using both reagents.

Organism		Expected Results	BD BBL Reagent	Commercial Reagent
<i>Enterobacter aerogenes</i>	ATCC 19615	–	–	–
<i>Acinetobacter anitratus</i>	ATCC 43608	–	–	–
<i>Klebsiella pneumoniae</i>	ATCC 33495	–	–	–
<i>Flavobacterium breve</i>	BD 8130	W+	W+	W+
<i>Providencia stuartii</i>	ATCC 33672	+	+	+
<i>Escherichia coli</i>	ATCC 25922	+	+	+
<i>Morganella morganii</i>	ATCC 25830	+	+	+
<i>Bacteroides thetaiotamicron</i>	ATCC 29741	+	+	+
<i>Weeksella virosa</i>	BD 5990	W+	W+	W+
<i>Klebsiella oxytoca</i>	BD 9214	+	+	+
<i>Vibrio parahaemolyticus</i>	ATCC 17802	W+	W+	W+

Key	
–	Negative
+	Positive
W+	Weak Positive

AVAILABILITY

Cat No. **Description**
261187 **BD BBL™** DMACA Indole Reagent Droppers, 50.

REFERENCES

1. Sutter, V.L., and W.T. Carter. 1972. Evaluation of media and reagents for indole-spot tests in anaerobic bacteriology. *Am. J. Clin. Pathol.* 58: 335-338.

2. Lowrance, B.L., P. Reich, and W.H. Traub. 1969. Evaluation of two spot-indole reagents. *Appl. Microbiol.* 17: 923-924.

3. Miller, J.M., and J.W. Wright. 1982. Spot indole test: evaluation of four reagents. *J. Clin. Microbiol.* 15: 589-592.

4. Pickett, M.J. 1989. Methods for identification of *Flavobacteria*. *J. Clin. Microbiol.* 27: 2309-2315.

5. Baron, E.J., and S.M. Finegold. 1990. Baily and Scott's diagnostic microbiology, 8th ed. The C.V. Mosby Company, St. Louis.

6. Balows, A., W.J. Hausler, Jr., K.L. Hermann, H.D. Isenberg, and H.J. Shadomy (ed.). 1991. Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.

7. MacFaddin, J.F. 2000. Biochemical tests for identification of medical bacteria. 3rd ed. Lippincott, Williams and Wilkins, Baltimore.

8. Vracko, R. and J.C. Sherris. 1963. Indole-spot test in bacteriology. *Am. J. Clin. Pathol.* 39: 429-432.

9. Data on file, BD Diagnostic systems.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

Australian Sponsor:

Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113
Australia

ATCC is a trademark of the American Type Culture Collection.

BD, BD Logo, and all other trademarks are property of Becton, Dickinson and Company.

© 2015 BD.