Part Number:	L0001480JAA	BALTS00191 Version 13.0 Template 4 Inserts	
Category and Description:	Package Insert, BD BBL™ DrySlide™ Indole	Rev from: 03	Rev to: 04
		Job Number:	839-18

Catalog Number: 231748

Blank (Sheet) Size: Length: 12" Width: 16.25"


Number of Pages: 10 Number of Sheets: 1

Page Size: Length: 12" Width: 3.25" Final Folded Size: 3" L X 3.25" W

Ink Colors: Number of Colors: 1 PMS #: Standard Black

Printed Two Sides: Yes: ☒ No: ☐

Style (see illustrations below): # 4

Vendor Printed: ☒ Online / In House Printed: ☐ Web Printed: ☐

See Specification control no. VSL-000148-0PR for material information.

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Company confidential. This document is the property of Becton, Dickinson and Company and is not to be used outside the company without written permission. Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Revised By: **REVISED BY**
By Sonia F.Thompson at 3:24 pm, Nov 05, 2018

Proofing Approved By: **PROOFING APPROVED BY**
By Terrence Means at 12:00 pm, Nov 06, 2018

Third Eye By: **THIRD EYE BY**
By Kim Stone at 12:45 pm, Nov 06, 2018

BD BBL™ DrySlide™ Indole

English: pages 1 – 2

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CE L0001480JAA(04)

IVD R_x Only 2018-11

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyts lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőtől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Inštrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

BD BBL™ DrySlide™ Indole is used for determining the indole reaction of bacteria.

SUMMARY AND EXPLANATION

The indole test is a qualitative procedure for determining the ability of bacteria to produce indole by reductive deamination of tryptophan.¹ Vracko and Sherris² showed that the indole test can be performed using filter paper moistened with *p*-dimethyl- aminobenzaldehyde (DMABA), an effective indicator of indole production.³⁻⁵

PRINCIPLES OF THE PROCEDURE

Indole, if present, reacts with *p*-dimethylaminobenzaldehyde to produce a pink-colored complex. The reaction is based on the pyrrole structure of indole, which permits a condensation reaction to occur between the CH₂ groups on the indole molecule and the aldehyde group of the indole reagent.¹

REAGENTS

BD BBL DrySlide Indole is a disposable slide with four filter paper reaction areas containing 5% *p*-dimethylaminobenzaldehyde (DMABA).

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Handle BD BBL DrySlide Indole only by the slide mount; do not touch the reaction area.

Follow established laboratory procedures in handling and disposing of infectious materials.

Storage Instructions: Store BD BBL DrySlide Indole at 2–8 °C. Store opened BD BBL DrySlide Indole pouches with unused slides at 15–30 °C.

Use BD BBL DrySlide Indole from opened, resealed pouches within 5 days. Discard unused slides after 5 days.

The expiration date applies to the product in its intact container when stored as directed.

Product Deterioration: Do not use a product if it fails to meet specifications for identity and performance.

SPECIMEN COLLECTION AND PREPARATION

Collect specimens or samples in sterile containers or with sterile swabs and transport immediately to the laboratory according to recommended guidelines.⁶⁻⁸

Process each specimen using procedures appropriate for that sample.⁶⁻⁸

Obtain the test organism from a primary isolation culture streaked on an appropriate solid medium (chocolate or tryptic soy blood agar).

Select well-isolated colonies. Mixed cultures require subculturing for purification. Insufficient growth indicates the need for further incubation.

Pick isolated colonies or a sweep of confluent growth from the pure culture.

PROCEDURE

Material Provided: BD BBL DrySlide Indole.

Materials Required But Not Provided: Inoculation loop (platinum or plastic) or wooden applicator stick and quality control organisms.

Test Procedure:

1. Equilibrate BD BBL DrySlide Indole to room temperature before opening the pouch.
2. Open the BD BBL DrySlide Indole pouch and remove the required number of slides.

After opening, the remaining BD BBL DrySlide Indole may be stored in the original pouch for up to 5 days at room temperature. Fold the top of the pouch over once and seal tightly with a self-adhesive sticker (provided). Discard unused slides 5 days after the package is opened.

Do not place used or partially used BD BBL DrySlide Indole (with organisms) back into the pouch.

3. Using an appropriate inoculation device, pick isolated colonies or a sweep of confluent growth from the culture to be tested.
4. Smear the specimen onto the reaction area of BD BBL DrySlide Indole.
5. Examine the reaction area for appearance of a pink color within 30 s.

User Quality Control:

Identity Specifications – Each pouch contains 3 slide mounts and 1 desiccant container. Each slide mount has 4 buff-yellow reaction areas.

Performance Response – Perform the spot indole test following the above procedure with the test organisms listed below.

Organism	ATCC®	Indole Reaction
<i>Bacteroides fragilis</i>	25285	–
<i>Escherichia coli</i>	25922	+
<i>Fusobacterium necrophorum</i>	25286	+
<i>Klebsiella pneumoniae</i>	13883	–

– buff-yellow + pink

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Indole-positive organisms change the color of the reaction area from buff-yellow to pink within 30 s. Indole-negative organisms produce no color change; the reaction area remains buff-yellow. A variable reaction may be seen in a small number of specimens due to marginal levels of tryptophan deaminase.

LIMITATIONS OF THE PROCEDURE

BD BBL DrySlide Indole may be used as an aid in organism identification. Additional biochemical testing is recommended for complete identification. For additional information, consult appropriate references.⁶⁻⁸

Indole-positive colonies have been reported to cause adjacent indole-negative colonies (within 5 mm) to appear falsely positive.^{2,5}

MacConkey agar has been reported to be unsuitable as the source of inoculum in spot test procedures.²⁻⁴

Some indole-positive strains of *Providencia* produce a color change different from the expected pink indole-positive reaction produced by other genera.^{2,3} An indole-positive reaction from these *Providencia* strains may appear brown/grey.

PERFORMANCE CHARACTERISTICS

Prior to release, all lots of BD BBL DrySlide Indole are tested for specific product specifications. A loopful of each of the following cultures is smeared onto the reaction area of samples of the lot: *Bacteroides fragilis* (ATCC 25285), *Fusobacterium necrophorum* (ATCC 25286), *Escherichia coli* (ATCC 25922), *Klebsiella oxytoca* (ATCC 8724), *Pseudomonas aeruginosa* (ATCC 27853) and *Providencia rettgeri* (ATCC 29944). The reaction areas of the slide are then examined for the appearance of a pink color within 30 s.

AVAILABILITY

Cat. No. Description

231748 BD BBL™ DrySlide™ Indole

REFERENCES

1. MacFadden, J.F. 1980. Biochemical tests for the identification of medical bacteria, 2nd ed. Williams and Wilkins, Baltimore.
2. Vracko, R., and J.C. Sherris. 1963. Indole spot-test in bacteriology. Am. J. Clin. Pathol. 39: 429-432.
3. Miller, J.M., and J.W. Wright. 1982. Spot indole test: evaluation of four reagents. J. Clin. Microbiol. 15: 589-592.
4. Lowrance, B.L., P. Reich, and W.H. Traub. 1969. Evaluation of two spot-indole reagents. Appl. Microbiol. 17: 923-924.
5. Peterson, W.C., D.C. Hale, and J.M. Matsin. 1982. An evaluation of the practicality of the spot-indole test for the identification of *Escherichia coli* in the clinical laboratory. Am. J. Clin. Pathol. 78: 755-758.
6. Leber, A.L (ed.) 2016. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology (ASM), Washington, D.C.
7. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., 2015. Manual of clinical microbiology, 11th ed. American Society for Microbiology (ASM), Washington, D.C.
8. Tille, P.M. & Forbes, B.A. 2014. Bailey & Scott's diagnostic microbiology, 13th ed. St. Louis, MO.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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R_x Only

This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."

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EC REP

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

Australian Sponsor:

Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113 Australia

ATCC is a trademark of the American Type Culture Collection.
DrySlide is a trademark of Difco Laboratories, a subsidiary of Becton, Dickinson and Company.

© 2018 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.