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231748 **BBL DrySlide** Indole, 25 x 3 lâminas.

BIBLIOGRAFIA

- MacFadden, J.F. 1980. Biochemical tests for the identification of medical bacteria, 2nd ed. Williams and Wilkins, Baltimore.
- Vracko, R., and J.C. Sherris. 1963. Indole spot-test in bacteriology. Am. J. Clin. Pathol. 39: 429-432.
- Miller, J.M., and J.W. Wright. 1982. Spot indole test: evaluation of four reagents. J. Clin. Microbiol. 15: 589-592.
- Lowrance, B.L., P. Reich, and W.H. Traub. 1969. Evaluation of two spot-indole reagents. Appl. Microbiol. 17: 923-924.
- Peterson, W.C., D.C. Hale, and J.M. Matsin. 1982. An evaluation of the practicality of the spot-indole test for the identification of *Escherichia coli* in the clinical laboratory. Am. J. Clin. Pathol. 78: 755-758.
- Isenberg, H.D. (ed.) 1994. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
- Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.). 1995. Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
- Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. Bailey & Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis.

Manufacturer / Výrobce / Producent / Fabrikant / Tootja / Valmistaja / Fabricant / Hersteller / Κατασκευαστής / Gyártó / Ditta produttrice / Gamintojas / Producent / Fabricante / Výrobca / Tillverkare / Производител / Producător / Üretici / Proizvođač / Производител / Аткарушы

Use by / Spotřebujte do / Anvendes før / Houdbaar tot / Kasutada enne / Viimeinkäyttöpäivä / A utiliser avant / Verwendbar bis / Ημερομηνία λήξης / Felhasználhatóság dátuma / Usare entro / Naudokite iki / Brukes før / Stosować do / Utilizar em / Použite do / Usar antes de / Använd före / Исполняйте до / A se utiliza până la / Son kulanma tarihi / Upotrebiti do / Исползовать до / дейи пайдаланура / Upotrijebiti do /
 YYYY-MM-DD / YYYY-MM (MM = end of month) /
 RRRR-MM-DD / RRRR-MM (MM = konec měsíce) /
 ÅÅÅÅ-MM-DD / ÅÅÅÅ-MM (MM = slutning af måned) /
 JJJJ-MM-DD / JJJJ-MM (MM = einde maand) /
 AAAA-KK-PP / AAAA-KK (KK = kuu lõpp) /
 VVVV-KK-PP / VVVV-KK (kuukauden loppuun mennessä) /
 AAAA-MM-JJ / AAAA-MM (MM = fin du mois) /
 JJJJ-MM-TT / JJJJ-MM (MM = Monatsende) /
 EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα) /
 ÉÉÉÉ-HH-NN / ÉÉÉÉ-HH (HH = hónap utolsó napja) /
 AAAA-MM-GG / AAAA-MM (MM = fine mese) /
 MMMM-MM-DD / MMMM-MM (MM = mensesio pabaiga) /
 ÅÅÅÅ-MM-DD / ÅÅÅÅ-MM (MM = slutten av måneden) /
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca) /
 AAAA-MM-DD / AAAA-MM (MM = fim do mês) /
 RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca) /
 aaaa-mm-dd / aaaa-mm (mm = fin del mes) /
 AAAA-MM-DD / AAAA-MM (MM = slutet på månaden) /
 ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = края на месеца) /
 AAAA-LL-ZZ / AAAA-LL (LL = sfârșitul lunii) /
 YYYY-AA-GG / YYYY-AA (AA = ayin sonu) /
 GGGG-MM-DD / GGGG-MM (MM = kraj meseca) /
 ГГГГ-ММ-ДД / ГГГГ-ММ (ММ = конец месяца) /
 ЖЖЖЖ-АА-КК / ЖЖЖЖ-АА (АА = айдың соңы) /
 GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)

Catalog number / Katalogové číslo / Katalognummer / Catalogusnummer / Kataloogi number / Tuotenumero / Numéro catalogue / Bestellnummer / Αριθμός καταλόγου / Katalógusszám / Numero di catalogo / Katalogo numeris / Numer katalogowy / Número do catálogo / Katalogové číslo / Número de catálogo / Каталоген номер / Număr de catalog / Katalog numarası / Kataloški broj / Номер по каталогу / Каталог нөмірі

Authorized Representative in the European Community / Autorizovaný zástupce pro Evropskou unii / Autoriseret repræsentant i EU / Erkend vertegenwoordiger in de Europese Unie / Volitatud esindaja Euroopa Nõukogus / Valtuutettu edustaja Euroopan yhteisössä / Représentant agréé pour la C.E.E. / Autorisierte EG-Vertretung / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Hivatalos képviselő az Európai Unióban / Rappresentante autorizzato nella Comunità europea / Jgaliotasis atstovas Europos Bendrijoje / Autorisert representant i EU / Autoryzowane przedstawicielstwo w Unii Europejskiej / Representante autorizado na União Europeia / Autorizovaný zástupca v Európskom spoločenstve / Representante autorizado en la Comunidad Europea / Auktoriserad representant i EU / Оторизирани представител в EU / Reprezentant autorizat în Uniunea Europeană / Avrupa Topluluğu Yetkilii Temsilcisi / Ovlašteni predstavnik u Evropskoj zajednici / Уполномоченный представитель в Европейском сообществе / Европа қауымдастығындағы уәкілетті өкіл / Autorizuirani predstavnik u EU

In Vitro Diagnostic Medical Device / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medisch hulpmiddel voor in vitro diagnose / In vitro diagnostika meditsiiniaparatuur / Lääkinnällinen in vitro -diagnostiikkalaite / Dispositif médical de diagnostic in vitro / Medizinisches In-vitro-Diagnostikum / In vitro διαγνωστική ιατρική συσκευή / In vitro diagnosztikai orvosi eszköz / Dispositivo medico diagnostico in vitro. / In vitro diagnostikos prietaisais / In vitro diagnostisk medisinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Medicinska pomôcka na diagnostiku in vitro / Dispositivo médico de diagnóstico in vitro / Medicinsk anordning för in vitro-diagnostik / Медицински уред за диагностика ин витро / Aparatură medicală de diagnosticare in vitro / In Vitro Diagnostik Tibbi Cihaz / Medicinski uređaj za in vitro dijagnostiku / Медицинский прибор для диагностики ин витро / Жасанды жағдайда жүргізетін медициналық диагностика аспабы / Medicinska pomagala za In Vitro Dijagnostiku

Temperature limitation / Teplotní omezení / Temperaturbegrænsning / Temperatuurlimiet / Temperatuuri piirang / Lämpötilarajoitus / Température limite / Zulässiger Temperaturbereich / Οριο θερμοκρασίας / Hőmérsékleti határ / Temperatura limite / Laikymo temperatūra / Temperaturbegrænsning / Ograniczenie temperatury / Limitação da temperatura / Ograničenje teploty / Limitación de temperatura / Temperaturbegrænsning / Температурни ограничения / Limitare de temperatură / Sıcaklık sınırlaması / Ograničenje temperature / Ограничение температуры / Температураны шектеу / Dozvoljena temperatura

Batch Code (Lot) / Kód (číslo) šarže / Batch code (Lot) / Chargennummer (lot) / Partii kood / Eräkoodi (LOT) / Code de lot (Lot) / Chargencode (Chargenbezeichnung) / Κωδικός παρτίδας (Παρτίδα) / Tétel száma (Lot) / Codice del lotto (partita) / Partijos numeris (Lot) / Batch-code (Serie) / Kod partii (seria) / Código do lote (Lote) / Kód série (šarža) / Código de lote (Lote) / Satskod (parti) / Код (Партида) / Număr lot (Lotul) / Parti Kodu (Lot) / Kod serije / Код партии (лот) / Топтама коды / Lot (kod)

Contains sufficient for <n> tests / Dostatečné množství pro <n> testů / Indeholder tilstrækkeligt til <n> test / Voldoende voor <n> tests / Küllaldane <n> testide jaoks / Sisältöön riittävä <n> testejä varten / Contenu suffisant pour <n> tests / Ausreichend für <n> Tests / Περιέχει επαρκή ποσότητα <n> εξετάσεις / <n> teszthez elegendő / Contenuto sufficiente per <n> test / Pakankamas kiekis atlikti <n> testų / Innholder tilstrekkelig for <n> tester / Zawiera ilość wystarczającą do <n> testów / Contêmo suficiente para <n> testes / Obsah vystačí na <n> testov / Contenido suficiente para <n> pruebas / Räckertill <n> antal tester / Съдържанието е достатъчно за <n> теста / Conține suficient pentru <n> teste / <n> testleri için yeterli miktarda içerir / Sadržaj dovoljan za <n> testova / Достаточно для <n> тестов(а) / <n> тесттері үшін жеткілікті / Sadržaj za (n) testova

Becton, Dickinson and Company
 7 Loveton Circle
 Sparks, MD 21152 USA
 800-638-8663
 www.bd.com/ds

Benex Limited
 Pottery Road, Dun Laoghaire
 Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection.

DrySlide is a trademark of Difco Laboratories, a subsidiary of Becton, Dickinson and Company.

BD, BD Logo and BBL are trademarks of Becton, Dickinson and Company. © 2013 BD.