BD Cepacia Medium BD OFPBL Agar

NAMJENA

BD Cepacia Medium (Hranjiva podloga BD Cepacia) i BD OFPBL Agar selektivne su diferencijalne podloge za izolaciju bakterije *Burkholderia cepacia* iz kliničkih uzoraka.

NAČELA I OBJAŠNJENJE POSTUPKA

Mikrobiološka metoda.

Burkholderia (prethodno poznata pod imenom *Pseudomonas*) *cepacia* opće je poznata kao bolnički patogen koji obično izaziva infekcije povezane s kontaminiranom opremom, sredstvima za dezinfekciju i lijekovima. Infekcije uključuju bakterijemiju, infekcije urinarnog i respiratornog trakta i ostale.¹ Uz to, organizam je važan patogen kod pacijenata s cističnom fibrozom (CF) koja se još naziva mukoviscidoza¹⁻³ te kod pacijenata s kroničnom granulomatoznom bolesti.⁴

Oporavak organizma na rutinskoj podlozi, poput krvnog agar ili agar MacConkey, ponekad je težak jer ostali organizmi poput *Staphylococcus* ili *Pseudomonas aeruginosa* često mogu prekomjerno narasti i prikriti njegovu prisutnost.^{1,2,5}

BD OFPBL (oksidacija/fermentacija-polimiksin-bacitracin-laktoza) **Agar** koristi se za izolaciju *B. cepacia*, a zabilježeno je da je manje selektivan od **BD Cepacia Medium**.^{1,6}

BD OFPBL Agar temelji se na O/F podlozi (podlozi za oksidaciju/fermentaciju). Indikator pH vrijednosti bromtimol plavo mijenja boju iz plave u žutu kad laktoza fermentira u kiseline. Fosfat se dodaje za stabilizaciju pH vrijednosti. Bacitracin i polimiksin B djeluju kao inhibitori za prateće bakterije.

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REFERENCE

1. Gilligan, P.H., G. Lum, P.A.R. Vandamme, and S. Whittier. 2003. *Burkholderia, Stenotrophomonas, Ralstonia, Brevundimonas, Comamonas, Delftia, Pandoraea, and Acidovorax*. In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Yolken (ed.). Manual of clinical microbiology, 8thed. American Society for Microbiology, Washington, D.C.
2. Gilligan, P.H. 1991. Microbiology of airway disease in patients with cystic fibrosis. Clin. Microbiol. Rev. 4: 35-51.
3. Gilligan, P.H., and D.V. Schidlow. 1984. The role of *Pseudomonas cepacia* in pulmonary disease of cystic fibrosis patients. Clin. Microbiol. Newsl. 6: 42-44.
4. O'Neil, K.M., et al. 1986. *Pseudomonas cepacia*: an emerging pathogen in chronic granulomatous disease. J. Pediatr. 108: 940-942.
5. Gilligan, P.H., et al. 1985. Isolation medium for the recovery of *Pseudomonas cepacia* from respiratory secretions of patients with cystic fibrosis. J. Clin. Microbiol. 25: 5-8.
6. Welch, D.F., et al. 1987. Selective and differential medium for recovery of *Pseudomonas cepacia* from the respiratory tracts of patients with cystic fibrosis. J. Clin. Microbiol. 25: 1730-1734.
7. Carson, L.A., et al. 1988. Comparative evaluation of selective media for isolation of *Pseudomonas cepacia* from cystic fibrosis patients and environmental sources. J. Clin. Microbiol. 26: 2096-2100.
8. Tablan, O.C., et al. 1987. Laboratory proficiency test results on use of selective media for isolating *Pseudomonas cepacia* from simulated sputum specimens of patients with cystic fibrosis. J. Clin. Microbiol. 25: 485-487.
9. Christenson, J.C., et al. 1989. Recovery of *Pseudomonas gladioli* from respiratory tract specimens of patients with cystic fibrosis. J. Clin. Microbiol. 27: 270-273.

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<http://www.bd.com>

<http://www.bd.com/europe/regulatory/>

ATCC is a trademark of the American Type Culture Collection

BD, BD logo and Stacker are trademarks of Becton, Dickinson and Company.

© 2014 BD

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Becton Dickinson GmbH