
Revisions

BALTSO0191 Version 10.0 Template 4
Inserts

Rev from	Rev to	JOB #
03	04	8464-16

NOTES:

- BD Catalog Number: 440910, 440911, 441951
- Blank (Sheet) Size: Length: 18" Width: 28.75"
- Number of Pages: 12 Number of Sheets: 1
- Page Size: Length: 18" Width: 4.75" Final Folded Size: 2.3 1/16" x 4.75"
- Ink Colors: No. of Colors: 1 PMS#: Blue 2755
- Printed two sides: Yes ☒ No ☐
- Style (see illustrations below): # 4

- Vendor Printed ☒ Online/In House Printed ☐ Web ☐
- See specification control no. N/A for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	REVISED BY By Tori Pagani at 3:43 pm, Jul 19, 2016	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	PROOFING APPROVED BY By Helaine Jones at 3:57 pm, Jul 19, 2016		
Checked By	THIRD EYE BY By Terrence Means at 4:09 pm, Jul 19, 2016		
Part Number: 8012507	Category and Description Package Insert, BD PhoenixSpec Nephelometer User's Guide	Sheet: 1 of 13	
		Scale: N/A	

BD PhoenixSpec™ Nephelometer

User's Guide

8012507(04)
2016-07

English: pages 1 – 2 Português: páginas 7 – 8
Français : pages 2 – 4 Español: páginas 8 – 10
Deutsch: Seiten 4 – 5 Japanese: ページ 10 – 11
Italiano: pagine 5 – 7

Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyts lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőjétől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o reprezentante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin.

INTENDED USE

The **BD PhoenixSpec™** Nephelometer is a portable device designed to measure turbidity of microbial suspensions equivalent to McFarland standards 0.10 to 4.50. The instrument may be used for measurement of inoculum density for the **BD BBL™ Crystal™** System and the **BD Phoenix™** System. The instrument is battery operated or can be used with an AC adapter.

SUMMARY AND EXPLANATION

The visual comparison of organism suspensions to turbidity standards is an accepted method of estimating organism densities. The most widely accepted standard is the McFarland standard.¹ A McFarland standard is prepared by adding barium chloride to aqueous sulfuric acid. The density of the resulting barium sulfate precipitate can be used to approximate the colony count of a prepared suspension; e.g., McFarland 1 is the equivalent of approximately 3.0 x 10⁸ colony forming units (CFU)/mL of *E. coli*, ATCC® 25922. Other standards have been used for density measurements, including titanium dioxide² and latex particle suspensions.³ The **BD PhoenixSpec** Nephelometer uses latex particle suspensions as calibration standards that do not require shaking before use.

The instrumented measurement of turbidity relies on the ability of particles to scatter light while in suspension. The measurement of this scattered light is referred to as nephelometry.⁴ In order to obtain an accurate density measurement, a reliable method of calibration must be used.

PRINCIPLES OF THE PROCEDURE

The **BD PhoenixSpec** Nephelometer is a solid state instrument that uses a tungsten-filament lamp as a light source, a 90° detector to monitor scattered light and a transmitted light detector. The instrument's microprocessor calculates the ratio of the signals for the 90° and transmitted light detectors. The ratiometric technique compensates for color and/or light absorbing materials, and for fluctuations in lamp intensity.

Before making measurements, the **BD PhoenixSpec** instrument is calibrated with the **BD PhoenixSpec** or **BD PhoenixSpec** AP Calibrator Kit. To determine the McFarland equivalent of a microbial suspension, a tube is placed in the reading chamber and the test button is pressed. The results are displayed on the LCD screen in McFarland units. On a daily basis, the 0.25 and 0.5 calibrators in the **BD PhoenixSpec** Calibrator Kit, are read as a sample to assure proper operation. On a daily basis, the 0.25, 0.5 and 2.0 calibrators in the **BD PhoenixSpec** AP Calibrator Kit, are read as a sample to assure proper operation.

SPECIFICATIONS

Range.....	McFarland 0.10–4.50	Ambient Operating	
Accuracy.....	0.5 ± 0.08 McFarland 2.0 ± 0.1 McFarland	Conditions.....	20 °C to 30 °C 0 to 90% RH at 30 °C 0 to 80% RH at 40 °C 0 to 70% RH at 50 °C
Repeatability.....	0.5 ± 0.04 McFarland 2.0 ± 0.1 McFarland		
Universal		Storage Temperature Range	
AC Adapter	Input:100-240VAC; 50-60 Hz, 0.5A Output:9.0VDC, 2.0A, 18W	Nephelometer:.....	-40 °C to 60 °C (-40 °F to 140 °F)
Battery	Four Size AA Alkaline Cells	Calibrators:.....	5 °C to 60 °C
Battery Life-		Sample Tube Size	BD™ "L" tube, 16 mm diameter, 75 mm height,
in-Use, # Tests.....	2300 Tests	Lamp Life.....	~ 400,000 readings

Warnings and Precautions

For *in vitro* Diagnostic Use.

1. The **BD PhoenixSpec** Nephelometer is designed to work with **BD "L"** tubes (16 x 75 mm) only. The minimum acceptable fill volume is 2.0 mL.
2. The **BD PhoenixSpec** Nephelometer is designed for use with **BD Phoenix** and **BD BBL Crystal** Systems only.

PROCEDURE

Materials Provided: **BD PhoenixSpec** Nephelometer, AC Adapter and 4 AA Alkaline batteries (not rechargeable). The **BD PhoenixSpec** Calibrator Kit (0.25, 0.5, 1.0 and 4.0 McFarland units) and **BD PhoenixSpec** AP Calibrator Kit (0.25, 0.5, 1.0, 2.0 and 4.0 McFarland units) are available separately - see **AVAILABILITY**.

Materials Required But Not Provided: Ancillary culture media and equipment required to prepare a bacterial suspension.

Instructions:

BATTERY INSTALLATION AND REPLACEMENT – Remove the instrument and batteries from the shipping box. Remove the battery compartment cover from the instrument bottom and install the batteries. Correct battery polarity is shown on the battery holder. Reinstall the battery compartment cover.

USING THE AC ADAPTER – Plug the adapter jack into the instrument's connector. The AC Adapter may be used with or without the batteries installed. The AC Adapter will not recharge the installed batteries. Position the equipment so the disconnecting device is easily accessible.

CALIBRATION PROCEDURE –

Notes:

- The instrument should be calibrated every three months, when the daily calibration verification check of the 0.25, 0.5 and/or 2.0 **BD PhoenixSpec** Calibrators read outside the 0.23–0.27, 0.45–0.55 and 1.8–2.2 McFarland range, respectively, or as dictated by the laboratory's experience or policy.
- Use only the **BD PhoenixSpec** or **BD PhoenixSpec** AP Calibrator Kit supplied for use with the **BD PhoenixSpec** Instrument.
- Keep the **BD PhoenixSpec** Calibrators in the supplied packaging to avoid scratching the tubes, as could happen in a metal test tube rack, for example. Scratches or dirt on the surface of the tubes can affect readings.
- Do not use the **BD PhoenixSpec** Calibrators beyond the expiration date which appears on the label.
- The instrument should be placed on a stationary, level surface.
- If an error code flashes on the LCD screen during the procedure, refer to the Troubleshooting section.
- The **BD PhoenixSpec** Calibrators do not require shaking before use.

1. Press the power button  to turn instrument ON.
2. Insert the 0.25 **BD PhoenixSpec** Calibrator into the reading chamber. Press the Calibration button. 
3. The screen will flash S1 (Standard 1).
4. Press the Read button.  The screen will count backwards from 10.
5. The screen will flash S2. Remove the 0.25 Calibrator and insert the 1.0 Calibrator. Press the Read button. 
6. The screen will count backwards from 10, then flash S3. Remove the 1.0 Calibrator, insert the 4.0 Calibrator and press the Read button. 
7. The screen will count backwards from 10, then go back to S1. Press the calibration button. 
8. Remove the 4.0 Calibrator. Calibration is now complete.

Note: If errors occur, refer to the TROUBLESHOOTING section.

TESTING PROCEDURES –

Notes:

- The instrument should be placed on a stationary, level surface.
- On a daily basis, the 0.25 and 0.5 **BD PhoenixSpec** Calibrators for the **BD PhoenixSpec** Calibrator Kit or the 0.25, 0.5 and 2.0 **BD PhoenixSpec** Calibrators for the **BD PhoenixSpec** AP Calibrator Kit must be read as samples and the readings must fall between 0.23–0.27 (for the 0.25), 0.45–0.55 (for the 0.5) and 1.8–2.2 (for the 2.0). If either reading falls outside its corresponding range, recalibrate the **BD PhoenixSpec** instrument.
- The **BD PhoenixSpec** Calibrators do not require shaking before use.

1. Press the power button to turn the instrument ON. 
2. Vortex the capped test sample and let bubbles disperse for 10 seconds.
3. Lift the lid, and insert the tube.
4. Press the Read button. 

5. Read the display after the lamp symbol turns off.
6. Remove the tube from the reading chamber.
7. Adjust the sample suspension with more organisms if the readout is lower than expected. Vortex the sample and reread.

Note: Refer to the BD Phoenix System User's Manual for the proper sample preparation for the BD Phoenix System if the reading exceeds the expected McFarland range or refer to the BD BBL Crystal system insert if the reading exceeds the expected McFarland range.

8. Repeat step 7 until the desired McFarland density is achieved.

9. When testing is complete, press the power button to turn the instrument OFF. 

RESULTS

The **BD PhoenixSpec** instrument displays the measurement in McFarland units. These units represent the optical density of the suspension and may be used to estimate the CFU/mL of the suspension tested.

LIMITATIONS OF THE PROCEDURE

1. The **BD PhoenixSpec** instrument must be calibrated with the **BD PhoenixSpec** Calibrators.
2. The **BD PhoenixSpec** Calibrators are to be used only for calibration of the **BD PhoenixSpec** Nephelometer and not as a visual, manual approximation to a barium sulfate McFarland standard.
3. Do not use the **BD PhoenixSpec** Nephelometer for measurements outside the range of 0.1–4.5 McFarland units.
4. Volumes below 2.0 mL must not be measured in the instrument.
5. **BD “L”** tubes - 16 mm diameter x 75 mm - must be used for samples.
6. Do not operate the instrument in direct sunlight.

PERFORMANCE CHARACTERISTICS

To determine the accuracy of the **BD PhoenixSpec** Nephelometer, colony counts were performed on suspensions of *E. coli* ATCC 25922 prepared in saline. Each test dilution was performed in replicates of six and the average of the plate counts is given in Table 1.

Table 1

McFarland	Expected CFU/mL x 10 ⁸	Adjusted* CFU/mL x 10 ⁸	% CV
0.25	0.75	0.9	14.1
0.5	1.5	1.7	13.2
1	3.0	3.2	7.4
2	6.0	6.1	4.2
3	9.0	9.3	3.2
4	12.0	12.6	2.7

*Because of difficulty in achieving the exact target McFarland density, these results were obtained by using the **BD PhoenixSpec** readings and adjusting the observed CFU/mL as if each dilution had been the target McFarland.

WARRANTY

The **BD PhoenixSpec** Nephelometer is warranted to be free from defects for one year from date of purchase.

MAINTENANCE

The general maintenance of the **BD PhoenixSpec** Instrument includes:

1. Wipe the exterior with a lint-free towel dampened with a mild, anti-bacterial solution.
2. Replace the batteries as needed.

TROUBLESHOOTING

Error messages

- E-2 Two calibrators that were read were too close together in value (i.e., the same calibrator was read twice during the calibration procedure).
- E-3 Low light error – check for obstructed light path.
- E-4 Memory failure. Press the on/off key and if the error reappears, contact service.
- E-5 A/D overrange. Check for obstructed light path.
- E-6 A/D underrange. Check for obstructed light path.
- E-7 Light leak. Ensure the tube is capped with a black cap. Press the tube down into the well to ensure it is fully inserted.
- E-8 Bad lamp or bad lamp circuit. Contact BD Technical Services.
- Cal? Do not use nephelometer. Contact BD Technical Services.

Low Battery Replace all four AA Alkaline batteries when the  icon is flashing at the lower left corner on the display.

Note: If errors E-4, E-5, or E-6 appear consistently, an internal hardware problem is indicated and BD Diagnostic Systems Product Support should be contacted.

AVAILABILITY

Cat. No.	Description
440910	BD PhoenixSpec™ Nephelometer
440911	BD PhoenixSpec™ Calibrator Kit (contains 0.25, 0.5, 1.0 and 4.0 calibrators)
441951	BD PhoenixSpec™ AP Calibrator Kit (contains 0.25, 0.5, 1.0, 2.0 and 4.0 calibrators)
440986	BD PhoenixSpec™ Universal AC Adapter

REFERENCES

1. McFarland, J. 1907. The nephelometer: an instrument for estimating the number of bacteria in suspensions used for calculating the opsonic index for vaccines. *JAMA* 49:1176-1178.
2. Roessler, W.G., and C.R. Brewer, 1967. Permanent turbidity standards. *Appl. Microbiol.* 15:1114-1121.
3. Pugh, T.L., and W. Heller, 1957. Density of polystyrene and polyvinyl toluene latex particles. *J. Colloid Sci.*, 12:173-180.
4. Mallette, M.F. 1969. XV. Evaluation of growth by physical and chemical means, p. 521-566. *In* J.R. Norris and D.W. Ribbons (ed.), *Methods in microbiology*, vol 1. Academic Press Inc., New York.

Technical Information: In the United States contact BD Technical Service and Support at 1.800.638.8663 or www.bd.com.

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1. McFarland, J. 1907. The nephelometer: an instrument for estimating the number of bacteria in suspensions used for calculating the opsonic index for vaccines. JAMA 49:1176-1178.
2. Roessler, W.G., and C.R. Brewer, 1967. Permanent turbidity standards. Appl. Microbiol. 15:1114-1121.
3. Pugh, T.L., and W. Heller, 1957. Density of polystyrene and polyvinyl toluene latex particles. J. Colloid Sci., 12:173-180.
4. Mallette, M.F. 1969. XV. Evaluation of growth by physical and chemical means, p. 521-566. In J.R. Norris and D.W. Ribbons (ed.), Methods in microbiology, vol 1. Academic Press Inc., New York.

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ATCC is a trademark of the American Type Culture Collection.
© 2016 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.

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Benex Limited Pottery Road, Dun Laoghaire Co. Dublin, Ireland