


Revisions

SO 0191-5


Rev from	Rev to	ECO #
0106	2010/07	5393-10

Notes:


- BD Cat. Number 295697
- Blank (Sheet) Size: Length: 11" Width: 8.5"
 Number of Pages: 2 Number of Sheets: 1
 Page Size: Length 11" Width 8.5" Final Folded Size: 4 1/4" x 1 3/8.
- Style (see illustrations below): # 1




#1




#2




#3




#4



#5




#6



#7

- See Specification Control Number N/A for Material Information
- Ink Colors: Printed two sides ☒ Yes ☐ No
 No. of Colors: 1 PMS# Black
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 8806871JAA		Category and Description Package Insert, Dubos Broth, Enriched	Sheet: 1 of 3	A
			Scale: N/A	

BD BBL™ Prepared Tubed Medium for Cultivation of Mycobacteria

Dubos Broth, Enriched



8806871JAA
2010/07

Pokyny vám poskytne miestni zástupce spoločnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіндік жергілікті BD өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

INTENDED USE

Dubos Broth, Enriched, is used in the cultivation of pure cultures of mycobacteria, particularly *Mycobacterium tuberculosis*.

SUMMARY AND EXPLANATION

Dubos Broth, Enriched, is a modified medium based on the formulation of Dubos et al.¹ This formulation differs from the original in that it has a strong buffering system and an acid pH.²

The particular value of Dubos Broth is that it provides dispersed growth, free of excessive clumps, which can be used to prepare a relatively uniform suspension of mycobacteria for use in bacterial studies.

It is also used as a subculture and enrichment medium for the rapid cultivation of *M. tuberculosis* and other mycobacterial species from treated clinical specimens and from direct inoculation of specimens that may yield pure cultures; e.g., spinal, pleural and peritoneal fluids.³

PRINCIPLES OF THE PROCEDURE

Dubos Broth, Enriched, contains enzymatic digest of casein and an amino acid, L-asparagine, as sources of nutrients. A variety of inorganic salts provide ions required for the metabolism of mycobacteria. Polysorbate 80, an oleic acid ester, supplies essential fatty acids for the replication of mycobacteria. Bovine albumin acts as a protective agent by binding free fatty acids that may be toxic to mycobacteria. The albumin is heat-treated to inactivate lipase, which may release fatty acids from the polysorbate 80. Phosphate buffers maintain the pH of the medium.

REAGENTS

Dubos Broth, Enriched

Approximate Formula* Per 900 mL Purified Water

Pancreatic Digest of Casein	0.5 g
L-Asparagine	2.0 g
Monopotassium Phosphate	1.0 g
Disodium Phosphate	2.5 g
Ferric Ammonium Citrate	0.05 g
Magnesium Sulfate	0.01 g
Polysorbate 80	0.2 g
Calcium Chloride	0.5 mg
Zinc Sulfate	0.1 mg
Copper Sulfate	0.1 mg
Bovine Albumin	100.0 mL

*Adjusted and/or supplemented as required to meet performance criteria.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Tubes with tight caps should be opened carefully to avoid injury due to breakage of glass.

Storage Instructions: On receipt, store tubes in the dark at 2 to 8°C. Avoid freezing and overheating. Do not open until ready to use. Minimize exposure to light. Tubed media stored as labelled until just prior to use may be inoculated up to the expiration date and incubated for up to 8 weeks. Allow the medium to warm to room temperature before inoculation.

Product Deterioration: Do not use tubes if they show evidence of microbial contamination, discoloration, precipitation, evaporation, or other signs of deterioration.

SPECIMEN COLLECTION AND HANDLING

This medium is not suitable for use directly with clinical specimens, except as stated above or as a "backup" enrichment medium in addition to isolation media. Consult appropriate texts for more information.³⁻⁷

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"⁸⁻¹¹ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids. Prior to discarding, sterilize specimen containers and other contaminated materials by autoclaving.

Biosafety Level 2 practices and procedures, containment equipment and facilities are required for non-aerosol-producing manipulations of clinical specimens such as preparation of acid-fast smears. All aerosol-generating activities must be conducted in a Class I or II biological safety cabinet. Biosafety Level 3 practices, containment equipment and facilities are required for laboratory activities in the propagation and manipulation of cultures of *M. tuberculosis* and *M. bovis*. Animal studies also require special procedures.¹⁰

PROCEDURE

Material Provided: Dubos Broth, Enriched

Materials Required But Not Provided: Ancillary culture media, reagents, quality control organisms and laboratory equipment as required for this procedure.

Test Procedure: Observe aseptic techniques.

The test procedures are those recommended by the Centers for Disease Control and Prevention (CDC) for primary isolation from specimens containing mycobacteria.⁴ N-acetyl-L-cysteine-sodium hydroxide (NALC-NaOH) solution is recommended as a gentle, but effective digesting and decontaminating agent. These reagents are provided in the BBL™ MycoPrep™ Specimen Digestion/Decontamination Kit. For detailed decontamination and culturing instructions, consult an appropriate text.³⁻⁷

Specimens that are less likely to be contaminated with other microorganisms (cerebrospinal fluid, pleural fluid, tissue biopsy, etc.) may be inoculated directly into the medium. Consult appropriate texts for recommended procedures.³⁻⁷

Incubate the tubes at 35 ± 2°C in a CO₂-enriched atmosphere. Keep the tube caps loosened for at least one week to permit circulation of CO₂, but tighten the caps thereafter to prevent dehydration. Loosen briefly once a week to replenish CO₂. Six to 8 weeks of incubation may be necessary for evidence of growth of many mycobacteria.

User Quality Control:

1. Examine the tubes for signs of deterioration as described under "Product Deterioration."
2. Check performance by inoculating a representative sample of each medium with pure cultures of stable control organisms that produce known, desired reactions. The following test strains are recommended:

TEST STRAIN	EXPECTED RESULTS
<i>Mycobacterium fortuitum</i> ATCC™ 6841	Growth
<i>Mycobacterium tuberculosis</i> ATCC 25177	Growth

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Growth in broth media is indicated by turbidity compared to an uninoculated control.

LIMITATIONS OF THE PROCEDURE

Enrichment broths should not be used as the sole isolation medium. They are to be used in conjunction with selective and nonselective plating media to increase the probability of isolating pathogens, especially when they may be present in small numbers.

For identification, organisms must be in pure culture. Morphological, biochemical and/or serological tests should be performed for final identification. Consult appropriate texts for detailed information and recommended procedures.³⁻⁷

PERFORMANCE CHARACTERISTICS

Since Dubos Broth provides dispersed growth, it is used to prepare a relatively uniform suspension of mycobacteria in bacterial studies. Dubos Broth was used to grow *Mycobacterium tuberculosis* to determine the enzymatic activity of alcohol dehydrogenase. After 7 days, the specific activity of the enzyme from cells grown in Dubos Broth was 2 to 3 times higher than from cells grown on Sauton medium.¹²

AVAILABILITY


Cat. No. Description


295697 BBL™ Dubos Broth, Enriched, Pkg of 10 size K tubes

REFERENCES

1. Dubos, R.J., F. Fenner, and C.H. Pierce. 1950. Properties of a culture of BCG grown in liquid media containing Tween 80 and the filtrate of heated serum. *Am. Rev. Tuberc.* 61:66.
2. MacFaddin, J.F. 1985. Media for isolation-cultivation-identification-maintenance of medical bacteria, vol. 1. Williams & Wilkins, Baltimore.
3. Cernoch, P.L., Enns, R.K., Saubolle, M.A. and R.J. Wallace. 1994. Cumitech 16A, Laboratory diagnosis of the mycobacterioses. Coordinating ed. A.S. Weissfeld, American Society for Microbiology, Washington, D.C.
4. Kent, P.T., and G.P. Kubica. 1985. Public health mycobacteriology: a guide for the level III laboratory. USDHHS. Centers for Disease Control, Atlanta.
5. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 1998. Bailey & Scott's diagnostic microbiology, 10th ed. Mosby, Inc., St. Louis.
6. Metchock B.G., Nolte, F.S., and R.J. Wallace. 1999. *Mycobacterium*, p. 399-437. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.
7. Master, R.N. 1992. Mycobacteriology, 3.01-3.16.4. In H.D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
8. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
9. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.

10. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
11. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
12. De Bruyn, J., A. Johannes, M. Weckx, and M.P. Beumer-Jochmans. 1980. Partial purification and characterization of an alcohol dehydrogenase of *Mycobacterium tuberculosis* var. *bovis* (BCG). J. of Gen. Microbiol. 124: 359-363.


 Manufacturer / Výrobce / Producent / Fabrikant / Tootja / Valmistaja / Fabricant / Hersteller / Κατασκευαστής / Gyártó / Ditta produttrice / Gamintojas / Producent / Fabricante / Výrobca / Tillverkare / Производитель / Producător / Üretici / Proizvođač / Производитель / Атқарушы

 Use by / Spotřebujte do / Anvendes før / Houdbaar tot / Kasutada enne / Viimeinkäyttöpäivä / A utiliser avant / Verwendbar bis / Ημερομηνία λήξης / Felhasználhatóság dátuma / Usare entro / Naudokite iki / Brukes før / Stosować do / Utilizar em / Použite do / Usar antes de / Använd före / Используйте до / A se utiliza până la / Son kullanma tarihi / Upotrebiti do / Исползовать до / дейін пайдалануға / Upotrijebiti do / YYYY-MM-DD / YYYY-MM (MM = end of month) / RRRR-MM-DD / RRRR-MM (MM = konec měsíce) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutning af måned) / JJJJ-MM-DD / JJJJ-MM (MM = einde maand) / AAAA-KK-PP / AAAA-KK (KK = kuu lõpp) / VVVV-KK-PP / VVVV-KK (kuukauden loppuun mennessä) / AAAA-MM-JJ / AAAA-MM (MM = fin du mois) / JJJJ-MM-TT / JJJJ-MM (MM = Monatsende) / EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα) / ÉÉÉÉ-HH-NN / ÉÉÉÉ-HH (HH = hónap utolsó napja) / AAAA-MM-GG / AAAA-MM (MM = fine mese) / MMMM-MM-DD / MMMM-MM (MM = mensesio pabaiga) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutten av måneden) / RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca) / AAAA-MM-DD / AAAA-MM (MM = fim do mês) / RRRR-MM-DD / RRRR-MM (MM = koniec miesiąca) / aaaa-mm-dd / aaaa-mm (mm = fin del mes) / ÁÁÁÁ-MM-DD / ÁÁÁÁ-MM (MM = slutet på månaden) / ГГГГ-MM-ДД / ГГГГ-MM (MM = края на месеца) / AAAA-LL-ZZ / AAAA-LL (LL = sfârșitul lunii) / YYYY-AA-GG / YYYY-AA (AA = ayın sonu) / GGGG-MM-DD / GGGG-MM (MM = kraj meseca) / ГГГГ-MM-ДД / ГГГГ-MM (MM = конец месяца) / ЖЖЖЖ-АА-КК / ЖЖЖЖ-АА (АА = айдың соңы) / GGGG-MM-DD / GGGG-MM (MM = kraj mjeseca)


REF Catalog number / Katalogové číslo / Katalognummer / Catalogusnummer / Kataloogi number / Tuotenumero / Numéro catalogue / Bestellnummer / Αριθμός καταλόγου / Katalógusszám / Numero di catalogo / Katalogo numeris / Numer katalogowy / Número do catálogo / Katalogové číslo / Número de catálogo / Каталоген номер / Număr de catalog / Katalog numarası / Kataloški broj / Номер по каталогу / Каталог нөмірі


EC REP Authorized Representative in the European Community / Autorizovaný zástupce pro Evropskou unii / Autoriseret repræsentant i EU / Erkend vertegenwoordiger in de Europese Unie / Volitatus esindaja Euroopa Nõukogus / Valituutettu edustaja Euroopan yhteisössä / Représentant agréé pour la C.E.E. / Autorisierte EG-Vertretung / Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Hivatalos képviselő az Európai Unióban / Rappresentante autorizzato nella Comunità europea / Igallotasis atstovas Europos Bendrijoje / Autorisert representant i EU / Autoryzowane przedstawicielstwo w Unii Europejskiej / Representante autorizado na União Europeia / Autorizovaný zástupca v Európskom spoločenstve / Representante autorizado en la Comunidad Europea / Auktoriserad representant i EU / Оторизирован представитель в EU / Reprezentant autorizat în Uniunea Europeană / Автура Топлүүлүгү Yetkili Temsilcisi / Ovlašćeni predstavnik u Evropskoj zajednici / Уполномоченный представитель в Европейском сообществе / Европа қауымдастығындағы уәкілетті өкіл / Autorizuirani predstavnik u EU

IVD In Vitro Diagnostic Medical Device / Lékařské zařízení určené pro diagnostiku in vitro / In vitro diagnostisk medicinsk anordning / Medisch hulpmiddel voor in vitro diagnose / In vitro diagnostika meditsiniaparatuur / Lääkinnällinen in vitro -diagnostiikkalaitte / Dispositif médical de diagnostic in vitro / Medizinisches In-vitro-Diagnostikum / In vitro διαγνωστική ιατρική συσκευή / In vitro diagnosztikai orvosi eszköz / Dispositivo medico diagnostico in vitro. / In vitro diagnostikos prietais / In vitro diagnostisk medisinsk utstyr / Urządzenie medyczne do diagnostyki in vitro / Dispositivo médico para diagnóstico in vitro / Medicinska pomůcka na diagnostiku in vitro / Dispositivo médico de diagnóstico in vitro / Medicinsk anordning för in vitro-diagnostik / Медицински уред за диагностика in vitro / Aparatură medicală de diagnosticare in vitro / In Vitro Diagnostik Tibbi Cihaz / Medicinski uređaj za in vitro dijagnostiku / Медицинский прибор для диагностики in vitro / Жасанды жағдайда жүргізетін медициналық диагностика аспабы / Medicinska pomagala za In Vitro Dijagnostiku

 Temperature limitation / Teplotní omezení / Temperaturbegrænsning / Temperatuurlimiet / Temperatuuri piirang / Lämpötilarajoitus / Température limite / Zulässiger Temperaturenbereich / Όριο θερμοκρασίας / Hölmérsékleti határ / Temperatura limite / Laikymo temperatūra / Temperaturbegrænsning / Ograniczenie temperatury / Limitação da temperatura / Ohraničenje teploty / Limitación de temperatura / Temperaturbegrænsning / Температурни ограничения / Limitare de temperatură / Sıcaklık sınırlaması / Ograničenje temperature / Ограничение температуры / Температураны шектеу / Dozvoljena temperatura

LOT Batch Code (Lot) / Kód (číslo) šarže / Batch kode (Lot) / Chargennummer (lot) / Partii kood / Eräkoodi (LOT) / Code de lot (Lot) / Chargencode (Chargenbezeichnung) / Κωδικός παρτίδας (Παρτίδα) / Tétel száma (Lot) / Codice del lotto (partita) / Partijos numeris (Lot) / Batch-kode (Serie) / Kod partii (seria) / Código do lote (Lote) / Kód série (šarža) / Código de lote (Lote) / Satskod (parti) / Код (Партия) / Număr lot (Lotul) / Parti Kodu (Lot) / Kod serije / Код партии (lot) / Топтама коды / Lot (kod)

 Consult Instructions for Use / Readpleeg gebruiksaanwijzing / Lugeda kasutusjuhendit / Tarkista käyttöohjeista / Consulter la notice d'emploi / Gebrauchsanweisung beachten / Συμβουλευτείτε τις οδηγίες χρήσης / Olvassa el a használati utasítást / Consultare le istruzioni per l'uso / Skaitykite naudojimo instrukcijas / Se i bruksanvisningen / Zobacz instrukcja użytkowania / Consulte as instruções de utilização / Pozri Pokyny na používanie / Consultar las instrucciones de uso / Se bruksanvisningen / Направте справка в инструкциите за употреба / Consultați instrucțiunile de utilizare / Kullanım Talimatları'na başvurun / Pogledajte uputstvo za upotrebu / См. руководство по эксплуатации / Пайдалану нұсқаулығымен танысып алыңыз / Koristi upute za upotrebu

 Becton, Dickinson and Company **EC REP** Benex Limited
7 Loveton Circle Rineanna House
Sparks, MD 21152 USA Shannon Free Zone
800-638-8663 Shannon, County Clare, Ireland
www.bd.com/ds

ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo, BBL and MycoPrep are trademarks of Becton, Dickinson and Company.
©2010 BD.