
Revisions

BALTSO0191 Version 8.0 Template 4
Inserts

Rev from	Rev to	JOB #
N/A	01	8038-15

NOTES:

- BD Catalog Number: 245015
- Blank (Sheet) Size: Length: 12" Width: 28.5"
- Number of Pages: 12 Number of Sheets: 1
- Page Size: Length: 12" Width: 4.75" Final Folded Size: 6" x 4.75"
- Ink Colors: No. of Colors: 1 PMS#: 2755 Blue
- Printed two sides: Yes ☒ No ☐
- Style (see illustrations below): # 4

- See specification control no. N/A for material information.
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

Label Design	REVISED BY By Nancy Carlsen at 2:30 pm, Dec 14, 2015	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION.	 BD Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	PROOFING APPROVED BY By Helaine Jones at 3:18 pm, Dec 14, 2015		
Checked By	THIRD EYE BY By Sarah Cork at 12:38 pm, Dec 16, 2015		
Part Number:	8809791JAA	Category and Description Package Insert CrystalSpec Nephelometer User's Guide	Sheet: 1 of 13 Scale: N/A A

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteys lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőjétől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielstwem BD. / Contacte o representante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinizle temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

The **CrystalSpec™** Nephelometer is a battery-powered, portable device designed to measure the turbidity of microbial suspensions equivalent to McFarland standards 0.5 through 4.0. The instrument can be used for any laboratory procedure that requires inoculum density adjustments within this range.

SUMMARY AND EXPLANATION

The visual comparison of organism suspensions to turbidity standards is an accepted method of estimating organism densities. The most widely accepted standard is the McFarland standard.¹ A McFarland standard is prepared by adding barium chloride to aqueous sulfuric acid. The density of the resulting barium sulfate precipitate can be used to approximate the colony count of a prepared suspension; e.g., McFarland 1 is the equivalent of approximately 3 x 10⁸ CFU/mL. Other standards have been used for density measurements, including titanium dioxide² and latex particle suspensions.³ The instrumented measurement of turbidity relies on the ability of particles to scatter light while in suspension. The measurement of this scattered light is referred to as nephelometry.⁴ In order to obtain an accurate density measurement, a reliable method of calibration must be used.

PRINCIPLES OF THE PROCEDURE

The **CrystalSpec** Nephelometer measures scattered light at a 90° angle of incidence. It is a solid state instrument that uses a light-emitting diode as a light source and a photodiode detector. The instrument was developed using bacterial suspensions at known concentrations (CFU/mL). Colony counts were correlated with McFarland units and the **CrystalSpec** instrument was programmed to display the measurement in McFarland units. The calibration standards supplied with the instrument have been developed to assure an accurate calibration. Before making a measurement on the **CrystalSpec** Nephelometer, the instrument is calibrated using the supplied standards. To determine the McFarland equivalent of a microbial suspension, a tube is placed in the reading chamber and the test button is pressed. The results are then displayed on the LCD screen in McFarland units.

SPECIFICATIONS:

Range.....	McFarland 0.5–4.0	Instrument Ambient	
Accuracy	± 0.2 McFarland Unit	Operating Conditions	20 to 30°C
Reproducibility	± 0.1 McFarland Unit	Storage Temperature Range	
Battery	9 V Alkaline	(Instrument and calibration tube).....	20 to 30°C
Battery Life.....	1 Year (typically)	Sample Tube Size.....	16–17 mm dia., flat bottomed glass tube

Precautions: *in vitro* Diagnostic
Read all instructions before use and follow precisely.

PROCEDURE

Materials Provided: **CrystalSpec** Nephelometer, including a **CrystalSpec** Calibration Blank (**B**), **CrystalSpec** Calibration Standard (**S**), Low Volume Sample Adapter, and 9-volt battery.
Materials Not Provided: Ancillary culture media and equipment required to prepare a bacterial suspension.

Instructions:

BATTERY INSTALLATION AND REPLACEMENT – Unpack the **CrystalSpec** Nephelometer and the 9 V battery. Invert the instrument and remove the battery cover plate. Connect the battery and replace the cover.
If replacing a discharged battery, press one of the buttons (either “Calibration” or “Test”) for 3 sec after discharged battery has been removed.
When a battery is initially installed, “**T**” will be displayed in the readout window while the internal evaluation is in progress. This “**T**” will also appear when a discharged battery is replaced with a new battery.
When installing a new battery, interrupting the contact before the battery clip is secure may cause the “**T**” display to not continue to show. The unit should be functioning properly, but if in doubt, re-install the battery per above.
The **CrystalSpec** instrument has an internal battery charge indicator. Along with sample readouts, a reading from “**B0**” to “**B5**” is displayed at the far right edge of the window to indicate the level of battery charge. A fully charged battery will display “**B5**”, a battery needing replacement will display “**B0**”.
When the battery charge is so low that the integrity of the readings is compromised, “**X**” will be displayed and the battery will be completely discharged by the instrument. The meter will not operate until the battery is replaced.
NOTE: *Anytime a battery is replaced or the battery is disconnected from the unit, the instrument must be re-calibrated.*
CALIBRATION PROCEDURES – The **CrystalSpec** instrument requires no warm-up period and is ready for calibration and immediate use. A **CrystalSpec** Calibration Blank (**B**) and **CrystalSpec** Calibration Standard (**S**) are supplied with the instrument; *use only these calibration standards for calibration.* Internal error checking is in place to assure proper calibration of the instrument.
The nephelometer has internal features that minimize the need for frequent calibrations. Calibration needs to be performed only under the following conditions: (a) Any time the battery is replaced or disconnected, (b) If the “**E**” message is displayed when taking a reading and (c) Whenever the daily reading using the Calibration Standard (**S**) fails to read 1.9, 2.0 or 2.1.

1. If the sample volume to be measured is between 2.0 mL and 4.0 mL, insert the Low Volume Sample Adapter provided prior to calibration. If the volume is greater than 4.0 mL, the Low Volume Sample Adapter is *not* needed.
2. For ease of calibrating, press the Calibration Button using a finger of your left hand and press the Test Button using a finger of your right hand.
3. The sequence of button press and release is important for the proper calibration and accurate functioning of the **CrystalSpec**.
4. If during calibration the expected displays do not show, after the display clears, re-start the calibration procedure with step 1.
5. If a calibration reading (step 8) of other than "2.0" is returned, remove the Calibration Standard and wipe the outside of the tube with a Kimwipe or clean cloth to remove any smudges, finger prints, etc. and re-calibrate.
6. Calibration between sample readings is not necessary; however, the **CrystalSpec** Calibration Standard (**S**) should be read as a sample on a daily basis. If the standard reads 1.9, 2.0 or 2.1, the instrument is calibrated and *does not need re-calibration*.
7. The calibration should be stable if ambient conditions of temperature and light remain relatively constant. If these conditions change, re-calibration is recommended.
8. Do not press the Test Button until the calibration procedure has been successfully completed (steps 1 - 8 are completed successfully) or the **CrystalSpec** will not be properly calibrated and readings will not be valid.

TESTING PROCEDURES – *Clear, flat bottomed glass tubes, 16–17 mm diameter must be used for samples.* Organisms may be suspended in water, saline or clear growth medium of yellow to pale brown color (e.g., Mueller-Hinton Broth, TSB, or BHI Broth). Use of other colored liquids as a suspending medium is not recommended. A minimum of 2.0 mL of liquid is required.

1. If the sample volume is between 2.0 mL and 4.0 mL, insert the Low Volume Sample Adapter. If the volume is greater than 4.0 mL, the Low Volume Sample Adapter is not needed.
2. Vortex the test sample and set aside briefly if bubbles are present. Gently tapping on the tube may aid in eliminating bubbles. Insert the test tube sample tube into the nephelometer. Make sure the tube is inserted as far as possible into the instrument.
3. Depress and release the test button [test tube graphic]. The result will be displayed along with the battery status readout.

EXAMPLE: For a sample of McFarland 1.0 and a battery that requires changing, the display will read:

"1.0 MCFARLAND B0".

4. Adjust the sample preparation with more organisms if the readout is lower than expected, or with sterile diluent if the reading is higher than expected. Vortex the sample and reread.
5. Repeat step 4 until the desired turbidity is achieved.
6. Repeat steps 2–5 (above) for additional samples.

RESULTS

The **CrystalSpec** Nephelometer displays the measurement in McFarland units. These readings represent the density of the suspension and can be used to estimate the CFU/mL of the suspension tested.

LIMITATIONS OF THE PROCEDURE

*The meter must be calibrated with the supplied **CrystalSpec** Calibration Blank (**B**) and **CrystalSpec** Calibration Standard (**S**). The use of other standards will result in miscalibration of the instrument.*

*The **CrystalSpec** Standard is to be used only for calibration of this instrument and not as a visual, i.e. manual approximation to a barium sulfate McFarland standard.*

Do not use the **CrystalSpec** Nephelometer for measurements outside the range of McFarland standards 0.5–4.0.

Volumes below 2.0 mL cannot be measured in the instrument.

Clear, flat bottomed glass tubes, 16–17 mm diameter must be used for samples.

For best results we recommend 9 V Alkaline batteries.

Do not use in direct sunlight.

PERFORMANCE CHARACTERISTICS

The reproducibility of the instrument was measured over time with standards at concentrations equivalent to McFarland 0.5, 1, 2, 3, and 4 turbidity standards. The error for any estimate was determined to be ± 0.1 McFarland unit.

To determine the accuracy of the **CrystalSpec** Nephelometer, colony counts were performed on suspensions of *E. coli* ATCC® 25922 prepared in saline. Each test was performed in replicates of five and the average of the plate counts is given in Table 1.

Table 1

McFarland	Expected CFU/mL x 10 ⁸	Adjusted* CFU/mL x 10 ⁸	%CV
0.5	1.5	1.7	11.6
1.0	3.0	3.3	10.7
2.0	6.0	6.5	14.8
3.0	9.0	10.4	8.2
4**0	12.0	13.0	9.9

* Because of difficulty in achieving the exact target McFarland density, these results were obtained by using the **CrystalSpec** readings and adjusting the observed CFU/mL as if each dilution had been the target McFarland.

**Four replicators were used to calculate the adjusted CFU/mL and % CV.

WARRANTY: The **CrystalSpec** Nephelometer is warranted to be free from defects for one year from the date of purchase.

MAINTENANCE AND TROUBLESHOOTING: No maintenance is required except for battery replacement when the battery indicator reads **"B0"**, or **"E"** is displayed, indicating a low battery.

If the **"E"** is displayed, the internal error checking has detected a problem due to improper use, discharged battery, or an electronic failure. If calibration and/or replacement of the battery does not correct the problem, (in the United States) contact Becton Dickinson, Technical Services, toll free (800) 638-8663.

AVAILABILITY

Cat. No.	Description
245015	BD CrystalSpec™ Calibration Blank and BD CrystalSpec™ Calibration Standard, carton of one each.

REFERENCES

1. McFarland, J. 1907. The nephelometer: an instrument for estimating the number of bacteria in suspensions used for calculating the opsonic index for vaccines. JAMA 49:1176-1178.

2. Roessler, W.G., and C.R. Brewer. 1967. Permanent turbidity standards. Appl. Microbiol. 15:1114-1121.

3. Pugh, T.L., and W. Heller. 1957. Density of polystyrene and polyvinyl toluene latex particles. J. Colloid Sci. 12:173-180.

4. Mallette, M.F. 1969. XV. Evaluation of growth by physical and chemical means, p. 521-566. In J.R. Norris and D.W. Ribbons (ed.). Methods in microbiology, vol. 1. Academic Press Inc., New York.

BD Diagnostics Technical Service and Support: contact your local BD representative or www.bd.com/ds.

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ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo and CrystalSpec are trademarks of Becton, Dickinson and Company. © 2015 BD
Benex Limited