

QUALITY CONTROL PROCEDURES

I INTRODUCTION

Trypticase™ Soy Agar is a general purpose medium which supports the growth of fastidious as well as nonfastidious microorganisms.

II PERFORMANCE TEST PROCEDURE

1. Liquefy the medium in the tubed deeps by heating in boiling water. Cool to 45–50 °C, add 1 mL of sterile defibrinated sheep blood to two tubes (for inoculation with *Streptococcus* strains) and pour into sterile Petri dishes. Mix well to evenly distribute the blood throughout the medium and allow to solidify for a minimum of 30 min.
2. Inoculate representative samples with the cultures listed below.
 - a. Using a 0.01 mL calibrated loop, inoculate the agar surfaces using 10⁻¹ dilutions of 18- to 24-h **Trypticase Soy Broth** cultures. Streak-inoculate the plates to assure the presence of well-isolated colonies.
 - b. Incubate plates or tubed slants (with loosened caps) at 35 ± 2 °C in an aerobic atmosphere. Blood plates should be incubated in the presence of carbon dioxide.
3. Examine plates or tubes after 18–24 and 42–48 h for amount of growth and pigmentation. Examine the blood agar plates for hemolysis.
4. Expected Results

Organisms	ATCC™	Recovery
Medium without the addition of blood.		
* <i>Shigella flexneri</i>	12022	Growth. Colonies medium to large, grayish-white, translucent, slightly convex and may be mucoid.
* <i>Escherichia coli</i>	25922	Growth
* <i>Staphylococcus aureus</i>	25923	Growth. Colonies medium to large, opaque, circular, entire with cream-yellow to gold pigment.
Medium with the addition of sterile defibrinated sheep blood (tubed deeps).		
* <i>Streptococcus pneumoniae</i>	6305	Growth. Colonies surrounded by zones of alpha hemolysis (green).
* <i>Streptococcus pyogenes</i>	19615	Growth. Colonies surrounded by zones of beta hemolysis (clear to hazy).

*Recommended organism strain for User Quality Control.

III ADDITIONAL QUALITY CONTROL

1. Examine tubes as described under "Product Deterioration."
2. Visually examine representative tubes to assure that any existing physical defects will not interfere with use.
3. Determine the pH potentiometrically at room temperature for adherence to the specification of 7.3 ± 0.2.
4. Incubate uninoculated representative tubes at 20–25 °C and 30–35 °C and examine after 7 days for microbial contamination.

PRODUCT INFORMATION

IV INTENDED USE

Trypticase Soy Agar is used for the isolation and cultivation of fastidious as well as nonfastidious microorganisms, including anaerobic as well as aerobic bacteria, although it is not the medium of choice for anaerobes.

V SUMMARY AND EXPLANATION

The nutritional composition of **Trypticase Soy Agar** has made it a popular medium for many years. It is the medium specified as Soybean-Casein Digest Agar Medium in *The United States Pharmacopeia* for the total aerobic microbial count portion of the microbial limit testing procedures.¹

The medium is used for a multitude of purposes including maintenance of stock cultures, plate counting, isolation of microorganisms from a variety of specimen types and as a base for media containing blood.²⁻⁴ It is included in the compendia of methods for the examination of water, wastewater and foods.^{5,6}

VI PRINCIPLES OF THE PROCEDURE

The combination of casein and soy peptones in **Trypticase Soy Agar** renders the medium highly nutritious by supplying organic nitrogen, particularly amino acids and longer-chained peptides. The sodium chloride maintains the osmotic equilibrium.

VII REAGENTS

Trypticase Soy Agar

Approximate Formula* Per Liter Purified Water

Pancreatic Digest of Casein.....	15.0 g
Papaic Digest of Soybean Meal.....	5.0 g
Sodium Chloride	5.0 g
Agar	15.0 g

*Adjusted and/or supplemented as required to meet performance criteria.

Warnings and Precautions: For *in vitro* Diagnostic Use.

Tubes with tight caps should be opened carefully to avoid injury due to breakage of glass.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"⁷⁻¹⁰ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids. After use, prepared tubes, specimen containers and other contaminated materials must be sterilized by autoclaving before discarding.

Storage Instructions: On receipt, store tubes in the dark at 2–25 °C. Avoid freezing and overheating. Do not open until ready to use. Minimize exposure to light. Tubed media stored as labeled until just prior to use may be inoculated up to the expiration date and incubated for the recommended incubation times. Allow the medium to warm to room temperature before inoculation.

Product Deterioration: Do not use tubes if they show evidence of microbial contamination, discoloration, drying or other signs of deterioration.

VIII SPECIMEN COLLECTION AND HANDLING

Specimens suitable for culture may be handled using various techniques. For detailed information, consult appropriate texts.^{3,11} Specimens should be obtained before antimicrobial agents have been administered. Provision must be made for prompt delivery to the laboratory.

IX PROCEDURE

Material Provided: **Trypticase** Soy Agar

Materials Required But Not Provided: Ancillary culture media, reagents, quality control organisms and laboratory equipment as required.

Test Procedure: Observe aseptic techniques.

Liquefy the medium in the tubed deeps by heating in boiling water. Cool to 45–50 °C, add blood, if desired, and pour into sterile Petri dishes.

For general use, streak the specimen as soon as possible after it is received in the laboratory. The streak plate is used primarily to isolate pure cultures from specimens containing mixed flora. Alternatively, if material is being cultured directly from a swab, roll the swab over a small area of the surface at the edge; then streak from this inoculated area. Since many pathogens require carbon dioxide on primary isolation, plates may be incubated in an atmosphere containing approximately 3–10% CO₂. Incubate plates at 35 ± 2 °C for 18–24 h.

Tubed slants of **Trypticase** Soy Agar are primarily used for the growth and maintenance of pure cultures. They should be inoculated with an inoculating loop and incubated under the same conditions as the plated medium.

User Quality Control: See "Quality Control Procedures."

Quality Control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

X RESULTS

After incubation, most plates will show an area of confluent growth. Because the streaking procedure is, in effect, a "dilution" technique, diminishing numbers of microorganisms are deposited on the streaked areas. Consequently, one or more of these areas should exhibit isolated colonies of the organisms contained in the specimen. In addition, growth of each organism may be semiquantitatively scored on the basis of growth in each of the streaked areas.

Hemolytic reactions should be noted for organisms inoculated on the medium containing blood.

The tubed slant containing pure cultures can be used for additional studies or as stock cultures as desired.

XI LIMITATIONS OF THE PROCEDURE

For identification, organisms must be in pure culture. Morphological, biochemical, and/or serological tests should be performed for final identification. Consult appropriate texts for detailed information and recommended procedures.^{3,11,12}

XII PERFORMANCE CHARACTERISTICS

Trypticase Soy Agar (TSA) with 5% Sheep Blood was used as a control in a study using broth-enhanced culture (Todd Hewitt) and Optical Immunoassay method for the diagnosis of b-hemolytic streptococcal infection. Five hundred two (502) specimens were tested. TSA with 5% Sheep Blood had a sensitivity and specificity of 92.5% and 99.4%, respectively.¹³ Nguyen et al. used **Trypticase** Soy Agar with 5% Sheep Blood as the 'gold standard' for the detection of group B *Streptococcus* from the lower genital tract of pregnant women.¹⁴ In another study, Rossmann et al. successfully reisolated *Lautropia mirabilis* on **Trypticase** Soy Agar with 5% Sheep Blood from the oral cavities of human immunodeficiency virus infected children.¹⁵ Of the 85 children evaluated in this study, 35 (41.4%) were positive for *L. mirabilis*. Isenberg et al. used **Trypticase** Soy Agar with 5% Sheep Blood as a control to evaluate the recovery of *Enterococcus* from a selective medium under study.¹⁶ Two hundred fifty (250) group D streptococcal strains isolated from clinical material and 8 strains obtained from the National Communicable Disease Center (Atlanta, Ga.) were used. Kantor et al. maintained stock cultures at room temperature using **Trypticase** Soy Agar slants covered with sterile mineral oil for a study on the identification of nonfermentative gram-negative bacteria in the clinical laboratory.¹⁷

XIII AVAILABILITY

Cat. No.	Description
221082	BBL™ Trypticase™ Soy Agar Deeps (Pour Tubes), 20 mL, Pkg. of 10 size A tubes
221086	BBL™ Trypticase™ Soy Agar Slants, Pkg. of 10 size K tubes
221087	BBL™ Trypticase™ Soy Agar Slants, Ctn. of 100 size K tubes

XIV REFERENCES

1. U.S. Pharmacopeial Convention, Inc. 2006. The U.S. pharmacopeia 29/The national formulary 24-2006. U.S. Pharmacopeial Convention, Inc., Rockville, Md.
2. MacFaddin, J.F. 1985. Media for isolation-cultivation-identification-maintenance of medical bacteria, vol. 1, Williams & Wilkins, Baltimore.
3. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2002. Bailey & Scott's diagnostic microbiology, 11th ed. Mosby, Inc., St. Louis.
4. Chapin, K.C., and P.R. Murray. 1999. Media, p. 1687-1707. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.
5. Clesceri, L.S., A.E. Greenberg, and A.D. Eaton (ed.) 1998. Standard methods for the examination of water and wastewater, 20th ed. American Public Health Association, Washington, D.C.
6. Downes, F.P. and K. Ito (ed.). 2001. Compendium of methods for the microbiological examination of foods, 4th ed. American Public Health Association, Washington, D.C.
7. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
8. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
9. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
10. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
11. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenover (ed.) 2003. Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
12. Holt, J.G., N.R. Krieg, P.H.A. Sneath, J.T. Staley, and S.T. Williams (ed.). 1994. Bergey's Manual™ of determinative bacteriology, 9th ed. Williams & Wilkins, Baltimore.
13. Fries, S.M. 1995. Diagnosis of group A streptococcal pharyngitis in a private clinic: comparative evaluation of an optical immunoassay method and culture. J. Ped. vol. 126, number 6.
14. Nguyen, T.M. et al. 1998. Detection of group B streptococcus: comparison of an optical immunoassay with direct plating and broth-enhanced culture methods. J. Matern. Fetal. Med. Jul-Aug; 7(4):172-176.
15. Rossmann, S.N. et al. 1998. Isolation of *Lautropia mirabilis* from oral cavities of human immunodeficiency virus infected children. J. Clin. Microbiol. 36:1756-1760.
16. Isenberg, H.D., D. Goldberg, and J. Sampson. 1970. Laboratory studies with a selective medium. Appl. Microbiol. Sept. 1970, p.433-436.
17. Kantor, L.T., D.K. Spyros, and R.B. Yee. 1975. Identification of nonfermentative gram-negative bacteria in the clinical laboratory. Amer. J. Med. Tech. vol. 41, no. 1.

Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

 Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

 Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo, BBL and Trypticase are trademarks of Becton, Dickinson and Company. ©2014 BD.

