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221730 **BD BBL GN Broth**, 8 mL, caixa com 100 tubos K

XIV BIBLIOGRAFIA

1. Hajna, A.A. 1955. A new specimen preservative for gram-negative organisms of the intestinal group. Public Health Lab. 13:59-62.
2. Hajna, A.A. 1955. A new enrichment broth medium for gram-negative organisms of the intestinal group. Public Health Lab. 13:83-89.
3. Croft, C.C., and M.J. Miller. 1956. Isolation of *Shigella* from rectal swabs with Hajna "GN" broth. Am. J. Clin. Pathol. 26:411-417.
4. Taylor, W.I., and D. Schelhart. 1967. Isolation of shigellae. IV. Comparison of plating media with stools. Am. J. Clin. Pathol. 48:356-362.
5. Downes, F.P. and K. Ito (ed.). 2001. Compendium of methods for the microbiological examination of foods, 4th ed. American Public Health Association, Washington, D.C.
6. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
7. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
8. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
9. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
10. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenover (ed.). 2003. Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.

11. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2002. Bailey & Scott's diagnostic microbiology, 11th ed. Mosby, Inc., St. Louis.
12. Ewing, W.H. 1986. Edwards and Ewing's identification of *Enterobacteriaceae*, 4th ed. Elsevier Science Publishing Co., Inc., New York.
13. MacFaddin, J.F. 1985. Media for isolation-cultivation-identification-maintenance of medical bacteria, vol. I. Williams & Wilkins, Baltimore.
14. Taylor, W.I. and D. Schelhart. 1968. Isolation of Shigellae. V. Comparison of enrichment broths with stools. Appl. Microbiol. 16:1383-1386.

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo, and all other trademarks are property of Becton, Dickinson and Company. © 2015 BD