BBL Trypticase Soy Agar

L007516 • Rev. 11 • Setembro 2014

PROCEDIMENTOS DE CONTROLO DE QUALIDADE

I INTRODUÇÃO

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XIV BIBLIOGRAFIA

1. U.S. Pharmacopeial Convention, Inc. 2006. The U.S. pharmacopeia 29/The national formulary 24-2006. U.S. Pharmacopeial Convention, Inc., Rockville, Md.
2. MacFaddin, J.F. 1985. Media for isolation-cultivation-identification-maintenance of medical bacteria, vol. 1, Williams & Wilkins, Baltimore.
3. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2002. Bailey & Scott's diagnostic microbiology, 11th ed. Mosby, Inc., St. Louis.
4. Chapin, K.C., and P.R. Murray. 1999. Media, p. 1687-1707. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.
5. Clesceri, L.S., A.E. Greenberg, and A.D. Eaton (ed.) 1998. Standard methods for the examination of water and wastewater, 20th ed. American Public Health Association, Washington, D.C.
6. Downes, F.P. and K. Ito (ed.). 2001. Compendium of methods for the microbiological examination of foods, 4th ed. American Public Health Association, Washington, D.C.
7. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
8. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
9. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
10. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
11. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenover (ed.) 2003. Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
12. Holt, J.G., N.R. Krieg, P.H.A. Sneath, J.T. Staley, and S.T. Williams (ed.). 1994. Bergey's Manual™ of determinative bacteriology, 9th ed. Williams & Wilkins, Baltimore.
13. Fries, S.M. 1995. Diagnosis of group A streptococcal pharyngitis in a private clinic: comparative evaluation of an optical immunoassay method and culture. J. Ped. vol. 126, number 6.
14. Nguyen, T.M. et al. 1998. Detection of group B streptococcus: comparison of an optical immunoassay with direct plating and broth-enhanced culture methods. J. Matern. Fetal. Med. Jul-Aug; 7(4):172-176.
15. Rossmann, S.N. et al. 1998. Isolation of *Lautropia mirabilis* from oral cavities of human immunodeficiency virus infected children. J. Clin. Microbiol. 36:1756-1760.
16. Isenberg, H.D., D. Goldberg, and J. Sampson. 1970. Laboratory studies with a selective medium. Appl. Microbiol. Sept. 1970, p.433-436.
17. Kantor, L.T., D.K. Spyros, and R.B. Yee. 1975. Identification of nonfermentative gram-negative bacteria in the clinical laboratory. Amer. J. Med. Tech. vol. 41, no. 1.

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection.

BD, BD Logo, BBL, GasPak and Trypticase are trademarks of Becton, Dickinson and Company. ©2014 BD.