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215102 **BBL BCYE Agar, Karton mit 100 Platten**

XIV LITERATUR

1. Feeley, J.C., R.J. Gibson, G.W. Gorman, N.C. Langford, J.K. Rasheed, D.C. Mackel, and W.B. Baine. 1979. Charcoal-yeast extract agar: primary isolation medium for *Legionella pneumophila*. *J. Clin. Microbiol.* 10:437-441.
2. Feeley, J.C., G.W. Gorman, R.E. Weaver, D.C. Mackel, and H.W. Smith. 1978. Primary isolation media for Legionnaires' disease bacterium. *J. Clin. Microbiol.* 8:320-325.
3. Pasculle, A.W., J.C. Feeley, R.J. Gibson, L.G. Cordes, R.L. Myerowitz, C.M. Patton, G.W. Gorman, C.L. Carmack, J.W. Ezzell, and J.N. Dowling. 1980. Pittsburgh pneumonia agent: direct isolation from human lung tissue. *J. Infect. Dis.* 141:727-732.
4. Edelstein, P.H. 1981. Improved semi-selective medium for isolation of *Legionella pneumophila* from contaminated clinical and environmental specimens. *J. Clin. Microbiol.* 14:298-303.
5. National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
6. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
7. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
8. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262, 17/10/2000, p. 0021-0045.*

9. Isenberg, H.D., FD. Schoenknecht, and A. von Graevenitz. 1979. Cumitech 9, Collection and processing of bacteriological specimens. Coordinating ed., S.J. Rubin. American Society for Microbiology, Washington, D.C.
10. Miller, J.M. and H.T. Holmes. 1999. Specimen collection, transport, and storage, p.33-63. *In* P.R. Murray, E.J. Baron, M.A. Tenover, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.
11. Winn, W.C., Jr. 1999. *Legionella*, p. 572-585. *In* P.R. Murray, E.J. Baron, M.A. Tenover, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C.
12. Weaver, R.E. 1978. Cultural and staining characteristics, p. 40-43. *In* G.L. Jones, and G.A. Herbert (ed.), "Legionnaires" the disease, the bacterium, and methodology. USDHEW, Center for Disease Control, Atlanta.
13. Koneman, E.W., S.D. Allen, W.M. Janda, P.C. Schreckenberger, and W.C. Winn, Jr. 1997. Color atlas and textbook of diagnostic microbiology, 5th ed. J.B. Lippincott Company, Philadelphia.
14. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2002. Bailey & Scott's diagnostic microbiology, 11th ed. Mosby, Inc., St. Louis.

 Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA
800-638-8663
www.bd.com/ds

 Benex Limited
Rineanna House
Shannon Free Zone
Shannon, County Clare, Ireland

ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo, BBL and Trypticase are trademarks of Becton, Dickinson and Company. ©2011 BD