BD MacConkey Agar with Sorbitol

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Bibliografia

1. Riley, L.W. et al. 1983: Hemorrhagic colitis associated with a rare *Escherichia coli* serotype. *New Engl. J. Med.* 308: 681-685.
2. Kaper, J.B., O'Brien, A.D. (eds.). 1998: *Escherichia coli* O157:H7 and Other Shiga Toxin-Producing *E. coli* strains. American Society for Microbiology, Washington, DC, USA.
3. Dorn, C.R., and E.J. Angrick. 1991. Serotype O157:H7 *Escherichia coli* from bovine and meat sources. *J. Clin. Microbiol.* 29: 1225-1231.
4. Wells, J.G. et al. 1983. Laboratory investigation of hemorrhagic colitis outbreaks associated with a rare *Escherichia coli* serotype. *J. Clin. Microbiol.* 18: 512-520.
5. Willshaw, G.A., et al. 1994: Vero cytotoxin-producing *Escherichia coli* O157 in beefburgers linked to an outbreak of diarrhea, haemorrhagic colitis and haemolytic uraemic syndrome in Britain. *Letters Appl. Microbiol.* 19: 304-307.
6. Ewing, W. H., and P. R. Edwards. 1954. Isolation and preliminary identification of *Escherichia coli* serotypes associated with cases of diarrhea of the newborn. *Public Health Lab.* 12:75-81.
7. March, S.B., and S. Ratman. 1986. Sorbitol-MacConkey medium for detection of *Escherichia coli* O157:H7 associated with hemorrhagic colitis. *J. Clin. Microbiol.* 23: 869-872.
8. Rappaport, F., and E. Henig. 1952. Media for the isolation and differentiation of pathogenic *Escherichia coli* (serotypes 0111 and 055). *J. Clin. Pathology.* 5:361-362.
9. Bopp, C. A., F. W. Brenner, P. I. Fields, J. G. Wells, and N. A. Stockbrine. 2003. *Escherichia, Shigella, and Salmonella*. In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). *Manual of clinical microbiology*, 8th ed. American Society for Microbiology, Washington, D.C.
10. Adams, S. 1991. Screening for verotoxin-producing *Escherichia coli*. *Clinical Lab Science* 4(1):19-20.
11. March, S. B., and S. Ratnam. 1986. Sorbitol-MacConkey medium for detection of *Escherichia coli* O157:H7 associated with hemorrhagic colitis. *J. Clin. Microbiol.* 23:869-872.

12. Meng, J., Feng, P., and M.P. Doyle. 2001. Pathogenic *Escherichia coli*. In: Downes, F.P., and K. Ito (ed.). Compendium of methods for the microbiological examination of foods., 4th edition. American Public Health Association, Washington. D.C.
13. Hitchins, A. D., P. Feng, W. D. Watkins, S. R. Rippey, and L. A. Chandler. 1995. *Escherichia coli* and the coliform bacteria. p. 4.01-4.29. In Bacteriological analytical manual, 8th ed. AOAC International, Gaithersburg, MD.
14. Mortlock, S. 1994. Recovery of *Escherichia coli* O157:H7 from mixed suspensions: evaluation and comparison of pre-coated immunomagnetic beads and direct plating. Brit. J. Biomed. Sci. 51: 207-214.
15. Ogden, I.D., Hepburn, N.F., and M. MacRae. 2001. The optimization of isolation media used in immunomagnetic separation methods for the detection of *Escherichia coli* O157 from foods. J. Appl. Microbiol. 91: 373-379.
16. Kist, M., et al. 2000. Infektionen des Darmes. In: Mauch, H., Lüttiken, R., and S. Gatermann (eds.): MiQ - Qualitätsstandards in der mikrobiologisch-infektiologischen Diagnostik, vol. 9. Urban & Fischer, Munich, Germany.

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<http://www.bd.com>

<http://www.bd.com/europe/regulatory/>

Dynabeads and Dynal are trademarks of Dynal

CHROMagar is a trademark of Dr. A. Rambach

ATCC is a trademark of the American Type Culture Collection

BD, BD logo, BBL and Trypticase are trademarks of Becton, Dickinson and Company

© 2013 BD

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Becton Dickinson GmbH