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XIV REFERANSLAR

1. Sabouraud, R. 1910. Les teignes, p. 553. Masson et Cie, Paris.
2. Ajello, L., L.K. Georg, W. Kaplan, and L. Kaufman. 1963. CDC Laboratory manual for medical mycology. PHS Publication No. 994, U.S. Government Printing Office, Washington, D.C.
3. Weitzman, I., J. Kane, and R.C. Summerbell. 1995. *Trichophyton*, *Microsporum*, *Epidermophyton*, and agents of superficial mycoses, p. 791-808. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
4. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
5. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol 17:53-80.
6. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
7. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from the risks related to exposure to the biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
8. Haley, L.D., H. Trandel, and M.B. Coyle. 1980. Cumitech 11, Practical methods for culture and identification of fungi in the clinical mycology laboratory. Coordinating ed., J.C. Sherris. American Society for Microbiology, Washington, D.C.
9. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 1998. Bailey & Scott's diagnostic microbiology, 10th ed. Mosby, Inc., St. Louis.
10. Merz, W.G., and G.D. Roberts. 1995. Detection and recovery of fungi from clinical specimens, p. 709-722. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
11. Kwon-Chung, K.J., and J.E. Bennett. 1992. Medical mycology. Lea & Febiger, Philadelphia.
12. Larone, D.H. 1995. Medically important fungi: a guide to identification, 3rd ed. American Society for Microbiology, Washington, D.C.
13. Koneman, E.W., S.D. Allen, W.M. Janda, P.C. Schreckenberger, and W.C. Winn, Jr. 1992. Mycology, p. 791-877. Color atlas and textbook of diagnostic microbiology, 4th ed., J.B. Lippincott Co., Philadelphia.

BD Diagnostics Teknik Desteęi: yerel BD temsilcinizle temasa gein veya www.bd.com/ds adresine bařvurun.

Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection.
BD, BD Logo and BBL are trademarks of Becton, Dickinson and Company. © 2015 BD.