**BBL Thioglycollate Medium,
Enriched with Vitamin K₁ and Hemin**
L007509 • Rev. 14 • Oktober 2015

KVALITETSKONTROLLPROSEDYRER (Valgfritt)

I INNLEDNING

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Kat. nr.	Beskrivelse
221742	BD BBL Thioglycollate Medium, Enriched with Vitamin K ₁ , and Hemin, 5 mL
221787	BD BBL Thioglycollate Medium, Enriched with Vitamin K ₁ , and Hemin, 8 mL
221788	BD BBL Thioglycollate Medium, Enriched with Vitamin K ₁ , and Hemin, 8 mL
297292	BD BBL Thioglycollate Medium, Enriched with Vitamin K ₁ , and Hemin, 10 mL

XIV REFERANSER

1. Brewer, J.H. 1940. A clear liquid medium for the "aerobic" cultivation of anaerobes. J. Bacteriol. 39:10.
2. Vera, H.D. 1944. Comparative study of materials suitable for the cultivation of clostridia. J. Bacteriol. 47:59-65.
3. Dowell, V.R., Jr., and T.M. Hawkins. 1987. Laboratory methods in anaerobic bacteriology. CDC laboratory manual. HHS Publication No. (CDC) 87-8272. Centers for Disease Control, Atlanta.
4. Dowell, V.R., Jr., G.L. Lombard, F.S. Thompson, and A.Y. Armfield. 1977. Media for isolation, characterization, and identification of obligately anaerobic bacteria. CDC laboratory manual. Center for Disease Control, Atlanta.
5. Forbes, B.A., and P.A. Granato. 1995. Processing specimens for bacteria, p. 265-281. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
6. Gibbons, R.J. and J.B. MacDonald. 1960. Hemin and vitamin K compounds as required factors for the cultivation of certain strains of *Bacteroides melaninogenicus*. J. Bacteriol. 80:164-170.
7. Wilkins, T.D., S.L. Chalgren, F. Jimenez-Ulate, C.R. Drake, Jr., and J.L. Johnson. 1976. Inhibition of *Bacteroides fragilis* on blood agar plates and reversal of inhibition by added hemin. J. Clin. Microbiol. 3:359-363.
8. Finegold, S.M., V.L. Sutter, H.R. Attebery, and J.E. Rosenblatt. 1974. Isolation of anaerobic bacteria, p. 365-375. In E.H. Lennette, G.H. Spaulding, and J.P. Truant (ed.), Manual of clinical microbiology, 2nd ed. American Society for Microbiology, Washington, D.C.
9. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
10. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
11. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
12. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
13. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry and M.A. Pfaller (ed.). 2007. Manual of clinical microbiology, 9th ed. American Society for Microbiology, Washington, D.C.
14. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey & Scott's diagnostic microbiology, 12th ed. Mosby, Inc., St. Louis.
15. Reischelderfer, C., and J.I. Mangels. 1994. Culture media for anaerobes, p. 2.3.1-2.3.8. In H.D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
16. Holt, J.G., N.R. Krieg, P.H.A. Sneath, J.T. Staley, and S.T. Williams (ed.). 1994. Bergey's Manual of determinative bacteriology, 9th ed. Williams & Wilkins, Baltimore.

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

ATCC is a trademark of the American Type Culture Collection.

BD, BD Logo, and all other trademarks are property of Becton, Dickinson and Company. © 2015 BD